


medi



medi Kidz range

Orthopaedic products from medi especially for children

medi. I feel better.





Orthopaedic products from medi especially for children

In this catalogue you will find a wide-ranging assortment of supports and braces. They have been designed especially to meet children's needs.

medi Kidz products promote healing after injuries such as bone fractures etc., and are used for a limited period of time. During the recovery period, they give children the support and safety they need to stay active.



Orthopaedics from medi

medi has been one of the leading manufacturers of medical devices for more than six decades.

Orthopaedic supports and braces promote the healing process and are good for the health of the entire musculoskeletal system.



6

Products for the upper limb

- 8 medi Collar soft Kidz
- 8 medi Clavicle support Kidz
- 9 medi Armschlinge Kidz
- 9 medi Shoulder sling Kidz
- 10 medi Epico ROMs Kidz
- 11 Manumed Kidz



12

Products for the lower limb

- 14 medi Classic Kidz
- 14 medi ROM Kidz
- 15 M.step Kidz
- 15 medi Walker boot Kidz



16

medi supports Perfection for children

- 18 Genumedi
- 18 Levamed
- 19 Manumed active





medi Kidz orthopaedic
products for the upper limb

medi Collar soft Kidz

Comfortable surgical collar

Indications

- Stiffness and pain of the neck and shoulders
- Osteoarthritis of the neck
- Wry neck muscle spasm (torticollis)
- Whiplash injury of the neck
- Instability of the neck
- Rheumatic pain

Mode of action

- Supports and stabilises the neck
- The warming effect helps treatment even further



Measuring for size: circumference of the neck

size	neck circumference	article number
1	18 – 22 cm (Height 5 cm)	R.222.011
2	22 – 28 cm (Height 6 cm)	R.222.012
3	28 – 33 cm (Height 7 cm)	R.222.013

medi Clavicle support Kidz

Clavicle support

Indications

- Post-traumatic and postoperative fractures
- Acromioclavicular dislocation

Contraindications

- Compression of the airways
- Compression of the neurovascular bundle

Mode of action

- Safe and reliable fixation of the shoulder girdle
- Promotes correct union of the fractured bone by holding the shoulder back



Measuring for size

size		article number
1	from 3 to 8 years of age	R.333.011
2	from 8 to 12 years of age	R.333.012

medi Armschlinge Kidz

Shoulder immobilisation brace

Indications

- Immobilisation after shoulder surgery
- Dislocation of the shoulder

Mode of action

- Relieves stresses on the shoulder, upper arm and elbow
- Adapts the shoulder and immobilises the arm



Measuring for size

size		article number
1	from 3 to 8 years of age	R.040.GK1
2	from 8 to 12 years of age	R.040.GK2

Note: children of about 12 of age and above can also use the medi Armschlinge.

size		article number
universal	from 12 years of age	R.040.G09

medi Shoulder sling Kidz

Shoulder joint support for immobilisation

Indications

- Shoulder contusions and sprains
- Shoulder dislocations
- Lesions of the rotator cuff
- Stable subcapital fractures of the humerus
- Scapular fractures

Mode of action

- Stable fixation of the shoulder and upper arm
- Additional stabilisation due to the chest strap



Measuring for size

size		article number
1	from 3 to 8 years of age	R.030.G29
2	from 8 to 12 years of age	R.030.G39

medi Epico ROMs Kidz

Elbow orthosis for mobilisation

Indications

- For conservative treatment of elbow dislocation without secondary injuries
- Orthotic care for fractures of the medial or lateral epicondyles
- Postoperative care after prosthesis implantation
- For functional follow-up treatment of stable fractures of the distal humerus
- For temporary immobilisation after arthrolysis of the elbow

Mode of action

- Orthosis for individually individual functional treatment
- Immobilises the elbow by fixing the upper and lower arm
- Flexion and extension stops in 10° steps (0° to 120°)



Measuring for size

model	arm length	article number
right	24 – 29 cm	R.031.3K9
left	24 – 29 cm	R.031.2K9

Note: The medi Epico ROMs can be used for arms 29 cm long and longer

model	arm length	article number
right	29 – 42 cm	R.032.302
left	29 – 42 cm	R.032.202

Manumed Kidz

Wrist support with stabilising splint

Indications

- Sprain / strain of the wrist
- Postoperative immobilisation of the wrist
- Inflammation and other disorders of the tendons in the wrist
- Inflammatory phases of rheumatic disorders

Mode of action

- Immobilisation and stabilisation of the wrist



Measuring for size: circumference of the wrist

size	wrist circumference	article number	
		left	right
1	10 – 11 cm	R.110.251	R.110.351
2	11 – 13 cm	R.110.252	R.110.352
3	13 – 14 cm	R.110.253	R.110.353
4	14 – 15 cm	R.110.254	R.110.354

Note: Manumed can be used for wrist circumferences of 15 cm and more.

colour	wrist circumference	article number	
		left	right
grey	15 – 18 cm	R.110.231	R.110.331
black	15 – 18 cm	R.110.201	R.110.301





medi Kidz orthopaedic
products for the lower limb

medi Classic Kidz

Knee immobilisation brace

Indications

- Postoperative immobilisation
- Preoperative pain conditions
- Status post luxation of the patella

Mode of action

- Immobilisation of the knee joint



Measuring for size: length of the orthosis

size		article number
1	21 cm	5.650.0K1
2	26 cm	5.650.0K2
3	31 cm	5.650.0K3
4	40 cm	5.650.0K4

medi ROM Kidz

Universal knee brace with flexion and extension stops and an immobilising function

Indications

- After an operation on the knee
- For immobilisation after knee injuries and for subsequent mobilisation
- For safe and reliable stabilisation of injured collateral ligaments or after surgery

Mode of action

- Collateral stabilisation of the knee joint
- Postoperative immobilisation
- Controlled mobilisation



Measuring for size: length of the orthosis

size		article number
Short	25 cm	G.180.0K1
Long	35 cm	G.180.0K2

M.step® Kidz

Ankle orthosis for stabilisation Indications

Indications

- Conservative and postoperative care of ligamentous injuries of the ankle joint (degree II and III)
- After sprains
- For chronic instability of the ankle joint

Mode of action

- Rigid half shells stabilise pronation and supination movements
- Autoadaptive cushion with integrated foam inlays for pressure relief



Measuring for size

size	article number
universal	G.092.009

medi Walker boot Kidz

Orthosis to immobilise the lower leg

Indications

- Postoperative treatment (rupture of the Achilles tendon, injuries of ligaments, soft-tissues and tendons, fracture of the fibula)
- Conservative treatment (fore- and midfoot and tarsus, Achilles tendon, fracture of the fibula)

Mode of action

- Fixes the ankle at 90°
- Rocker sole promotes the natural gait
- Simple wound inspection and hygiene



Measuring for size: shoe size

size	EU	article number
S	18 – 22	G.190.0K1
M	23 – 25	G.190.0K2
L	26 – 28	G.190.0K3
XL	29 – 32	G.190.0K4

Note: The medi Walker boot can be used for shoe sizes of 32 and more.

EU	article number
32 – 38	G.190.002





medi supports Perfection for children

Love of detail creates supports that provide relief for young and old patients alike.

Genumedi®


Knee support for soft tissue compression

Indications

- Mild instability
- Joint effusion and swelling
- Irritative conditions (chronic, post-traumatic, postoperative)
- Osteoarthritis and arthritis
- Patellofemoral Pain syndrome


Mode of action

- Integrated silicone pads for a massage effect
 - > Reduce swelling
- Proprioceptive compression
 - > Guides and stabilises the joints



Measuring for size

colour	thigh circumference	calf circumference	size	article number
silver	37 – 40 cm	25 – 28 cm	I	K.100.201
black	37 – 40 cm	25 – 28 cm	I	K.100.501
sand	37 – 40 cm	25 – 28 cm	I	K.100.801



Levamed

Support for compression of the soft tissues at the ankle joint

Indications

- Mild instability
- Sprains, contusions
- Swelling and joint effusion (osteoarthritis / arthritis)
- Irritative conditions (chronic, postoperative, post-traumatic)

Mode of action

- Integrated silicone pads for a massage effect
 - > Reduce swelling (oedema, haematomas)
- Proprioceptive compression
 - > Guides and stabilises the joints



Measuring for size: ankle circumference

colour	ankle circumference	size	article number
silver	17 – 20 cm	I	K.010.201
black	17 – 20 cm	I	K.010.501
sand	17 – 20 cm	I	K.010.801



Manumed active

Soft wrist support with stabilising splint

Indications

- Wrist injuries, e.g. sprains, strains
- Osteoarthritis of the wrist
- Instability
- Tendovaginitis
- Irritative conditions (post-traumatic, post-operative and rheumatic)

Mode of action

- Compressive fabric for massage effect and improved blood circulation
- The volar splint stabilises and reduces stress on the wrist
- The circular strap provides additional stability
- Pinpoint pressure relief due to the anatomically designed strap



Measuring for size: circumference wrist

colour	circumference wrist	size	article number	
			left	right
silver	13 – 15 cm	I	K.412.201	K.413.201
sand	13 – 15 cm	I	K.412.801	K.413.801



Headquarters
medi GmbH & Co KG
Medicusstraße 1
95448 Bayreuth
Germany
T +49 921 912-0
F +49 921 912-781
export@medi.de
www.medi.de/en

