



Respiratory Therapy

MINIPLUS

Comfortable Compressor Nebulizer

The Mini-Plus is a compressor nebulizer designed for easy portability, especially for people on the move. Chronic respiratory problems are unpredictable and it is important for people to have their treatment devices close at hand. The compact size, light weight and convenient carrying handle of the Mini-Plus facilitates portability and allows each person to keep their Mini-Plus accessible at all times.

Small and Easy-to-carry!



Respiratory Therapy



Features:

- It's easy to carry and ready to go right out of the box.
- Updated sleek design makes it easy to pack and store.
- Integrated handle makes it convenient for travel and mobility.
- Childproof cap prevents getting lost or being swallowed by a small child.
- Motor-piston drive for reliable, consistent performance.



Delivers aerosolized particles of 0.5-5 microns.
 Dispenses 3 cc's of medication in 10 minutes or less.
 Reliable output for efficient treatment.

Convenient Tube Holder design.

Flow outlet fitting, air filter and on/off switch are positioned around the periphery on the front panel and recessed for extra protection.

One button operation for easy set up within seconds.

Accessories Ordering Information:



Mouthpiece set : Anti-drool "T" mouthpiece with 6" reservoir tube and 7' supply tube.

Mouthpiece set and 2 masks : Anti-drool "T" mouthpiece with 6" reservoir tube, 7' supply tube and one each of adult & pediatric mask

Specifications:

Electrical Rating	100VAC,60Hz,1.2A	100VAC,50Hz,1.2A
(Note: Refer to the rating label on the product)	120VAC,60Hz,1.2A	230VAC,50Hz,0.6A
	230VAC,60Hz,0.6A	230V-240VAC, 50Hz, 0.6A
Maximum Compressor Pressure	≥ 40 psi	≥ 40 psi
Operation Flow Rate	≥ 7.0lpm	≥ 5.5 lpm
Nebulization Rate	0.3 ml/min	
Particle size	0.5 ~ 5 microns	
Dimensions (W x D x H)	14.0 x 18.8 x 10.0 cm / 5.5" x 7.4" x 3.9"	
Weight	1.8 kg / 4.0 lb	



APEX MEDICAL CORP.

9, Min Sheng St. Tu-Cheng,
 Taipei County, 236, Taiwan, R.O.C.
 TEL: 886-2-22683100
 FAX: 886-2-22686525
 E-mail:marketing@apexmedical.com.tw
 http://www.apexbrand.com

Distributed by:

