




# **NAGI™** Stent

*for pancreatic pseudocyst drainage*



■ Fluoroscopic image




*"An endoscopic approach that uses this new FCSEMS(Nagi) is feasible in the treatment of PFCs of both pancreatic pseudocysts and walled-off pancreatic necrosis"*

Natsuyo Yamamoto, MD, PhD [Vol.77 No.5 : GIE 2013]

▲ Placement of a transnasal drainage tube for irrigation and direct endoscopic necrosectomy through the stent

# Taewoong Niti-S™ Nagi™ Stent for pancreatic pseudocyst drainage




- **Wide and acute flare ends** prevent the risk of migration. Also, the smooth edges reduce the possibility of stent related luminal damages
- **Large diameter** of the stent (up to 16mm) enough to insert the scope for walled-off pancreatic necrosis
- **Retrieval String** at proximal end helps for repositioning or easy removal
- **Radiopaque markers** on the stent provide clear visualization under the fluoroscopy  
[\*3(three) at both ends & 2(two) in the middle]


## Ordering Information

### Endoscopic Approach

Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
BS1001FW	10	1	10.5	180
BS1002FW		2		
BS1003FW		3		
BS1201FW	12	1		
BS1202FW		2		
BS1203FW		3		
BS1401FW	14	1		
BS1402FW		2		
BS1403FW		3		
BS1601FW	16	1		
BS1602FW		2		
BS1603FW		3		



A large amount of pus came out through the stent



Conventional upper GI endoscope advances into the pseudocyst through the stent. Then performing of necrosectomy is using the snare forceps

► J Hepatobiliary Pancreat Sci.2012;DOI 10.1007

## Released Articles

\* **Preliminary report on a new, fully covered, metal stent designed for the treatment of pancreatic fluid collections**

Gastrointestinal Endoscopy 2013 Hiroyuki Isayama, MD, PhD, Natsuyo Yamamoto, MD, PhD, et al.

\* **Treatment of infected pancreatic pseudocysts using a novel, dedicated covered self-expandable metal stent (CSEMS) with an effective antimigration system**

Endoscopy 2012 I.Tarantino, M. Traina, L.Barresi, M.D.Pisa et al.

\* **New fully-covered self-expandable metal stent for endoscopic ultrasonography-guided intervention in infectious walled-off pancreatic necrosis**

J Hepatobiliary Pancreat Sci.2012 Takao Itoi.D. Nageshwar Reddy. Ichiro Yasuda



**Aquilant Interventional**

### UK Distributors

Aquilant Interventional  
Aquilant House, Units B1-B2  
Bond Close, Kingsland Business Park  
Basingstoke, Hampshire, RG24 8PZ

T: +44 (0) 845 6024 007  
F: +44 (0) 1256 365486  
E: [contactus@aquilantinterventional.com](mailto:contactus@aquilantinterventional.com)  
W: [www.aquilantinterventional.com](http://www.aquilantinterventional.com)

