



Devon™ Needle Counters

The safer way to count, contain and dispose of contaminated scalpel blades and suture needles.




COVIDIEN

positive results for life™

Counting on Safety

Added safety during handling and disposal





Blade Remover

Fast easy blade removal system with magnetic strips to ensure blade remains within counter.

The range of Devon™ Needle Counters offers an extensive choice of different capacity sizes and options that include configurable interiors using foam block, foam strips, magnets, or adhesive strips to secure, contain and accurately count needles and scalpel blades.

- Rigid plastic box – puncture resistant, closed on all sides
- Hinged containers – may be taken apart to save space
- Choose magnetic, foam or combination of both for 10 to 100 needles
- Halves can be detached for use in two separate locations
- Adhesive backing secures foam block to any surface
- Outer adhesive tabs secure container halves to any work surface



“peel off”
adhesive backing
for securing



High density foam
block ensures
needle retention

Boxlocks

The unique Boxlocks feature on all Devon™ needle counters allows you to lock after each count but re-open should a recount be required.



Devon™ Foam block needle counters – high density printed foam block



1. 31142352

Foam block with numbered magnet
maximum 100 count and blade removal



2. 31142386

Foam block with single magnet and
blade removal



3. 31142394

Foam block for small counts with
single magnet and blade removal



4. 31142360

Foam block with adhesive strip and
blade removal



5. 31142311

Foam block with single black magnet



6. 31142493

Foam block with single black magnet



7. 31142469

Double foam block (maximum 80 count)



8. 31142402

Foam block with adhesive strip



9. 31142501

Foam block with adhesive strip

Devon™ Magnetic needle counters – numbered magnetic strips



10. 31142113

Double black magnet with blade removal



11. 31142154

Large double black magnet



12. 31142147

Double black magnet



13. 31142196

Double black magnet



14. 31142139

Clear top, single black magnet



15. 31142410

Clear top, single black magnet



16. 31142238

Double white magnet with blade removal

Devon™ Foam strip needle counters



17. 31142451
Double foam strips with blade remover



18. 31142295
Single foam strips with magnet
and blade remover



19. 31142188
Large foam strips with single
black magnet



20. 31142303
Double foam strip



21. 31142212
Single foam strips with clear top

Devon™ Double blade removers with size 3 and size 4 blade removal systems



22. 31181459
Single foam block with magnet



23. 31181467
Double magnet



24. 31313425
Single foam strips with single magnet

Ordering Information

Devon™ Foam block needle counters – high density printed foam block

	Code	Description	Count / Capacity	Sales unit
1.	31142352	Foam block with numbered magnet maximum 100 count	100/100	8 box/64 case
2.	31142386	Foam block with single magnet and blade removal	20/40	12 box/96 case
3.	31142394	Foam block for small counts with single magnet	10/20	12 box/96 case
4.	31142360	Foam block with adhesive strip	20/40	12 box/96 case
5.	31142311	Foam block with single black magnet	40/70	8 box/64 case
6.	31142493	Foam block with single black magnet	20/40	12 box/96 case
7.	31142469	Double foam block (maximum 80 count)	80/80	8 box/64 case
8.	31142402	Foam block with adhesive strip	40/70	8 box/64 case
9.	31142501	Foam block with adhesive strip	20/40	12 box/96 case

Devon™ Magnetic needle counters – numbered magnetic strips

	Code	Description	Count / Capacity	Sales unit
10.	31142113	Double black magnet with blade removal	30/30	12 box/96 case
11.	31142154	Large double black magnet	60/60	8 box/64 case
12.	31142147	Double black magnet	30/30	12 box/96 case
13.	31142196	Double black magnet	20/20	21 box/168 case
14.	31142139	Clear top, single black magnet	30/30	8 box/64 case
15.	31142410	Clear top, single black magnet	10/10	21 box/168 case
16.	31142238	Double white magnet with blade removal	30/30	12 box/96 case

Devon™ Foam strip needle counters

	Code	Description	Count / Capacity	Sales unit
17.	31142451	Double foam strips with blade remover	30/30	12 box/96 case
18.	31142295	Single foam strips with magnet and blade remover	15/30	12 box/96 case
19.	31142188	Large foam strips with single black magnet	30/60	8 box/64 case
20.	31142303	Double foam strip	30/30	12 box/96 case
21.	31142212	Single foam strips with clear top	15/15	12 box/96 case

Devon™ Double blade removers with size 3 and size 4 blade removal systems

	Code	Description	Count / Capacity	Sales unit
22.	31181459	Single foam block with magnet	30/30	12 box/96 case
23.	31181467	Double magnet	25/25	12 box/96 case
24.	31313425	Single foam strips with single magnet	25/25	12 box/96 case



COVIDIEN

positive results for life™

COVIDIEN, COVIDIEN with logo, Covidien logo and *positive results for life* are U.S. and internationally registered trademarks of Covidien AG. Other brands are trademarks of a Covidien company. © 2011 Covidien.

MS_2011_MS_0249

COVIDIEN (UK) LIMITED
154 FAREHAM ROAD
GOSPORT HAMPSHIRE
PO13 0AS
UNITED KINGDOM (UK)

+44 (0)203 027 1757 [T]
+44 (0)203 027 1758 [F]

WWW.COVIDIEN.COM