

## Infant Flow LP Nasal CPAP System

- New Low Pressure dual jet generator to help reduce patient's work of breathing by up to 70%
- Delivers Nasal CPAP at a consistent pressure
- Softer and ergonomically designed prongs and masks for superior fixation
- Option between Headgear and Bonnet for ease of fit

### LP Generator

- Dual jet generator incorporating fluidic flip technology
- Low momentum impinging jets reduce WOB during inspiratory and expiratory cycle
- Infant Flow LP generator utilizes 80% less driving pressure
- Support cradle stabilizes generator and minimizes torque during infant movement



### LP Prongs

- Soft, comfort-fit prongs provide an effective seal and minimizes leaks
- High profile design reduces pressure on the infant's septum minimizing skin necrosis
- Available in 5 sizes to accommodate 450gm to 5kg infants



### Aquilant

Unit B1-B2 Bond Close, Kingsland Business Park  
Basingstoke, Hampshire RG24 8PZ  
T: +44 (0) 1256 365525 | F: +44 (0) 1256 365486  
E: [contactus@aquilantservices.com](mailto:contactus@aquilantservices.com) | [www.aquilant.net](http://www.aquilant.net)

## Infant Flow LP Nasal CPAP System

### LP Masks

- Flexible bellow allows mask to find natural position over infant's nose minimizing leaks and maintaining consistent nCPAP level
- Anatomically sized with deep nasal cavity to accommodate various nasal shapes and reduce pressure points
- Available in 5 colour coded sizes – Extra Small to Extra Large



### LP Headgear

- Intuitive comfort wrap design with adjustable straps for a custom fit helps minimize infant discomfort
- Fixation straps provide comfort and flexibility and keeps device away from infant's eyes
- Comfort wrap design secures generator at proper angle minimizing leaks
- Available in 6 sizes – Extra Small to Extra Large



### LP Bonnets

- Soft cotton offers a little stretch to conform to the infant's head and provide warmth
- Open top offers scalp access for IV therapy, phototherapy and ultrasound scans
- Available in 10 colour coded sizes to provide a comfortable fit



### Aquilant

Unit B1-B2 Bond Close, Kingsland Business Park  
Basingstoke, Hampshire RG24 8PZ  
T: +44 (0) 1256 365525 | F: +44 (0) 1256 365486  
E: [contactus@aquilantservices.com](mailto:contactus@aquilantservices.com) | [www.aquilant.net](http://www.aquilant.net)