

iVascular®

therapies for living

oceanus 14 pro

PTA balloon dilatation catheter

To reach the **most challenging** lesions
with the **highest visibility**

Best flexibility
and trackability*

Smallest balloon
diameters

Enhanced
visibility*

Lowest
profiles



Proprietary hydrophilic
durable coating



www.ivascular.global

oceanus 14 pro features

* Radiopaque markers	
	Ø ≤ 2.0mm, Length ≤ 20mm: 1 tungsten marker
	Ø ≤ 2.0mm, Length > 20mm: 2 tungsten markers
	Ø ≥ 2.5mm: 2 Pt/Ir markers

- > Over the wire catheter (OTW)
- > Usable Catheter length: 100 cm & 150 cm
- > Guidewire compatibility: 0.014"
- > Tip profile: 0.017"
- > Crossing profile: 0.021" to 0.033"
- > Nominal pressure: 7 atm
- > Rated Burst Pressure (RBP): 16 atm
- > Average Burst Pressure (ABP): 22 atm

Product with CE mark, certified by Notified Body 0318

Usable catheter length	Balloon Ø (mm)	Balloon Length (mm)									
		10	15	20	40	60	80	120	150	200	
100 cm	1.25	BP PCI4N 100 125 010	BP PCI4N 100 125 015	BP PCI4N 100 125 020	-	-	-	-	-	-	-
	1.5	BP PCI4N 100 150 010	BP PCI4N 100 150 015	BP PCI4N 100 150 020	BP PCI4N 100 150 040	BP PCI4N 100 150 060	BP PCI4N 100 150 080	-	-	-	-
	2.0	BP PCI4N 100 200 010	BP PCI4N 100 200 015	BP PCI4N 100 200 020	BP PCI4N 100 200 040	BP PCI4N 100 200 060	BP PCI4N 100 200 080	BP PCI4N 100 200 120	BP PCI4N 100 200 150	BP PCI4N 100 200 200	-
	2.5	-	-	-	BP PCI4N 100 250 040	BP PCI4N 100 250 060	BP PCI4N 100 250 080	BP PCI4N 100 250 120	BP PCI4N 100 250 150	BP PCI4N 100 250 200	-
	3.0	-	-	-	BP PCI4N 100 300 040	BP PCI4N 100 300 060	BP PCI4N 100 300 080	BP PCI4N 100 300 120	BP PCI4N 100 300 150	BP PCI4N 100 300 200	-
	3.5	-	-	-	BP PCI4N 100 350 040	BP PCI4N 100 350 060	BP PCI4N 100 350 080	BP PCI4N 100 350 120	BP PCI4N 100 350 150	BP PCI4N 100 350 200	-
	4.0	-	-	-	BP PCI4N 100 400 040	BP PCI4N 100 400 060	BP PCI4N 100 400 080	BP PCI4N 100 400 120	-	-	-
150 cm	1.25	BP PCI4N 150 125 010	BP PCI4N 150 125 015	BP PCI4N 150 125 020	-	-	-	-	-	-	-
	1.5	BP PCI4N 150 150 010	BP PCI4N 150 150 015	BP PCI4N 150 150 020	BP PCI4N 150 150 040	BP PCI4N 150 150 060	BP PCI4N 150 150 080	-	-	-	-
	2.0	BP PCI4N 150 200 010	BP PCI4N 150 200 015	BP PCI4N 150 200 020	BP PCI4N 150 200 040	BP PCI4N 150 200 060	BP PCI4N 150 200 080	BP PCI4N 150 200 120	BP PCI4N 150 200 150	BP PCI4N 150 200 200	-
	2.5	-	-	-	BP PCI4N 150 250 040	BP PCI4N 150 250 060	BP PCI4N 150 250 080	BP PCI4N 150 250 120	BP PCI4N 150 250 150	BP PCI4N 150 250 200	-
	3.0	-	-	-	BP PCI4N 150 300 040	BP PCI4N 150 300 060	BP PCI4N 150 300 080	BP PCI4N 150 300 120	BP PCI4N 150 300 150	BP PCI4N 150 300 200	-
	3.5	-	-	-	BP PCI4N 150 350 040	BP PCI4N 150 350 060	BP PCI4N 150 350 080	BP PCI4N 150 350 120	BP PCI4N 150 350 150	BP PCI4N 150 350 200	-
	4.0	-	-	-	BP PCI4N 150 400 040	BP PCI4N 150 400 060	BP PCI4N 150 400 080	BP PCI4N 150 400 120	-	-	-

● Introducer compatibility: 4F

Distributed by:



Aquilant

Aquilant House, Unit B1-B2 Bond Close,
Kingsland Business Park, Basingstoke,
Hampshire, RG248PZ, United Kingdom

T: +44 (0) 1256 365 454
F: +44 (0) 1256 365 486

E: contactus@aquilantservices.com
W: www.aquilantservices.com



Manufactured by:

Life Vascular Devices Biotech, S.L.
www.ivascular.global
info@ivascular.global