


Steerable ureteral stent set, **optimatic, multilength (OptiPur)**

Steerable ureteral stent with double pigtail for variable lengths and simple application system - facilitates the manipulation (turning, pushing, pulling) within the urinary tract

assembled system





ureteral stent set, open distal end, with guide wire

Catheter Ø	Length A	Length B	GW-Ø	length guid wire	Order-Nr.
CH	cm	cm	inch	cm	
6	24-32	20-28	.035	150	3044-2060
7	24-32	20-28	.035	150	3044-3060
8	24-32	20-28	.038	150	3044-4060

Box: 10 PCS

ureteral stent set, open distal end, without guide wire

Catheter Ø	Length A	Length B	Adapted to guide wire	Order-Nr.
CH	cm	cm	inch	
6	24-32	20-28	.035	3044-2065*
7	24-32	20-28	.035	3044-3065
8	24-32	20-28	.035	3044-4065*

Box: 10 PCS