

At Siemens Healthineers, our purpose is to enable healthcare providers to increase value by empowering them on their journey toward expanding precision medicine, transforming care delivery, and improving patient experience, all made possible by digitalizing healthcare.

An estimated 5 million patients globally benefit every day from our innovative technologies and services in the areas of diagnostic and therapeutic imaging, laboratory diagnostics, and molecular medicine, as well as digital health and enterprise services.

We are a leading medical technology company with over 120 years of experience and 18,000 patents globally. Through the dedication of more than 50,000 colleagues in 75 countries, we will continue to innovate and shape the future of healthcare.

Acute Care, ADVIA, ADVIA Centaur, CLINITEK Advantus, CLINITEK Status, DCA Vantage, epoc, RAPIDComm, RAPIDLab, RAPIDPoint, Stratus, Xprecia Stride, and all associated marks are trademarks of Siemens Healthcare Diagnostics Inc., or its affiliates. All other trademarks and brands are the property of their respective owners.

Product availability may vary from country to country and is subject to varying regulatory requirements. Please contact your local representative for availability.

Siemens Healthineers Headquarters
Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen, Germany
Phone: +49 9131 84-0
siemens-healthineers.com

Published by
Siemens Healthcare Diagnostics Inc.
Point of Care Diagnostics
2 Edgewater Drive
Norwood, MA 02062-4637
USA
Phone: +1 781-269-3000

Point of Care Informatics

Open Connectivity Device Interface List

siemens-healthineers.com/poc-informatics

Vendor	Device	Connection	Application
77 Elektronika	DocUReader 2 Pro	Bidirectional Download operator/patient Configuration Positive patient identification	Urinalysis
Abaxis	Piccolo Xpress	Unidirectional	Clinical chemistry
Abbott Diabetes Care	FreeStyle Precision Pro (W-LAN)	Bidirectional Download operator/patient Configuration Positive patient identification	Glucose Ketone
	PCx Precision	Unidirectional	Glucose
	Precision Xceed Pro	Bidirectional Download operator/patient Configuration Positive patient identification	Glucose Ketone
Abbott Point of Care	Afinion AS 100 Analyzer	Bidirectional Configuration Download operator	Multitest
	Afinion 2	Bidirectional Configuration Download operator	Multitest
	Alere Reader	Unidirectional	Multitest
	Alere i	Unidirectional	Multitest
	CELL-DYN Emerald	Unidirectional	Hematology
	Cholestech LDX System	Unidirectional	Multitest
	Pima	Unidirectional	HIV
	i-STAT 1	Bidirectional (via i-STAT-DE-DMI software) Download operator/patient Configuration Positive patient identification	Multitest
	i-STAT Alinity	Bidirectional Download operator/patient Configuration Positive patient identification	Multitest
	i-STAT/DE/DMI	Bidirectional Download operator/patient Configuration Positive patient identification	Supplier software
Acon Laboratories	Mission U120	Unidirectional	Urinalysis
A.MENARINI Diagnostics	Zenit RA	Unidirectional	Immunoassay
Analyticon	CombiScan 100	Unidirectional	Urinalysis

Vendor	Device	Connection	Application
Analyticon	Urilyzer 100	Bidirectional Download operator	Urinalysis
	Urilyzer 100 Pro	Bidirectional Configuration Download operator/patient Positive patient identification	Urinalysis
Arkray	ADAMS HA-8160	Unidirectional	HbA1c
	Spotchem EZ SP-4430	Unidirectional	Clinical chemistry (no QC!)
Beckman Coulter	Coulter Ac-T diff	Unidirectional	Hematology (no QC!)
	Coulter Ac-T diff2	Unidirectional	Hematology
	Coulter Ac-T 5diff	Unidirectional	Hematology
	Coulter Ac-T 8diff	Unidirectional	Hematology
	Coulter Ac-T 10diff	Unidirectional	Hematology
	DxH 500	Unidirectional	Hematology
Becton Dickinson	BD Veritor Plus	Unidirectional	Multitest
biosurfit	spinit	Unidirectional	Multitest
Boule Medical	Medonic M-series M16/M20	Unidirectional	Hematology
	Medonic M-series M16/M20M	Unidirectional	Hematology
	Medonic M-series M16M/ M20M	Unidirectional	Hematology
	Medonic M-series M16C/M20C	Unidirectional	Hematology
	Medonic M-series M16S/M20S	Unidirectional	Hematology
	Swelab Alfa Basic	Unidirectional	Hematology
	Swelab Alfa Standard	Unidirectional	Hematology
	Swelab Alfa Cap Piercer	Unidirectional	Hematology
	Swelab Alfa Plus	Unidirectional	Hematology
Swelab Alfa AutoSampler	Unidirectional	Hematology	
CARE diagnostica	eco Solo II	Bidirectional Configuration Download operator/patient	Glucose
Cepheid	GeneXpert PCR	Unidirectional	Molecular diagnostics
Chirus	DXpress Reader	Unidirectional	Multitest
concile	Q100	Unidirectional	Immunoassay
Dräger	Jaundice Meter JM-105	Unidirectional	Bilirubin
Drucker Diagnostics	QBC STAR	Unidirectional	Hematology
Enigma Diagnostics Ltd	Enigma MiniLab	Unidirectional	Multitest
Eurolyser Diagnostica	Eurolyser cube	Unidirectional	Multitest
	Eurolyser solo	Unidirectional	Multitest
	Eurolyser smart	Unidirectional	Multitest
Fuji	FUJI DRI-CHEM 4000	Unidirectional	Clinical chemistry
	FUJI DRI-CHEM 7000	Unidirectional	Clinical chemistry
	FUJI DRI-CHEM NX500	Unidirectional	Clinical chemistry
GenePOC Diagnostics	Revogene	Bidirectional Download operator Configuration	Multitest

Vendor	Device	Connection	Application
Haemonetics	TEG 6s	Bidirectional Positive patient identification	Hematology
Helena Laboratories	Cascade Abrazo	Bidirectional Configuration	Coagulation
	Cascade POC	Bidirectional Download operator	Coagulation
HemoCue	Glucose 201 DM	Bidirectional Download operator/patient Configuration Positive patient identification	Glucose
	Glucose 201 DM RT	Bidirectional Download operator/patient Configuration Positive patient identification	Glucose
	HbA1c 501	Unidirectional	Hemoglobin
	Hb 201 DM	Bidirectional Download operator/patient Configuration Positive patient identification	Hemoglobin
	WBC Diff	Unidirectional	Hematology
	HemoSonics	Quantra	Bidirectional Download operator/patient
Hologic	Rapid fFN	Unidirectional	Immunoassay
	PeriLynx	Unidirectional	Immunoassay
HORIBA	ABX Micros 45	Unidirectional	Hematology
	ABX Micros 60	Unidirectional	Hematology
	ABX Micros ES 60	Unidirectional	Hematology
	ABX Micros CRP	Unidirectional	Multitest
	ABX Micros CRP 200 (+ ABX e-Sat)	Unidirectional	Multitest
	ABX Microsemi CRP	Unidirectional	Multitest
	ABX Pentra 60	Unidirectional	Hematology
	ABX Pentra 60 C+	Unidirectional	Hematology
Yumizen H500	Unidirectional	Hematology	
LifeHealth	IRMA TRUPOINT	Bidirectional Positive patient identification	Blood gas
Lite-On Technology Corp.	skyla HB1	Bidirectional Download operator	Clinical chemistry
	skyla Hi	Bidirectional Download operator	Clinical chemistry
Macherey-Nagel	Uryxxon Relax ^{NEW}	Bidirectional Download operator Positive patient identification	Urinalysis
Medtronic	ACT Plus	Unidirectional	Coagulation
	HMS Plus	Unidirectional	Coagulation
Mitsubishi	PATHFAST	Unidirectional	Immunoassay
Nihon Koden	Celltac Alpha ^{NEW} (MEK 6500/6510)	Unidirectional	Hematology

Vendor	Device	Connection	Application
Nihon Kodon	Celltac Chemi ^{NEW} (CHM-4100K)	Unidirectional	Hematology
Norma Diagnostika	Norma Icon 3/5	Unidirectional	Hematology
Nova Biomedical	NovaNet*	Bidirectional Download operator/patient Positive patient identification	Supplier software
	Stat Profile pHox/pHOx Plus/Ultra	Bidirectional Configuration Positive patient identification	Blood gas
	Stat Profile CCX	Bidirectional Configuration Positive patient identification	Blood gas
	Stat Profile Prime	Unidirectional	Blood gas/multitest
	Stat Profile Prime Plus ^{NEW}	Bidirectional Download operator/patient Positive patient identification (via NovaNet)	Blood gas
	StatStrip Glucose	Bidirectional Download operator/patient Configuration Positive patient identification	Glucose
	StatStrip Glucose/Ketone	Bidirectional Download operator/patient Configuration Positive patient identification	Glucose Ketone
	StatStrip HH ^{NEW}	Bidirectional Download operator/patient Positive patient identification (via NovaNet)	Hemoglobin
	StatStrip Lactate	Bidirectional Download operator/patient Configuration Positive patient identification	Lactate
	StatSensor Creatinine	Bidirectional Download operator/patient Configuration Positive patient identification	Creatinine
OPTI Medical	OPTI CCA-TS2	Bidirectional Download operator/patient	Blood gas
Orion Diagnostica	QuikRead 101	Unidirectional	CRP
	QuikRead go	Bidirectional Download operator/patient Configuration	CRP
Philips	BiliChek	Unidirectional	Bilirubin
Quidel	Sofia	Unidirectional	Immunoassay
	Triage MeterPlus	Bidirectional Download operator	Multitest
	Triage MeterPro	Bidirectional Download operator	Multitest
Radiometer	ABL5xx series	Bidirectional Positive patient identification	Blood gas

Vendor	Device	Connection	Application
Radiometer	ABL6xx series	Bidirectional Positive patient identification	Blood gas
	ABL7xx series	Bidirectional Positive patient identification	Blood gas
	ABL8xx series	Bidirectional Positive patient identification	Blood gas
	ABL77	Bidirectional Positive patient identification	Blood gas
	ABL80 FLEX series	Bidirectional Positive patient identification	Blood gas
	ABL90 FLEX/FLEX Plus	Bidirectional Positive patient identification	Blood gas
	AQT90 FLEX	Bidirectional Positive patient identification	Multitest
	AQURE	Bidirectional Positive patient identification Download operator	Supplier software
	NPT7	Bidirectional Positive patient identification	Blood gas
	RADIANCE	Bidirectional Positive patient identification Download operator	Supplier software
Response Biomedical	RAMP Reader	Unidirectional	Cardiac markers
	RAMP 200	Bidirectional Download operator Configuration	Multitest
Roche Diagnostics	ACCU-CHEK Inform	Bidirectional Download operator/patient Configuration Positive patient identification	Glucose
	ACCU-CHEK Inform II	Bidirectional Download operator/patient Configuration Positive patient identification	Glucose
	CARDIAC reader	Bidirectional Configuration	Cardiac markers
	cobas b101	Bidirectional Download operator/patient Configuration Positive patient identification	Multitest
	cobas b 123	Bidirectional Download operator/patient Positive patient identification	Blood gas
	cobas b 221/121	Unidirectional	Blood gas
	cobas h 232	Bidirectional Download operator/patient Positive patient identification Configuration	Cardiac markers
	cobas c 111	Unidirectional	Clinical chemistry
	cobas liat	Bidirectional Download operator	Multitest

Vendor	Device	Connection	Application	
Roche Diagnostics	cobas 6000	Unidirectional	Clinical chemistry	
	cobas IT 1000	Bidirectional Download operator/patient Positive patient identification	Supplier software	
	CoaguChek Pro II	Bidirectional Download operator/patient Configuration Positive patient identification	Coagulation	
	CoaguChek XS	Unidirectional	Coagulation	
	CoaguChek XS Plus	Bidirectional Download operator/patient Configuration Positive patient identification	Coagulation	
	CoaguChek XS Pro	Bidirectional Download operator/patient Configuration Positive patient identification	Coagulation	
	OMNI 3,4,5,6,7,8,9	Bidirectional Positive patient identification	Blood gas	
	OMNI Classic	Bidirectional Positive patient identification	Blood gas	
	OMNI S	Bidirectional Positive patient identification	Blood gas	
	Reflotron Plus	Unidirectional	Clinical chemistry	
	Reflotron Sprint	Unidirectional	Clinical chemistry	
	Urisys 1100	Bidirectional (with firmware version 5.x or higher) Download operator	Urinalysis	
	Samsung	LABGEO ^{HC10}	Bidirectional Download operator	Hematology
		LABGEO ^{IB10}	Bidirectional Download operator	Immunoassay
LABGEO ^{PT10}		Bidirectional Download operator	Clinical chemistry	
SD Biosensor	Standard F200	Unidirectional	Multitest	
	Standard F2400	Unidirectional	Multitest	
Siemens Healthcare Diagnostics	ADVIA [®] 1650 Clinical Chemistry System	Unidirectional	Clinical chemistry	
	ADVIA [®] 360 Hematology System	Unidirectional	Hematology	
	ADVIA Centaur [®] CP/XP Immunoassay System	Unidirectional	Immunoassay	
	CLINITEK Advantus [®] Urine Chemistry Analyzer	Unidirectional	Urinalysis	
	CLINITEK Status [®] Analyzer	Unidirectional	Urinalysis	
	CLINITEK Status [®] + Urine Chemistry Analyzer	Bidirectional Download operator Configuration	Urinalysis	

Vendor	Device	Connection	Application	
Siemens Healthcare Diagnostics	DCA Vantage [®] Analyzer	Bidirectional Download operator	HbA1c	
	epoc [®] Blood Analysis System	Bidirectional Download operator Configuration Positive patient identification	Blood gas	
	RAPIDComm [®] Data Management System	Bidirectional Download operator Positive patient identification	Supplier software	
	RAPIDLab [®] 8xx series Blood Gas Systems	Bidirectional Download operator Configuration Positive patient identification	Blood gas	
	Atellica [®] UAS 1500 ^{NEW} Automated Urinalysis System (Clinitek Novus [®] Analyzer)	Unidirectional	Urinalysis	
	Atellica [®] UAS 800 ^{NEW} Automated Urinalysis System (Clinitek Novus [®] Analyzer)	Unidirectional	Urinalysis	
	RAPIDLab [®] 248/348/348EX Blood Gas Systems	Bidirectional	Blood gas	
	RAPIDLab 1200 series Blood Gas Systems	Bidirectional Download operator Configuration Positive patient identification	Blood gas	
	RAPIDPoint [®] 340/350 Blood Gas Systems	Bidirectional Download operator Configuration Positive patient identification	Blood gas	
	RAPIDPoint 400, 405 Blood Gas Systems	Bidirectional Download operator Configuration Positive patient identification	Blood gas	
	RAPIDPoint 500 Blood Gas Systems	Bidirectional Download operator Configuration Positive patient identification	Blood gas	
	Stratus [®] CS Acute Care [™] Diagnostic System	Unidirectional	Cardiac markers	
	Stratus CS 200 Acute Care Diagnostic System	Bidirectional Download operator	Cardiac markers	
	Xprecia Stride [™] Coagulation Analyzer	Bidirectional Download operator Configuration	Coagulation	
	Sphere Medical	Proxima	Unidirectional	Blood gas
	Stat-Dx	DiagCORE	Unidirectional	Multitest
Sysmex	Alifax Roller	Unidirectional	Blood gas	
	CA-560	Unidirectional	Coagulation	
	CA-600	Unidirectional	Coagulation	
Sysmex	K-1000	Unidirectional	Hematology	
	K-4500	Unidirectional	Hematology	

Vendor	Device	Connection	Application
Sysmex	KX-21	Unidirectional	Hematology
	poch-100i	Unidirectional	Hematology
	XN-350, XN-450, XN-550	Unidirectional (via SUIT Manager)	Hematology
	XP-300	Unidirectional	Hematology
	XS-800i	Unidirectional (via SUIT Manager)	Hematology
	XS-1000i	Unidirectional (via SUIT Manager)	Hematology
Thermo Fisher Scientific	B·R·A·H·M·S Direct Reader	Bidirectional Configuration	Multitest
Veda.Lab	Easy Reader Plus	Bidirectional Download operator	Multitest
Werfen	ACL TOP	Unidirectional	Coagulation
	AVOXimeter 1000e	Unidirectional	Blood gas
	AVOXimeter 4000	Unidirectional	Blood gas
	BIO-FLASH	Unidirectional	Immunoassay
	GEM Premier 3000	Bidirectional Positive patient identification	Blood gas
	GEM Premier 3500	Bidirectional Positive patient identification	Blood gas
	GEM Premier 4000/5000	Bidirectional Download operator Positive patient identification	Blood gas
	GEM OPL	Unidirectional (via GEMweb Plus)	CO-oximetry
	GEM PCL	Unidirectional	Coagulation
	GEM PCL Plus	Unidirectional	Coagulation
	GEMweb Plus	Bidirectional Download operator Positive patient identification	Supplier software
	HEMOCHRON Jr. Signature	Unidirectional	Coagulation
	HEMOCHRON Jr. Signature + (Plus)	Unidirectional	Coagulation
	HEMOCHRON Signature Elite	Bidirectional Download operator Configuration	Coagulation
	HEMOCHRON Response	Unidirectional	Coagulation
	IL Synthesis 35	Bidirectional Positive patient identification	Blood gas
	IL Synthesis 45	Bidirectional Positive patient identification	Blood gas
	ProTime InRhythm	Unidirectional	Coagulation
	ROTEM delta	Bidirectional Positive patient identification	Coagulation
	ROTEM gamma	Bidirectional Positive patient identification	Coagulation
ROTEM sigma	Bidirectional Positive patient identification	Coagulation	