



150 Cassette Printer

Thermo Scientific Shandon PrintMate Cassette Printer Specifications

Description	Order Number
Shandon PrintMate 150 Cassette Printer, 100-240 volts, 50/60 Hz, 100 VA, 342 BTU <i>Included with Product</i> Two-hopper PrintMate, USB cable, power cable, two hoppers, one roll of tape (installed), basic software CD, operator guide	A79500001
Shandon PrintMate 450 Cassette Printer, 100-240 volts, 50/60 Hz, 100 VA, 342 BTU <i>Included with Product</i> Six-hopper PrintMate, USB cable, power cable, six hoppers, one roll of tape (installed), basic software CD, operator guide	A79500002
Shandon PrintMate 900 Cassette Printer, 100-240 volts, 50/60 Hz, 100 VA, 342 BTU <i>Included with Product</i> Twelve-hopper PrintMate, USB cable, power cable, twelve hoppers, one roll of tape (installed), basic software CD, operator guide	A79500003

Thermo Scientific Shandon PrintMate Cassette Printing Systems

Description	Order Number
Shandon PrintMate 150 Cassette Printing System, 100-240 volts, 50/60 Hz, 150 VA, 513 BTU <i>Included with Product</i> Two-hopper PrintMate, USB cable, power cable, two hoppers, one roll of tape (installed), basic software CD, operator guide, PC with monitor, digital imager	A79500004
Shandon PrintMate 450 Cassette Printing System, 100-240 volts, 50/60 Hz, 150 VA, 513 BTU <i>Included with Product</i> Six-hopper PrintMate, USB cable, power cable, six hoppers, one roll of tape (installed), basic software CD, operator guide, PC with monitor, digital imager	A79500005
Shandon PrintMate 900 Cassette Printing System, 100-240 volts, 50/60 Hz, 150 VA, 513 BTU <i>Included with Product</i> Twelve-hopper PrintMate, USB cable, power cable, twelve hoppers, one roll of tape (installed), basic software CD, operator guide, PC with monitor, digital imager	A79500006

Thermo Scientific Shandon PrintMate Accessories

Description	Order Number
Collection system	A79510014
Collection trays (box of 4)	A79530063
Hopper	A79510004
Hot foil tape, box of 5 (roll of tape prints approx. 20,000 cassettes)	A79510013
Left exit chute	A79520007
Right exit chute	A79520008
Interface software package	A79510019
Protocol software package	A79510020

Dimensions for PrintMate 900 (without cassette tubes): (D x W x H): 15.75" x 19" x 9" (40 x 48 x 23cm)

Dimensions for PrintMate 900 (with cassette tubes): (D x W x H): 15.75" x 19" x 35.4" (40 x 48 x 90cm)

Weight: 26.5lbs (14kg)

© 2007 Thermo Fisher Scientific Inc. All rights reserved. LabWriter is a trademark of Raymond A. Lamb Limited. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

United States
1-800-522-7270
(269) 372-2674 fax

United Kingdom
0800 0189396 +44 (0) 1928 562501
+44 (0) 1928 562627 fax

www.thermo.com/pathology

M41000 R11/07

Thermo Scientific Shandon PrintMate™ Automated Cassette Printing System



Thermo Scientific Shandon PrintMate 450 Cassette Printer

- Flexibility
- Connectivity
- Reliability
- Speed
- Security

The Thermo Scientific Shandon PrintMate Saves Time, Money, and Reduces Specimen Mix-ups in the Lab

Flexibility

The Shandon PrintMate can grow with your laboratory workload. From two to six or six to twelve hoppers, your new labeler will never be obsolete by your changing workload.

Small footprint—fits easily into any workflow pattern

The Shandon PrintMate was designed with you in mind. It combines ruggedness with the capacity to work reliably as a batch printer in high volume laboratories.

- Small enough to be used at a single grossing workstation.
- Intelligent enough to be shared between two workstations while maintaining specimen integrity.

Connectivity

- Compatible with all LIS interfaces to Thermo Scientific printers.
- Compatible with all Thermo Scientific Shandon MicroWriter™ software.
- Connects via USB or RJ45, ethernet connection.

Collection system

The dual track collection system collects up to 100 cassettes in organizer trays. The trays are easily stacked and simplify matching cassettes to specimen containers at the grossing workstation.

Print-on-demand

The LabWriter™ printing software allows you to easily enter data for large batches or single cassettes using standard or special numbering configurations.

- Uses bar codes from requisitions to automatically print cassettes—eliminating transcription errors.
- Collection chute enables you to place one Shandon PrintMate between two technicians without risk of cassette mix-up.
- Easily prints additional cassettes.
- Print one case at a time, eliminating the risk of cassette mix-up.
- Fits easily into the grossing workstation.
- Provides hard copy record of printed cassettes.
- Complies with lean workflow requirements.

CAP Requirements

As requirements for cassette labeling change, Shandon PrintMate can help your laboratory meet the newer regulations and maintain your laboratory accreditation. The latest changes, “requiring a second unique identifier” can be easily applied to the marking face of the cassette.



Reliability—No Jams

The patented dispense system is able to reliably dispense a wide range of standard cassettes year after year without jamming or misfeeding. The Shandon PrintMate is the only labeler to let you choose your favorite cassettes and use them reliably.



History in the Marketplace, A Name you can Trust...

Thermo Scientific is a name you can trust. For over twenty years, Thermo Scientific instruments have been manufactured to help improve sample identification and result integrity.

We continue to use the same basic printing method that has been reliable and long lasting. Our new printing system builds on this success with newer, faster printhead technology that you trust and rely on for accurate, readable sample identification year after year.



Getting Accurate Results Doesn't Have to Slow you Down!

Whatever workflow your lab demands, the choice is yours...
...you select what is best suited for your laboratory workflow.

Thermo Scientific Shandon PrintMate 900 Cassette Printer

- **Thermo Scientific Shandon PrintMate 150**
Holds up to 150 cassettes
- **Thermo Scientific Shandon PrintMate 450**
Holds up to 450 cassettes
- **Thermo Scientific Shandon PrintMate 900**
Holds up to 900 cassettes



Specimen Security—Reduced Errors

The high quality print of the Shandon PrintMate means that every cassette is clearly marked. The higher level of specimen identification, accuracy and sample integrity contributes to improved patient care.



High-Quality Printing System

Prints easy to read, permanent black type. The print is perfect for printing linear and 2D bar codes as well as graphics and logos.

- No ink cartridges
- No cleaning cartridges
- No special coatings
- No special curing processes