

*Direct to the  
right choice*

## Diagnostics Electrodes



**COVIDIEN**

*positive results for life™*

### Experts on electrodes

Our complete range of high-quality electrodes offers options for all applications. We can help you select the right electrode for your specific purposes.



# Covidien recommends...

## For stress test/Holter ECG – H34LG

50x45 mm

### Key highlights

- Liquid gel
- Suitable for long- and short-term diagnostics
- Extra strong adhesive
- Snap adapter – A34 attaches electrode to 4 mm plug cables



### Further details

- Liquid gel results in very low skin impedance and enhances safety due to short recovery time of the ECG electrode
- Pressure sensitive adhesive is extra strong and skin friendly – particularly suitable for stress or Holter testing
- Foam backing is water resistant and washable
- High quality Ag/AgCl sensor
- Press stud connection

| Product | Order code | PC/Pouch | PC/Box | PC/Carton |
|---------|------------|----------|--------|-----------|
| H34LG   | 31.1345.21 | 30       | 300    | 3600      |

# Covidien recommends...

## For resting ECG – Q-Trace™ 5400

22x31 mm

### Key highlights

- Universal size
- Radiolucent
- Full surface hydrogel – no adhesive
- 10 electrodes on a plate for 1 test



### Further details

- Disposable tab style, short-term diagnostic electrode
- Universal size, suitable for adults and infants
- Radiolucent – no need to change electrodes for x-ray testing, cost effective
- Full surface hydrogel – conductive adhesive hydrogel covering all the contact surface
- No adhesive means electrode is non allergic and repositionable
- 10 electrodes on a plate for 1 test
- Clear tape (polyester or polyethylene) backing
- High quality Ag/AgCl sensor
- Suitable for crocodile connectors; can be attached to 4 mm plug cables using A34 adapter

| Product       | Order code | PC/Pouch | PC/Box | PC/Carton |
|---------------|------------|----------|--------|-----------|
| Q-Trace™ 5400 | 31433538   | 100      | 1000   | 4000      |



| Product     | Order code | PC/Pouch | PC/Box | PC/Carton |
|-------------|------------|----------|--------|-----------|
| A34 adapter | 22660A     | 10       | -      | -         |

# Quick guide

## Overview of Diagnostic Electrodes

| Product            | Order code | Size      | Recommended use              | Gel        |           |      | Material   |            |       |     | X-ray | MRI | Connector |      |  | Other features |
|--------------------|------------|-----------|------------------------------|------------|-----------|------|------------|------------|-------|-----|-------|-----|-----------|------|--|----------------|
|                    |            |           |                              | Liquid gel | Hydro gel | Foam | Clear tape | Micro pore | Cloth | Tab |       |     | 4 mm plug | Stud |  |                |
| H1354LG            | 31.1414.21 | Ø 40 mm   | Resting ECG                  | ●          |           | ●    |            |            |       |     |       |     | ●         |      | Strong, skin friendly adhesive                             |                |
| H664LG             | 31.1664.21 | Ø 55 mm   | Stress test/ Holter          | ●          |           |      |            |            |       |     |       |     | ●         |      | Skin friendly adhesive                                     |                |
| H9145G             | 31.1914.21 | 57x34 mm  | Resting/ Stress test/ Holter |            | ●         | ●    |            |            |       |     |       |     | ●         |      |  |                |
| H914TSG            | 31.1914821 | 57x34 mm  | Resting/ Stress test         |            | ●         | ●    |            |            |       | ●   | ●     |     | ●         |      |  |                |
| H34LG              | 31.1345.21 | 50x45 mm  | Stress test/ Holter          | ●          |           | ●    |            |            |       |     |       |     |           | ●    | Extra strong, skin friendly adhesive                       |                |
| H135LG             | 31.1413.21 | Ø 40 mm   | Stress test/ Holter          | ●          |           | ●    |            |            |       |     |       |     |           | ●    | Strong, skin friendly adhesive                             |                |
| H66LG              | 31.1663.21 | Ø 55 mm   | Stress test/ Holter          | ●          |           | ●    |            |            |       |     |       |     |           | ●    | Strong, skin friendly adhesive                             |                |
| H98LG              | 31.1985.21 | Ø 55 mm   | Stress test/ Holter          | ●          |           |      |            | ●          |       |     |       |     |           | ●    |  |                |
| H99LG              | 31.1995.21 | Ø 55 mm   | Stress test/ Holter          | ●          |           |      |            |            | ●     |     |       |     |           | ●    |  |                |
| H715G              | 31.1715.21 | 58x51 mm  | Stress test/ Holter          |            | ●         | ●    |            |            |       |     |       |     |           | ●    | Strong, skin friendly adhesive<br>Backing has cable fixing |                |
| Q-Trace™ 5400      | 31433538   | 22x31 mm  | Resting ECG                  |            | ●†        |      | ●          |            |       | ●   | ●     | ●   |           |      |  |                |
| Q-Trace™ Gold 5500 | 30807732   | 23x25 mm  | Resting ECG                  |            | ●†        |      | ●          |            |       | ●   | ●     | ●   |           |      |  |                |
| A34 adapter        | 22660A     |           | Resting/ Stress test         |            |           |      |            |            |       |     |       |     | ●         | ●    |  |                |
| Nu-Prep™*          | 30806726   | 114g tube | Universal                    |            |           |      |            |            |       | ●   | ●     |     |           |      |  |                |

● Covidien recommended electrodes

□ Additional electrode options

† Full surface hydrogel, no adhesive – non allergic, and allows repositioning and repeated use of an electrode on one patient.



FOR MORE INFORMATION AND SUPPORT, PLEASE CONTACT YOUR LOCAL COVIDIEN REPRESENTATIVE

COVIDIEN (UK) COMMERCIAL LTD.  
4500 PARKWAY  
WHITELEY, FAREHAM  
HAMPSHIRE, PO15 7NY  
UNITED KINGDOM

+44(0)203 027 1757 (T)      WWW.COVIDIEN.COM  
+44(0)203 027 1758 (F)  
UKSALES@COVIDIEN.COM

COVIDIEN, COVIDIEN with logo, Covidien logo and *positive results for life* are U.S. and internationally registered trademarks of Covidien AG. TM\* Trademark of its respective owner. Other brands are trademarks of a Covidien company. ©2013 Covidien.

MS\_2011\_SW\_0283