





# Biliary stent




# Taewoong Niti-S™ S Biliary stent

## for benign and malignant biliary stricture



### Feature

- Silicone coated on both inner and outer surface :
  - to prevent the risk of tumor ingrowth
  - to help smooth bile flow
- Retrieval string facilitates safe and smooth removal



### Radiopaque marker : 3(three) at both ends & 2(two) in the middle

### Ordering Information

#### Endoscopic Approach

#### Percutaneous Approach

#### ► Fully covered

Code	Stent		Delivery		
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)	
BS0604F	6	4	8.5	180	
BS0605F		5			
BS0606F		6			
BS0607F		7			
BS0608F		8			
BS0609F		9			
BS0610F		10			
BS0612F		12			
BS0804F		8			4
BS0805F					5
BS0806F					6
BS0807F					7
BS0808F	8				
BS0809F	9				
BS0810F	10				
BS0812F	12				
BS1004F	10	4			
BS1005F		5			
BS1006F		6			
BS1007F		7			
BS1008F		8			
BS1009F		9			
BS1010F		10			
BS1012F		12			

#### ► Fully covered

Code	Stent		Delivery		
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)	
TS0604F	6	4	8.5	50	
TS0605F		5			
TS0606F		6			
TS0607F		7			
TS0608F		8			
TS0609F		9			
TS0610F		10			
TS0612F		12			
TS0804F		8			4
TS0805F					5
TS0806F					6
TS0807F					7
TS0808F	8				
TS0809F	9				
TS0810F	10				
TS0812F	12				
TS1004F	10	4			
TS1005F		5			
TS1006F		6			
TS1007F		7			
TS1008F		8			
TS1009F		9			
TS1010F		10			
TS1012F		12			

\* Available both bare and uncovered stent

Coding : BS\*\*\_B: both bare / B\*\*\_: Uncovered

### Released Article

\* Polyurethane-Covered Self-Expandable Nitinol Stent for Malignant Biliary Obstruction : Preliminary Results Cardio Vascular and Int  
Young-Min Han, Seung-Bae Hwang, Soo-Teik Lee, Jeong-Min Lee, Gyung-Ho Chung

#### UK Distributors



**Aquilant Interventional**

Aquilant Interventional

Aquilant House, Units B1-B2

Bond Close, Kingsland Business Park

Basingstoke, Hampshire, RG24 8PZ

T: +44 (0) 845 6024 007

F: +44 (0) 1256 365486

E: [contactus@aquilantinterventional.com](mailto:contactus@aquilantinterventional.com)

W: [www.aquilantinterventional.com](http://www.aquilantinterventional.com)

