






# S Esophageal stent



Distal release



Proximal release



# Taewoong Niti-S™ S Esophageal stent for benign and malignant esophageal strictures

## ● Features

- **Fixed cell with braided construction**
  - High flexibility and optimal radial force
  - Both head ends (8mm larger than trunk) help to minimize migration
- **Silicone covering and soft round ends**
  - Reduce tissue ingrowth and hyperplasia reaction
- **Visible green suture for easy removal**
- **Radiopaque marker : 4(four) at both ends & 2(two) in the middle**

## ● Ordering Information

Distal Release					
Product	Code	Stent		Delivery	
		Diameter (mm)	Length (cm)	Profile (Fr)	Usable length (cm)
Fully covered	ES**__F	16, 18	6, 8, 10, 12, 14, 15	16	70
		20, 22		20	
		24, 28		22	
Both ends bare	ES**__B	16, 18	6, 8, 10, 12, 14, 15	16	
		20, 22		20	
		24, 28		22	
Fully covered / Removal string on both ends	ES**__F2	16, 18	6, 8, 10, 12, 14, 15	16	
		20, 22		20	
		24, 28		22	

\* Available in USA  
(The sizes marked in blue)

Proximal Release					
Product	Code	Stent		Delivery	
		Diameter (mm)	Length (cm)	Profile (Fr)	Usable length (cm)
Fully covered	ESP**__F	16, 18	6, 8, 10, 12, 14, 15	16	70
		20, 22		20	
		24, 28		22	
Both ends bare	ESP**__B	16, 18	6, 8, 10, 12, 14, 15	16	
		20, 22		20	
		24, 28		22	
Fully covered / Removal string on both ends	ESP**__F2	16, 18	6, 8, 10, 12, 14, 15	16	
		20, 22		20	
		24		22	

## Released Article

\* **Fully covered, retrievable self-expanding metal stents(Niti-s) in palliation of malignant dysphagia: Long-term results of a prospective study**  
Gastroenterology, 2011; 46: 875-880 Sung Jun Choi, Jin Hong Kim, MD



**Aquilant Interventional**

### UK Distributors

Aquilant Interventional  
Aquilant House, Units B1-B2  
Bond Close, Kingsland Business Park  
Basingstoke, Hampshire, RG24 8PZ

T: +44 (0) 845 6024 007  
F: +44 (0) 1256 365486  
E: [contactus@aquilantinterventional.com](mailto:contactus@aquilantinterventional.com)  
W: [www.aquilantinterventional.com](http://www.aquilantinterventional.com)

