





S Enteral colonic stent (Covered)

for benign and malignant colorectal obstruction



Feature

- **Fixed cell with braided construction**
 - High flexibility and optimal radial force
- **Silicone covering and soft round ends**
 - Reduce tissue ingrowth and hyperplasia reaction
- **Visible green suture for easy removal**
- **Radiopaque marker** : 4(four) at both ends & 2(two) in the middle



- ◀ **A.** Complete disunion of the ileoanal anastomosis
B. The 5-week endoscopic control showing complete repair of the anastomotic disunion

[GASTROINTESTINAL ENDOSCOPY 2009;69:1282-7]

Ordering Information

TTS (Through The Scope)

► Fully covered

Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
CST1806F	18	6	10.5	220
CST1808F		8		
CST1810F		10		
CST1812F		12		
CST1814F		14		
CST1815F		15		
CST2006F	20	6		
CST2008F		8		
CST2010F		10		
CST2012F		12		
CST2014F		14		
CST2015F		15		

OTW (Over The Wire)

► Fully covered

Code	Stent		Delivery			
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)		
CS1806F	18	6	16	70		
CS1808F		8				
CS1810F		10				
CS1812F		12				
CS1814F		14				
CS1815F		15				
CS2006F	20	6	20			
CS2008F		8				
CS2010F		10				
CS2012F		12				
CS2014F		14				
CS2015F		15				
CS2206F	22	6			22	
CS2208F		8				
CS2210F		10				
CS2212F		12				
CS2214F		14				
CS2215F		15				
CS2406F	24	6	22			
CS2408F		8				
CS2410F		10				
CS2412F		12				
CS2414F		14				
CS2415F		15				
CS2606F	26	6				22
CS2608F		8				
CS2610F		10				
CS2612F		12				
CS2614F		14				
CS2615F		15				
CS2806F	28	6			22	
CS2808F		8				
CS2810F		10				
CS2812F		12				
CS2814F		14				
CS2815F		15				

Released Article

* *From iatrogenic digestive perforation to complete anastomotic disunion: endoscopic stenting as a new concept of stent-guided regeneration and re-epithelialization* GASTROINTESTINAL ENDOSCOPY Vol. 69, No. 7: 2009 Laila Amrani

Ordering Information

TTS (Through The Scope)

► Both ends bare

Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
CST1806B	18	6	10.5	220
CST1808B		8		
CST1810B		10		
CST1812B		12		
CST1814B		14		
CST1815B		15		
CST2006B	20	6		
CST2008B		8		
CST2010B		10		
CST2012B		12		
CST2014B		14		
CST2015B		15		

OTW (Over The Wire)

► Both ends bare

Code	Stent		Delivery			
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)		
CS1806B	18	6	16	70		
CS1808B		8				
CS1810B		10				
CS1812B		12				
CS1814B		14				
CS1815B		15				
CS2006B	20	6	20			
CS2008B		8				
CS2010B		10				
CS2012B		12				
CS2014B		14				
CS2015B		15				
CS2206B	22	6			22	
CS2208B		8				
CS2210B		10				
CS2212B		12				
CS2214B		14				
CS2215B		15				
CS2406B	24	6	22			
CS2408B		8				
CS2410B		10				
CS2412B		12				
CS2414B		14				
CS2415B		15				
CS2606B	26	6				22
CS2608B		8				
CS2610B		10				
CS2612B		12				
CS2614B		14				
CS2615B		15				
CS2806B	28	6			22	
CS2808B		8				
CS2810B		10				
CS2812B		12				
CS2814B		14				
CS2815B		15				

Ordering Information

TTS (Through The Scope)

► Distal end bare

Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
CST1806H	18	6	10.5	220
CST1808H		8		
CST1810H		10		
CST1812H		12		
CST1814H		14		
CST1815H		15		
CST2006H	20	6		
CST2008H		8		
CST2010H		10		
CST2012H		12		
CST2014H		14		
CST2015H		15		

OTW (Over The Wire)

► Distal end bare

Code	Stent		Delivery				
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)			
CS1806H	18	6	16	70			
CS1808H		8					
CS1810H		10					
CS1812H		12					
CS1814H		14					
CS1815H		15					
CS2006H	20	6	20				
CS2008H		8					
CS2010H		10					
CS2012H		12					
CS2014H		14					
CS2015H		15					
CS2206H	22	6			22		
CS2208H		8					
CS2210H		10					
CS2212H		12					
CS2214H		14					
CS2215H		15					
CS2406H	24	6				22	
CS2408H		8					
CS2410H		10					
CS2412H		12					
CS2414H		14					
CS2415H		15					
CS2606H	26	6					22
CS2608H		8					
CS2610H		10					
CS2612H		12					
CS2614H		14					
CS2615H		15					
CS2806H	28	6	22				
CS2808H		8					
CS2810H		10					
CS2812H		12					
CS2814H		14					
CS2815H		15					



Aquilant Interventional

UK Distributors

Aquilant Interventional
Aquilant House, Units B1-B2
Bond Close, Kingsland Business Park
Basingstoke, Hampshire, RG24 8PZ

T: +44 (0) 845 6024 007
F: +44 (0) 1256 365486
E: contactus@aquilantinterventional.com
W: www.aquilantinterventional.com