


Kendall SCD™ Comfort Sleeve

Sizing guide



Thigh length - Kendall SCD™ Comfort Sleeve



Measure **upper thigh** circumference at the gluteal furrow.

Once measurement has been obtained, use the table to identify the correct sleeve size and product code.

Thigh Circumference 1	Description	Size	Order Code	Case Qty
Up to 40.6 cm	Thigh Length	X-Small	74010	5 pairs
40.6 cm to 55.9 cm	Thigh Length	Small	74011	5 pairs
55.9 cm to 71.1 cm	Thigh Length	Medium	74012	5 pairs
71.1 cm to 91.4 cm	Thigh Length	Large	74013	3 pairs

Thigh length - Kendall SCD™ Tear-Away Sleeve




Measure **upper thigh** circumference at the gluteal furrow.

Once measurement has been obtained, use the table to identify the correct sleeve size and product code.

Thigh Circumference 1	Description	Size	Order Code	Case Qty
40.6 cm to 55.9 cm	Thigh Length	Small	74041	5 pairs
55.9 cm to 71.1 cm	Thigh Length	Medium	74042	5 pairs
71.1 cm to 91.4 cm	Thigh Length	Large	74043	3 pairs

Knee length - Kendall SCD™ Comfort Sleeve



Measure **calf** circumference at the greatest portion.

Once measurement has been obtained, use the table to identify the correct sleeve size and product code.

Calf Circumference 2	Description	Size	Order Code	Case Qty
Up to 35.6 cm	Knee Length	Small	74021	5 pairs
35.6 cm to 53.3 cm	Knee Length	Medium	74022	5 pairs
53.3 cm to 66.0 cm	Knee Length	Large	74023	5 pairs

Kendall SCD™ Sequential Compression Tubing

Description	Order Code	Case Qty
Kendall SCD Tubing Set	9528	1 pair
Kendall SCD Tubing Set (Ambulatory)	9595	1 pair