



Kendall SCD™ Sequential Compression Comfort Sleeves


Proven clinically. Now proven comfortable.



Your partner in VTE prevention.

The Kendall SCD™ comfort sleeves have been designed to provide an increased level of comfort, and have been designed to increase patient compliance with mechanical prophylaxis.¹

Factors that Lead to Patient Discomfort



Sweat and heat

Cool, venting inner fabric helps to wick sweat away from the skin, evaporating through the soft edges. Through this process, cool air is also transferred through the vent holes.



Itchiness

The sleeve material's fibres have diameters similar to silk and cashmere, helping to deliver softness and smoothness.



Pressure

Improved bladder geometry helps to evenly distribute pressure.



Skin irritation

The sleeve's pillowcase design may help address skin irritation by allowing the layers to glide freely among each other as patient shifts in bed.



Reference:

1. Conway D et al. A controlled study evaluating a newly designed sequential compression comfort sleeve in hospitalized patients. Presented at the Academy of Medical-Surgical Nurses (AMSNS), 22 October, 2010, in Las Vegas, USA.

All Kendall™ Sequential Compression Comfort Sleeves are indicated to help prevent DVT and PE in obese and morbidly obese patients.

Ordering information

Kendall SCD™ Sequential Compression Comfort Sleeve - Knee Length

Item Code	Size	Calf Circumference	Quantity/Case
74021	Small	Up to 35.6 cm	5 pairs
74022	Medium	35.6 cm to 53.3 cm	5 pairs
74023	Large	53.3 cm to 66.0 cm	5 pairs

Kendall SCD™ Sequential Compression Comfort Sleeve - Thigh Length

Item Code	Size	Thigh Circumference	Quantity/Case
74010	X-Small	Up to 40.6 cm	5 pairs
74011	Small	40.6 cm to 55.9 cm	5 pairs
74012	Medium	55.9 cm to 71.1 cm	5 pairs
74013	Large	71.1 cm to 91.4 cm	3 pairs

Kendall SCD™ Sequential Compression Tear-Away Sleeve - Thigh Length

Item Code	Size	Thigh Circumference	Quantity/Case
74041	Small	40.6 cm to 55.9 cm	5 pairs
74042	Medium	55.9 cm to 71.1 cm	5 pairs
74043	Large	71.1 cm to 91.4 cm	3 pairs

IMPORTANT: Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

© 2018 Cardinal Health. All rights reserved. CARDINAL HEALTH, the Cardinal Health LOGO, ESSENTIAL TO CARE and KENDALL SCD are trademarks or registered trademarks of Cardinal Health. All other marks are the property of their respective owners. 2PATR18-797953