



**Aquilant**  
A HC21 Group Company

**HC21**  
Healthcare

**SecurAcath**

## Patient Information leaflet

**Your line/ drain is secured with a SecurAcath device.**

About SecurAcath:

- The SecurAcath has small blunt, flexible securement feet that sit beneath the skin holding the catheter/drain in place
- You should not feel any discomfort or pain from the SecurAcath
- The SecurAcath is designed to stay in place for the entire duration of your therapy and does not need to be changed
- Your dressing will be changed and the site cleaned once a week
- You can shower as directed by your clinicians who may provide you with a shower sleeve, you must keep the dressing dry
- The SecurAcath will be removed when the line/drain is re-moved
- **The SecurAcath device is MRI compatible**

If your SecurAcath becomes painful, the site becomes red or you feel generally unwell/feverish please contact your care provider as soon as possible.

For more information about SecurAcath visit [www.securacath.com](http://www.securacath.com)

**Distributed By Aquilant**  
Unit B1-B2 Bond Close,  
Kingsland Business Park, Basingstoke,  
Hampshire, RG24 8PZ, United Kingdom  
T: +44 (0) 1256 306 506  
E: [contactus@aquilant-services.com](mailto:contactus@aquilant-services.com)

Interrad Medical  
181 Cheshire Lane,  
Suite 100, Plymouth,  
MN 55441 USA  
+1.763.225.6699  
[www.securacath.com](http://www.securacath.com)

