




Nurse Information Leaflet

Your patients line/ drain is secured with SecurAcath.

- SecurAcath has small blunt, flexible securement feet that sit beneath the skin holding the catheter/drain in place
- SecurAcath stays in place for the 'Life of the Line'
- Change film dressing and clean site as per local policy
- You do not need to apply any additional securement devices
- The SecurAcath device is MRI conditional


For more information about SecurAcath visit www.securacath.com

SecurAcath® app!



**For the
Life of
the Line™**



Use QR codes to access videos or search for SecurAcath in  **YouTube**



Dressing change
and cleaning

Removal; Fold
Method



Removal; Fold
or Split Method

Minimise pain
on Removal



SecurAcath® app!

