







CLINICIAN-INSPIRED DESIGN. ENHANCED PATIENT CARE.

Non-DEHP Shiley™ neonatal and pediatric tracheostomy tubes, available with TaperGuard™ cuff*



The new Shiley™ neonatal and pediatric tracheostomy tubes feature non-DEHP material and expanded size offerings, as well as several new design enhancements.

| | DESIGN ENHANCEMENT | PREVIOUS VERSION | |
|--|--|--|---|
|  | Flange Softer**, clear material on neck flange for enhanced visualization of the underlying skin and stoma site. Softer material may increase patient comfort. | Opaque, white flange |  |
|  | Cuff Low-profile (low-volume), low-pressure air-filled cuff to provide sealing, TaperGuard™ cuff technology | Traditional barrel-shaped cuff |  |
|  | Material Made with materials that do not contain DEHP | Made with materials containing DEHP |  |
| | Flange Eyelet Raised outer eyelet portions for easier access when inserting the trach tube holder/tape tie through the flange eyelet | Relatively flat profile | |
| | 15 mm Connector Increased length on both cuffed and uncuffed devices | Approximately 2.5 mm shorter than new design | |
| | Sizes Additional smaller sizes as well as a range of cuffed neonatal tubes | Not available in size 2.5 | |

* Shiley™ neonatal and pediatric tracheostomy tubes with TaperGuard™ cuff are not available in all markets. Shiley™ neonatal and pediatric tracheostomy tubes with TaperGuard™ cuff is CE marked and commercially available in European Union countries. Please refer to local product labeling.

** Compared to legacy Shiley™ product

Specifications

On average, there is a 0.2-mm difference in the outer diameter of the new Shiley™ tracheostomy tubes for neonatal and pediatric patients compared with the previous version, as shown in the table below. The new product offers a consistent outer diameter for each product size across all configurations.

Shiley™ neonatal tracheostomy tubes

| NEW | | | PREVIOUS | |
|--|-------------|--------|-------------|--------|
| Description | Product SKU | O.D. | Product SKU | O.D. |
| Neonatal tracheostomy tube, cuffless | 2.5NEF | 4.2 mm | N/A | N/A |
| | 3.0NEF | 4.8 mm | 3.0NEO | 4.5 mm |
| | 3.5NEF | 5.4 mm | 3.5NEO | 5.2 mm |
| | 4.0NEF | 6.0 mm | 4.0NEO | 5.9 mm |
| | 4.5NEF | 6.7 mm | 4.5NEO | 6.5 mm |
| Neonatal tracheostomy tube with TaperGuard™ cuff | 2.5NCF | 4.2 mm | N/A | N/A |
| | 3.0NCF | 4.8 mm | N/A | N/A |
| | 3.5NCF | 5.4 mm | N/A | N/A |
| | 4.0NCF | 6.0 mm | N/A | N/A |
| | 4.5NCF | 6.7 mm | N/A | N/A |

Shiley™ pediatric tracheostomy tubes

| NEW | | | PREVIOUS | |
|---|-------------|--------|-------------|--------|
| Description | Product SKU | O.D. | Product SKU | O.D. |
| Pediatric tracheostomy tube, cuffless | 2.5PEF | 4.2 mm | N/A | N/A |
| | 3.0PEF | 4.8 mm | 3.0PED | 4.5 mm |
| | 3.5PEF | 5.4 mm | 3.5PED | 5.2 mm |
| | 4.0PEF | 6.0 mm | 4.0PED | 5.9 mm |
| | 4.5PEF | 6.7 mm | 4.5PED | 6.5 mm |
| | 5.0PEF | 7.3 mm | 5.0PED | 7.1 mm |
| | 5.5PEF | 7.9 mm | 5.5PED | 7.7 mm |
| Extra long pediatric tracheostomy tube, cuffless | 5.0PELF | 7.3 mm | 5.0PDL | 7.1 mm |
| | 5.5PELF | 7.9 mm | 5.5PDL | 7.7 mm |
| | 6.0PELF | 8.5 mm | 6.0PDL | 8.3 mm |
| | 6.5PELF | 9.0 mm | 6.5PDL | 9.0 mm |
| Pediatric tracheostomy tube with TaperGuard™ cuff* | 2.5PCF | 4.2 mm | N/A | N/A |
| | 3.0PCF | 4.8 mm | N/A | N/A |
| | 3.5PCF | 5.4 mm | N/A | N/A |
| | 4.0PCF | 6.0 mm | 4.0PDC | 5.9 mm |
| | 4.5PCF | 6.7 mm | 4.5PDC | 6.5 mm |
| | 5.0PCF | 7.3 mm | 5.0PDC | 7.1 mm |
| Extra long pediatric tracheostomy tube with TaperGuard™ cuff* | 5.5PCF | 7.9 mm | 5.5PDC | 7.7 mm |
| | 5.0PLCF | 7.3 mm | 5.0PLC | 7.1 mm |
| | 5.5PLCF | 7.9 mm | 5.5PLC | 7.7 mm |
| | 6.0PLCF | 8.5 mm | 6.0PLC | 8.3 mm |
| | 6.5PLCF | 9.0 mm | 6.5PLC | 9.0 mm |

*Shiley™ neonatal and pediatric tracheostomy tubes with TaperGuard™ cuff are not available in all markets. Shiley™ neonatal and pediatric tracheostomy tubes with TaperGuard™ cuff is CE marked and commercially available in European Union countries. Please refer to local product labeling.

Photo credit: iStock

IMPORTANT: Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

© 2018 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. 18-weu-shiley-trach-neo-ped-sales-sheet-2233773

Medtronic