



CE
0120

Easy to handle • Reduce tissue charring • Make ideal ablation zone

star RF ELECTRODE



star RF Electrode _Uni-faced

The star RF electrodes of **Uni-faced** type are intended to create the half moon ablation zone by generating RF energy at square expose range. Uni-faced type is able to avoid damage of normal tissue and heat sink effect from the contiguous vessel.



Half Moon Ablation Zone

Useful to avoid damage to normal tissue and organs



Half moon ablation zone

star RF Electrode _Uni-faced SPECIFICATION

Unit : mm

17-07s05U	17		05	17-20s05U	17		05
17-07s07U	17		07	17-20s07U	17		07
17-07s10U	17	70	10	17-20s10U	17	200	10
17-07s15U	17		15	17-20s15U	17		15
17-10s05U	17		05	17-25s05U	17		05
17-10s07U	17		07	17-25s07U	17		07
17-10s10U	17	100	10	17-25s10U	17	250	10
17-10s15U	17		15	17-25s15U	17		15
17-15s05U	17		05	17-30s05U	17		05
17-15s07U	17		07	17-30s07U	17		07
17-15s10U	17	150	10	17-30s10U	17	300	10
17-15s15U	17		15	17-30s15U	17		15