



Seating

R82[®]



User height
750–1500 mm



Max. load
40 kg



Environment
Indoor & outdoor

Stingray[®]

Stingray is a super stylish buggy with tool free adjustment of postural support and unique swivel action for forward/rearward facing options.

The Stingray is available in 2 sizes.

Highlights about Stingray



Foldable

The Stingray is folded quickly and fits into almost any car. Fold the seat and frame units separately or altogether and pack it into your car - soon you will be on your way.



180° turnable seat unit



Tilt in space and back recline

Perfect seating position due to adjustable seat tilt and back angle adjustment. The seat unit tilts 45° and the back reclines down, allowing your child to rest while out and about.



Safe transportation in vehicles

The Stingray complies with ISO 7176-19 for safe transportation in motor vehicles. The Stingray is approved for transport in vehicles of a user weight of up to 40 kgs.



Standard

Tilt in space, back recline, adjustable seat depth, adjustable back height, height and angle adjustable foot support, turnable seat, frame one size.

Options

Sizes

- Size 1
- Size 2

Front wheels

- 7" solid or air, with swivel and directional stabilizer
- 12" solid or air, no swivel

Rear wheels

- 12" solid or air, no swivel

Foot brake

- Individual
- Central

Hand brake

- Hand brakes
- No hand brakes

Foot rest

- Long bar
- Short bar

Cushion colours



Blue



Red



Green



Pink



Grey



Yellow



Orange












Standard Stingray with 12" no swivel front wheels



Standard Stingray with 7" turnable front wheels, with directional stabilizer

Measurements

		1	2
	Effective seat width	300 mm	350 mm
	Effective seat depth	180-300 mm	250-400 mm
	Backrest height	410-550 mm	500-650 mm
	Seat surface height at front edge**	520 mm	520 mm
	Max. load	40 kg	40 kg
	Max. load at transportation	40 kg	40 kg
	Backrest angle	170°	170°
	Seat plane angle	45°	45°
	Length*	740 mm	740 mm
	Overall width	590 mm	590 mm
	Folded length	900 mm	900 mm
	Folded width	590 mm	590 mm
	Folded height	490 mm	490 mm
	Total mass	19 kg	20 kg

*Note: The configuration of the product. **Measured to adaptorplate. Further measurements at our website

Nice to know

The big range of accessories, flexible adjustment of the foot support and hip and side supports on the Stingray, gives the opportunity of making a comfortable and supportive seating solution for the individual user.

The Stingray has supreme driving properties, fantastic suspension and with the safe and stable brakes, security always has pride of place.



Accessories



Folding top
8740110, 8740120



Head rest cushion
8740810, 8740820



Sculpted head support
872083x



Anatomic head rest
8720862, 8720863




Two-part head support
8720870



Grip handle
872040x



Hip supports
87402x0



Pommel
872050x



Tray
87200x0, 8776350



Transport fittings
8720902



Side supports
872030x



Foam back plate
877632x

Covers for knobs on side supports
8720360

For more accessories - visit our website



Hip belt
85427-x



Oxygen bottle holder
8776295



Bunting bag
874013x



Vent tray
8776290



Shopping basket
8720600



Head rest fitting
8720850



Thermo bag
87460-xxx



Backpack
87043183-x



Rain cover
872014x



Rain cape with hood
87465-xxx

