




Argyle™
Suction Catheters



CardinalHealth
Essential to care™












Argyle™ Suction Catheters

A comprehensive family of suction catheters for a range of procedures with choice of tip designs and optional connectors for different suction sources. Constructed with a frosted surface and to an effective level of hardness to resist kinking.



EASIER IDENTIFICATION
Universal size colour coding for tips and connectors across the range.
 BS EN ISO 8836; 2009

Catheter size – colour code system

CH/FR	Colour		Tube ID	Diameter
05CH/FR	Grey		1.7 mm	○
06CH/FR	Light Green		2.0 mm	○
07CH/FR	Violet		2.3 mm	○
08CH/FR	Blue		2.7 mm	○
10CH/FR	Black		3.3 mm	○
12CH/FR	White		4.0 mm	○
14CH/FR	Green		4.7 mm	○
16CH/FR	Orange		5.3 mm	○
18CH/FR	Red		6.0 mm	○
20CH/FR	Yellow		6.7 mm	○

Suction Catheters - GENTLE FLOW for Atraumatic Suction*

Tracheobronchial suction catheter

- Atraumatic* tube for bronchial suction
- With one central lumen and 4 elliptical side holes to protect the mucosa from invagination and minimize suction trauma
- Uniformly distributed flow
- Available with 2 different kinds of connector: Funnel or vacuum suction control system
- Colour coded per size
- Medical grade PVC, particularly soft and thermosensitive with a SATIN SLIP surface for minimal trauma during insertion and extraction



Tip 2168



Moulded Funnel 2168

GENTLE FLOW with Funnel Connector

CODE	DESCRIPTION	QTY
2058	GENTLE FLOW with Funnel Connector – ● 5Fr x 39 cm – Grey	100
2068	GENTLE FLOW with Funnel Connector – ● 6Fr x 39 cm – Clear Green	100
2088	GENTLE FLOW with Funnel Connector – ● 8Fr x 39 cm – Blue	100
2108	GENTLE FLOW with Funnel Connector – ● 10Fr x 50 cm – Black	100
2128	GENTLE FLOW with Funnel Connector – ○ 12Fr x 50 cm – White	100
2148	GENTLE FLOW with Funnel Connector – ● 14Fr x 50 cm – Dark Green	100
2168	GENTLE FLOW with Funnel Connector – ● 16Fr x 50 cm – Orange	100
2188	GENTLE FLOW with Funnel Connector – ● 18Fr x 61 cm – Red	100
2108L	GENTLE FLOW with Funnel Connector – ● 10Fr x 61 cm – Black	100
2128L	GENTLE FLOW with Funnel Connector – ○ 12Fr x 61 cm – White	100
2148L	GENTLE FLOW with Funnel Connector – ● 14Fr x 61 cm – Dark Green	100
2168L	GENTLE FLOW with Funnel Connector – ● 16Fr x 61 cm – Orange	100
2188L	GENTLE FLOW with Funnel Connector – ● 18Fr x 61 cm – Red	100



Tip 2181



Vacuum Control 2181

GENTLE FLOW with Vacuum Control

CODE	DESCRIPTION	QTY
2051	GENTLE FLOW with Vacuum Control – ● 5Fr x 39 cm – Grey	100
2061	GENTLE FLOW with Vacuum Control – ● 6Fr x 39 cm – Clear Green	100
2081	GENTLE FLOW with Vacuum Control – ● 8Fr x 39 cm – Blue	100
2101	GENTLE FLOW with Vacuum Control – ● 10Fr x 50 cm – Black	100
2121	GENTLE FLOW with Vacuum Control – ○ 12Fr x 50 cm – White	100
2141	GENTLE FLOW with Vacuum Control – ● 14Fr x 50 cm – Dark Green	100
2161	GENTLE FLOW with Vacuum Control – ● 16Fr x 50 cm – Orange	100
2181	GENTLE FLOW with Vacuum Control – ● 18Fr x 50 cm – Red	100
2101L	GENTLE FLOW with Vacuum Control – ● 10Fr x 61 cm – Black	100
2121L	GENTLE FLOW with Vacuum Control – ○ 12Fr x 61 cm – White	100
2141L	GENTLE FLOW with Vacuum Control – ● 14Fr x 61 cm – Dark Green	100
2161L	GENTLE FLOW with Vacuum Control – ● 16Fr x 61 cm – Orange	100
2181L	GENTLE FLOW with Vacuum Control – ● 18Fr x 61 cm – Red	100

* it provides a smooth tip minimizing trauma

Suction Catheters – AERO FLOW for Atraumatic* Suction

Tracheobronchial suction catheter with Atraumatic* tip

- Available with the Vacuum suction control system and funnel connector
- SLIDE-TEX™ external finish
- Colour coded
- Also available with DIRECTROL™ undulated tip to facilitate placement into the bronchial cavities or Coude Tip
- Medical grade PVC
- Tube for bronchial suction with “ring” tip
- With one central hole and four side holes to protect the mucosa from invagination and minimize suction trauma
- Uniformly distributed flow

AERO FLOW with Vacuum Control

CODE	DESCRIPTION	QTY
8888250092	AERO FLOW with Vacuum Control – ● 5Fr x 25 cm – Lavender	50
8888250118	AERO FLOW with Vacuum Control – ● 8Fr x 40 cm – Blue	50
8888250126	AERO FLOW with Vacuum Control – ● 10Fr x 52 cm – Black	50
8888250134	AERO FLOW with Vacuum Control – ○ 12Fr x 52 cm – White	50
8888250142	AERO FLOW with Vacuum Control – ● 14Fr x 52 cm – Green	50
1180250159	AERO FLOW with Vacuum Control – ● 16Fr x 52 cm – Orange	50



Tip 8888250134



Vacuum Control 8888250134

AERO FLOW with Integrated Funnel Connector

CODE	DESCRIPTION	QTY
1180250415	AERO FLOW with Integrated Funnel Connector – ● 8Fr x 48 cm – Blue	50
1180250423	AERO FLOW with Integrated Funnel Connector – ● 10Fr x 48 cm – Black	50
1180250431	AERO FLOW with Integrated Funnel Connector – ○ 12Fr x 48 cm – White	50
1180250449	AERO FLOW with Integrated Funnel Connector – ● 14Fr x 48 cm – Green	50
1180250464	AERO FLOW with Integrated Funnel Connector – ● 16Fr x 48 cm – Orange	50
1180250472	AERO FLOW with Integrated Funnel Connector – ● 16Fr x 63 cm – Orange	50




Tip 1180250464



Funnel End 1180250464

AERO FLOW with Vacuum Control and Coude Tip

CODE	DESCRIPTION	QTY
1180251736	AERO FLOW with Vacuum Control and Coude Tip – ● 10Fr x 51 cm – Black	50
1180251744	AERO FLOW with Vacuum Control and Coude Tip – ○ 12Fr x 51 cm – White	50
1180251751	AERO FLOW with Vacuum Control and Coude Tip – ● 14Fr x 51 cm – Green	50
1180251769	AERO FLOW with Vacuum Control and Coude Tip – ● Black	50



Close-up of Tip




Coude Tip 1180251744




Vacuum Control 1180251744

* it provides a smooth tip minimizing trauma




AERO FLOW with Funnel Connector and Coude Tip

CODE	DESCRIPTION	QTY
1180250811	AERO FLOW with Funnel Connector and Coude Tip – ● 8Fr x 40 cm – Blue	50
1180250829	AERO FLOW with Funnel Connector and Coude Tip – ● 10Fr x 52 cm – Black	50
1180250837	AERO FLOW with Funnel Connector and Coude Tip – ○ 12Fr x 52 cm – White	50
1180250845	AERO FLOW with Funnel Connector and Coude Tip – ● 14Fr x 52 cm – Green	50
1180250860	AERO FLOW with Funnel Connector and Coude Tip – ● 14Fr x 65 cm – Green	50
1180250852	AERO FLOW with Funnel Connector – Coude Tip – ● 16Fr x 52 cm – Orange	50



AERO FLOW with Vacuum Control and Undulated Tip

CODE	DESCRIPTION	QTY
1180250514	AERO FLOW Suction Catheter with Vacuum Control and Undulated Tip – ○ 12Fr x 53 cm – White	50
1180250522	AERO FLOW Suction Catheter with Vacuum Control and Undulated Tip – ● 14Fr x 53 cm – Green	50



AERO FLOW with Funnel Connector and Undulated Tip

CODE	DESCRIPTION	QTY
1180250001	AERO FLOW Suction Catheter with Integrated Funnel and Undulated Tip – ● 10Fr x 53 cm – Black	50
1180250019	AERO FLOW Suction Catheter with Integrated Funnel and Undulated Tip – ○ 12Fr x 53 cm – White	50
1180250027	AERO FLOW Suction Catheter with Integrated Funnel and Undulated Tip – ● 14Fr x 53 cm – Green	50

Argyle™ AERO JET Suction Catheter for Bronchial Lavage

Tracheobronchial suction and Lavage Catheter

- Two-way bronchial suction and lavage tube
- AERO FLOW Atraumatic* "ring" tip*
- With one central hole and four side holes to protect the mucosa from invagination and minimize suction trauma
- With a second lumen for bronchial lavage with a Lumen connector
- Available with the Vacuum suction control system
- Radiopaque
- Medical grade PVC

Dual Lumen Bronchial Suction Catheter with Integrated Funnel Connector

CODE	DESCRIPTION	QTY
1180250936	ARGYLE™ AERO-JET Dual Lumen Bronchial Suction Catheter with Integrated Funnel Connector – ○ 12Fr x 47 cm – White	50
1180250944	ARGYLE™ AERO-JET Dual Lumen Bronchial Suction Catheter with Integrated Funnel Connector – ● 14Fr x 47 cm – Green	50
1180250951	ARGYLE™ AERO-JET Dual Lumen Bronchial Suction Catheter with Integrated Funnel Connector – ● 16Fr x 47 cm – Orange	50



* it provides a smooth tip minimizing trauma

Tracheobronchial Suction Catheters

- Standard suction catheter for bronchial suction
- With 2 side eyes
- Available with 2 different kinds of connector: integrated Funnel, sure-grip funnel connector or vacuum breaker connector
- Colour coded per size



Argyle™ suction catheter with vacuum breaker

CODE	DESCRIPTION	QTY
1180851055	Standard suction catheter with vacuum breaker and 2 side eyes – ● 5Fr x 30 cm – Violet	50
1180851063	Standard suction catheter with vacuum breaker and 2 side eyes – ● 6Fr x 44 cm – light green	100
1180851071	Standard suction catheter with vacuum breaker and 2 side eyes – ● 6Fr x 30 cm – light green	50
1180851089	Standard suction catheter with vacuum breaker and 2 side eyes – ● 8Fr x 45 cm – Blue	100
1180851097	Standard suction catheter with vacuum breaker and 2 side eyes – ● 8Fr x 30 cm – Blue	50
1180851105	Standard suction catheter with vacuum breaker and 2 side eyes – ● 10Fr x 45 cm – Black	100
1180851121	Standard suction catheter with vacuum breaker and 2 side eyes – ○ 12Fr x 50 cm – White	100
1180851147	Standard suction catheter with vacuum breaker and 2 side eyes – ● 14Fr x 50 cm – Green	100
1180851162	Standard suction catheter with vacuum breaker and 2 side eyes – ● 16Fr x 50 cm – Orange	100
1180851188	Standard suction catheter with vacuum breaker and 2 side eyes – ● 18Fr x 50 cm – Red	100



Argyle™ suction catheter with molded funnel connector

CODE	DESCRIPTION	QTY
1180853085	Standard suction catheter with funnel connector and 2 side eyes – ● 8Fr x 50 cm – Blue	100
1180853101	Standard suction catheter with funnel connector and 2 side eyes – ● 10Fr x 50 cm – Black	100
1180853127	Standard suction catheter with funnel connector and 2 side eyes – ○ 12Fr x 50 cm – White	100
1180853143	Standard suction catheter with funnel connector and 2 side eyes – ● 14Fr x 50 cm – Green	100
1180853168	Standard suction catheter with funnel connector and 2 side eyes – ● 16Fr x 50 cm – Orange	100
1180853184	Standard suction catheter with funnel connector and 2 side eyes – ● 18Fr x 50 cm – Red	100



Argyle™ suction catheter with integrated funnel connector

CODE	DESCRIPTION	QTY
1180830067	Standard suction catheter with funnel connector and 1 side eye – ● 6Fr x 45 cm – Green	100
1180830083	Standard suction catheter with funnel connector and 1 side eye – ● 8Fr x 45 cm – Blue	100
1180830109	Standard suction catheter with funnel connector and 1 side eye – ● 10Fr x 50 cm – Black	100
1180830125	Standard suction catheter with funnel connector and 1 side eye – ○ 12Fr x 50 cm – White	100
1180830141	Standard suction catheter with funnel connector and 1 side eye – ● 14Fr x 50 cm – Green	100

Suction Catheters

For oral, tracheal and endobronchial suction procedures. To be used by connection to the suction source with Argyle™ suction connecting tubing

- Appropriate degree of hardness provides kink resistance
- Metric type has accurate positioning depth marks at each cm starting 4 cm from the distal end
- Soft rounded, open tip or bevelled tip or bevelled open tip
- Two small relieving eyes at tip / two lateral eyes with smooth edges
- Vacuum control connector – Transparent for easy observation of secretions
- Marker on funnel connector indicating direction of angled tip
- Angled tip for endobronchial suction procedures, providing easier access to the right or left bronchus
- Colour coded connector
- Effective length: 50 cm

Suction Catheter with Vacuum Control Rounded tip – 50 cm

CODE	DESCRIPTION	LENGTH	QTY
019 5201 1	Suction Catheter, Vacuum Control, Rounded Tip – ● 5Fr – Grey	50 cm	100
019 5202 1	Suction Catheter, Vacuum Control, Rounded Tip – ● 6Fr – Light Green	50 cm	100
019 5203 1	Suction Catheter, Vacuum Control, Rounded Tip – ● 7Fr – Violet	50 cm	100
019 5204 1	Suction Catheter, Vacuum Control, Rounded Tip – ● 8Fr – Blue	50 cm	100
019 5205 1	Suction Catheter, Vacuum Control, Rounded Tip – ● 10Fr – Black	50 cm	100
019 5206 1	Suction Catheter, Vacuum Control, Rounded Tip – ○ 12Fr – White	50 cm	100
019 5207 1	Suction Catheter, Vacuum Control, Rounded Tip – ● 14Fr – Green	50 cm	100
019 5208 1	Suction Catheter, Vacuum Control, Rounded Tip – ● 16Fr – Orange	50 cm	100
019 5209 1	Suction Catheter, Vacuum Control, Rounded Tip – ● 18Fr – Red	50 cm	100



Rounded Tip



Vacuum control

Suction Catheter with Molded Funnel Connector and Bevelled tip – 50 cm

CODE	DESCRIPTION	LENGTH	QTY
019 5101 1	Suction Catheter, Moulded Funnel, Bevelled Tip – ● 6Fr – Light Green	50 cm	100
019 5102 1	Suction Catheter, Moulded Funnel, Bevelled Tip – ● 8Fr – Blue	50 cm	100
019 5103 1	Suction Catheter, Moulded Funnel, Bevelled Tip – ● 10Fr – Black	50 cm	100
019 5104 1	Suction Catheter, Moulded Funnel, Bevelled Tip – ○ 12Fr – White	50 cm	100
019 5105 1	Suction Catheter, Moulded Funnel, Bevelled Tip – ● 14Fr – Green	50 cm	100
019 5106 1	Suction Catheter, Moulded Funnel, Bevelled Tip – ● 16Fr – Orange	50 cm	100
019 5107 1	Suction Catheter, Moulded Funnel, Bevelled Tip – ● 18Fr – Red	50 cm	100
019 5108 1	Suction Catheter, Moulded Funnel, Bevelled Tip – ● 20Fr – Yellow	50 cm	75



Bevelled Tip



Molded Funnel Connector



Rounded Tip



Molded Funnel Connector



Coude tip



Molded Funnel Connector



Bevelled tip



Vacuum Control

Suction Catheter with Molded Funnel Connector and Rounded tip – 50 cm

CODE	DESCRIPTION	LENGTH	QTY
019 5141 1	Suction Catheter, Moulded Funnel, Rounded Tip – ● 6Fr – Light Green	50 cm	100
019 5142 1	Suction Catheter, Moulded Funnel, Rounded Tip – ● 8Fr – Blue	50 cm	100
019 5143 1	Suction Catheter, Moulded Funnel, Rounded Tip – ● 10Fr – Black	50 cm	100
019 5144 1	Suction Catheter, Moulded Funnel, Rounded Tip – ○ 12Fr – White	50 cm	100
019 5145 1	Suction Catheter, Moulded Funnel, Rounded Tip – ● 14Fr – Green	50 cm	100
019 5146 1	Suction Catheter, Moulded Funnel, Rounded Tip – ● 16Fr – Orange	50 cm	100
019 5147 1	Suction Catheter, Moulded Funnel, Rounded Tip – ● 18Fr – Red	50 cm	100
019 5148 1	Suction Catheter, Moulded Funnel, Rounded Tip – ● 20Fr – Yellow	50 cm	75

Suction Catheter with Molded Funnel Connector and Coude tip – 50 cm

CODE	DESCRIPTION	LENGTH	QTY
019 5121 1	Suction Catheter, Moulded Funnel, Coude Tip – ● 10Fr – Black	50 cm	100
019 5122 1	Suction Catheter, Moulded Funnel, Coude Tip – ○ 12Fr – White	50 cm	100
019 5123 1	Suction Catheter, Moulded Funnel, Coude Tip – ● 14Fr – Green	50 cm	100
019 5124 1	Suction Catheter, Moulded Funnel, Coude Tip – ● 16Fr – Orange	50 cm	100

Suction Catheter with Vacuum Control and Bevelled Tip – 50 cm

CODE	DESCRIPTION	LENGTH	QTY
019 5301 1	Suction Catheter, Vacuum Control, Bevelled Tip – ● 6Fr – Light Green	50 cm	100
019 5302 1	Suction Catheter, Vacuum Control, Bevelled Tip – ● 8Fr – Blue	50 cm	100
019 5303 1	Suction Catheter, Vacuum Control, Bevelled Tip – ● 10Fr – Black	50 cm	100
019 5304 1	Suction Catheter, Vacuum Control, Bevelled Tip – ○ 12Fr – White	50 cm	100
019 5305 1	Suction Catheter, Vacuum Control, Bevelled Tip – ● 14Fr – Green	50 cm	100
019 5306 1	Suction Catheter, Vacuum Control, Bevelled Tip – ● 16Fr – Orange	50 cm	100
019 5307 1	Suction Catheter, Vacuum Control, Bevelled Tip – ● 18Fr – Red	50 cm	100

Suction Catheters – Graduation Marks

Available with molded funnel connector or vacuum control rounded tip. Length 50 cm or 40 cm

Suction Catheter with Vacuum Control, Graduation Mark & Rounded tip

CODE	DESCRIPTION	LENGTH	QTY
019 5310 1	Suction Catheter with Vacuum Control, Graduation Mark & Rounded tip – 5Fr	50 cm	100
019 5311 1	Suction Catheter with Vacuum Control, Graduation Mark & Rounded tip – 6Fr	50 cm	100
019 5312 1	Suction Catheter with Vacuum Control, Graduation Mark & Rounded tip – 7Fr	50 cm	100
019 5313 1	Suction Catheter with Vacuum Control, Graduation Mark & Rounded tip – 8Fr	50 cm	100
019 5314 1	Suction Catheter with Vacuum Control, Graduation Mark & Rounded – 10Fr	50 cm	100
019 5315 1	Suction Catheter with Vacuum Control, Graduation Mark & Rounded – 12Fr	50 cm	100
019 5316 1	Suction Catheter with Vacuum Control, Graduation Mark & Rounded – 14Fr	50 cm	100
019 5317 1	Suction Catheter with Vacuum Control, Graduation Mark & Rounded – 16Fr	50 cm	100
019 5318 1	Suction Catheter with Vacuum Control, Graduation Mark & Rounded – 18Fr	50 cm	100



Rounded Tip



Graduation marks



Vacuum control

Suction Catheter with Molded Funnel, Graduation Mark & Rounded tip – 50 cm

CODE	DESCRIPTION	LENGTH	QTY
019 5131 1	Suction Catheter with Molded Funnel, Graduation Mark & Rounded tip – 6Fr	50 cm	100
019 5132 1	Suction Catheter with Molded Funnel, Graduation Mark & Rounded tip – 8Fr	50 cm	100
019 5133 1	Suction Catheter with Molded Funnel, Graduation Mark & Rounded tip – 10Fr	50 cm	100
019 5134 1	Suction Catheter with Molded Funnel, Graduation Mark & Rounded tip – 12Fr	50 cm	100
019 5135 1	Suction Catheter with Molded Funnel, Graduation Mark & Rounded tip – 14Fr	50 cm	100
019 5136 1	Suction Catheter with Molded Funnel, Graduation Mark & Rounded tip – 16Fr	50 cm	100
019 5137 1	Suction Catheter with Molded Funnel, Graduation Mark & Rounded tip – 18Fr	50 cm	100
019 5138 1	Suction Catheter with Molded Funnel, Graduation Mark & Rounded tip – 20Fr	50 cm	75



Rounded Tip



Graduation marks



Molded Funnel Connector



Rounded Tip



Graduation marks



Molded Funnel Connector

Suction Catheter with Molded Funnel, Graduation Mark & Rounded tip – 40 cm

CODE	DESCRIPTION	LENGTH	QTY
019 5161 1	Suction Catheter with Molded Funnel, Graduation Mark & Rounded tip – 6Fr	40 cm	100
019 5162 1	Suction Catheter with Molded Funnel, Graduation Mark & Rounded tip – 8Fr	40 cm	100
019 5163 1	Suction Catheter with Molded Funnel, Graduation Mark & Rounded tip – 10Fr	40 cm	100
019 5164 1	Suction Catheter with Molded Funnel, Graduation Mark & Rounded tip – 12Fr	40 cm	100
019 5165 1	Suction Catheter with Molded Funnel, Graduation Mark & Rounded tip – 14Fr	40 cm	100
019 5166 1	Suction Catheter with Molded Funnel, Graduation Mark & Rounded tip – 16Fr	40 cm	100
019 5167 1	Suction Catheter with Molded Funnel, Graduation Mark & Rounded tip – 18Fr	40 cm	100
019 5168 1	Suction Catheter with Molded Funnel, Graduation Mark & Rounded tip – 20Fr	40 cm	75

Argyle™ Specimen Trap

Argyle™ specimen traps feature one-cc graduations for accurate measurement of samples.

- Screw-on cap
- Removable male connector
- Full 40 ml or 70 ml capacity with 1 ml graduations
- Patient label
- Sterile
- Latex-free

Argyle™ Specimen Trap

CODE	DESCRIPTION	QTY
8884724504	Argyle™ Specimen Trap to be used with Bronchoscope – 70 ml	50
8884724500	Argyle™ Specimen Trap DeLee Neo Safe Bubble – 40 ml	50



8884724000

8884724504

Argyle™ SAF-T-MARK Paediatric Suction Catheter

Paediatric Tracheobronchial suction catheter

- Vacuum Control version has one central hole, two side eyes and 10 mm metric markings
- In cuffed suction control version, it has one central hole and one side hole, radiopaque markers and a latex cuff
- Medical grade PVC

Argyle™ SAF-T-MARK Paediatric Suction Catheter

CODE	DESCRIPTION	QTY
1180251801	Paediatric Suction Catheter with Vacuum Control Connector and 1 cm Depth Marks – ● 5Fr x 24 cm – Lavender	50
1180251819	Paediatric Suction Catheter with Vacuum Control Connector and 1 cm Depth Marks – ● 6Fr x 24 cm – Green	50
1180251827	Paediatric Suction Catheter with Vacuum Control Connector and 1 cm Depth Marks – ● 8Fr x 29 cm – Blue	50
8888256016	Paediatric Suction Catheter with Cuffed Suction Control – ● 5Fr x 32 cm – Mauve	50
8888256024	Paediatric Suction Catheter with Cuffed Suction Control – ● 8Fr x 32 cm – Green	50



Bevelled tip



Vacuum Control



Bevelled tip 8888256016



Vacuum Control 8888256016

Argyle™ Lukens Mucus Collection

The Argyle™ Lukens specimen container comes with a whistle tip suction catheter with a tapered connector and a smooth Slide-Tex finish that makes it easy to insert and withdraw.

- 14 Fr whistle tip suction catheter with Slide-Tex finish
- 20 ml specimen trap graduated at 2 ml intervals
- Extra trap cap
- Suction tube with integral tapered connector
- Patient label
- Sterile
- Latex-free

Argyle™ Lukens Collection

CODE	DESCRIPTION	QTY
8888258608	Argyle™ Lukens Mucus Collection – 14Fr Catheter – 20 ml	50



Tip 8888258608



Tapered Connector 8888258608

QUICK REFERENCE GUIDE

ARGYLE™ SUCTION CATHETERS

Suction Catheters - GENTLE FLOW for Atraumatic* Suction

GENTLE FLOW with Funnel Connector

CODE	DESCRIPTION	QTY
2058	GENTLE FLOW with Funnel Connector – 5Fr x 39 cm – Grey	100
2068	GENTLE FLOW with Funnel Connector – 6Fr x 39 cm – Clear Green	100
2088	GENTLE FLOW with Funnel Connector – 8Fr x 39 cm – Blue	100
2108	GENTLE FLOW with Funnel Connector – 10Fr x 50 cm – Black	100
2128	GENTLE FLOW with Funnel Connector – 12Fr x 50 cm – White	100
2148	GENTLE FLOW with Funnel Connector – 14Fr x 50 cm – Dark Green	100
2168	GENTLE FLOW with Funnel Connector – 16Fr x 50 cm – Orange	100
2188	GENTLE FLOW with Funnel Connector – 18Fr x 61 cm – Red	100
2108L	GENTLE FLOW with Funnel Connector – 10Fr x 61 cm – Black	100
2128L	GENTLE FLOW with Funnel Connector – 12Fr x 61 cm – White	100
2148L	GENTLE FLOW with Funnel Connector – 14Fr x 61 cm – Dark Green	100
2168L	GENTLE FLOW with Funnel Connector – 16Fr x 61 cm – Orange	100
2188L	GENTLE FLOW with Funnel Connector – 18Fr x 61 cm – Red	100

GENTLE FLOW with Vacuum Control

CODE	DESCRIPTION	QTY
2051	GENTLE FLOW with Vacuum Control – 5Fr x 39 cm – Grey	100
2061	GENTLE FLOW with Vacuum Control – 6Fr x 39 cm – Clear Green	100
2081	GENTLE FLOW with Vacuum Control – 8Fr x 39 cm – Blue	100
2101	GENTLE FLOW with Vacuum Control – 10Fr x 50 cm – Black	100
2121	GENTLE FLOW with Vacuum Control – 12Fr x 50 cm – White	100
2141	GENTLE FLOW with Vacuum Control – 14Fr x 50 cm – Dark Green	100
2161	GENTLE FLOW with Vacuum Control – 16Fr x 50 cm – Orange	100
2181	GENTLE FLOW with Vacuum Control – 18Fr x 50 cm – Red	100
2101L	GENTLE FLOW with Vacuum Control – 10Fr x 61 cm – Black	100
2121L	GENTLE FLOW with Vacuum Control – 12Fr x 61 cm – White	100
2141L	GENTLE FLOW with Vacuum Control – 14Fr x 61 cm – Dark Green	100
2161L	GENTLE FLOW with Vacuum Control – 16Fr x 61 cm – Orange	100
2181L	GENTLE FLOW with Vacuum Control – 18Fr x 61 cm – Red	100

Suction Catheters – AERO FLOW for Atraumatic Suction*

AERO FLOW with Vacuum Control

CODE	DESCRIPTION	QTY
8888250092	AERO FLOW with Vacuum Control – 5Fr x 25 cm – Lavender	50
8888250118	AERO FLOW with Vacuum Control – 8Fr x 40 cm – Blue	50
8888250126	AERO FLOW with Vacuum Control – 10Fr x 52 cm – Black	50
8888250134	AERO FLOW with Vacuum Control – 12Fr x 52 cm – White	50
8888250142	AERO FLOW with Vacuum Control – 14Fr x 52 cm – Green	50
1180250159	AERO FLOW with Vacuum Control – 16Fr x 52 cm – Orange	50

AERO FLOW with Integrated Funnel Connector

CODE	DESCRIPTION	QTY
1180250415	AERO FLOW with Integrated Funnel Connector – 8Fr x 48 cm – Blue	50
1180250423	AERO FLOW with Integrated Funnel Connector – 10Fr x 48 cm – Black	50
1180250431	AERO FLOW with Integrated Funnel Connector – 12Fr x 48 cm – White	50
1180250449	AERO FLOW with Integrated Funnel Connector – 14Fr x 48 cm – Green	50
1180250464	AERO FLOW with Integrated Funnel Connector – 16Fr x 48 cm – Orange	50
1180250472	AERO FLOW with Integrated Funnel Connector – 16Fr x 63 cm – Orange	50

AERO FLOW with Vacuum Control and Coude Tip

CODE	DESCRIPTION	QTY
1180251736	AERO FLOW with Vacuum Control and Coude Tip – 10Fr x 51 cm – Black	50
1180251744	AERO FLOW with Vacuum Control and Coude Tip – 12Fr x 51 cm – White	50
1180251751	AERO FLOW with Vacuum Control and Coude Tip – 14Fr x 51 cm – Green	50
1180251769	AERO FLOW with Vacuum Control and Coude Tip – Black	50

AERO FLOW with Funnel Connector and Coude Tip

CODE	DESCRIPTION	QTY
1180250811	AERO FLOW with Funnel Connector and Coude Tip – 8Fr x 40 cm – Blue	50
1180250829	AERO FLOW with Funnel Connector and Coude Tip – 10Fr x 52 cm – Black	50
1180250837	AERO FLOW with Funnel Connector and Coude Tip – 12Fr x 52 cm – White	50
1180250845	AERO FLOW with Funnel Connector and Coude Tip – 14Fr x 52 cm – Green	50
1180250860	AERO FLOW with Funnel Connector and Coude Tip – 14Fr x 65 cm – Green	50
1180250852	AERO FLOW with Funnel Connector – Coude Tip – 16Fr x 52 cm – Orange	50

AERO FLOW with Vacuum Control and Undulated Tip

CODE	DESCRIPTION	QTY
1180250514	Atraumatic Suction Catheter* with Vacuum Control and Undulated Tip – 12Fr x 53 cm – White	50
1180250522	Atraumatic Suction Catheter* with Vacuum Control and Undulated Tip – 14Fr x 53 cm – Green	50

AERO FLOW with Funnel Connector and Undulated Tip

CODE	DESCRIPTION	QTY
1180250001	Atraumatic Suction Catheter* with Integrated Funnel and Undulated Tip – 10Fr x 53 cm – Black	50
1180250019	Atraumatic Suction Catheter* with Integrated Funnel and Undulated Tip – 12Fr x 53 cm – White	50
1180250027	Atraumatic Suction Catheter* with Integrated Funnel and Undulated Tip – 14Fr x 53 cm – Green	50

QUICK REFERENCE GUIDE

Argyle™ AERO JET Suction Catheter for Bronchial Lavage

Dual Lumen Bronchial Suction Catheter with Integrated Funnel Connector

CODE	DESCRIPTION	QTY
1180250936	ARGYLE™ AERO-JET Dual Lumen Bronchial Suction Catheter with Integrated Funnel Connector – 12Fr x 47 cm – White	50
1180250944	ARGYLE™ AERO-JET Dual Lumen Bronchial Suction Catheter with Integrated Funnel Connector – 14Fr x 47 cm – Green	50
1180250951	ARGYLE™ AERO-JET Dual Lumen Bronchial Suction Catheter with Integrated Funnel Connector – 16Fr x 47 cm – Orange	50

Tracheobronchial Suction Catheter

Argyle™ suction catheter with vacuum breaker

CODE	DESCRIPTION	QTY
1180851055	Standard suction catheter with vacuum breaker and 2 side eyes – 5Fr x 30 cm – Violet	50
1180851063	Standard suction catheter with vacuum breaker and 2 side eyes – 6Fr x 44 cm – light green	100
1180851071	Standard suction catheter with vacuum breaker and 2 side eyes – 6Fr x 30 cm – light green	50
1180851089	Standard suction catheter with vacuum breaker and 2 side eyes – 8Fr x 45 cm – Blue	100
1180851097	Standard suction catheter with vacuum breaker and 2 side eyes – 8Fr x 30 cm – Blue	50
1180851105	Standard suction catheter with vacuum breaker and 2 side eyes – 10Fr x 45 cm – Black	100
1180851121	Standard suction catheter with vacuum breaker and 2 side eyes – 12Fr x 50 cm – White	100
1180851147	Standard suction catheter with vacuum breaker and 2 side eyes – 14Fr x 50 cm – Green	100
1180851162	Standard suction catheter with vacuum breaker and 2 side eyes – 16Fr x 50 cm – Orange	100
1180851188	Standard suction catheter with vacuum breaker and 2 side eyes – 18Fr x 50 cm – Red	100

Argyle™ suction catheter with molded funnel connector

CODE	DESCRIPTION	QTY
1180853085	Standard suction catheter with funnel connector and 2 side eyes – 8Fr x 50 cm – Blue	100
1180853101	Standard suction catheter with funnel connector and 2 side eyes – 10Fr x 50 cm – Black	100
1180853127	Standard suction catheter with funnel connector and 2 side eyes – 12Fr x 50 cm – White	100
1180853143	Standard suction catheter with funnel connector and 2 side eyes – 14Fr x 50 cm – Green	100
1180853168	Standard suction catheter with funnel connector and 2 side eyes – 16Fr x 50 cm – Orange	100
1180853184	Standard suction catheter with funnel connector and 2 side eyes – 18Fr x 50 cm – Red	100

Argyle™ suction catheter with integrated funnel connector

CODE	DESCRIPTION	QTY
1180830067	Standard suction catheter with funnel connector and 1 side eye – 6Fr x 45 cm – Green	100
1180830083	Standard suction catheter with funnel connector and 1 side eye – 8Fr x 45 cm – Blue	100
1180830109	Standard suction catheter with funnel connector and 1 side eye – 10Fr x 50 cm – Black	100
1180830125	Standard suction catheter with funnel connector and 1 side eye – 12Fr x 50 cm – White	100
1180830141	Standard suction catheter with funnel connector and 1 side eye – 14Fr x 50 cm – Green	100

Suction Catheters

Suction Catheter with Vacuum Control Rounded tip – 50 cm

CODE	DESCRIPTION	LENGTH	QTY
019 5201 1	Suction Catheter, Vacuum Control, Rounded Tip – 5Fr – Grey	50 cm	100
019 5202 1	Suction Catheter, Vacuum Control, Rounded Tip – 6Fr – Light Green	50 cm	100
019 5203 1	Suction Catheter, Vacuum Control, Rounded Tip – 7Fr – Violet	50 cm	100
019 5204 1	Suction Catheter, Vacuum Control, Rounded Tip – 8Fr – Blue	50 cm	100
019 5205 1	Suction Catheter, Vacuum Control, Rounded Tip – 10Fr – Black	50 cm	100
019 5206 1	Suction Catheter, Vacuum Control, Rounded Tip – 12Fr – White	50 cm	100
019 5207 1	Suction Catheter, Vacuum Control, Rounded Tip – 14Fr – Green	50 cm	100
019 5208 1	Suction Catheter, Vacuum Control, Rounded Tip – 16Fr – Orange	50 cm	100
019 5209 1	Suction Catheter, Vacuum Control, Rounded Tip – 18Fr – Red	50 cm	100

Suction Catheter with Molded Funnel Connector and Bevelled tip – 50 cm

CODE	DESCRIPTION	LENGTH	QTY
019 5101 1	Suction Catheter, Moulded Funnel, Bevelled Tip – 6Fr – Light Green	50 cm	100
019 5102 1	Suction Catheter, Moulded Funnel, Bevelled Tip – 8Fr – Blue	50 cm	100
019 5103 1	Suction Catheter, Moulded Funnel, Bevelled Tip – 10Fr – Black	50 cm	100
019 5104 1	Suction Catheter, Moulded Funnel, Bevelled Tip – 12Fr – White	50 cm	100
019 5105 1	Suction Catheter, Moulded Funnel, Bevelled Tip – 14Fr – Green	50 cm	100
019 5106 1	Suction Catheter, Moulded Funnel, Bevelled Tip – 16Fr – Orange	50 cm	100
019 5107 1	Suction Catheter, Moulded Funnel, Bevelled Tip – 18Fr – Red	50 cm	100
019 5108 1	Suction Catheter, Moulded Funnel, Bevelled Tip – 20Fr – Yellow	50 cm	75

Suction Catheter with Molded Funnel Connector and Rounded tip – 50 cm

CODE	DESCRIPTION	LENGTH	QTY
019 5141 1	Suction Catheter, Moulded Funnel, Rounded Tip – 6Fr – Light Green	50 cm	100
019 5142 1	Suction Catheter, Moulded Funnel, Rounded Tip – 8Fr – Blue	50 cm	100
019 5143 1	Suction Catheter, Moulded Funnel, Rounded Tip – 10Fr – Black	50 cm	100
019 5144 1	Suction Catheter, Moulded Funnel, Rounded Tip – 12Fr – White	50 cm	100
019 5145 1	Suction Catheter, Moulded Funnel, Rounded Tip – 14Fr – Green	50 cm	100
019 5146 1	Suction Catheter, Moulded Funnel, Rounded Tip – 16Fr – Orange	50 cm	100
019 5147 1	Suction Catheter, Moulded Funnel, Rounded Tip – 18Fr – Red	50 cm	100
019 5148 1	Suction Catheter, Moulded Funnel, Rounded Tip – 20Fr – Yellow	50 cm	75

Suction Catheter with Molded Funnel Connector and Coude tip – 50 cm

CODE	DESCRIPTION	LENGTH	QTY
019 5121 1	Suction Catheter, Moulded Funnel, Coude Tip – 10Fr – Black	50 cm	100
019 5122 1	Suction Catheter, Moulded Funnel, Coude Tip – 12Fr – White	50 cm	100
019 5123 1	Suction Catheter, Moulded Funnel, Coude Tip – 14Fr – Green	50 cm	100
019 5124 1	Suction Catheter, Moulded Funnel, Coude Tip – 16Fr – Orange	50 cm	100

Suction Catheter with Vacuum Control and Bevelled Tip – 50 cm

CODE	DESCRIPTION	LENGTH	QTY
019 5301 1	Suction Catheter, Vacuum Control, Bevelled Tip – 6Fr – Light Green	50 cm	100
019 5302 1	Suction Catheter, Vacuum Control, Bevelled Tip – 8Fr – Blue	50 cm	100
019 5303 1	Suction Catheter, Vacuum Control, Bevelled Tip – 10Fr – Black	50 cm	100
019 5304 1	Suction Catheter, Vacuum Control, Bevelled Tip – 12Fr – White	50 cm	100
019 5305 1	Suction Catheter, Vacuum Control, Bevelled Tip – 14Fr – Green	50 cm	100
019 5306 1	Suction Catheter, Vacuum Control, Bevelled Tip – 16Fr – Orange	50 cm	100
019 5307 1	Suction Catheter, Vacuum Control, Bevelled Tip – 18Fr – Red	50 cm	100

QUICK REFERENCE GUIDE

Suction Catheters – Graduation Marks

Suction Catheter with Vacuum Control, Graduation Mark & Rounded tip

CODE	DESCRIPTION	LENGTH	QTY
019 5310 1	Suction Catheter with Vacuum Control, Graduation Mark & Rounded tip – 5Fr	50 cm	100
019 5311 1	Suction Catheter with Vacuum Control, Graduation Mark & Rounded tip – 6Fr	50 cm	100
019 5312 1	Suction Catheter with Vacuum Control, Graduation Mark & Rounded tip – 7Fr	50 cm	100
019 5313 1	Suction Catheter with Vacuum Control, Graduation Mark & Rounded tip – 8Fr	50 cm	100
019 5314 1	Suction Catheter with Vacuum Control, Graduation Mark & Rounded – 10Fr	50 cm	100
019 5315 1	Suction Catheter with Vacuum Control, Graduation Mark & Rounded – 12Fr	50 cm	100
019 5316 1	Suction Catheter with Vacuum Control, Graduation Mark & Rounded – 14Fr	50 cm	100
019 5317 1	Suction Catheter with Vacuum Control, Graduation Mark & Rounded – 16Fr	50 cm	100
019 5318 1	Suction Catheter with Vacuum Control, Graduation Mark & Rounded – 18Fr	50 cm	100

Suction Catheter with Molded Funnel, Graduation Mark & Rounded tip – 50 cm

CODE	DESCRIPTION	LENGTH	QTY
019 5131 1	Suction Catheter with Molded Funnel, Graduation Mark & Rounded tip – 6Fr	50 cm	100
019 5132 1	Suction Catheter with Molded Funnel, Graduation Mark & Rounded tip – 8Fr	50 cm	100
019 5133 1	Suction Catheter with Molded Funnel, Graduation Mark & Rounded tip – 10Fr	50 cm	100
019 5134 1	Suction Catheter with Molded Funnel, Graduation Mark & Rounded tip – 12Fr	50 cm	100
019 5135 1	Suction Catheter with Molded Funnel, Graduation Mark & Rounded tip – 14Fr	50 cm	100
019 5136 1	Suction Catheter with Molded Funnel, Graduation Mark & Rounded tip – 16Fr	50 cm	100
019 5137 1	Suction Catheter with Molded Funnel, Graduation Mark & Rounded tip – 18Fr	50 cm	100
019 5138 1	Suction Catheter with Molded Funnel, Graduation Mark & Rounded tip – 20Fr	50 cm	75

Suction Catheter with Molded Funnel, Graduation Mark & Rounded tip – 40 cm

CODE	DESCRIPTION	LENGTH	QTY
019 5161 1	Suction Catheter with Molded Funnel, Graduation Mark & Rounded tip – 6Fr	40 cm	100
019 5162 1	Suction Catheter with Molded Funnel, Graduation Mark & Rounded tip – 8Fr	40 cm	100
019 5163 1	Suction Catheter with Molded Funnel, Graduation Mark & Rounded tip – 10Fr	40 cm	100
019 5164 1	Suction Catheter with Molded Funnel, Graduation Mark & Rounded tip – 12Fr	40 cm	100
019 5165 1	Suction Catheter with Molded Funnel, Graduation Mark & Rounded tip – 14Fr	40 cm	100
019 5166 1	Suction Catheter with Molded Funnel, Graduation Mark & Rounded tip – 16Fr	40 cm	100
019 5167 1	Suction Catheter with Molded Funnel, Graduation Mark & Rounded tip – 18Fr	40 cm	100
019 5168 1	Suction Catheter with Molded Funnel, Graduation Mark & Rounded tip – 20Fr	40 cm	75

Argyle™ Specimen Trap

CODE	DESCRIPTION	QTY
8884724000	Argyle™ Specimen Trap with Screw Top Lid without Catheter – 40 ml	25
8884724504	Argyle™ Specimen Trap to be used with Bronchoscope – 70 ml	50
8884724500	Argyle™ Specimen Trap DeLee Neo Safe Bubble – 40 ml	50

Argyle™ SAF-T-MARK Paediatric Suction Catheter

CODE	DESCRIPTION	QTY
1180251801	Paediatric Suction Catheter with Vacuum Control Connector and 1 cm Depth Marks – ● 5Fr x 24 cm – Lavender	50
1180251819	Paediatric Suction Catheter with Vacuum Control Connector and 1 cm Depth Marks – ● 6Fr x 24 cm – Green	50
1180251827	Paediatric Suction Catheter with Vacuum Control Connector and 1 cm Depth Marks – ● 8Fr x 29 cm – Blue	50
8888256016	Paediatric Suction Catheter with Cuffed Suction Control – ● 5Fr x 32 cm – Mauve	50
8888256024	Paediatric Suction Catheter with Cuffed Suction Control – ● 8Fr x 32 cm – Green	50

Argyle™ Lukens Mucus Collection

CODE	DESCRIPTION	QTY
8888258608	Argyle™ Lukens Mucus Collection – 14Fr Catheter – 20 ml	50



CardinalHealth
Essential to care™

For healthcare professionals only.

© 2020 Cardinal Health. All Rights Reserved. CARDINAL HEALTH, the Cardinal Health LOGO, ESSENTIAL TO CARE and ARGYLE are trademarks of Cardinal Health and may be registered in the US and/or in other countries. All other marks are the property of their respective owners. Important information: Prior to use, refer to the instructions for use supplied with this device for indications, contraindications, side effects, suggested procedure, warnings and precautions. For healthcare professionals only. As part of its continuous product development policy, Cardinal Health reserves the right to change product specifications without prior notification. Please contact your Cardinal Health representative for additional product availability information. 2SRV19-1071559 01/2020