



Argyle™  
Surgical Suction  
Connecting Tubes

# Argyle™ Surgical Suction Connecting Tubes

A comprehensive range of surgical suction connecting tubes available with molded and funnel connectors.



## Argyle™ Bubble Tubing

Can be easily cut to exact length and is available in conductive and non-conductive variations.

- PVC formulation
- Low memory prevents twisting and kinking
- Various lumen sizes
- Bubble at every 36" (91 cm) can be cut to form a tapered or funnel connector
- Dispensed in a single roll for ease of use

### Argyle™ Bubble Tubing

CODE	DESCRIPTION	QTY
8888280214	Argyle™ Bubble Tubing 5 mm x 30.5 m	1
8888290213	Argyle™ Bubble Tubing 5 mm x 30.5 m – Conductive	1
8888280602	Argyle™ Bubble tubing 7 mm x 30 m – Non Conductive	1
8888280412	Argyle™ Bubble Tubing 6 mm x 30.5 m	1
1180839035	Argyle™ Bubble Tubing 6 mm x 50 m	1
8888280610	Argyle™ Bubble Tubing 7 mm x 30.5 m	1
8888290619	Argyle™ Bubble Tubing 7 mm x 30.5 m – Conductive	1
8888280206	Argyle™ Bubble tubing 5 mm x 30 m – Non Conductive	1



8888290213



8888230201

## Argyle™ Oxygen Bubble Tubing

- PVC formulation
- Low memory prevents twisting and kinking
- Bubble at every 36" (91 cm) can be cut to form a tapered or funnel connector
- Dispensed in a single roll for ease of use

### Argyle™ Oxygen Bubble Tubing

CODE	DESCRIPTION	QTY
8888230201	Argyle™ Oxygen Bubble Tubing 3 mm x 30 m – Green	1

## Argyle™ Suction Tubing

Argyle™ surgical suction connecting tubing is available with female molded connectors (male connector also included). Argyle™ connecting tubing has the following features:

- Wide range of sizes
- Clear PVC construction
- Striated walls to avoid tube collapse
- Minimal coil memory prevents kinking
- Sterile
- Latex free



Funnel connection Suction tubing

### Argyle™ Suction Tubing with Funnel Connectors

CODE	DESCRIPTION	QTY
8888284505	Argyle™ Connecting Tube Funnel Connectors 5 mm x 50 cm	50
8888284513	Argyle™ Connecting Tube Funnel Connectors 5 mm x 1.8 m	50
8888284539	Argyle™ Connecting Tube Funnel Connectors 5 mm x 3.1 m	50
8888284547	Argyle™ Connecting Tube Funnel Connectors 5 mm x 3.7 m	20
8888284554	Argyle™ Connecting Tube Funnel Connectors 5 mm x 6.1 m	10
8888284604	Argyle™ Connecting Tube Funnel Connectors 6 mm x 1.8 m	50
8888284612	Argyle™ Connecting Tube Funnel Connectors 6 mm x 3.1 m	50
8888284620	Argyle™ Connecting Tube Funnel Connectors 6 mm x 3.7 m	20
8888284638	Argyle™ Connecting Tube Funnel Connectors 6 mm x 6.1 m	10
8888284703	Argyle™ Connecting Tube Funnel Connectors 7 mm x 1.8 m	50
8888284737	Argyle™ Connecting Tube Funnel Connectors 7 mm x 3.7 m	50
8888284745	Argyle™ Connecting Tube Funnel Connectors 7 mm x 6.1 m	10



8888284521

### Argyle™ Suction Tubing with Molded to Tapered Connectors

CODE	DESCRIPTION	QTY
8888284521	Argyle™ Connecting Tube Funnel / Tapered Connectors 5 mm x 1.8 m	50
8888284711	Argyle™ Connecting Tube Funnel / Tapered Connectors 7 mm x 1.8 m	50



### Argyle™ Suction Tubing with Molded Connectors

CODE	DESCRIPTION	QTY
8888301507	Argyle™ Connecting Tube w/Molded Connectors F/F with Male adaptor* 5 mm x 50 cm	50
8888301523	Argyle™ Connecting Tube w/Molded Connectors F/F with Male adaptor* 5 mm x 3.1 m	50
8888301531	Argyle™ Connecting Tube w/Molded Connectors F/F with Male adaptor* 5 mm x 3.7 m	20
8888301606	Argyle™ Connecting Tube w/Molded Connectors F/F with Male adaptor* 6 mm x 1.8 m	50
8888301614	Argyle™ Connecting Tube w/Molded Connectors F/F with Male adaptor* 6 mm x 3.1 m	50
8888301622	Argyle™ Connecting Tube w/Molded Connectors F/F with Male adaptor* 6 mm x 3.7 m	20
8888301705	Argyle™ Connecting Tube w/Molded Connectors F/F with Male adaptor* 7 mm x 1.8 m	50
8888301713	Argyle™ Connecting Tube w/Molded Connectors F/F with Male adaptor* 7 mm x 3.1 m	50
8888301721	Argyle™ Connecting Tube w/Molded Connectors F/F with Male adaptor* 7 mm x 3.7 m	20



## Medi-Vac™ Suction Tubing

- Medical grade PVC  
Provides easy handling and drapability with low coil memory
- Thick tubing walls  
Offer safety and collapse-resistance at high vacuum pressures

- Maxi-Grip™ Connectors  
Provides a secure fit on a wide range of surgical instruments and equipment
- Rigid male connector
- Single use & sterile

### Medi-Vac™ Suction Tubing

CODE	DESCRIPTION	QTY
NEX510	5 mm x 3.05 m Sterile Tubing	30 ea/cs
NEX56A	5 mm x 1.8 m Sterile Tubing	50 ea/cs
NEX610	6 mm x 3.05 m Sterile Tubing	25 ea/cs
NEX66A	6 mm x 1.8 m Sterile Tubing	45 ea/cs
NEX710	7 mm x 3.05 m Sterile Tubing	20 ea/cs
NEX712	7 mm x 3.60 m Sterile Tubing	20 ea/cs
NEX76A	7 mm x 1.8 m Sterile Tubing	40 ea/cs

Sure-Grip female molded connectors provides exceptional grip



# QUICK REFERENCE GUIDE

## Argyle™ Bubble Tubing PAGE 20

CODE	DESCRIPTION	QTY
8888280214	Argyle™ Bubble Tubing 5 mm x 30.5 m	1
8888290213	Argyle™ Bubble Tubing 5 mm x 30.5 m – Conductive	1
8888280602	Argyle™ Bubble tubing 5 mm x 30 m – Non Conductive	1
8888280412	Argyle™ Bubble Tubing 6 mm x 30.5 m	1
1180839035	Argyle™ Bubble Tubing 6 mm x 50 m	1
8888280610	Argyle™ Bubble Tubing 7 mm x 30.5 m	1
8888290619	Argyle™ Bubble Tubing 7 mm x 30.5 m – Conductive	1
8888280206	Argyle™ Bubble tubing 7 mm x 30 m – Non Conductive	1

## Argyle™ Oxygen Bubble Tubing PAGE 21

CODE	DESCRIPTION	QTY
8888230201	Argyle™ Oxygen Bubble Tubing 3 mm x 30 m – Green	1

## Argyle™ Suction Tubing PAGE 21-22

### Argyle™ Suction Tubing with Funnel Connectors

CODE	DESCRIPTION	QTY
8888284505	Argyle™ Connecting Tube Funnel Connectors 5 mm x 50 cm	50
8888284513	Argyle™ Connecting Tube Funnel Connectors 5 mm x 1.8 m	50
8888284539	Argyle™ Connecting Tube Funnel Connectors 5 mm x 3.1 m	50
8888284547	Argyle™ Connecting Tube Funnel Connectors 5 mm x 3.7 m	20
8888284554	Argyle™ Connecting Tube Funnel Connectors 5 mm x 6.1 m	10
8888284604	Argyle™ Connecting Tube Funnel Connectors 6 mm x 1.8 m	50
8888284612	Argyle™ Connecting Tube Funnel Connectors 6 mm x 3.1 m	50
8888284620	Argyle™ Connecting Tube Funnel Connectors 6 mm x 3.7 m	20
8888284638	Argyle™ Connecting Tube Funnel Connectors 6 mm x 6.1 m	10
8888284703	Argyle™ Connecting Tube Funnel Connectors 7 mm x 1.8 m	50
8888284737	Argyle™ Connecting Tube Funnel Connectors 7 mm x 3.7 m	50
8888284745	Argyle™ Connecting Tube Funnel Connectors 7 mm x 6.1 m	10

### Argyle™ Suction Tubing with Molded to Tapered Connectors

CODE	DESCRIPTION	QTY
8888284521	Argyle™ Connecting Tube Funnel / Tapered Connectors 5 mm x 1.8 m	50
8888284711	Argyle™ Connecting Tube Funnel / Tapered Connectors 7 mm x 1.8 m	50

## Argyle™ Suction Tubing with Molded Connectors

CODE	DESCRIPTION	QTY
8888301507	Argyle™ Connecting Tube w/Molded Connectors F/F with Male adaptor* 5 mm x 50 cm	50
8888301523	Argyle™ Connecting Tube w/Molded Connectors F/F with Male adaptor* 5 mm x 3.1 m	50
8888301531	Argyle™ Connecting Tube w/Molded Connectors F/F with Male adaptor* 5 mm x 3.7 m	20
8888301606	Argyle™ Connecting Tube w/Molded Connectors F/F with Male adaptor* 6 mm x 1.8 m	50
8888301614	Argyle™ Connecting Tube w/Molded Connectors F/F with Male adaptor* 6 mm x 3.1 m	50
8888301622	Argyle™ Connecting Tube w/Molded Connectors F/F with Male adaptor* 6 mm x 3.7 m	20
8888301705	Argyle™ Connecting Tube w/Molded Connectors F/F with Male adaptor* 7 mm x 1.8 m	50
8888301713	Argyle™ Connecting Tube w/Molded Connectors F/F with Male adaptor* 7 mm x 3.1 m	50
8888301721	Argyle™ Connecting Tube w/Molded Connectors F/F with Male adaptor* 7 mm x 3.7 m	20

## Medi-Vac™ Suction Tubing PAGE 22

### Medi-Vac™ Suction Tubing

CODE	DESCRIPTION	QTY
NEX510	5 mm x 3.05 m Sterile Tubing	30 ea/cs
NEX56A	5 mm x 1.8 m Sterile Tubing	50 ea/cs
NEX610	6 mm x 3.05 m Sterile Tubing	25 ea/cs
NEX66A	6 mm x 1.8 m Sterile Tubing	45 ea/cs
NEX710	7 mm x 3.05 m Sterile Tubing	20 ea/cs
NEX712	7 mm x 3.60 m Sterile Tubing	20 ea/cs
NEX76A	7 mm x 1.8 m Sterile Tubing	40 ea/cs



## ... and this was just a portion of our Complete Surgical Suction Portfolio

---

Fluid management systems and surgical suction instruments  
devices that meet the needs of every Operating Room



### **Medi-Vac™ Suction and Fluid Collection Products**

Fluid management products, all in one place. Our Flex Advantage™ single-use systems deliver industry-leading performance.



### **Argyle™ Surgical Suction Yankauers**

Flexible and rigid Yankauers with a choice of tip design and the option of pre-connected suction tubes.



### **Argyle™ Suction Connecting Tubes**

Surgical suction connecting tubes available with molded and funnel connectors.



### **Argyle™ Suction Catheters**

Atraumatic suction catheters with choice of tip designs and optional connectors for different suction sources.



### **Accessories, Connectors and Adaptors**

Accessories and adaptors for the Medi-Vac™ fluid management system.



**CardinalHealth**  
*Essential to care™*

For healthcare professionals only.

© 2020 Cardinal Health. All Rights Reserved. CARDINAL HEALTH, the Cardinal Health LOGO, ESSENTIAL TO CARE, MEDI-VAC, FLEX ADVANTAGE and ARGYLE are trademarks of Cardinal Health and may be registered in the US and/or in other countries. All other marks are the property of their respective owners. Important information: Prior to use, refer to the instructions for use supplied with this device for indications, contraindications, side effects, suggested procedure, warnings and precautions. For healthcare professionals only. As part of its continuous product development policy, Cardinal Health reserves the right to change product specifications without prior notification. Please contact your Cardinal Health representative for additional product availability information. 2SRV19-1071567 01/2020