

Your **Power** for Health


greiner bio-one



VACUETTE® QUICKSHIELD Complete PLUS

See the 'flash' upon
successful venipuncture


VACUETTE® QUICKSHIELD Complete PLUS

Disposable Safety Holder with VACUETTE® VISIO PLUS Multi-Drawing Needle

VACUETTE® QUICKSHIELD Complete PLUS combines the safety of QUICKSHIELD Safety Tube Holder with rotational shield and the innovation of VISIO PLUS Multi-Drawing Needle.

- The safety shield is attached directly to the holder and can be rotated to the user-preferred position. The safety feature is activated by pressing on a solid surface or by thumb activation. Hands are always behind the needle and the needle is locked into the shield after activation.
- The multi-drawing needle has a clear translucent plastic hub/view window which provides visual “flash” confirming vein penetration.
- QUICKSHIELD Complete PLUS is individually wrapped and sterilized to ensure the ultimate hygienic product delivery.
- The unit is fully assembled to save time.

Unique clear translucent hub lets you see ‘red’ upon successful venipuncture



Ordering Information

Item Number	Description	Inner	Outer
450235	VACUETTE® QUICKSHIELD Complete PLUS (Safety Tube Holder and VISIO PLUS 21G x 1" Multi-Drawing Needle Combination)	30 pcs.	300 pcs.
450239	VACUETTE® QUICKSHIELD Complete PLUS (Safety Tube Holder and VISIO PLUS 21G x 1.5" Multi-Drawing Needle Combination)	30 pcs.	300 pcs.
450226	VACUETTE® QUICKSHIELD Complete PLUS (Safety Tube Holder and VISIO PLUS 22G x 1" Multi-Drawing Needle Combination)	30 pcs.	300 pcs.
450228	VACUETTE® QUICKSHIELD Complete PLUS (Safety Tube Holder and VISIO PLUS 22G x 1.5" Multi-Drawing Needle Combination)	30 pcs.	300 pcs.

Greiner Bio-One, Inc.
4238 Capital Drive
Monroe, NC 28110
USA

Phone 888.286.3883
Fax 800.726.0052
office@us.gbo.com
www.gbo.com