



Argyle™
Surgical Suction
Connecting Tubes

Argyle™ Surgical Suction Connecting Tubes

A comprehensive range of surgical suction connecting tubes available with molded and funnel connectors.



Argyle™ Bubble Tubing

Can be easily cut to exact length and is available in conductive and non-conductive variations.

- PVC formulation
- Low memory prevents twisting and kinking
- Various lumen sizes
- Bubble at every 36" (91 cm) can be cut to form a tapered or funnel connector
- Dispensed in a single roll for ease of use

Argyle™ Bubble Tubing

| CODE | DESCRIPTION | QTY |
|------------|--|-----|
| 8888280214 | Argyle™ Bubble Tubing 5 mm x 30.5 m | 1 |
| 8888290213 | Argyle™ Bubble Tubing 5 mm x 30.5 m – Conductive | 1 |
| 8888280602 | Argyle™ Bubble tubing 7 mm x 30 m – Non Conductive | 1 |
| 8888280412 | Argyle™ Bubble Tubing 6 mm x 30.5 m | 1 |
| 1180839035 | Argyle™ Bubble Tubing 6 mm x 50 m | 1 |
| 8888280610 | Argyle™ Bubble Tubing 7 mm x 30.5 m | 1 |
| 8888290619 | Argyle™ Bubble Tubing 7 mm x 30.5 m – Conductive | 1 |
| 8888280206 | Argyle™ Bubble tubing 5 mm x 30 m – Non Conductive | 1 |





8888230201

Argyle™ Oxygen Bubble Tubing

- PVC formulation
- Low memory prevents twisting and kinking
- Bubble at every 36" (91 cm) can be cut to form a tapered or funnel connector
- Dispensed in a single roll for ease of use

Argyle™ Oxygen Bubble Tubing

| CODE | DESCRIPTION | QTY |
|------------|--|-----|
| 8888230201 | Argyle™ Oxygen Bubble Tubing 3 mm x 30 m – Green | 1 |

Argyle™ Suction Tubing

Argyle™ surgical suction connecting tubing is available with female molded connectors (male connector also included). Argyle™ connecting tubing has the following features:

- Wide range of sizes
- Clear PVC construction
- Striated walls to avoid tube collapse
- Minimal coil memory prevents kinking
- Sterile
- Latex free



Funnel connection Suction tubing

Argyle™ Suction Tubing with Funnel Connectors

| CODE | DESCRIPTION | QTY |
|------------|--|-----|
| 8888284505 | Argyle™ Connecting Tube Funnel Connectors 5 mm x 50 cm | 50 |
| 8888284513 | Argyle™ Connecting Tube Funnel Connectors 5 mm x 1.8 m | 50 |
| 8888284539 | Argyle™ Connecting Tube Funnel Connectors 5 mm x 3.1 m | 50 |
| 8888284547 | Argyle™ Connecting Tube Funnel Connectors 5 mm x 3.7 m | 20 |
| 8888284554 | Argyle™ Connecting Tube Funnel Connectors 5 mm x 6.1 m | 10 |
| 8888284604 | Argyle™ Connecting Tube Funnel Connectors 6 mm x 1.8 m | 50 |
| 8888284612 | Argyle™ Connecting Tube Funnel Connectors 6 mm x 3.1 m | 50 |
| 8888284620 | Argyle™ Connecting Tube Funnel Connectors 6 mm x 3.7 m | 20 |
| 8888284638 | Argyle™ Connecting Tube Funnel Connectors 6 mm x 6.1 m | 10 |
| 8888284703 | Argyle™ Connecting Tube Funnel Connectors 7 mm x 1.8 m | 50 |
| 8888284737 | Argyle™ Connecting Tube Funnel Connectors 7 mm x 3.7 m | 50 |
| 8888284745 | Argyle™ Connecting Tube Funnel Connectors 7 mm x 6.1 m | 10 |



8888284521

Argyle™ Suction Tubing with Molded to Tapered Connectors

| CODE | DESCRIPTION | QTY |
|------------|--|-----|
| 8888284521 | Argyle™ Connecting Tube Funnel / Tapered Connectors 5 mm x 1.8 m | 50 |
| 8888284711 | Argyle™ Connecting Tube Funnel / Tapered Connectors 7 mm x 1.8 m | 50 |



Argyle™ Suction Tubing with Molded Connectors

| CODE | DESCRIPTION | QTY |
|------------|---|-----|
| 8888301507 | Argyle™ Connecting Tube w/Molded Connectors F/F with Male adaptor* 5 mm x 50 cm | 50 |
| 8888301523 | Argyle™ Connecting Tube w/Molded Connectors F/F with Male adaptor* 5 mm x 3.1 m | 50 |
| 8888301531 | Argyle™ Connecting Tube w/Molded Connectors F/F with Male adaptor* 5 mm x 3.7 m | 20 |
| 8888301606 | Argyle™ Connecting Tube w/Molded Connectors F/F with Male adaptor* 6 mm x 1.8 m | 50 |
| 8888301614 | Argyle™ Connecting Tube w/Molded Connectors F/F with Male adaptor* 6 mm x 3.1 m | 50 |
| 8888301622 | Argyle™ Connecting Tube w/Molded Connectors F/F with Male adaptor* 6 mm x 3.7 m | 20 |
| 8888301705 | Argyle™ Connecting Tube w/Molded Connectors F/F with Male adaptor* 7 mm x 1.8 m | 50 |
| 8888301713 | Argyle™ Connecting Tube w/Molded Connectors F/F with Male adaptor* 7 mm x 3.1 m | 50 |
| 8888301721 | Argyle™ Connecting Tube w/Molded Connectors F/F with Male adaptor* 7 mm x 3.7 m | 20 |



Medi-Vac™ Suction Tubing

- Medical grade PVC
Provides easy handling and drapability with low coil memory
- Thick tubing walls
Offer safety and collapse-resistance at high vacuum pressures

- Maxi-Grip™ Connectors
Provides a secure fit on a wide range of surgical instruments and equipment
- Rigid male connector
- Single use & sterile

Medi-Vac™ Suction Tubing

| CODE | DESCRIPTION | QTY |
|--------|------------------------------|----------|
| NEX510 | 5 mm x 3.05 m Sterile Tubing | 30 ea/cs |
| NEX56A | 5 mm x 1.8 m Sterile Tubing | 50 ea/cs |
| NEX610 | 6 mm x 3.05 m Sterile Tubing | 25 ea/cs |
| NEX66A | 6 mm x 1.8 m Sterile Tubing | 45 ea/cs |
| NEX710 | 7 mm x 3.05 m Sterile Tubing | 20 ea/cs |
| NEX712 | 7 mm x 3.60 m Sterile Tubing | 20 ea/cs |
| NEX76A | 7 mm x 1.8 m Sterile Tubing | 40 ea/cs |

Sure-Grip female molded connectors provides exceptional grip



QUICK REFERENCE GUIDE

Argyle™ Bubble Tubing PAGE 20

| CODE | DESCRIPTION | QTY |
|------------|--|-----|
| 8888280214 | Argyle™ Bubble Tubing 5 mm x 30.5 m | 1 |
| 8888290213 | Argyle™ Bubble Tubing 5 mm x 30.5 m – Conductive | 1 |
| 8888280602 | Argyle™ Bubble tubing 5 mm x 30 m – Non Conductive | 1 |
| 8888280412 | Argyle™ Bubble Tubing 6 mm x 30.5 m | 1 |
| 1180839035 | Argyle™ Bubble Tubing 6 mm x 50 m | 1 |
| 8888280610 | Argyle™ Bubble Tubing 7 mm x 30.5 m | 1 |
| 8888290619 | Argyle™ Bubble Tubing 7 mm x 30.5 m – Conductive | 1 |
| 8888280206 | Argyle™ Bubble tubing 7 mm x 30 m – Non Conductive | 1 |

Argyle™ Oxygen Bubble Tubing PAGE 21

| CODE | DESCRIPTION | QTY |
|------------|--|-----|
| 8888230201 | Argyle™ Oxygen Bubble Tubing 3 mm x 30 m – Green | 1 |

Argyle™ Suction Tubing PAGE 21-22

Argyle™ Suction Tubing with Funnel Connectors

| CODE | DESCRIPTION | QTY |
|------------|--|-----|
| 8888284505 | Argyle™ Connecting Tube Funnel Connectors 5 mm x 50 cm | 50 |
| 8888284513 | Argyle™ Connecting Tube Funnel Connectors 5 mm x 1.8 m | 50 |
| 8888284539 | Argyle™ Connecting Tube Funnel Connectors 5 mm x 3.1 m | 50 |
| 8888284547 | Argyle™ Connecting Tube Funnel Connectors 5 mm x 3.7 m | 20 |
| 8888284554 | Argyle™ Connecting Tube Funnel Connectors 5 mm x 6.1 m | 10 |
| 8888284604 | Argyle™ Connecting Tube Funnel Connectors 6 mm x 1.8 m | 50 |
| 8888284612 | Argyle™ Connecting Tube Funnel Connectors 6 mm x 3.1 m | 50 |
| 8888284620 | Argyle™ Connecting Tube Funnel Connectors 6 mm x 3.7 m | 20 |
| 8888284638 | Argyle™ Connecting Tube Funnel Connectors 6 mm x 6.1 m | 10 |
| 8888284703 | Argyle™ Connecting Tube Funnel Connectors 7 mm x 1.8 m | 50 |
| 8888284737 | Argyle™ Connecting Tube Funnel Connectors 7 mm x 3.7 m | 50 |
| 8888284745 | Argyle™ Connecting Tube Funnel Connectors 7 mm x 6.1 m | 10 |

Argyle™ Suction Tubing with Molded to Tapered Connectors

| CODE | DESCRIPTION | QTY |
|------------|--|-----|
| 8888284521 | Argyle™ Connecting Tube Funnel / Tapered Connectors 5 mm x 1.8 m | 50 |
| 8888284711 | Argyle™ Connecting Tube Funnel / Tapered Connectors 7 mm x 1.8 m | 50 |

Argyle™ Suction Tubing with Molded Connectors

| CODE | DESCRIPTION | QTY |
|------------|---|-----|
| 8888301507 | Argyle™ Connecting Tube w/Molded Connectors F/F with Male adaptor* 5 mm x 50 cm | 50 |
| 8888301523 | Argyle™ Connecting Tube w/Molded Connectors F/F with Male adaptor* 5 mm x 3.1 m | 50 |
| 8888301531 | Argyle™ Connecting Tube w/Molded Connectors F/F with Male adaptor* 5 mm x 3.7 m | 20 |
| 8888301606 | Argyle™ Connecting Tube w/Molded Connectors F/F with Male adaptor* 6 mm x 1.8 m | 50 |
| 8888301614 | Argyle™ Connecting Tube w/Molded Connectors F/F with Male adaptor* 6 mm x 3.1 m | 50 |
| 8888301622 | Argyle™ Connecting Tube w/Molded Connectors F/F with Male adaptor* 6 mm x 3.7 m | 20 |
| 8888301705 | Argyle™ Connecting Tube w/Molded Connectors F/F with Male adaptor* 7 mm x 1.8 m | 50 |
| 8888301713 | Argyle™ Connecting Tube w/Molded Connectors F/F with Male adaptor* 7 mm x 3.1 m | 50 |
| 8888301721 | Argyle™ Connecting Tube w/Molded Connectors F/F with Male adaptor* 7 mm x 3.7 m | 20 |

Medi-Vac™ Suction Tubing PAGE 22

Medi-Vac™ Suction Tubing

| CODE | DESCRIPTION | QTY |
|--------|------------------------------|----------|
| NEX510 | 5 mm x 3.05 m Sterile Tubing | 30 ea/cs |
| NEX56A | 5 mm x 1.8 m Sterile Tubing | 50 ea/cs |
| NEX610 | 6 mm x 3.05 m Sterile Tubing | 25 ea/cs |
| NEX66A | 6 mm x 1.8 m Sterile Tubing | 45 ea/cs |
| NEX710 | 7 mm x 3.05 m Sterile Tubing | 20 ea/cs |
| NEX712 | 7 mm x 3.60 m Sterile Tubing | 20 ea/cs |
| NEX76A | 7 mm x 1.8 m Sterile Tubing | 40 ea/cs |



... and this was just a portion of our Complete Surgical Suction Portfolio

Fluid management systems and surgical suction instruments
devices that meet the needs of every Operating Room



Medi-Vac™ Suction and Fluid Collection Products

Fluid management products, all in one place. Our Flex Advantage™ single-use systems deliver industry-leading performance.



Argyle™ Surgical Suction Yankauers

Flexible and rigid Yankauers with a choice of tip design and the option of pre-connected suction tubes.



Argyle™ Suction Connecting Tubes

Surgical suction connecting tubes available with molded and funnel connectors.



Argyle™ Suction Catheters

Atraumatic suction catheters with choice of tip designs and optional connectors for different suction sources.



Accessories, Connectors and Adaptors

Accessories and adaptors for the Medi-Vac™ fluid management system.



CardinalHealth
Essential to care™

For healthcare professionals only.

© 2020 Cardinal Health. All Rights Reserved. CARDINAL HEALTH, the Cardinal Health LOGO, ESSENTIAL TO CARE, MEDI-VAC, FLEX ADVANTAGE and ARGYLE are trademarks of Cardinal Health and may be registered in the US and/or in other countries. All other marks are the property of their respective owners. Important information: Prior to use, refer to the instructions for use supplied with this device for indications, contraindications, side effects, suggested procedure, warnings and precautions. For healthcare professionals only. As part of its continuous product development policy, Cardinal Health reserves the right to change product specifications without prior notification. Please contact your Cardinal Health representative for additional product availability information. 2SRV19-1071567 01/2020