

# UVENTA™

## Ureteral Stent

*for malignant ureteral obstruction*

**A quantum leap in stent technology for  
a new evolution in palliative therapy**



#### HEAD OFFICE

14, Gojeong-ro, Wolgot-myeon, Gimpo-si, Gyeonggi-do, Korea 10022

#### SALES OFFICE

Bldg. C 4F, 158 Haneulmaeul-ro, Ilsandong-gu, Goyang-si, Gyeonggi-do, Korea 10355  
E: lea@stent.net T: +82 31 904 6153 F: +82 31 904 6157

[www.taewoongmedical.com](http://www.taewoongmedical.com)

**TaeWoong**  
MEDICAL 

# UVENTA™ Ureteral Stent




## Indication

- Malignant ureteral obstruction

## Features

- Triple layers (Mesh–Membrane–Mesh)
- Self-expandable
- Optimal radial force
- High conformability
- Large diameter for re-intervention
- High visibility under fluoroscopic view
- Pre-mounted on the delivery system for easy deployment
- Repositioning suture



## Ordering Information

Specification of UVENTA™ Ureteral Stent

Code	Stent		Delivery		Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)		Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
UC0606FB	6	6	9	70	UC0806FB	8	6	9	70
UC0607FB		7			UC0807FB		7		
UC0608FB		8			UC0808FB		8		
UC0609FB		9			UC0809FB		9		
UC0610FB		10			UC0810FB		10		
UC0611FB		11			UC0811FB		11		
UC0612FB	12	UC0812FB	12		10				
UC0706FB	7	6	9			UC0814FB	14		
UC0707FB		7				UC0816FB	16		
UC0708FB		8				UC0818FB	18		
UC0709FB		9				UC1006FB	6	9	
UC0710FB		10				UC1007FB	7		
UC0711FB		11		UC1008FB	8				
UC0712FB	12	UC1009FB	9						
UC0714FB	14	UC1010FB	10						
UC0716FB	16	UC1011FB	11	10					
UC0718FB	18	UC1012FB	12						
		UC1014FB	14		10				
		UC1016FB	16						
		UC1018FB	18						

## Released Article

- \* Comparison of Efficacy and Safety Between a Segmental Thermo-expandable Metal Alloy Spiral Stent (Memokath 051) and a Self-Expandable Covered Metallic Stent (UVENTA™) in the Management of Ureteral Obstructions [Journal of Laparoendoscopic & Advanced Surgical Techniques / 2014]
- \* Multicenter Experience of the Newly Designed Covered Metallic Ureteral Stent for Malignant Ureteral Occlusion: Comparison With Double J Stent Insertion [Cardiovasc Intervent Radiol / 08 August 2013]
- \* Efficacy and Safety of a Novel, Double-Layered, Coated, Self-Expandable Metallic Mesh Stent (UVENTA™) in Malignant Ureteral Obstructions [Journal of Endourology / Vol.27, No.7, July 2013]
- \* Palliative Care of Malignant Ureteral Obstruction with Polytetrafluoroethylene Membrane-Covered Self-Expandable Metallic Stents: Initial Experience [Korean Journal of Urology / 2012;53:625-631]