

Breathing circuits product catalog



6

Anesthesia

24

ICU

34

PICU/NICU

38

Anesthetic and
intensive care tubes

42

Catheter
mounts


49

Anesthetic and intensive
care accessories

Contents

Anesthesia	6	Intensive care adult (ICU)	24
Anesthetic systems	7	Intensive care systems	25
Adult anesthesia circuits	9	Adult standard circuits	25
Adult polyvinyl chloride (PVC) – smoothbore breathing systems	9	Polyethylene (PE) – corrugated breathing systems	26
Adult polyvinyl chloride (PVC) – smoothbore breathing systems with accessories	10	Adult standard circuit with 1 water trap	26
Adult polyvinyl chloride (PVC) – smoothbore breathing systems with swivel Y-piece	10	Polyvinyl chloride (PVC) – smoothbore breathing systems	26
Adult polypropylene (PP) – extendable breathing systems	11	Polypropylene (PP) – extendable breathing systems	27
Adult polyethylene (PE) – corrugated breathing systems	12	Adult standard circuit with 2 water traps	27
Pediatric anesthesia circuits	13	Polyvinyl chloride (PVC) – smoothbore breathing systems	27
Polyvinyl chloride (PVC) smoothbore breathing systems	13	Polypropylene (PP) – extendable breathing systems	28
Neonatal anesthesia circuits	14	Polyethylene (PE) – corrugated breathing systems	28
Polyvinyl chloride (PVC) smoothbore breathing systems	14	Adult insulated smoothbore breathing systems	29
Mapleson and Bain breathing systems adult, pediatric, neonatal	15	Adult PVC heated circuits for active humidifiers	30
Adult Mapleson C breathing systems	16	Adult PVC heated circuits for Fisher and Paykel™ MR730 active humidifier	31
Pediatric Mapleson C breathing systems	17	Adult heated circuits for Fisher and Paykel™ MR850 active humidifier	31
Infant Mapleson C breathing systems	18	Adult heated circuits for Kendall™ Aerodyne 2000 active humidifier	32
Mapleson F – pediatric and infant breathing systems	19	Adapter cables	33
Coaxial breathing circuits	20	Intensive care pediatric/neonatal (PICU/NICU)	34
3-in-1 coaxial breathing circuits	21	Pediatric standard circuits	35
Accessories	22	Pediatric circuits with water traps	35
Anesthesia face masks	22	Pediatric heated circuits	36
Guedel airways	23	Pediatric single heated circuits used with Fisher and Paykel™ MR730 and/or MR850 active humidifiers	36
		Neonatal standard circuits	37
		Neonatal circuits with water traps	37
		Neonatal single heated circuits	37
		Neonatal single heated circuits used with Fisher and Paykel™ MR730 and MR850 active humidifiers	37

Anesthetic and intensive care tubes	38	Anesthetic and intensive care accessories	49
Adult tubings	39	Connectors	50
Adult PVC tubings	39	Angled patient connectors	50
Adult PVC tubings with accessories	39	Y-pieces	50
Adult PP extendable tubings	39	Straight, short connectors	50
Adult PP extendable tubings with accessories	40	Straight, long connectors.....	50
Adult PE tubings	40	L adapters	51
Pediatric tubings	41	T adapters	52
Pediatric PVC tubings	41	Caps	53
Neonatal tubings	41	Water traps	53
Neonatal PVC tubings	41	Valves	53
Catheter mounts	42	Mouthpieces	53
Catheter mounts	43	Breathing bags and accessories	54
PVC catheter mounts	43	With 22 female flexible connector	54
PVC catheter mounts with double swivel elbow connector	43	With 22 female rigid connector	54
PVC catheter mounts with 120-degree elbow connector	44	Accessories	54
PVC catheter mounts with 90-degree elbow connector	45	Gas sampling lines	55
PVC catheter mounts with CO ₂ sampling port on		PVC pressure monitoring lines	55
90-degree elbow connector	45		
PVC straight catheter mounts	46		
PVC Carlens catheter mounts	46		
PP catheter mounts (extendable)	47		
Extendable PP catheter mounts with double-swivel elbow			
connector	47		
Extendable PP catheter mounts with 120-degree elbow			
connector	47		
Extendable PP catheter mounts with 90-degree elbow			
connector	47		
Extendable PP catheter mounts with CO ₂ sampling port on			
90-degree elbow connector	48		
Extendable PP straight catheter mounts	48		



Medtronic is a global leader in medical technology, services, and solutions.

We're on a mission to create and deliver innovative healthcare solutions, developed in collaboration with healthcare providers, to enhance the quality of life for patients and improve experiences for customers.

For more than 70 years, our brands have defined their categories in the healthcare field. Today we support the lifesaving work of medical professionals in more than 150 countries.

6

Anesthesia

24

ICU

34

PICU/NICU

38

Anesthetic and
intensive care tubes

42

Catheter
mounts

49

Anesthetic and intensive
care accessories

DAR™ breathing systems

Under the DAR™ brand, Medtronic offers a full range of breathing circuits and other respiratory care accessories for the OR, recovery rooms, intensive care units, and the home care environments.

DAR™ products include different types of single use systems to meet the requirements of medical staff and patients. Some system types are firsts in the industry and exclusive to the DAR™ brand, such as the heated-wire circuits with integrated heating element or the 3-in-1 DAR™ coaxial anesthesia circuit.

DAR™ breathing circuits are available in three different diameters:

- Adult (ø 22 mm)
- Pediatric (ø 15 mm)
- Neonatal (ø 10 mm)

Quality solutions fitted to your needs

Our respiratory disposables combine innovation and efficiency with exceptional reliability.

DAR™ quality system procedures involve an extensive series of controls at every stage of production, meeting all recognized quality assurance standards.

The devices are manufactured in a stringently controlled clean room environment.

Beyond the frequent controls made on raw materials, semi-finished, and finished products, all products undergo integrity testing at a pressure of 100 cm H₂O to detect any leaks or weak points.

Products are CE marked in compliance with European Medical Device Regulation (MDR) 2017/745/EU and manufactured according to ISO 13485.

DAR™ products are a natural extension to other Medtronic product lines, such as Shiley™ endotracheal tubes, Shiley™ tracheostomy tubes, inhalation therapy products – making Medtronic the complete solution to all your respiratory care requirements.



6

Anesthesia

24

ICU

34

PICU/NICU

38

Anesthetic and
intensive care tubes

42

Catheter
mounts

49

Anesthetic and intensive
care accessories

Anesthesia



6

Anesthesia

24

ICU

34

PICU/NICU

38

Anesthetic and
intensive care tubes

42

Catheter
mounts

49

Anesthetic and intensive
care accessories

Anesthetic systems

Medtronic offers a range of anesthetic breathing systems that meet all your day-to-day requirements in the anesthetic room, the operating room, and the recovery unit.

Circuits include basic circuits suitable for use with most ventilators or circle absorber systems and several versions of the most popular Mapleson configurations.



Circuits are available in adult, pediatric, and neonatal sizes.

We can also provide custom-made breathing systems.

All products are latex free.

Materials include:

- PVC smoothbore (high mechanical resistance and low compliance)
- Extendable polypropylene (minimal storage space required)
- Polyethylene (lightweight and cost-effective)



Polyvinyl chloride (PVC) - smoothbore breathing systems

DAR™ smoothbore breathing circuits ensure optimal ventilatory support by offering the following benefits:

- Totally clear tube makes constant monitoring of the circuit possible.
- Smooth inner surface allows an easy drainage of condensed water, reducing risk of bacterial growth.
- Reinforced external helix keeps the dimension of the tube diameter constant, even when bent.
- Extremely flexible and crush proof, reducing the risk of interruptions in ventilation
- Very low compliance
- Soft, leak safe, flexible connectors ensures a quality seal
- Single patient use

Polypropylene (PP) - extendable breathing systems

- Can extend up to three times their packed length, overcoming storage problems; even circuits of considerable length require very little storage room.
- "Cemory" feature: the system will maintain the position set by the operator for as long as required, minimizing strain on the ET tube.
- Available in adult size only (22 mm diameter)

Polyethylene (PE) - corrugated breathing systems

- Lightweight and cost-effective
- Connectors directly bonded to the tube to ensure gas tight connections and prevent leakage
- Available in 50 m rolls
- Available in blue color for scavenging applications
- Supplied clean or sterile



6

Anesthesia

24

ICU

34

PICU/NICU

38

Anesthetic and
intensive care tubes

42

Catheter
mounts


49

Anesthetic and intensive
care accessories

Adult anesthesia circuits

Adult polyvinyl chloride (PVC) - smoothbore breathing systems

Code	Description	Sterile	Notes
300/6042	120 cm	✓	
300/6043	150 cm	✓	
300/6044	180 cm	✓	
300/6037	180 cm, cap on Y-piece ^A	✓	
300/13783	180 cm, elbow on Y-piece	✓	
300/13774	150 cm, elbow with capno port on Y-piece ^B	✓	
300/6099	200 cm	✓	
300/6049	220 cm	✓	
300/6449	300 cm	✓	
300/14147	150 cm, detachable Y-piece	✓	
302S14192	280 cm, 2 mechanical filter large	✓	
301/13404	200 cm, ports on Y-piece	✓	
301/13786	180 cm, ports and cap on Y-piece, D angled flex connector for Puritan Bennett™ 7200 ventilators	✓	
301/6130	180 cm, ports and cap on Y-piece	✓	
301/6378	150 cm, ports on detachable Y-piece	✓	
301/6326	150 cm, ports on Y-piece	✓	
301/14105	180 cm, ports on detachable Y-piece	✓	



Box quantity: 10

6

Anesthesia

24

ICU

34

PICU/NICU

38

Anesthetic and
intensive care tubes

42


Catheter
mounts

49

Anesthetic and intensive
care accessories

Adult polyvinyl chloride (PVC) - smoothbore breathing systems with accessories

300/6365




Code	Description	Sterile	Notes
300/6365	PVC breathing system 120 cm, extra limb 60 cm	✓	2 L breathing bag
Also available			
300S14224	150 cm, extra limb 160 cm	✓	2 L breathing bag
300S14385	150 cm, extra limb 100 cm	✓	2 L breathing bag with hook
300S14269	180 cm, extra limb 120 cm, elbow with capno	✓	2 L breathing bag

Box quantity: 10

Adult polyvinyl chloride (PVC) - smoothbore breathing systems with swivel Y-piece

300R14421




Code	Description	Sterile	Notes
300R14421	PVC breathing system 200 cm, extra limb 200 cm with 2 L breathing bag, straight connector 22M-22M/19F		For multi-patient use
Also available			
300R14422	300 cm, extra limb 300 cm with 2 L breathing bag, straight connector 22M-22M/19F		For multi-patient use

Box quantity: 10

10


Adult polypropylene (PP) - extendable breathing systems

300/13316



Code	Description	Sterile	Notes
300/13316	PP breathing system, 60/180 cm		
Also available			
300/13532	60/180 cm	✓	
300/13620	60/180 cm, 20 cm PVC catheter mount	✓	Swivel mount
300/13315	60/180 cm, elbow with capno on Y-piece		
300/13321	60/180 cm, angled elbow on Y-piece		
300/13811	60/180 cm, detachable Y-piece with angled elbow		
Box quantity: 10			

300/13330




Code	Description	Sterile	Notes
300/13330	PP extendable breathing system 60/180 cm, extra limb 45/120 cm, elbow with capno on Y-piece		2 L breathing bag
Also available			
300/13324	60/180 cm, extra limb 45/120 cm, elbow on Y-piece		2 L breathing bag
300/13329	60/180 cm, extra limb 45/120 cm		2 L breathing bag
300/13815	60/180 cm, extra limb 60/180 cm, elbow with capno on Y-piece	✓	2 L breathing bag
*300P16000	60/180 cm, extra limb 60/180cm, 90 degree elbow with capno on Y-piece, with electrostatic filter	✓	3 L breathing bag
Box quantity: 10			

*Note: This product is only available in Latin America.


Adult polyethylene (PE) - corrugated breathing systems

300NP14344




Code	Description	Sterile	Notes
300NP14344	PE breathing system, 150 cm, without ports, cap on Y-piece		For multi-patient use
Also available			
300NP14345	190 cm		For multi-patient use
300NP14348	190 cm, elbow with capno and cap on Y-piece		For multi-patient use
Box quantity: 10			

301NP14324



Code	Description	Sterile	Notes
301NP14324	PE breathing system 150 cm, ports on Y-piece, elbow with capno and cap on Y-piece		For multi-patient use
Box quantity: 10			

301NP14330




Code	Description	Sterile	Notes
301NP14330	PE breathing system 150 cm, 2 L breathing bag, extra limb 92 cm, elbow with capno and cap on Y-piece		For multi-patient use
Also available			
301NP14331	190 cm, extra limb 92 cm, 2 L breathing bag		For multi-patient use
301NP14332	240 cm, extra limb 92 cm, 2 L breathing bag		For multi-patient use
301NP14333	150 cm, extra limb 92 cm, 2 L breathing bag, electrostatic filter small		For multi-patient use
Box quantity: 10			

Pediatric anesthesia circuits


Polyvinyl chloride (PVC) smoothbore breathing systems

305/6748




Code	Description	Sterile	Notes
305/6748	100 cm	✓	
Also available			
305/6712	120 cm	✓	
305/6713	150 cm	✓	
305/6714	180 cm	✓	
305/6967	220 cm	✓	
305/6722	300 cm	✓	
Box quantity: 10			

305/8804



Code	Description	Sterile	Notes
305/8804	180 cm, ports on Y-piece	✓	
Box quantity: 10			

305/8773




Code	Description	Sterile	Notes
305/8773	120 cm, 1 L breathing bag, extra limb 120 cm	✓	
Box quantity: 10			

Neonatal anesthesia circuits

Polyvinyl chloride (PVC) smoothbore breathing systems


305/6980



Code	Description	Sterile	Notes
305/6980	Neonatal PVC smoothbore circuit 120 cm, no ports on Y-piece	✓	22 flex ventilator connections
Also available			
305/8310	150 cm	✓	22 flex ventilator connections

Box quantity: 10

305/6852



Code	Description	Sterile	Notes
305/6852	120 cm	✓	22 flex ventilator connections
305/8173	150 cm	✓	22 flex ventilator connections
305/6745	180 cm	✓	22 flex ventilator connections

Box quantity: 10

14

Mapleson and Bain breathing systems: adults, pediatric, neonatal

Our range of compact disposable Mapleson systems offer a choice of open anesthetic breathing systems.

APL valve

The adjustable pressure limiting valve maintains a constant pressure in the breathing system. Fresh gas flow, along with manual compression of the breathing bag, deliver an increasing pressure. The excess gas will be vented through the APL valve to the scavenging system without allowing entry of the outside air into the breathing system during inspiration.

The APL valve also inhibits overinflation of the reservoir bag allowing a lower resistance to expiration and preventing trauma to the airways.

Unique geometric design without spring helps to eliminate valve blockage in the event of mechanical failure and allows fine control of respiratory pressures. It is also fitted with a shrouded magnetic pressure relief valve to prevent barotrauma.

Mapleson C breathing systems

Lightweight and compact, DAR™ Mapleson C systems allow manual administration of gas mixtures and are indicated for application during anesthesia, resuscitation, patient transport, tube replacement, and bagging. They represent a possible alternative to the self-inflating resuscitation bags in ICU, recovery, and emergency departments.

Mapleson D - Bain breathing systems


Adult coaxial anesthesia systems with 90° angled elbow at patient side are ideal for use on both spontaneous breathing as well as mechanically ventilated patients.

Mapleson F - pediatric and infant breathing systems

The Mapleson F systems have a low dead space as the fresh gas flow enters directly at the removable 90° elbow connector. For that reason, the Mapleson F are the preferred Anesthetic systems for pediatric and infant patients, allowing manual administration of gas mixtures in ICU, recovery, and emergency departments.

Adult Mapleson C breathing systems

309/14059




Code	Description	Sterile	Notes
309/14059	Adult Mapleson C with lateral gas inlet, 180 cm tube, 2 L breathing bag, elbow and APL valve with overpressure relief valve	✓	

Also available

309/6612	Adult Mapleson C with lateral gas inlet, 180 cm tube, 2 L breathing bag, elbow and APL valve	✓	APL valve without overpressure limiting valve
309/13362	Adult Mapleson C with lateral gas inlet, 180 cm tube, 2 L breathing bag, straight connector with APL valve	✓	APL valve without overpressure limiting valve
309R14254	Adult Mapleson C with lateral gas inlet, 180 cm tube, 2 L breathing bag, straight connector with cap and APL valve with overpressure limiting valve		For multi-patient use

Box quantity: 10

309/6310



Code	Description	Sterile	Notes
309/6310	Adult Mapleson C with vertical gas inlet, 180 cm tube, 2 L breathing bag, elbow and APL valve	✓	APL valve without overpressure limiting valve


Also available

309S14217	Adult Mapleson C with vertical gas inlet, 180 cm tube, 2 L breathing bag, elbow and APL valve, cap	✓	APL valve without overpressure limiting valve
309/6498	Adult Mapleson C with vertical gas inlet, 300 cm tube, 2 L breathing bag, elbow and APL valve	✓	APL valve without overpressure limiting valve

Box quantity: 10

Pediatric Mapleson C breathing systems

309/6353




Code	Description	Sterile	Notes
309/6353	Pediatric Mapleson C with lateral gas inlet, 180 cm tube, 1 L breathing bag, elbow and APL valve	✓	APL valve with overpressure limiting valve and one-way valve on elbow

Also available			
309S14219	Pediatric Mapleson C with lateral gas inlet, 180 cm tube, 1 L breathing bag, elbow and APL valve	✓	APL valve without overpressure limiting valve

Box quantity: 10

309/6563




Code	Description	Sterile	Notes
309/6563	Pediatric Mapleson C with vertical gas inlet, 180 cm tube, 1 L breathing bag, elbow and APL valve	✓	APL valve without overpressure limiting valve

Box quantity: 10

Infant Mapleson C breathing systems

309/6351




Code	Description	Sterile	Notes
309/6351	Infant Mapleson C with lateral gas inlet, 180 cm tube, 1/2 L breathing bag, elbow and APL valve	✓	APL valve with overpressure limiting valve and one-way valve on elbow

Also available

309S14220	Infant Mapleson C with lateral gas inlet, 180 cm tube, 1/2 L breathing bag, elbow and APL valve, cap	✓	APL valve with overpressure limiting valve and one-way valve on elbow
-----------	--	---	---

Box quantity: 10


309/6355



Code	Description	Sterile	Notes
309/6355	Infant Mapleson C with vertical gas inlet, 180 cm tube, 1/2 L breathing bag, elbow and APL valve	✓	APL valve without overpressure limiting valve

Box quantity: 10

Mapleson F - pediatric and infant breathing systems



Code	Description	Sterile	Notes
286/7814	50 cm Mapleson F, Jackson Rees modification, with 1/2 L open tail loose breathing bag, scavenging tube 180 cm		
Also available			
286R25691	50 cm Mapleson F, Jackson Rees modification, with 1/2 L open tail breathing bag, scavenging tube 180 cm		For multi-patient use
Box quantity: 10			

Straight connector with APL valve, sterile

616/9048

Box quantity: 50

Elbow connector with APL valve, sterile

616/5256


Box quantity: 25

Coaxial breathing circuits

Our 3-in-1 coaxial circuit is a circuit where inspiratory and CO₂ sampling lines are enclosed in the larger bore expiratory line. The fresh gases passing through the inner tube are warmed by the expired gases passing through the outer tube.


Main features

- The circuit is less bulky and easier to handle as it permits the elimination of one parallel limb.
- A line dedicated to capnography sampling is included in the outer tubing, thus allowing CO₂ sampling without additional tubing.
- Also available as the coaxial breathing circuit (without CO₂ sampling line)
- The coaxial breathing circuit can also be used with circle absorbent.
- Coaxial breathing circuit can be used in ICU, provided that an HME or a filter/HME is in place at the patient end.



3-in-1 coaxial breathing circuits

285/25408



Code	Description	Sterile	Notes
285/25408	130 cm, with pressure port	✓	
Also available			
285/25417	200 cm, with pressure port	✓	
285P25417	200 cm, with pressure port		
285/25446	260 cm, with pressure port	✓	

Box quantity: 10

Coaxial Breathing Circuits

Code	Description	Sterile	Notes
285/25557	200 cm, with pressure port	✓	

Box quantity: 10


Accessories

Anesthesia face masks

- Standard and inflatable configurations provide all-around solutions for all patients
- Available for neonates to large adults
- Transparent body for easier visualization of vital signs
- Color coded for ease of identification
- Extra soft cushion for better anatomic fit
- Sterile or clean versions available
- Latex free
- Individually packed, in boxes of 10



Standard anesthesia face masks

Large adult size 5 - blue ring	Medium adult size 4 - yellow ring	Small adult size 3 - white ring	Pediatric size 2 - green ring	Infant size 1	Neonate size 0
 22F	 22F	 22F	 22F	 15M	 15M
8800DAR	8801DAR	8802	8803	8804	8805
Box quantity: 10	Box quantity: 10	Box quantity: 10	Box quantity: 10	Box quantity: 10	Box quantity: 10

Inflatable anesthesia face masks

Large adult size 5 - blue ring	Medium adult size 4 - yellow ring	Small adult size 3 - white ring	Large adult size 5 - blue ring sterile	Medium adult size 4 - yellow ring sterile	Small adult size 3 - white ring sterile
 22F	 22F	 22F	 22F	 22F	 22F
8806	8807	8808	8809	8810	8811
Box quantity: 10	Box quantity: 10	Box quantity: 10	Box quantity: 10	Box quantity: 10	Box quantity: 10

Guedel airways

Available in sterile and non-sterile clean pack options.

Single-use Guedel airways

Code	Length cm	Color Code	Sterile
287/7839	3		✓
287/7838	4		✓
287/7636	5		✓
287/7635	6		✓
287/7634	7		✓
287/7633	8		✓
287/7632	9		✓
287/7631	10		✓
287/7843	11		✓
287P7839SP	3		
287P7838SP	4		
287P7636SP	5		
287P7635SP	6		
287P7634SP	7		
287P7633SP	8		
287P7632SP	9		
287P7631SP	10		
287P7843SP	11		
Box quantity: 10			



Intensive care (ICU) adult



6

Anesthesia

24

ICU

34

PICU/NICU

38

Anesthetic and
intensive care tubes

42

Catheter
mounts

49

Anesthetic and intensive
care accessories

Intensive care systems

Medtronic offers a range of breathing systems that meet all your day-to-day requirements for intensive care units.

Circuits include basic systems suitable for use with heat and moisture exchangers as well as several versions of heated and non-heated circuits for use with most available active humidification systems. Circuits are available in adult, pediatric, and neonatal sizes.

We can also provide custom-made breathing systems.


All products are latex free.

Materials include:

- PVC smoothbore (high mechanical resistance and low compliance)
- Extendable polypropylene (minimal storage space required)
- Polyethylene (lightweight and cost-effective)

Adult standard circuits


Polyvinyl chloride (PVC) - smoothbore breathing systems



Code	Description	Sterile	Notes
301/6326	PVC circuit, 150 cm, ports on Y-piece	✓	22 flex ventilator connections
Also available			
301/6378	150 cm, detachable Y-piece	✓	22 flex ventilator connections
301/13404	200 cm	✓	22 flex ventilator connections
301/6130	180 cm, cap on Y-piece	✓	22 flex ventilator connections
301/14105	180 cm, detachable Y-piece	✓	22 flex ventilator connections
Box quantity: 10			

Polyethylene (PE) - corrugated breathing systems

301NP14324




Code	Description	Sterile	Notes
301NP14324	PE breathing system 150 cm, ports on Y-piece, elbow with capno and cap on Y-piece		For multi-patient use

Box quantity: 10

Adult standard circuit with 1 water trap

Polyvinyl chloride (PVC) - smoothbore breathing systems

303/6112



Code	Description	Sterile	Notes
303/6112	PVC Breathing system without ports, on Y-piece inspiratory line 180 cm, expiratory line 180 cm with water trap	✓	22 flex ventilator connections

Also available

304/6486	160 cm, ports on Y-piece, extra limb 60 cm	✓	22 flex ventilator connections
304/13567	210 cm, ports on Y-piece	✓	22 flex ventilator connections

Box quantity: 10

Polypropylene (PP) - extendable breathing systems

303/13609

60/180 cm
22M/15F
22F
45/120 cm
30/70 cm
22F

Code	Description	Sterile	Notes
303/13609	PP extendable breathing system with ports, on Y-piece inspiratory line 60/180 cm, expiratory line 90/240 cm with water trap, extra limb 30/70 cm		22F ventilator connections

Box quantity: 10

Adult standard circuit with 2 water traps

Polyvinyl chloride (PVC) - smoothbore breathing systems

304/13686

80 cm
22M/15F
80 cm
80 cm
22 flex
100 cm
22 flex

Code	Description	Sterile	Notes
304/13686	PVC breathing system, ports on Y-piece, inspiratory and expiratory line 160 cm, extra limb 100 cm	✓	22 flex ventilator connections


Also available

304/13693	160 cm, extra limb 60 cm, detachable Y-piece	✓	22 flex ventilator connections
304/6400	160 cm, cap on Y-piece	✓	22 flex ventilator connections
304S14365	210 cm, detachable water traps, elbow on Y-piece	✓	22 flex ventilator connections
303/6117	160 cm, no ports on Y-piece, detachable water traps, extra limb 40 cm	✓	22 flex ventilator connections
303/6376	200 cm, no ports on Y-piece, extra limb 60 cm, straight connector 22M-22M/19F	✓	22 flex ventilator connections

Box quantity: 10

Polypropylene (PP) - extendable breathing systems

304/13320



Code	Description	Sterile	Notes
304/13320	PP extendable breathing system, ports on Y-piece, 60/140 cm, extra limb 30/70 cm		22F ventilator connections


Also available

304/13614	90/240 cm, extra limb 45/120 cm, 22M-22 flex connectors set		22F ventilator connections
304/13607	60/140 cm, extra limb 45/120 cm, 3 L latex free bag		22F ventilator connections

Box quantity: 10

Polyethylene (PE) - corrugated breathing systems

304/13565



Code	Description	Sterile	Notes
304/13565	PE breathing system, ports and cap on Y piece, 158 cm, extra limb 79 cm, straight connectors 22M-22M/19F		22F ventilator connections

Box quantity: 10


28

Adult insulated smoothbore breathing systems


To maintain a constant temperature of the gas inside the circuit, a range of insulated PVC smoothbore circuits is available in the adult diameter only. In these systems, the external spiral of the tube is coated with a further layer of plastic material which forms an insulating air cavity.

PVC insulated smoothbore tube

1. Smoothbore internal surface
2. Armouring spiral
3. External wall
4. Insulating air space



304/6085



Code	Description	Sterile	Notes
304/6085	PVC thermo insulated breathing system, ports on Y-piece, 160 cm, extra limb 60 cm	✓	22 flex ventilator connections
Also available			
303/5140	160 cm, no ports on Y-piece, 60 cm extra limb	✓	22 flex ventilator connections
Box quantity: 10			

Adult PVC heated circuits for active humidifiers

Whenever an active humidifier is used, large quantities of condensation in the circuit increases the likelihood of contamination. Previously, this problem has been overcome by fitting the circuit with a heating element consisting of an electrical wire that passes inside the circuit itself. The DAR™ PVC breathing circuits with integrated heating element employ an innovative and exclusive solution for overcoming the problem. In addition to the advantages offered by PVC systems, these circuits have the heating element embedded in the outer armouring helix of the limb. Circuits are available for most common active humidifiers.



Main features

- Uniformly heated limb, without “cold spots,” preventing condensation buildup inside the circuit (the main cause of circuit contamination).
- Embedded heating element increases safety, for both operators and patients.
- Available for adult, pediatric, and neonatal
- The embedded wire can be easily connected to the power supply through an adapter cable.

Adult heated circuits used with the Kendall™ Aerodyne 2000 active humidifier

The Aerodyne 2000 active humidifier is an ideal partner for today's ventilators, such as the Puritan Bennett™ 840 ventilator, and uses a closed sterile water system, an effective way to minimize infection risk and increase safety.


Aerodyne 2000 is versatile system that can be used for several application and ventilator modes, and therefore is a cost-effective solution together with DAR™ heated wire circuit range.



Adult PVC heated circuits for Fisher and Paykel™ MR730 active humidifier

Single heated circuits

304/13031




Code	Description	Sterile	Notes
304/13031	Heated breathing system, ports on detachable Y-piece, inspiratory heated line 120 cm, expiratory line 160 cm, extra limb 30 cm	✓	22 flex ventilator connections
Also available			
304/13421	Heated line 120 cm, expiratory line 160 cm, extra limb 60 cm	✓	22 flex ventilator connections
304/13358	Inspiratory heated line 150 cm, expiratory line 160 cm, extra limb 30 cm cap on Y-piece	✓	22 flex ventilator connections

Box quantity: 10

Adult heated circuits for Fisher and Paykel™ MR850 active humidifier

Single heated circuits

304S14292



Code	Description	Sterile	Notes
304S14292	Heated breathing system, ports and cap on detachable Y-piece, inspiratory heated line 120 cm, expiratory line 160 cm, extra limb 30 cm	✓	22 flex ventilator connections
Also available			
304S14291	Inspiratory heated line 150 cm, expiratory line 160 cm, extra limb 60 cm	✓	22 flex ventilator connections

Box quantity: 10

31

Dual heated wire circuit

304S14300

Code	Description	Sterile	Notes
304S14300	Dual heated breathing system, detachable Y-piece with ports and elbow with cap, inspiratory heated line 150 cm, expiratory heated line 180 cm, extra limb 40 cm	✓	22 flex ventilator connections

Box quantity: 10

Adult heated circuits for Kendall™ Aerodyne 2000 active humidifier

Single heated circuits

304S14262


Code	Description	Sterile	Notes
304S14262	Heated breathing system, ports on detachable Y-piece, inspiratory heated line 140 cm, expiratory line 160 cm, extra limb 60 cm	✓	22 flex ventilator connections

Box quantity: 10

32

Dual heated circuit

304S14263



Code	Description	Sterile	Notes
304S14263	Dual heated breathing system, ports on detachable Y-piece, 140 cm, extra limb 80 cm	✓	22 flex ventilator connections

Box quantity: 10

Adapter cables



111/1146



111/1149



111/1130

Code	Description	Sterile	Notes
111/1146	Adapter cable for single heated circuit Fisher and Paykel™ MR850 active humidifier		
111/1149	Adapter cable for dual heated circuit Fisher and Paykel™ MR850 active humidifier		
111/1130	Adapter cable for inspiratory line for Kendall™ Aerodyne 2000 active humidifier		

Box quantity: 1

Intensive care pediatric/neonatal (PICU/NICU)



6

Anesthesia

24

ICU

34

PICU/NICU

38

Anesthetic and
intensive care tubes

42


Catheter
mounts

49

Anesthetic and intensive
care accessories

Pediatric standard circuits

305/8804




Code	Description	Sterile	Notes
305/8804	Pediatric basic circuit, ports on Y-piece, 180 cm	✓	22 flex ventilator connections

Box quantity: 10

Pediatric circuits with water traps

306/8787



Code	Description	Sterile	Notes
306/8787	Pediatric breathing system 160 cm, extra limb 80 cm	✓	22 flex ventilator connections

Also available

306/6948	160 cm, extra limb 60 cm	✓	22 flex ventilator connections
306/6899	300 cm, no ports on Y-piece, no extra limb	✓	22 flex ventilator connections
306/6705	120 cm, one water trap, no extra limb	✓	22 flex ventilator connections

Box quantity: 10

Pediatric heated circuits

Pediatric single heated circuits used with Fisher and Paykel™ MR730 and/or MR850 active humidifiers

306/6843

22 flex 60 cm 22 flex

22M/15F 22F

60 cm 60 cm 22 flex

22M/15F

135 cm

Embedded heated wire

Code	Description	Sterile	Notes
306/6843	Single heated breathing system, ports on detachable Y-piece, inspiratory heated line 135 cm, expiratory line 120 cm, extra limb 60 cm	✓	22 flex ventilator connections, compatible with Fisher and Paykel™ MR730 and MR850 humidifiers

Also available


306S8987	Ports on detachable Y-piece, inspiratory heated line 135 cm, expiratory line 160 cm, extra limb 50 cm, straight connectors	✓	22 flex ventilator connections, compatible with Fisher and Paykel™ MR730 and MR850 humidifiers
----------	--	---	--

Box quantity: 10

36

Neonatal standard circuits

305/6852



Code	Description	Sterile	Notes
305/6852	Neonatal PVC smoothbore circuit 120 cm, ports on Y-piece	✓	22 flex ventilator connections
Also available			
305/8173	150 cm	✓	22 flex ventilator connections
305/6745	180 cm	✓	22 flex ventilator connections
Box quantity: 10			

Neonatal circuits with water traps

Code	Description	Sterile	Notes
307/8448	160 cm, ports and monitoring line 180 cm on Y-piece, extra limb 60 cm	✓	22 flex ventilator connections
307S8957	120 cm, no ports on Y-piece, extra limb 60 cm	✓	22 flex ventilator connections
307/6747	160 cm, ports on Y-piece, extra limb 60 cm	✓	22 flex ventilator connections
Box quantity: 10			

Neonatal single heated circuits

Neonatal single heated circuits used with Fisher and Paykel™* MR730 and MR850 active humidifiers

Code	Description	Sterile	Notes
307/8447	Y-piece with ports and pressure line 180 cm, inspiratory heated line 135 cm with 40 cm additional limb, expiratory line 160 cm, extra limb 60 cm, straight connector set	✓	22 flex ventilator connections
Box quantity: 10			

6

Anesthesia

24

ICU

34

PICU/NICU

38

Anesthetic and
intensive care tubes

42

Catheter
mounts

49

Anesthetic and intensive
care accessories

Anesthetic and intensive care tubes



6

Anesthesia

24

ICU

34

PICU/NICU

38

Anesthetic and
intensive care tubes

42

Catheter
mounts


49

Anesthetic and intensive
care accessories

Adult tubings

Adult PVC tubings

285/5058




Code	Description	Sterile	Notes
285/5058	40 cm	✓	
285/5059	60 cm	✓	
285/5060	80 cm	✓	
285/5061	100 cm	✓	
285/5062	120 cm	✓	
285/5063	150 cm	✓	
285/5064	180 cm	✓	
285/7032	200 cm	✓	

Box quantity: 10

Adult PVC tubings with accessories

285/7515



Code	Description	Sterile	Notes
285/7515	150 cm, 22 flex-22M connections, 3 L latex bag on 22M connector	✓	
Also available			
285/7283	100 cm, 22 flex-22 flex connections, 2 L breathing bag on 22M connector	✓	
285S25640	150 cm, 22 flex-22M connections, 2 L breathing bag on 22M connector	✓	

Box quantity: 10

Adult PP extendable tubings

Code	Description	Sterile	Notes
285/7770	60/180 cm	✓	

Box quantity: 25

6

Anesthesia

24

ICU

34

PICU/NICU

38

Anesthetic and
intensive care tubes

42


Catheter
mounts

49

Anesthetic and intensive
care accessories

Adult PP extendable tubings with accessories

285/7743



2 L

Latex free

22M/19F

45/120 cm


22F

22F

Code	Description	Sterile	Notes
285/7743	45/120 cm, 22F-22F connections, 2 L latex-free breathing bag		

Box quantity: 25

285/7728



Adult water trap

60/180 cm

22F

22F

60/180 cm

Code	Description	Sterile	Notes
285/7728	60/180 cm - adult water trap - 60/180 cm, 22F-22F connections	✓	

Box quantity: 10

Adult PE tubings

290/7264



22F

22F

Code	Description	Sterile	Notes
290/7264	50 mt with cuffs every 39.5 cm		Packed in box of 1 roll (50 mt)
290/5035	50 mt with cuffs every 16.5 cm		Packed in box of 1 roll (50 mt)


Box quantity: 10

40

Pediatric tubings

Pediatric PVC tubings

286/5066



Code	Description	Sterile	Notes
286/5066	60 cm, 22 flex-22 flex connections	✓	
286/5068	100 cm, 22 flex-22 flex connections	✓	
286/5069	120 cm, 22 flex-22 flex connections	✓	
286/5070	150 cm, 22 flex-22 flex connections	✓	
286/5071	180 cm, 22 flex-22 flex connections	✓	
286/7040	180 cm, 15 flex-22 flex connections		
286/7438	200 cm, 22 flex-15 flex connections	✓	
286/25081	300 cm, 15 flex-22 flex connections	✓	

Box quantity: 10

Neonatal tubings

Neonatal PVC tubings

286/25222



Code	Description	Sterile	Notes
286/25222	120 cm, 15M-15F connections		

Box quantity: 10

6

Anesthesia

24

ICU

34

PICU/NICU

38

Anesthetic and
intensive care tubes

42

Catheter
mounts

49

Anesthetic and intensive
care accessories

Catheter mounts



6

Anesthesia

24

ICU

34

PICU/NICU

38

Anesthetic and
intensive care tubes

42

Catheter
mounts

49

Anesthetic and intensive
care accessories

Catheter mounts

Catheter mounts are used to connect a breathing system to an endotracheal or a tracheostomy tube and help provide optimal gas delivery, administer drugs if needed, and perform endotracheal suctioning. They allow adjustments to the breathing systems while keeping the endotracheal tube in place. Our catheter mounts include a wide choice of smoothbore PVC, PP, and PE configurations.


Main features

- Lightweight, atoxic, transparent catheter mounts with compliance less than 2%
- Spiral armouring avoids occlusion under torsion or flexion
- Low dead space
- Smoothbore internal surface prevents fluid buildup reducing risks of bacterial growth
- PP expandable material allows the operator to set and fix the catheter mount in several positions and saves storage space
- Bonded connectors
- 15M or 22F adapters available
- Individually packed, sterile, in boxes of 25

PVC catheter mounts

PVC catheter mounts with double swivel elbow connector

331/5389



10 cm

15M


Double-swivel elbow with cap for suctioning and bronchoscopy

22M/15F

Code	Description	Sterile	Notes
331/5389	10 cm long, 15M adapter – 22M/15F patient connections	✓	
Also available			
331/5390	15 cm long, 15M adapter	✓	
331/5391	10 cm long, 22F adapter	✓	
331/5373	20 cm long, 15M adapter	✓	
331/5385	15 cm long, 22F adapter	✓	
331/5395	20 cm long, 22F adapter	✓	
331/5626	5 cm long, 15M adapter	✓	

Box quantity: 25


43

331/5654

Code	Description	Sterile	Notes
------	-------------	---------	-------

331/5654	15 cm long, 15M adapter - 22M/15F patient connections, with fixed cap	✓	
----------	---	---	--

Box quantity: 25


331/5629

Code	Description	Sterile	Notes
------	-------------	---------	-------

331/5629	20 cm long, 22F adapter - 22M/15F patient connections, with fixed cap	✓	
----------	---	---	--

Box quantity: 25

PVC catheter mounts with 120-degree elbow connector

331/5305

Code	Description	Sterile	Notes
------	-------------	---------	-------

331/5305	10 cm long, 15M adapter - 22M/15F patient connections	✓	
----------	---	---	--

Also available

331/5306	15 cm long, 15M adapter	✓	
----------	-------------------------	---	--

331/5112	10 cm long, 22F adapter	✓	
----------	-------------------------	---	--

331/5113	15 cm long, 22F adapter	✓	
----------	-------------------------	---	--

331/5353	5 cm long, 15M adapter	✓	
----------	------------------------	---	--

Box quantity: 25

6

Anesthesia

24

ICU

34

PICU/NICU

38

Anesthetic and
intensive care tubes

42


Catheter
mounts

49

Anesthetic and intensive
care accessories

PVC catheter mounts with 90-degree elbow connector


330/5307



Code	Description	Sterile	Notes
330/5307	10 cm long, 15M adapter - 22M/15F patient connections	✓	
Also available			
330/5109	10 cm long, 22F adapter	✓	
330/5110	15 cm long, 22F adapter	✓	
330/5308	15 cm long, 15M adapter	✓	
Box quantity: 25			

PVC catheter mounts with CO₂ sampling port on 90-degree elbow connector

330/5321




Code	Description	Sterile	Notes
330/5321	10 cm long, 15M adapter - 22M/15F patient connections	✓	
Also available			
330/5322	15 cm long, 15M adapter	✓	
330/5322TC	15 cm long, 15M adapter tethered cap	✓	
330/5645	10 cm long, 22F adapter	✓	
330/5646	15 cm long, 22F adapter	✓	



Box quantity: 25

PVC straight catheter mounts


332/5114



Code	Description	Sterile	Notes
332/5114	10 cm long, 15M-15F	✓	
Also available			
332/5115	15 cm long, 15M-15F	✓	
332/5309	10 cm long, 22F-15F	✓	
332/5310	15 cm long, 22F-15F	✓	
332/5316	5 cm long, 22F-15F	✓	
332/5317	5 cm long, 15M-15F	✓	
332/5332	20 cm long, 15M-15F	✓	

Box quantity: 25

332S5685



Code	Description	Sterile	Notes
332S5685	20 cm long, 15M adapter-22M/15F patient connections	✓	
Also available			
332/5339	30 cm long, 15M adapter	✓	

Box quantity: 25

PVC Carlens catheter mounts

This catheter mount with double-channeled tube provides independent lung ventilation. It also enables the Carlens or White connection endobronchial tubes, according to specific need, through an obturator in one of two ways. The obturator is included in the packaging.

6

Anesthesia

24

ICU

34

PICU/NICU

38

Anesthetic and
intensive care tubes

42

Catheter
mounts

49


Anesthetic and intensive
care accessories

PP catheter mounts (extendable)

Extendable PP catheter mounts with double-swivel elbow connector

331/5661

Double-swivel elbow with cap for suctioning and bronchoscopy




Code	Description	Sterile	Notes
331/5661	7/16 cm long, 15M adapter-22M/15F patient connections	✓	
Also available			
331/5665	7/16 cm long, 15M adapter, polyethylene cap on connector	✓	
331/5669	7/16 cm long, 22F adapter	✓	

Box quantity: 25

Extendable PP catheter mounts with 120-degree elbow connector

331/5662

Cap for suctioning




Code	Description	Sterile	Notes
331/5662	7/16 cm long, 15M adapter-22M/15F patient connections	✓	
Also available			
331/5670	7/16 cm long, 22F adapter	✓	

Box quantity: 25

Extendable PP catheter mounts with 90-degree elbow connector

330/5337



Code	Description	Sterile	Notes
330/5337	7/16 cm long, 15M adapter - 22M/15F patient connections	✓	
Also available			
330/5678	7/16 cm long, 22F adapter	✓	

Box quantity: 25

6

Anesthesia

24

ICU

34

PICU/NICU

38

Anesthetic and
intensive care tubes

42


Catheter
mounts

49

Anesthetic and intensive
care accessories

Extendable PP catheter mounts with CO₂ sampling port on 90-degree elbow connector

330/5664



Code	Description	Sterile	Notes
330/5664	7/16 cm long, 15M adapter - 22M/15F patient connections	✓	
Also available			
330/5667	7/16 cm long, 22F adapter	✓	

Box quantity: 25

Extendable PP straight catheter mounts

332/5663



Code	Description	Sterile	Notes
332/5663	7/16 cm long, 15M-15F	✓	
Also available			
332/5666	7/16 cm long, 22F-15F	✓	

Box quantity: 25

48

Anesthetic and intensive care accessories



6

Anesthesia

24

ICU

34

PICU/NICU

38

Anesthetic and
intensive care tubes

42

Catheter
mounts


49

Anesthetic and intensive
care accessories

Connectors

Angled patient connectors

Double swivel, port for suctioning and bronchoscopy with rubber cap, sterile




610/5257

Box quantity : 10

Y-Pieces


Adult, sterile



600/5087

Box quantity: 10

Adult, with pressure and temperature ports, sterile




601/5089

Box quantity: 10

Straight, short connectors


Sterile



608/5261

Box quantity: 10


Sterile



609/5244

Box quantity: 10

Sterile



609/5248

Box quantity: 10


Sterile



609/9262

Box quantity: 10


Sterile



609/5252

Box quantity: 10

Sterile




609/5176

Box quantity: 10

Straight, long connectors

With pressure port, sterile



608/5093A

Box quantity: 10


With metered-dose inhaler adapter, nonsterile



608/5105

Box quantity: 10


Sterile



608/5094

Box quantity: 10


With pressure port, sterile



608/5095A

Box quantity: 10

Sterile



608/5092

Box quantity: 10

6

Anesthesia

24

ICU

34

PICU/NICU

38

Anesthetic and intensive care tubes

42


Catheter mounts

49

Anesthetic and intensive care accessories

L adapters


Sterile



620/5173

Box quantity: 10


With male luer lock, sterile



620/9019

Box quantity: 10

With male luer lock, sterile



620/9019TC

Box quantity: 10


**For Fisher and Paykel™
MR730, sterile**



620S9286

Box quantity: 10


Flex, sterile



620/9121

Box quantity: 10


With APL valve, sterile



616/5256

Box quantity: 25


Nonsterile, autoclavable



620/5183

Box quantity: 2 bags × 25 pcs


Nonsterile, autoclavable



620/5217

Box quantity: 2 bags × 25 pcs

**With temperature port,
nonsterile, autoclavable**



620/5180B

Box quantity: 2 bags × 25 pcs

6

Anesthesia

24

ICU

34

PICU/NICU

38

Anesthetic and
intensive care tubes

42


Catheter
mounts

49

Anesthetic and intensive
care accessories

T adapters

Sterile




22F 22M 22M

615/5100

Box quantity: 10

Sterile




22M 22F 22M/15F

615/5245

Box quantity: 10

Sterile




22F 22M/15F 22M

615/5254

Box quantity: 10

Sterile




22M 22M 22M/15F

615/5103

Box quantity: 10

With pressure port, sterile




22M 22M 22M/15F

615/5103A

Box quantity: 10

With temperature port, sterile




22M 22M 22M/15F

615/5103B

Box quantity: 10

Sterile




22M/15F 22F 22M

615/9234

Box quantity: 10

Caps

For 15F connections, sterile	Sterile	Sterile	Bulk packed, nonsterile
			
608/9196	608/9068	608/9138	608B9099
Box quantity: 10	Box quantity: 10	Box quantity: 10	Box quantity: 300

Water traps

Adult water trap with 22M-22M connections, sterile	PEEP valve with 22M connection, nonsterile	Suction control valve, sterile
		
652/5406	609/9034	291/7617
Box quantity: 50	Box quantity: 10	Box quantity: 10

Mouthpieces

Mouthpiece on T-piece with 22M connections, nonsterile	Flex mouthpiece with 33F round shaped connection, nonsterile
	
615/5288	500P30002
Box quantity: 10	Box quantity: 25

Breathing bags and accessories

A complete range of reservoir breathing bags available in four sizes, all meet ISO standards

- 22mm I.D. connector
- Single patient use
- Latex free
- Individually packed, nonsterile, in boxes of 10

With 22 female flexible connector

3 L, nonsterile	2 L, nonsterile	1 L, nonsterile	0.5 L, nonsterile
			
10911	10910	10909	10908
Box quantity: 10	Box quantity: 10	Box quantity: 10	Box quantity: 10

With 22 female rigid connector

3 L, nonsterile	2 L, nonsterile	1 L, nonsterile	0.5 L, nonsterile
			
10907	10906	10905	10904
Box quantity: 10	Box quantity: 10	Box quantity: 10	Box quantity: 10

Accessories

Circuit holder, nonsterile	Bag hook, nonsterile
	
111P1114	111/1153
Box quantity: 50	Box quantity: 25

Gas sampling lines

Code	Description	Sterile	Notes
291/7565	PVC line 200 cm, Ø 1,2 × 2,5 mm, male luer lock connections	✓	
291/7731	PVC line 300 cm, Ø 1,2 × 2,5 mm, male luer lock connections	✓	
291/7733	PVC line 300 cm, Ø 1,2 × 2,5 mm, male-female luer lock connections	✓	
Box quantity: 25			
291/7801	PVC line 200 cm, Ø 1,2 × 2,5 mm, male luer lock-universal connector	✓	
Box quantity: 50			

PVC pressure monitoring lines

Code	Description	Tube diameter	Sterile	Notes
291/25539	200 cm with angled-male luer lock connections	5 × 7,5		
291/7365	200 cm	6,36 × 9,54	✓	
Box quantity: 10				

IMPORTANT: Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

©2024 Medtronic. Medtronic, Medtronic logo, and Engineering the extraordinary are trademarks of Medtronic. TM* Third-party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.
emea-re-2400063-dar-patient-circuits-catalogue-en-gb-emea-13916595

[medtronic.co.uk/acm](https://www.medtronic.co.uk/acm)

Medtronic