

Custom Fit Measuring Form

For quotes please email: quotes@hc21.group
To order, please email Customer Services on: info@hc21.group

Date: _____ Patient Name: _____

Measured By: _____ Tel: _____ Email: _____

Delivery Address: _____

Invoice Address: _____



☐ JOBST® Confidence	
Compression Class (RAL) ☐ CCL 1 (18-21mmHg) ☐ CCL 2 (23-32mmHg) ☐ CCL 3 (34-46mmHg)	Quantity Left Right
Style and Options ☐ AD Knee high ☐ Slipform ☐ Ankle functional zone ☐ Hallux valgus (bunion relief zone)* ☐ Decor-line ☐ Initials (max 2 letters)	
☐ AG Thigh high ☒ Slipform ☐ Ankle functional zone ☐ Knee functional zone ☐ Decor-line ☐ Initials (max 2 letters) ☐ Hallux valgus (bunion relief zone)*	
Colour ☐ Beige ☐ Black ☐ Caramel ☐ Red heather ☐ Anthracite heather ☐ Jeans heather Note: the functional zone for red, anthracite and jeans heather coloured garments is more visible	Top band options Knee high A-D ☐ 5cm SoftFit Thigh high A-G ☒ 5cm silicone top band
Remarks *Only available with slant foot.	

Please note: Compression garments carry a 6-month warranty – based on 2 pairs worn alternately. Garments should be replaced annually to maintain optimal treatment and fit.

Please note: All Measurements in cms.

Left **Right**

Length **Circumference** **Circumference** **Length**

a-K _____
Required for thigh high garment

a-G _____ G _____ a-G _____

a-F _____ F _____ a-F _____

a-E _____ E _____ a-E _____

a-D _____ D _____ a-D _____
(Knee high)

a-C _____ C _____ a-C _____

a-B₁ _____ B₁ _____ a-B₁ _____

a-B _____ B _____ a-B _____

_____ Y _____

_____ A _____

a*

*a = floor / base of measuring board

Foot Length - Left	Foot Length - Right
For open toe _____ cm	For open toe _____ cm
For closed toe _____ cm (longest toe)	For closed toe _____ cm (longest toe)

The diagram illustrates the anatomical planes for the foot. On the left, a foot is shown from a lateral perspective, with a vertical double-headed arrow labeled 'lateral' indicating the side view, and a vertical double-headed arrow labeled 'medial' indicating the inner view. On the right, a foot is shown from a medial perspective, with a vertical double-headed arrow labeled 'medial' indicating the inner view, and a vertical double-headed arrow labeled 'lateral' indicating the side view.

NOTE: Please turn page for garments options and application video

JOBST® Confidence Knee High

Compression Class / Option



Basic style

CCL 1 (18-21mmHg)
CCL 2 (23-32mmHg)
CCL 3 (34-46mmHg)

Options

Closed toe
SoftFit
Slipform
Ankle functional zone

JOBST® Confidence Thigh High

Compression Class / Option



Basic style

CCL 1 (18-21mmHg)
CCL 2 (23-32mmHg)
CCL 3 (34-46mmHg)

Options

Closed toe
5cm silicone band
Slipform
Ankle functional zone
Knee functional zone



Application of JOBST®
Confidence video



Please note: Compression garments carry a 6-month warranty – based on 2 pairs worn alternately. Garments should be replaced annually to maintain optimal treatment and fit.

Options



Ankle Functional Zone

Circumferential functional zone allows greater movement and flexibility at the ankle



Knee Functional Zone

Circumferential functional zone allows greater movement and flexibility at the knee



Slipform

Slipform provides a higher shape to the outside of the garment, at the top, conforming to the limb shape



Deco Line

A decorative seam line can be added to the front of the garment on request



Hallux Vallgus (bunion relief zone)*

Creates a looser knitting pattern around the ball of the foot to provide relief to painful bunions



SoftFit Top Band

Soft and secure top band that has been uniquely knitted into the garment top band for a more comfortable secure fit



Initials

Patient initials can be knitted into the garment below the JOBST logo for a personal touch