

JOBST® ELVAREX® PLUS

Foot Cap Measuring Form


For quotes please email: quotes@hc21.group
To order, please email Customer Services on: info@hc21.group



Date: _____ Patient Name: _____
Measured By: _____ Tel: _____ Email: _____
Delivery Address: _____
Invoice Address: _____



PLEASE NOTE: All measurements in cms



Note: Circumferential measurements are required for all 5 toes however length measurement not required for little toe unless 5th toe required


Note: All foot caps will be manufactured with **4 toes** unless noted in this remarks box that 5 toes are required.

Left Foot:

- ☐ 5 Toes
☐ Length of little toe

Right Foot:

- ☐ 5 Toes
☐ Length of little toe



JOBST Elvarex Plus

Compression Class:

- ☐ CCL 1 (18-21mmHg)
☐ CCL 2 (23-32mmHg)
☐ CCL 3 (34-46mmHg)*

* Note: this will be JOBST Elvarex (seamed)

Colour:

- ☐ Beige ☐ Black
☐ Dark blue ☐ Grey
☐ Dark brown
☐ Cranberry

Style:

- ☒ AA Footcap
☒ 4 Toes

Quantity:

Left
Right

Options:

- ☒ 1.5cm slanted little toe opening †

Remarks:

† This is standard. If different size is required, please note in remarks

JOBST Elvarex Soft Seamless

Compression Class:

- ☐ CCL 1 (18-21mmHg)
☐ CCL 2 (23-32mmHg)

Colour:

- ☐ Beige
☐ Black

Style:

- ☒ AA Footcap
☒ 4 Toes

Quantity:


Left
Right

Options: ☒ 1.5cm slanted little toe opening

JOBST Elvarex Plus / JOBST Elvarex Soft Seamless Custom-Fit Information



JOBST Elvarex Plus Foot Cap with Toes



Compression Class/Option

Basic style – foot cap with toes in beige*

CCL 1 – 18-21mmHg


CCL 2 – 23-32mmHg

Options

Non-standard colour

*1.5cm slanted little toe

JOBST Elvarex Foot Cap with Toes



Basic style – foot cap with toes in beige*


CCL 3 – 34-46mmHg

Options

Non-standard colour

*1.5cm slanted little toe

JOBST Elvarex Soft Seamless Foot Cap with Toes



Basic style – foot cap with toes in beige*

CCL 1 – 18-21mmHg

CCL 2 – 23-32mmHg

Options

Non-standard colour

*1.5cm slanted little toe opening



Please note: Compression garments carry a 6-month warranty – based on 2 pairs worn alternately. Garments should be replaced annually to maintain optimal treatment and fit.

*The HC21 Group is fully compliant with the GDPR and applies the appropriate security to the processing of all personal data. By completing this order form, you are confirming that you are aware of your obligation to obtain consent from the patient on the processing of their data for the production of their JOBST® compression garment in accordance with the Data Protection Act 2018, and the EU General Data Protection Regulation [GDPR]. For more information about our Privacy Policy, please contact gdp@hc21.group

