

SEE CLEARLY. EVERY TIME.

The next generation McGRATH™ MAC video laryngoscope** features enhanced optics and is durable for routine use.

CLINICIAN INSPIRED

Based on your feedback, we designed the McGRATH™ MAC video laryngoscope because your first attempt should be your best.

SEE THE DIFFERENCE

The McGRATH™ MAC video laryngoscope provides a simple and convenient solution for intubation. The latest design offers:

- Enhanced optics
- Increased durability
- Intelligent battery management

To learn more about the next generation McGRATH™ MAC video laryngoscope, visit

medtronic.co.uk/mcgrath-mac-next-gen



THE EVOLUTION CONTINUES

2×
light spread*

3×
brighter*

4×
resolution*

Next Generation
McGRATH™ MAC
video laryngoscope



McGRATH™
MAC video
laryngoscope



ENHANCED DESIGN

LEGACY DESIGN

Visualization

Light color	Natural, warmer light	Bright white light
Light brightness	3× brighter light source ¹	Standard light source
White balance	Automatically adjusts	Fixed
Camera detail	Enhanced resolution	Standard resolution

Durability

McGRATH™ MAC Handle	Reinforced screen hinge	Standard screen hinge
IP rating	IPX8	IPX7

Battery management

Auto-off feature	Manual power on/off feature
------------------	-----------------------------

Ordering information

SKU number	301-000-000	300-000-000
------------	-------------	-------------

1. Based on internal test report #RE001500188, V&V plan.

* As compared to the previous version of the McGRATH™ MAC video laryngoscope

** May not be available for sale in all markets

Medtronic

IMPORTANT: Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

© 2020 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic.

™** Third party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.
19-weu-mcgrath-new-gen-comparison-3781099

medtronic.co.uk/rms

HC21

csuk@hc21.group
01256 306506
healthcare21.eu