


Niti-S Biliary Covered Stent [Giobor]

for malignant biliary strictures



Features

- **Partially covered design**
 - (1) The long covered portion (70%)
 - (2) The bare portion (30%)
- **Anti-migration flared end**
- **Radiopaque markers:** Three (3) at both ends & Two (2) at the end of the covered part

Ordering Information

| Endoscopic | | | | |
|------------|---------------|-------------|--------------|--------------------|
| Code | Stent | | Delivery | |
| | Diameter (mm) | Length (cm) | Profile (Fr) | Usable Length (cm) |
| BS0806FP | 8 | 6 | 8.5 | 180 |
| BS0808FP | | 8 | | |
| BS0810FP | | 10 | | |
| BS0812FP | | 12 | | |
| BS1006FP | 10 | 6 | 8.5 | 180 |
| BS1008FP | | 8 | | |
| BS1010FP | | 10 | | |
| BS1012FP | | 12 | | |

Aquilant
A HCB Company

info@hc21.group
00 353 214 875 055
healthcare21.eu