


Niti-S Biliary Uncovered Stent [S-Type]

for malignant biliary strictures



Features

- **Fixed cell with braided construction**
 - Flexible and resistant to fracture
- **Radiopaque marker: 3 (three) at both ends & 2 (two) in the middle**

Ordering Information

Percutaneous				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
T0604	6	4	7	50
T0605		5		
T0606		6		
T0607		7		
T0608		8		
T0609		9		
T0610		10		
T0612		12		
T0804	8	4	7	50
T0805		5		
T0806		6		
T0807		7		
T0808		8		
T0809		9		
T0810		10		
T0812		12		
T1004	10	4	7	50
T1005		5		
T1006		6		
T1007		7		
T1008		8		
T1009		9		
T1010		10		
T1012		12		

Endoscopic				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
B0604	6	4	7	180
B0605		5		
B0606		6		
B0607		7		
B0608		8		
B0609		9		
B0610		10		
B0612		12		
B0804	8	4	7	180
B0805		5		
B0806		6		
B0807		7		
B0808		8		
B0809		9		
B0810		10		
B0812		12		
B1004	10	4	7	180
B1005		5		
B1006		6		
B1007		7		
B1008		8		
B1009		9		
B1010		10		
B1012		12		

Endoscopic, Short-wire				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
BM0604	6	4	8.5	180
BM0605		5		
BM0606		6		
BM0607		7		
BM0608		8		
BM0609		9		
BM0610		10		
BM0612		12		
BM0804	8	4	8.5	180
BM0805		5		
BM0806		6		
BM0807		7		
BM0808		8		
BM0809		9		
BM0810		10		
BM0812		12		

Endoscopic, Short-wire				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
BM1004	10	4	8.5	180
BM1005		5		
BM1006		6		
BM1007		7		
BM1008		8		
BM1009		9		
BM1010		10		
BM1012		12		

Aquilant
A HCGI Company

csuk@hc21.group
01256 306506
healthcare21.eu