

Cocoon Patient Warming Blankets and Gowns

One solution for pre, peri and
post operative



Pre Operative Patient Warming Gowns



Intraoperative & Underbody Blankets



Bariatric & Extended Blankets



Neonatal & Paediatric Underbody Blankets



Neonatal & Paediatric Overbody Blankets

Speciality Blankets



Cardiac Blankets



Cocoon Convective Warming Machine CWS 4000/CWS 5000

Cocoon Warming Blankets Codes



Pre Operative Patient Warming Gowns

Latex Free




Easily convert to use as a body blanket during and post procedure

To meet the new NICE CG65 guidelines Cocoon have developed a range of pre-operative patient warming gowns that easily converts for use during and post procedure negating the need to open further blankets.

Keeping a patient warm before, during and after surgery is a cornerstone of patient care. Patient temperature management is not just a question of keeping a patient comfortable while they are in the hospital. If not properly managed, a patient's core temperature may drop below normal levels, which can lead to 'inadvertent perioperative hypothermia' (IPH).


Model - **CLM 096**
Pre Warming Gown XL
Size 135 x 100cm
Hose ports - TBC

- Comfort and warmth for the patient before operation
- Can be converted into adult lower body blanket after prewarming
- Can be used for post-operative care




Model - **CLM 097**
Pre Warming Gown L
Size 100 x 135cm
Hose ports - TBC

- Comfort and warmth for the patient before operation
- Can be converted into adult lower body blanket after prewarming
- Can be used for post-operative care



Model - **CLM 098**
Pre Warming Gown M
Size 130 x 92cm
Hose ports - TBC

- Comfort and warmth for the patient before operation
- Can be converted into adult lower body blanket after prewarming
- Can be used for post-operative care



Model - **CLM 099**
Pre Warming Gown S
Size 98 x 60cm
Hose ports - TBC


- Comfort and warmth for the patient before operation
- Can be converted into paediatric blanket after pre-warming
- Can be used for post-operative care



Intraoperative & Underbody Blankets

Latex Free



Model - CLM 0101 Full Body Size 100 x 195cm Hose ports - 1	• Ample dimensions to accommodate all patients • Universal inlet port ties ensure the hose remains firmly in the blanket	 FULL BODY 100 x 195cm CLM 0101
Model - CLM 0102 Lower Body Size 100 x 120cm Hose ports - 1	• Generous amount of tape to secure blanket over abdomen • Universal inlet port ties ensure the hose remains firmly in the blanket	 LOWER BODY 100 x 120cm CLM 0102
Model - CLM 0103 Upper Body-Arms in Size 100 x 124cm Hose ports - 2	• Generous amount of tape to secure blanket over abdomen • Two hose inlets for operator convenience • Universal inlet port ties ensure the hose remains firmly in the blanket	 UPPER BODY, ARMS IN 100 x 124cm CLM 0103
Model - CLM 0104 Upper Body-Arms out Size 57 x 196cm Hose ports - 2	• Two hose inlets for operator convenience • Arm board ties to secure blanket once in position • Clear face shield to retain warmth around the head • Universal inlet port ties ensure the hose remains firmly in the blanket • Ample tape (latex free) to secure blanket to torso	 UPPER BODY, ARMS OUT 57 x 196cm CLM 0104
Model - CLM 0127 Adult Full Access Under Body Size 206 x 87cm Hose ports - 4	• Under Body, full, Unrestricted surgical access to the patient • Fluid drainage holes in the blanket allow fluids to be drained out of the blanket and not pool on the blanket surface • Arms slits for easier patient positioning • 4 universal inlet ports at each corner allow flexibility for placement of hose insertion	 ADULT FULL ACCESS UNDER BODY BLANKET 206 x 87cm CLM 0127
Model - CLM 0138 Large Adult Full Access Under Body Size 235 x 200cm Hose ports - 4	• Under Body, full, Unrestricted surgical access to the patient • Fluid drainage holes in the blanket allow fluids to be drained out of the blanket and not pool on the blanket surface • Arms slits for easier patient positioning • 4 universal inlet ports at each corner allow flexibility for placement of hose insertion	 Large Adult Full Access Under Body Size 235 x 200cm CLM 0138
Model - CLM 0122 Adult Tube Size 100 x 217cm Hose ports - 4	• Full surgical access as patient lies between two tubular sections • Warm air is directed both across and down the length of the patient • Universal inlet port ties ensure the hose remains firmly in the blanket	 ADULT, TUBE 100 x 217cm CLM 0122

Bariatric & Extended Blankets

Latex Free



Cocoon patient warming provides effective temperature management with our full range of warming blankets. Our unique microscopic perforations provide for even air distribution without the risk of residue.


Model - **CLM 0126**
Bariatric Half Body
Size 116 x 150cm
Hose ports - 1

- New larger dimensions to ensure coverage on all patients
- Universal inlet port ties ensure the hose remains firmly in the blanket
- Larger size blanket to suite bariatric patients




Model - **CLM 0128**
Full Torso Bariatric Blanket
Size 147 x 195cm
Hose ports - 1

- Ample Dimensions to accommodate all patients
- Universal inlet port ties ensure the hose remains firmly in the blanket




Model - **CLM 0129**
Extended Lower Body Blanket
Size 98 x 160cm
Hose ports - 1

- Generous amount of tape to secure blanket over abdomen
- Universal inlet port ties ensure the hose remains firmly in the blanket




Model - **CLM 0130**
Extended Full Torso Blanket
Size 100 x 200cm
Hose ports - 1

- Ample Dimensions to accommodate all patients
- Universal inlet port ties ensure the hose remains firmly in the blanket
- More inflated area for extra warmth to the patient




Model - **CLM 0131**
Extended Upper Torso - Arms In Blanket
Size 137.5 x 97cm
Hose ports - 2

- Generous amount of tape to secure blank
- Two hose inlets for operator convenience
- Universal inlet port ties ensure the hose remains firmly in the blanket



Model - **CLM 0132**
Extended Upper Torso - Arms Out Blanket
Size 76 x 197.5cm
Hose ports - 2

- Two hose inlets for operator convenience
- Arm board ties to secure blanket once in position
- Clear face shield to retain warmth around the head
- Universal inlet port ties ensure the hose remains firmly in the blanket



Neonatal & Paediatric Underbody Blankets

Latex Free



Cocoon patient warming provides effective temperature management for the infants and children in any medical environment. The microscopic perforations keep the infants and children safe.


Model - **CLM 0105**

Paediatric - Neonatal

Size 87 x 89 cm

Hose ports - 2

- Full surgical access as patient lies on top of quilted area
- Two hose inlets for operator convenience
- Universal inlet port ties ensure the hose remains firmly in the blanket
- 2 sheets of clear plastic (for additional warmth, if required)



Model - **CLM 0107**

Paediatric - Child

Size 130 x 92cm

Hose ports - 2

- Full surgical access as patient lies on top of quilted area
- Two hose inlets for operator convenience
- Universal inlet port ties ensure the hose remains firmly in the blanket
- 2 sheets of clear plastic (for additional warmth, if required)




Model - **CLM 0115**

Large Paediatric Under Body

Size 100 x 120cmcm

Hose ports - 1

- Full surgical access as patient lies on top of quilted area
- Universal inlet port ties ensure the hose remains firmly in the blanket



Model - **CLM 0116**

Small Paediatric Under Body

Size 59 x 87cm

Hose ports - 1

- Full surgical access as patient lies on top of quilted area
- Universal inlet port ties ensure the hose remains firmly in the blanket



Neonatal & Paediatric Overbody Blankets

Latex Free



Cocoon patient warming provides effective temperature management for the infants and children in any medical environment. The microscopic perforations keep the infants and children safe.

Model - **CLM 0117**
Paediatric Lower Body
 Size 59 x 87cm
 Hose ports - 1

- Over body blanket
- Universal inlet port ties ensure the hose remains firmly in the blanket



Model - **CLM0109**
Paediatric Full Body
 Size 120 x 100cm
 Hose ports - 1

- Ample dimensions to accommodate all paediatric patients
- Universal inlet port ties ensure the hose remains firmly in the blanket




Model - **CLM 0106**
Paediatric - Tube
 Size 117 x 78cm
 Hose ports - 1

- Full surgical access as patient lies between two tubular sections
- Warm air is directed both across and down the length of the patient
- 1 sheet of clear plastic (for additional warmth, if required)
- Universal inlet port ties ensure the hose remains firmly in the blanket



Model - **CLM 0112**
Multi Purpose Air Sock
 Size 100 x 12cm
 Hose ports - 1

- Full surgical access as patient lies next to the tubular Air Sock
- By placing a covering sheet over the Air Sock warm air is retained and directed both across and down the length of the patient
- Universal inlet port ties ensure the hose remains firmly in the blanket



Speciality Blankets


Latex Free



Cocoon provides speciality blankets to suit every procedure.

Model - **CLM 0118**
Abdominal Access
Size 100 x 195cm
Hose ports - 3


- Cut out section over the chest and abdomen allowing for easy access
- Ample dimensions to accommodate all patients and positions
- Three hose inlets for operator convenience
- Arm board ties to secure blanket once in position
- Ample tape (latex free) around fenestration to secure blanket to torso
- Universal inlet port ties ensure the hose remains firmly in the blanket



ABDOMINAL ACCESS
190 x 180cm
CLM 0118

Model - **CLM 0124**
Surgical Access
- Arms In
Size 100 x 195cm
Hose ports - 1


- Larger cut out section over the chest and abdomen allowing for easy access
- Ample dimensions to accommodate all patients and positions
- Ample tape (latex free) around fenestration to secure blanket to torso
- Universal inlet port ties ensure the hose remains firmly in the blanket



SURGICAL ACCESS
100 x 195cm
CLM0124

Model - **CLM 0125**
Lithotomy and
Laparoscopy
Size 190 x 180cm
Hose ports - 2


- Larger cut out section over the chest and abdomen allowing for easy access
- Ample dimensions to accommodate all patients and positions
- Ample tape (latex free) around fenestration to secure blanket to torso
- Universal inlet port ties ensure the hose remains firmly in the blanket



LITHOTOMY & LAPAROSCOPY
190 x 180cm
CLM 0125

Model - **CLM0133**
Under Body, Spinal
Blanket
Size 83 x 184cm
Hose ports -

- For use in spinal surgery
- Provides unrestricted full patient access
- Drainage hole provides drain fluid on the surface of the blanket
- Universal tie system provides easy fixing of the blanket the table
- Providing comfort and warmth to patient



Under Body, Spinal Blanket
83 x 184cm
1 Hose Entry Port
1 Face cover, 12 Ties
CLM 0133

Cardiac Blankets


Latex Free



Cocoon provides speciality blankets to suit every procedure.

Model - **CLM 0108**
Cardiac Access
Size 100 x 195cm
Hose ports - 1


- 40cm slit from the neck area down the length of the blanket over the chest and abdomen allowing for easy access
- Universal inlet port ties ensure the hose remains firmly in the blanket



CARDIAC ACCESS
100 x 195cm
CLM 0108

Model - **CLM 0121**
Multi Access
Size 100 x 195cm
Hose ports - 3

- Flexible design to allow surgical access by folding back sections of the blanket
- Universal inlet port ties ensure the hose remains firmly in the blanket



MULTI ACCESS
100 x 195cm
CLM 0121

Model - **CLM0202**
Sterile Cardiac Blanket
Size 100 x 120cm
Hose ports - 1

- Lower Body, Sterile Cardiac Blanket
- Generous amount of adhesive tape to secure blanket
- Can be repositioned and will still adhere
- Universal inlet port ties ensure the hose remains firmly in the blanket



STERILE CARDIAC BLANKET
100 x 120cm
CLM 0202

Model - **CLM 0138**
Large Adult Full Access Under Body
Size 235 x 200cm
Hose ports - 4

- Under Body, full, Unrestricted surgical access to the patient
- Fluid drainage holes in the blanket allow fluids to be drained out of the blanket and not pool on the blanket surface
- Arms slits for easier patient positioning
- 4 universal inlet ports at each corner allow flexibility for placement of hose insertion



Large Adult Full Access Under Body
Size 235 x 200cm
CLM 0138

Cocoon Convective Warming Machine CWS 4000/CWS 5000

10 Key Points

- Dual, self-resetting, 53 degree celcius thermostats assure patient safety
- Heating element incorporating thermal switch and thermal fuse provides system protection in case of air flow failure
- Hepa filter ensures clean air delivery to patient
- Thoughtfully designed air intake does not draw air from the floor area maximising filter life



- Filter status indicator takes the guesswork out of filter maintenance
- Microprocessor control delivers excellent safety and performance
- Built-in computer communication port facilitates easy testing and service
- Retractable metal bed rails offer a strong, convenient, fold-away bed mounting
- Flexible mounting options of IV pole or stand to suit your needs
- Protective hose cover delivers thermal insulation and ease of cleaning

Cocoon Warming Blankets Codes

Product Code	NHS SC Code *	Product Description	Size	UOI
CLM096	VCB745	Patient PRE WARM Gown XL with Extra Fabric	100 x 135cm	10
CLM097	VCB746	Patient PRE WARM Gown L	100 x 135cm	10
CLM098	VCB747	Patient RE WARM Gown M	92 x 130cm	10
CLM099	VCB748	Patient PRE WARM Gown S	60 x 98cm	10
CLM0101	VCB749	Adult Full Body Blanket	100 x 195cm	30
CLM0102	VCB750	Lower Body Blanket	100 x 120cm	30
CLM0103	VCB751	Upper Body-Arms In Blanket	100 x 124cm	30
CLM0104	VCB752	Upper Body-Arms Out Blanket	57 x 196cm	30
CLM0105	VCB738	Paediatric-Neonatal Underbody Blanket	87 x 89cm	10
CLM0106	VCB739	Paediatric-Tube Underbody Blanket	117 x 78cm	10
CLM0107	VCB740	Paediatric-Child Underbody Blanket	130 x 92cm	10
CLM0108	VCB755	Cardiac Access Blanket	100 x 195cm	10
CLM0109	VCB761	Paediatric Full Body Blanket	120 x 100cm	10
CLM0112	VCB753	Multi Purpose Air Sock	100 x 12cm	30
CLM0115	VCB741	Large Paediatric Underbody Surgical Access Blanket	100 x 120cm	10
CLM0116	VCB742	Small Paediatric Underbody Surgical Access Blanket	59 x 87cm	10
CLM0117	VCB743	Paediatric Over Body Blanket	59 x 87cm	10
CLM0118	VCB756	Abdominal Access Blanket	190 x 180cm	10
CLM0121	VCB757	Multi Access Blanket	100 x 195cm	10
CLM0122	VCB737	Adult Full Access Underbody Tube	100 x 217cm	10
CLM0124	VCB758	Surgical Access Blanket	100 x 195cm	10
CLM0125	VCB759	Lithotomy and Laparoscopy	190 x 180cm	10
CLM0126	VCB762	Bariatric Lower Body Blanket	150 x 116cm	10
CLM0127	VCB760	Adult Underbody Full Surgical Access	206 x 87cm	10
CLM0128	VCB763	Full Body Bariatric Blanket	147 x 195cm	10
CLM0129	VCB764	Extended Lower Body Blannket	98 x 160cm	10
CLM0130	VCB765	Extended Full Torso Blannket	100 x 200cm	10
CLM0131	VCB766	Extended Upper Torso Arms-in Blannket	137.5 x 7cm	10
CLM0132	VCB767	Extended Upper Torso Arms-out Blannket	76 x 97.5cm	10
CLM0133	VCB768	Underbody Spinal Blanket	83 x 184 cm	10
CLM0138	VCB770	Large Adult Full Access Underbody Blanket	235 x 200cm	10
CLM0202	VCB772	STERILE Cardiac Access Blanket	100 x 120cm	10

*NHS Codes only applicable to UK



csuk@hc21.group
01256 306506
healthcare21.eu