

ordering information

One of the most complete coil offers on the market

Prestige Plus Complex coils	
Reference	Coil sizes
PRES0202CXPPL	2mm x 2cm
PRES0204CXPPL	2mm x 4cm
PRES0206CXPPL	2mm x 6cm
PRES0307CXPPL	3mm x 7cm
PRES0310CXPPL	3mm x 10cm
PRES0407CPKPL	4mm x 7cm
PRES0410CPKPL	4mm x 10cm
PRES0413CPKPL	4mm x 13cm
PRES0510CPKPL	5mm x 10cm
PRES0517CPKPL	5mm x 17cm
PRES0620CPKPL	6mm x 20cm
PRES0720CPKPL	7mm x 20cm
PRES0830CPKPL	8mm x 30cm
PRES0930CPKPL	9mm x 30cm
PRES1030CPKPL	10mm x 30cm
PRES1040CPKPL	10mm x 40cm
PRES1240COM18	12mm x 40cm
PRES1447COM18	14mm x 47cm
PRES1655COM18	16mm x 55cm
PRES1860COM18	18mm x 60cm
PRES2065COM18	20mm x 65cm
PRES2265COM18	22mm x 65cm
PRES2465COM18	24mm x 65cm

Prestige Plus Helical coils	
Reference	Coil sizes
PRES0204HXPPL	2mm x 4cm
PRES0206HXPPL	2mm x 6cm
PRES0208HXPPL	2mm x 8cm
PRES0310HXPPL	3mm x 10cm
PRES0315HXPPL	3mm x 15cm
PRES0410HPKPL	4mm x 10cm
PRES0420HPKPL	4mm x 20cm
PRES0515HPKPL	5mm x 15cm
PRES0520HPKPL	5mm x 20cm
PRES0620HPKPL	6mm x 20cm
PRES0720HPKPL	7mm x 20cm
PRES0830HPKPL	8mm x 30cm
PRES0930HPKPL	9mm x 30cm
PRES1040HPKPL	10mm x 40cm

ULTRA Detachment Controller	
Model No.	Box content
ULTRA	1 detachment controller

Distributed in the UK by:

Aquilant

A HCS Company

csuk@hc21.group

01256 306506

healthcare21.eu

*Data on file

The Prestige Coil System is intended for use in the peripheral vasculature to endovascularly obstruct or occlude blood flow in vascular abnormalities of the peripheral vessels. Class III CE0297 in compliance with Medical Device Directive (MDD 93/42/EEC amended by 2007/47/EC). Manufactured by BALT USA LLC, 29 Parker, Irvine, CA 92618. Carefully read the instructions for use before use. First CE marking: 2021. For a list of indications, warnings, precautions and contraindications, refer to IFU-060. This document with associated pictures is non-contractual and is solely dedicated to healthcare professionals and BALT's distributors (BALT's supplier's distributors). The products commercialized by BALT shall exclusively be used in accordance with the instructions for use included in the boxes. The content of this document, in particular data, information, trademarks and logos is BALT S.A.S and affiliates' sole property. Consequently, all representation and/or reproduction, whether in part or in full, is forbidden and would be considered a violation of BALT S.A.S and affiliates' copyrights and other intellectual proprietary rights ©2022 BALT S.A.S and affiliates all rights reserved. This document with associated pictures is non-contractual and is solely dedicated to healthcare professionals and BALT S.A.S and affiliates' distributors. The products commercialized by BALT S.A.S and affiliates shall exclusively be used in accordance with the package inserts included in the boxes. DC057GB (07/2024)



Balt
10, rue de la Croix Vigneron,
95160 Montmorency, France
Tél. : +33 (0)1 39 89 46 41
www.baltgroup.com




prestige plus

volume and versatility



prestige plus

Intended for use in the peripheral vasculature to endovascularly obstruct or occlude blood flow in vascular abnormalities of the peripheral vessels.



Prestige design

Coil softness

premium coil softness through variations of filar and primary wind for deliverability and conformability

Complete offering

with a full range of Complex & Helical configurations to provide options in a variety of clinical scenarios

Pusher design

with spiral cut hypotube & progressively softer body coil ensuring smooth **pushability** & microcatheter **stability**

Prestige detachment

Instantaneous detachment

in under 1 second allows for rapid treatment & efficient case times

Reliable detachment system



entire size range
deliverable through
.0165"-.027" catheters*

new pusher
length of **185cm**

peripheral treatment