

fabian® Ventilators


Product Catalog



fabian® Ventilators

Product Catalog

- 4 Devices
- 6 Software Options
- 8 Accessories/Replacement Parts
 - 8 Flow Sensors
 - 9 Exhalation Valve Parts
 - 9 Patient Circuit Spare Parts
 - 10 Dual-Limb Ventilation
 - 10 Test Lungs
 - 11 Humidifiers
 - 12 Air Compressors
 - 13 Air and O₂ Sets
 - 13 Hose Sets
 - 13 Water Traps
 - 13 Adapters
 - 14 Carts
 - 14 Cart Spares
 - 14 Cart Accessories
 - 15 Power Cords
- 16 Aerogen® Accessories
- 18 Infant Flow® LP Nasal CPAP system
- 20 User Manuals



fabian® Ventilators

Devices

113001

fabian HFO

Neonatal & Paediatric HFO ventilator

Backup time: Conventional mode 2.5 hours
HFO mode 1 hour

111001.01

fabian HFO Light (without HFO module)

Neonatal & Paediatric battery

Backup time: Conventional mode 2.5 hours



121001

fabian Therapy Evolution

NIV ventilator

Backup time: 3 hours

Mode: Time cycled, pressure controlled

Ventilation modes: NCPAP, DUOPAP, O₂ Therapy

Included capability: Data trend 5 days



122001

fabian +nCPAP Evolution

Neonatal & Paediatric invasive and non-invasive unit

Backup time: 3 hours

Mode: Time cycled, pressure controlled

Ventilation modes: IPPV, SIPPV, SIMV, SIMV+PSV, PSV, CPAP, NCPAP, DUOPAP, O₂ Therapy

Included capability: Volume Limit (VL) & Volume Guarantee (VG), data trend 5 days

User interface: Colour touchscreen

Mode: Time cycled, pressure controlled

Ventilation modes: HFO, IPPV/IMV, SIPPV, SIMV, SIMV+PSV, PSV, CPAP, NCPAP, DUOPAP, O₂ Therapy,

Included capability: Volume Limit (VL) & Volume Guarantee (VG), data trend 5 days

Gas connection type (air and oxygen): NIST F

Power supply: 100-240 V, 50/60 Hz

All fabian battery backups are nickel-metal hydride battery

All platforms include software:

- 150212 – SW module NCPAP and DUOPAP, 1 piece
- 1502111 – SW module O₂ high-flow therapy, 1 piece
- 150226 – SW module data trend 5 days, 1 piece

fabian HFO, HFO Light, and fabian +nCPAP evolution include:

- 150221 – SW module VG & VL, 1 piece
- 100001.02 Accessory kit “invasive” includes:
 - 151100.01 – Flow sensor cable, 1 piece, reusable
 - 151120 – Flow sensor MPU neonatal, 1 piece
 - 151900 – Test lung , 1 piece
- 153011.02 – Infant Flow connecting tube, 2 pieces

fabian Therapy evolution includes accessory:

- 153001 – Infant Flow connecting tube, 2 pieces

Please choose a power cord when placing order

- 497180.XX* – Power cord, 1 piece

Extended Warranties

	1 year	2 years	3 years
fabian HFO ventilator	111191	111192	111193
fabian Therapy evolution	121191	121192	121193
fabian +nCPAP evolution	122191	122192	122193

fabian® Ventilators

Software Options

PRICO (Predictive Intelligent Control of Oxygenation)

150216

AT

PRICO-NIV Closed-Loop Control FiO₂ _SpO₂
4th Edition – includes ODU cable 3849 and SW key 150238

150210

AT

PRICO-NIV Closed-Loop Control FiO₂
3rd Edition – includes Masimo Module USB 152401 and SW key 150238

150217

AL

AN

PRICO-CNV Closed-Loop Control FiO₂ _SpO₂
4th Edition – includes ODU cable 3849 and SW key 150239

150220

AL

AN

PRICO-CNV Closed-Loop Control FiO₂
3rd Edition – includes Masimo Module USB 152401 and SW key 150239

150218

AI

PRICO-HFO Closed-Loop Control FiO₂ _SpO₂
4th Edition – includes ODU cable 3849 and SW key 150240

150230

AI

PRICO-HFO Closed-Loop Control FiO₂
3rd Edition – includes Masimo Module USB 152401 and SW key 150240

Device Abbreviations			
AI	fabian HFO	AL	fabian HFO Light
AN	fabian +nCPAP evolution	AT	fabian Therapy evolution

150213

AT

NIV Trigger for fabian Therapy evolution
Available for devices with SN > AT-001920

- Option for trigger in NCPAP/DuoPAP modes with single-limb ventilation and apnea detection in nCPAP mode
- Flow triggering (Infant Flow LP generators) in DuoPAP
- Pressure triggering (for Infant Flow LP Generator) in DuoPAP
- Allows S-DuoPAP functionalities

- Includes:**
- 150215 – software licence key, 1 piece
 - 151100.01 – flow sensor cable, reusable, 1 piece
 - 151111 – flow sensor, SPU, neonatal, 1 box of 10

150234

AI

Lung Recruitment Software for HFO-Mode

- Function to open the lung by increasing the pressure level cyclically

150223

AI

AL

AN

NIV Trigger for fabian +nCPAP & HFO

- Option for NCPAP/DuoPAP modes in single-limb ventilation
- Apnea detection in nCPAP mode
- Flow triggering (Infant Flow LP generators) in DuoPAP mode
- Pressure triggering in DuoPAP mode
- *Only for software version >=5.0*

150236

AI

AL

FOT Software Option License Key

- Forced Oscillation Technique (FOT) to measures the respiratory reactance in nearly all ventilation modes

150225

AI

AL

AN

AT

PDMS Licence Key

- For VueLink, AcuLink and Intellibridge (communication interface)

150214

AT

NIV Trigger Upgrade Kit
Available for devices with SN > AT-001920

- Option for trigger in NCPAP/DuoPAP modes with single-limb ventilation and apnea detection in nCPAP mode
- Flow triggering (Infant Flow LP generators) in DuoPAP
- Allows S-DuoPAP functionalities

- Includes:**
- 150215 – software licence key, 1 piece
 - 122222 – connector PCB for AT device, 1 piece
 - 121200.22 – main board rev. 2.2.1, 1 piece
 - 151100.01 – flow sensor cable, reusable, 1 piece
 - 151111 – flow sensor, SPU, neonatal, 1 box of 10

fabian® Ventilators

Accessories

Flow Sensors



151100.01 AI AL AN AT

Flow Sensor Cable, Reusable
Compatible with all listed flow sensors in this pricelist
Cable length: 2m Weight: 82.1g



151111 AI AL AN AT

Flow Sensor, SPU, Neonatal
Compatible with flow sensor cable (151100.01)
Weight: 10.5g 10 pieces per box



151120 AI AL AN AT

Flow Sensor, MPU, Neonatal
Compatible with flow sensor cable (151100.01)
Autoclavable at 134°C, 30 cycles maximum
Individually packed
Refer to Instructions for Use for complete instructions for autoclaving



151121 AI AL AN AT

Flow Sensor, MPU, Paediatric
Compatible with flow sensor cable (151100.01)
Autoclavable at 134°C, 30 cycles maximum
Individually packed

Device Abbreviations

AI	fabian HFO	AL	fabian HFO Light
AN	fabian +nCPAP evolution	AT	fabian Therapy evolution

Exhalation Valve Parts

151501.01 AI AL AN

Exhalation Valve Membrane, MPU
Autoclavable at 134° C, 30 cycles maximum



151520.01 AI AL AN

Membrane Holder, 22OD/15ID, MPU, PC
Autoclavable at 134°C for 30 cycles maximum
Made of polycarbonate (PC)



Patient Circuit Spares

155200 AI AL AN AT

Connector 15OD/9-11OD, SPU
■ 1 piece per box



156206 AL

Straight Connector 22/10, MPU
■ 10 mm male to 22 mm female
■ Autoclavable at 134°C
■ Material: Polysulfone



156209 AI AL AN

Straight Connector 15M/22F, MPU
■ Autoclavable at 134°C



Dual-Limb Ventilation

154308.02 AI AL AN

Neo Tubing Set, Dual Heated, SPU

Made by ALTERA

- 10 sets per box, includes:
- Proximal pressure line, 1 piece
 - Incubator extension tube, 1 piece
 - Infant Flow connecting tube, 1 piece
 - Adapter to convert to NIV ventilation, 1 piece
 - Autofeed AL-0005 humidifier chamber, 1 piece



154902 AI AL AN AT

iNO Kit for fabian HFO

- 10 sets per box
- With Luer-lock adapters, one-way valve and 22F-22F straight connector
 - For connection to NOxBOX system, e.g., single-limb ventilation

153011.02 AI AL AN AT

Infant Flow Tube, 60cm, SPU

- 10 pieces per box

Test Lungs

151900 AI AL AN

Test Lung, 1 Bellows, Neonatal

Autoclavable at 134° C
Material: Silicone



151901 AI AL AN

Test Lung, 2 Bellows, Paediatric

Test lung for neonatal ventilators, with 2 bellows, to simulate one lung ventilation and ET tube leakage



Device Abbreviations			
AI	fabian HFO	AL	fabian HFO Light
AN	fabian +nCPAP evolution	AT	fabian Therapy evolution

Humidifiers



152120.01EU AI AL AN AT

Humidifier, WILAmated, AIRcon G2, 220V/50H

- EU power cord, type C, 1 piece
- 497090.XX – angled cords (See page 15)

152120.01UK AI AL AN AT

Humidifier WILAmated, AIRcon G2, 220V/50Hz

- UK power cord, type G, 1 piece
- 152120.01EU – additional power cords

Each colour LCD respiratory gas humidifier with dual servo control for airway temperature also includes:

- User manual, 1 piece
- Fixation for bracket, 1 piece
- Heater wire adapter (i+e, Wilamed type), 1 piece
- Short instructions, 1 piece
- Temperature probe, 1 piece



Air Compressors

490052.EU AI AL AN AT

Compressor Silent 230V/50Hz

Device on wheels

All air compressors include:

- 2 m tube, 1 piece
- 497180.XX – power cord, 1 piece

See page 15



Compressor with wheels

490058 AI AL AN AT

Compressor Silent 230V/60Hz*

490059.EU AI AL AN AT

Compressor Silent 230V/50Hz*

490069.US AI AL AN AT

Compressor Silent 115V/60Hz*



Compressor without wheels

*For use with Smart Trolley part #493000

Device Abbreviations

AI	fabian HFO	AL	fabian HFO Light
AN	fabian +nCPAP evolution	AT	fabian Therapy evolution

Air and O₂ Sets

499090 AI AL AN AT

Filter Set for Medical Grade Air-NIST

To be mounted on the trolley or on the wall, screws are included



Hose Sets

499230.BS AI AL AN AT

Air/O₂ Hose Set, 3m, BL/WH, NIST-BS5682

- Air hose, BK-WH, 3 m, NIST-BS5682, 1 piece
- O₂ hose, WH, 3 m, NIST-BS5682, 1 piece

499230.DE AI AL AN AT

Air/O₂ Hose set, 3m, NIST-DIN

- Air hose, BK-WH, 3 m, NIST-DIN, 1 piece
- O₂ hose, WH, 3 m, NIST-DIN, 1 piece

499200.DE

Air/O₂ Hose Set 5M BL/WH, DIN120 - NIST

- Air hose, BK-WH, 5m, NIST-DIN, 1 piece
- O₂ hose, WH, 5m, NIST-DIN, 1 piece



Water Traps

499091.EU AI AL AN AT

Air Water Trap NIST 90°

Water trap with female NIST adapter on ventilator side and male NIST adapter on input side



499301.US AI AL AN AT

Air Water Trap NIST 90° Male #A3/DISS

Canada and Japan



Adapters

499391.US AI AL AN AT

Air Adapter, NIST-f / DISS-m, 90°

Canada and Japan



499392.US AI AL AN AT

O₂ Adapter, NIST-f / DISS-m, 90°

Canada and Japan



Carts

400001 AI AL

Cart for fabian HFO*



400002 AN AT

Cart for fabian EVO*



*Mounting plate is not included and needs to be ordered separately.

Cart Spares

401101 AI AL

HFO Mounting Plate

Material: Stainless steel, powder-coated



Includes:

- DIN rails (square tube 25/10), 2 pieces
- Fixing material

401202 AN AT

EVO Mounting Plate

Material: Stainless steel, powder-coated



Includes:

- DIN rails (square tube 25/10), 2 pieces
- Fixing material

Cart Accessories

400011.01

IV Pole with DIN Clamp

- Short IV-pole length 1.2 m
- Adjustable height
- Universal clamp for 25 x 10 mm standard bar and 25 mm round pole



410001

Circuit Support Arm with DIN Rail Clamp

- One knob for fixation and adjustment
- DIN rail clamp bracket



400013

Basket

AMS cart & light cart for trolley #401100 or #401200

400014

Bottle Holder

AMS cart & light cart for trolley #401100 or #401200

Device Abbreviations

AI	fabian HFO	AL	fabian HFO Light
AN	fabian +nCPAP evolution	AT	fabian Therapy evolution

Power Cords

90°Angle

Straight

Australia (AUS)

497090.AU

AI AL AN AT

- C13 / AS3112
- 2.5 meters



497180.AU

AI AL AN AT

- C13 / AS3112
- 2.5 meters



Switzerland (CH)

497090.CH

AI AL AN AT

- C13 / T12
- 2.5 meters



497180.CH

AI AL AN AT

- C13 / T12
- 2.5 meters



Europe / SCHUKO (EU)

497090.EU

AI AL AN AT

- CEE 7/7
- 2.5 meters



497180.EU

AI AL AN AT

- CEE 7/7
- 2.5 meters



United Kingdom (UK)

497090.UK

AI AL AN AT

- C13 / BS1363
- 2.5 meters



497180.UK

AI AL AN AT

- C13 / BS1363
- 2.5 meters



United States (USA) and Japan (JP)

497090.JP

AI AL AN AT

- C13 / T12
- 2 meters



497180.01US

AI AL AN AT

- Type-B - NEMA



fabian® Ventilators

Aerogen® accessories*

AG-UC1000-XX

Aerogen® USB Controller

- Set includes 1 piece



AG-AS3300-XX

Aerogen Solo® Convenience Kit - Nebuliser & Adult T-piece

- Set includes 5 pieces



AG-152222-XX

Aerogen Solo® Nebuliser

- Set includes 5 pieces



AG-AS3010-XX

Aerogen Solo® Adult T-piece - 22mm

- Set includes 10 pieces



AG-AS3035-XX

Aerogen Solo® Paediatric T-piece - 12mm

- Set includes 10 pieces



Regional Suffix Guide (XX)

The SKU suffix (XX) designates the approved country/market version of the product. Orders should be placed using the suffix corresponding to the destination country.

For example:

Code	Market	Code	Market
UK	UK & Ireland	CH	China
EU	Europe	IC	India
LA	Latin America	BR	Brazil
AP	Malaysia	SA	South Africa
JP	Japan	AU	Australia/ New Zealand

* Product availability may vary by market. Please contact your local representative for details.



Infant Flow® LP


Nasal CPAP system†

7772020LP-40

Infant Flow LP Generator Circuit Kit

Includes generator, sizing guide, and small, medium, and large nasal prongs

- Single-patient use
- Set includes 20 kits



Infant Flow LP Generator Circuit Kit

For use with F&P MR 730 and 850 Humidifier with RT132 circuit

- 7772011-40 ▪ Set includes 10 kits
- 7772021-40 ▪ Set includes 20 kits

China PN suffix: -15

7772000LP

Infant Flow Generator without Prongs

- Single-patient use
- Set includes 10 packs

Infant Flow LP Nasal Prongs

Each package includes 10 pieces.



China PN suffix: -15

Infant Flow LP Nasal Mask

Each package includes 10 pieces.



China PN suffix: -15

Infant Flow LP Headgear

Each package includes 10 pieces.

777040XS-40	Extra small	17–21 cm
777040S-40	Small	21–26 cm
777040SM-40	Small/med	24–28 cm
777040M-40	Medium	26–32 cm
777040L-40	Large	32–37 cm
777040XL-40	Extra large	37–42 cm



Infant Flow LP Bonnets

Each package includes 10 pieces.

777010-40	Size 000	18–20 cm
777012-40	Size 00	20–22 cm
777014-40	Size 0	22–24 cm
777016-40	Size 1	24–26 cm
777018-40	Size 2	26–28 cm
777020-40	Size 3	28–30 cm
777022-40	Size 4	30–32 cm
777024-40	Size 5	32–34 cm
777026-40	Size 6	34–36 cm
777028-40	Size 7	36–38 cm



†Product availability may vary by market. Please contact your local representative for details.

†Product availability may vary by market. Please contact your local representative for details.



fabian HFO

113003.AR	Arabic	113003.DE	German	113003.PT	Portuguese
113003.BG	Bulgarian	113003.GR	Greek	113003.PTBR	Brazilian Portuguese
113003.HR	Croatian	113003.hu	Hungarian	113003.RO	Romanian
113003.CZ	Czech	113003.IT	Italian	113003.RU	Russian
113003.DK	Danish	113003.JP	Japanese	113003.RS	Serbian
113003.NL	Dutch	113003.KO	Korean	113003.SK	Slovakian
113003.EN	English	113003.LV	Latvian	113003.ES	Spanish
113003.EE	Estonian	113003.LT	Lithuanian	113003.SE	Swedish
113003.FI	Finnish	113003.NO	Norwegian	113003.TR	Turkish
113003.FR	French	113003.PL	Polish		

fabian +nCPAP Evolution

122003.AR	Arabic	122003.FR	French	122003.PT	Portuguese
122003.BG	Bulgarian	122003.DE	German	122003.PTBR	Brazilian Portuguese
122003.CN	Chinese	122003.GR	Greek	122003.RO	Romanian
122003.HR	Croatian	122003.HU	Hungarian	122003.RU	Russian
122003.CZ	Czech	122003.IT	Italian	122003.SE	Slovakian
122003.DK	Danish	122003.JP	Japanese	122003.ES	Spanish
122003.NL	Dutch	122003.LT	Lithuanian	122003.SE	Swedish
122003.EN	English	122003.LV	Latvian	122003.TR	Turkish
122003.EE	Estonian	122003.NO	Norwegian	122003.UA	Ukrainian
122003.FI	Finnish	122003.PL	Polish		

fabian Therapy Evolution

121003.AR	Arabic	121003.FR	French	121003.PT	Portuguese
121003.BG	Bulgarian	121003.DE	German	121003.PTBR	Brazilian Portuguese
121003.CN	Chinese	121003.GR	Greek	121003.RO	Romanian
121003.HR	Croatian	121003.HU	Hungarian	121003.RU	Russian
121003.CZ	Czech	121003.IT	Italian	121003.SE	Slovakian
121003.DK	Danish	121003.JP	Japanese	121003.ES	Spanish
121003.NL	Dutch	121003.LT	Lithuanian	121003.SE	Swedish
121003.EN	English	121003.LV	Latvian	121003.TR	Turkish
121003.EE	Estonian	121003.NO	Norwegian	121003.UA	Ukrainian
121003.FI	Finnish	121003.PL	Polish		



For more information about [fabian ventilators](#),
please visit qr.zoll.com/fabianHFO or call (800) 348-9011.

Aquilant
A HC21 Company

csuk@hc21.group
01256 306506
healthcare21.eu

ZOLL Medical Corporation

269 Mill Road | Chelmsford, MA 01824 | 978-421-9655 | 800-804-4356 | zoll.com

Not all products are available in all markets. Availability is subject to regulatory approvals and local market conditions.

Copyright © 2026 ZOLL Medical Corporation. All rights reserved. fabian, iFlow, Infant Flow, and ZOLL are trademarks or registered trademarks of ZOLL Medical Corporation in the United States and/or other countries. All other trademarks are the property of their respective owners. MCN VP 2503 0041-05

ZOLL®