



Medtronic

McGRATH™ video laryngoscopes

Comparison guide




McGRATH™ MAC video laryngoscope



McGRATH™ MAC+ video laryngoscope

4.5x

higher screen resolution



45%

larger screen size

	McGRATH™ MAC video laryngoscope	McGRATH™ MAC+ video laryngoscope
Specification		
Weight	0.180kg	0.180kg
Useful life	4 years or 2000 decontamination cycles	4 years or 2000 decontamination cycles
Visualization		
Screen resolution	320x240	720x480 (4.5x increase)
Screen Size	2.4 inch	3 inch (45% increase)
Durability		
Handle	Reinforced screen hinge	Reinforced screen hinge
IP rating	IPX8	IPX8
Battery management		
Battery type	Disposable and rechargeable batteries with auto-off feature	Rechargeable battery with auto-off feature
Disposable battery life	250 minutes	N/A
Rechargeable battery life	10,000 minutes	10,000 minutes
Run time on full charge	120 min (with rechargeable battery)	60 min (with rechargeable battery)
Ordering information and SKU numbers		
Handle	301-000-000	302-000-000
3.6V disposable battery	340-000-000	N/A
Rechargeable battery	341-100-00	341-100-00
Relay	N/A	302-100-000
Connected screen	N/A	302-200-000
MAC + Features		
Streaming and file management	N/A	Yes
Video recording	N/A	Yes
Touch screen	N/A	Yes
Still image capture	N/A	Yes
File storage (video and still image)	N/A	Yes
Wireless connectivity	N/A	Yes

For trained personnel only. For specific indications and instructions for use, please refer to the product manual.

medtronic.com/mcgrathmacplus