
# ComVi Pyloric/Duodenal Stent [Flare-Type]

for intrinsic and/or extrinsic malignant pyloric/duodenal stricture



## Features

- **Triple layered construction:** A biocompatible PTFE membrane tube is held between the inner and outer mesh
- An uncovered proximal flare end to reduce migration
- The PTFE membrane prevents the risk of tissue invasions

- **Large diameters (up to 26mm)** are loaded into a 10.5Fr delivery system

- **Radiopaque marker:** 3 (three) at the Flared end and 2 (two) on both ends of the body

TaeWoong MEDICAL

# Ordering Information

| TTS (Through The Scope) | | | | |
|---|---|---|---|---|
| Distal End Bare Type | | | | |
| Code | Stent | | Delivery | |
| | Diameter (mm) | Length (cm) | Profile (Fr) | Usable Length (cm) |
| PCT1806P | 18 | 6 | 10.5 | 180 |
| PCT1808P | | 8 | | |
| PCT1810P | | 10 | | |
| PCT1812P | | 12 | | |
| PCT2006P | 20 | 6 | 10.5 | 180 |
| PCT2008P | | 8 | | |
| PCT2010P | | 10 | | |
| PCT2012P | | 12 | | |
| PCT2206P | 22 | 6 | 10.5 | 180 |
| PCT2208P | | 8 | | |
| PCT2210P | | 10 | | |
| PCT2212P | | 12 | | |

| TTS (Through The Scope) | | | | |
|---|---|---|---|---|
| Distal End Bare Type | | | | |
| Code | Stent | | Delivery | |
| | Diameter (mm) | Length (cm) | Profile (Fr) | Usable Length (cm) |
| PCT2406P | 24 | 6 | 10.5 | 180 |
| PCT2408P | | 8 | | |
| PCT2410P | | 10 | | |
| PCT2412P | | 12 | | |
| PCT2606P | 26 | 6 | 10.5 | 180 |
| PCT2608P | | 8 | | |
| PCT2610P | | 10 | | |
| PCT2612P | | 12 | | |

| TTS (Through The Scope) | | | | |
|---|---|---|---|---|
| Distal End Bare Type | | | | |
| Code | Stent | | Delivery | |
| | Diameter (mm) | Length (cm) | Profile (Fr) | Usable Length (cm) |
| PCT1806P-12 | 18 | 6 | 10.5 | 120 |
| PCT1808P-12 | | 8 | | |
| PCT1810P-12 | | 10 | | |
| PCT1812P-12 | | 12 | | |
| PCT2006P-12 | 20 | 6 | 10.5 | 120 |
| PCT2008P-12 | | 8 | | |
| PCT2010P-12 | | 10 | | |
| PCT2012P-12 | | 12 | | |
| PCT2206P-12 | 22 | 6 | 10.5 | 120 |
| PCT2208P-12 | | 8 | | |
| PCT2210P-12 | | 10 | | |
| PCT2212P-12 | | 12 | | |
| PCT2406P-12 | 24 | 6 | 10.5 | 120 |
| PCT2408P-12 | | 8 | | |
| PCT2410P-12 | | 10 | | |
| PCT2412P-12 | | 12 | | |
| PCT2606P-12 | 26 | 6 | 10.5 | 120 |
| PCT2608P-12 | | 8 | | |
| PCT2610P-12 | | 10 | | |
| PCT2612P-12 | | 12 | | |

| TTS (Through The Scope) | | | | |
|---|---|---|---|---|
| Distal End Bare Type | | | | |
| Code | Stent | | Delivery | |
| | Diameter (mm) | Length (cm) | Profile (Fr) | Usable Length (cm) |
| PCT1806P-22 | 18 | 6 | 10.5 | 220 |
| PCT1808P-22 | | 8 | | |
| PCT1810P-22 | | 10 | | |
| PCT1812P-22 | | 12 | | |
| PCT2006P-22 | 20 | 6 | 10.5 | 220 |
| PCT2008P-22 | | 8 | | |
| PCT2010P-22 | | 10 | | |
| PCT2012P-22 | | 12 | | |
| PCT2206P-22 | 22 | 6 | 10.5 | 220 |
| PCT2208P-22 | | 8 | | |
| PCT2210P-22 | | 10 | | |
| PCT2212P-22 | | 12 | | |
| PCT2406P-22 | 24 | 6 | 10.5 | 220 |
| PCT2408P-22 | | 8 | | |
| PCT2410P-22 | | 10 | | |
| PCT2412P-22 | | 12 | | |
| PCT2606P-22 | 26 | 6 | 10.5 | 220 |
| PCT2608P-22 | | 8 | | |
| PCT2610P-22 | | 10 | | |
| PCT2612P-22 | | 12 | | |

Distributed in the UK by:

Aquilant
A HC2I Company
csuk@hc21.group
01256 306506
healthcare21.eu

ⓒ 2024, Taewoong Medical Co., Ltd. All rights reserved. TW-MBPD-CF (Rev.03) Printed in Korea. May 16. 2024