JOBST® FarrowWrap® 4000

A-D Below Knee

For quotes please email: **quotes@hc21.group**
To order please email Customer Services on: **info@hc21.group**

Date:_____ Patient Name:_____

Measured By:_____ Tel:_____ Email:_____

Delivery Address:_____

HC21 Healthcare

Measuring Guide Video

Measuring Instructions

Legpiece


1. Measure ankle **B** and calf **C** circumferences
2. Measure length posteriorly (back of leg) following the contour from **B** to **D** (two finger widths below the crease of the knee)
3. Determine size **(XS-L)** & length **(Regular or Tall)** using chart below

NOTE: 1 x JOBST® FarrowHybrid (Liners) wide is provided. If additional or alternative sized liners are required, **please order separately.**

JOBST FarrowHybrid (Liners)

1. Measure **A1** circumference
2. Record **shoe size**
3. Determine size **(S-XL)**, width **(Regular or Wide)**

NOTE: Do not apply JOBST® FarrowHybrid (Liners) if calf circumference exceeds **60cm**. When a large JOBST® FarrowWrap® 4000 is ordered, an additional non-compressive liner will be provided to accommodate a calf circumference greater than **60cm.** Consider a JOBST® FarrowWrap® footpiece in this circumstance.



NOTE: 3cm allowance on either side of measurement - **B & C**

	X Small	Small	Medium	Large
JOBST® FarrowWrap® 4000 Legpiece				
C Calf	25-38cm	32-48cm	37-56cm	43-67cm
B Ankle	17-23cm	19-28cm	25-37cm	29-43cm
B-D Regular	30-33cm	30-33cm	30-33cm	30-33cm
B-D Tall	34-37cm	34-37cm	34-37cm	34-37cm

☐ JOBST® FarrowWrap® 4000 Legpiece

Size ☐ X Small ☐ Small ☐ Medium ☐ Large
Length **B-D** ☐ Regular (30-33cm) ☐ Tall (34-37cm)
Colour ☐ Tan ☐ Black
Product Code ☐☐☐☐☐☐—☐☐
Quantity Required ☐☐☐
Repeat prescription required every ☐☐ months

☐ JOBST® FarrowHybrid® (Liners)

Size ☐ Small (2-5) ☐ Medium (6-8) ☐ Large (9-11) ☐ X Large (12-13) (shoe size)
Width **A1** ☐ Regular (20-27cm) ☐ Wide (28-35cm)
Colour ☐ Tan ☐ Black (Available in Wide variant only)
Product Code ☐☐☐☐☐☐—☐☐
Quantity Required ☐☐☐
Repeat prescription required every ☐☐ months

DEFY GRAVITY

JOBST

NOTE: Please turn page for product codes required for ordering

JOBST® FarrowWrap® 4000 Order Codes

JOBST® FarrowWrap® 4000 Legpiece (30-40mmHg)

For patients with **mild to moderate** oedema, with or without a leg ulcer



Product Code	Size	Length	Colour
76661-00	Extra small	Regular	Beige
76661-08	Extra small	Regular	Black
76661-04	Extra small	Tall	Beige
76661-12	Extra small	Tall	Black
76661-01	Small	Regular	Beige
76661-09	Small	Regular	Black
76661-05	Small	Tall	Beige
76661-13	Small	Tall	Black
76661-02	Medium	Regular	Beige
76661-10	Medium	Regular	Black
76661-06	Medium	Tall	Beige
76661-14	Medium	Tall	Black
76661-03	Large	Regular	Beige
76661-11	Large	Regular	Black
76661-07	Large	Tall	Beige
76661-15	Large	Tall	Black

JOBST® FarrowHybrid® (20-30mmHg)

For the treatment of **mild to moderate** oedema in the foot and ankle



Product Code	Size	Colour
76664-00	Small - Regular (shoe size 2-5)	Tan
76664-01	Medium - Regular (shoe size 6-8)	Tan
76664-02	Large - Regular (shoe size 9-11)	Tan
76664-03	Extra Large - Regular (shoe size 12-13)	Tan
76664-12	Small - Wide (shoe size 2-5)	Tan
76664-13	Medium - Wide (shoe size 6-8)	Tan
76664-14	Large - Wide (shoe size 9-11)	Tan
76664-15	Extra Large - Wide (shoe size 12-13)	Tan
76664-16	Small - Wide (shoe size 2-5)	Black
76664-17	Medium - Wide (shoe size 6-8)	Black
76664-18	Large - Wide (shoe size 9-11)	Black
76664-19	Extra Large - Wide (shoe size 12-13)	Black

Do not apply JOBST® FarrowHybrid® if calf circumference exceeds **60cm**, consider a JOBST® FarrowWrap® footpiece in this circumstance

Application Videos

Scan to access HC21 Compression Products available on PCRS

JOBST Code Finder

Scan the QR code or search for **codefinder.jobst.com** on your mobile device.

HC21
info@hc21.group
00 353 214 875 055
healthcare21.eu

essity JOBST®, an Essity brand

*The HC21 Group is fully compliant with the GDPR and applies the appropriate security to the processing of all personal data. By completing this order form, you are confirming that you are aware of your obligation to obtain consent from the patient on the processing of their data for the production of their JOBST® compression garment in accordance with the Data Protection Act 2018, and the EU General Data Protection Regulation [GDPR]. For more information about our Privacy Policy, please contact **gdpr@hc21.group**