

COMPRESSION THERAPY CATALOGUE

DEFY GRAVITY

JOBST

Find ordering and prescribing information fast with JOBST Code Finder.

No more searching through a catalogue, use our web app on your mobile device to find the information you need.

CODEFINDER.JOBST.COM



About JOBST® Compression Garments

Starting in 1950 with its founder, Conrad Jobst, who suffered from venous disease, JOBST® has always had the strongest commitment to improving the quality of life for those with venous disease or oedema-related conditions, such as lymphoedema.

Essity develop and manufacture medical products used in the care and management of venous and lymphatic diseases including: varicose veins, deep vein thrombosis, chronic venous insufficiency, venous ulcers, peripheral oedemas and lymphoedema. We manufacture and distribute a wide range of innovative, high quality products ranging from compression garments to bandages.

Our commitment to quality is reflected in our superior product design, vigorous product testing, and our state-of-the-art manufacturing procedures and facilities. Essity is committed to total satisfaction for its clinicians and patients, which is why each JOBST® product comes with a six-month guarantee against garment failure due to manufacturing or material defect. In addition, we strive to deliver all JOBST® garments within five working days from receipt of order.

Contents

Compression Bandages / Foam Pad

Comprilan®	4
Tensopress®	4
JOBST® Foam pad	4
JOBST® Compri2	5
JOBST® Compri2 lite	5
JOBST® Comprifore	6-7

Ready-to-Wear Compression Garments

JOBST® Classic	8-10
JOBST® Bella™ Lite	11-13
JOBST® ForMen Ambition	14-15
JOBST® ForMen Explore	16-17
JOBST® Opaque	18-39
JOBST® UltraSheer	40-55
JOBST® UlcerCare	56-57

Wrap Compresion Systems

JOBST® FarrowWrap®	58-67
JOBST® Farrow® Toe Cap	68
JOBST® Farrow® Glove	69

Custom-Fit Compression Garments

JOBST® Elvarex®	70-74
JOBST® Elvarex® Plus	75-76
JOBST® Elvarex® Soft	77-79
JOBST® Elvarex® Soft Seamless	80-81
JOBST® Custom-Fit Styles and Options	82-85
JOBST® Confidence	86-87

Others

JOBST® Relax	88-89
JOBST® Surgical Vest	90
Bellisse®	91-92
JOBST® JoViPads	93-96

Please note: we are unable to include all styles, options and colours for certain products therefore if you cannot find the exact product you require, please contact our Customer Service team for further assistance. Additional information can also be found on our website at **www.jobst.co.uk**

Comprilan®

100% Cotton Short-Stretch Bandage

Available on Drug Tariff

Comprilan can be used for phlebological and lymphological indications to support the muscle pump action and improve the haemodynamics of venous return. This short-stretch bandage is durable and extremely resistant.

- Washable and re-usable
- Air permeable
- Bandage ends do not fray



Code	Size	Pack Qty	Case Qty	PIP Code	NHS SC Code
01026-00	6cm x 5m	1 x 5	40 x 5	204-1895	EBA046
01027-00	8cm x 5m	1 x 5	34 x 5	204-1903	EBA047
01028-00	10cm x 5m	1 x 5	28 x 5	204-1911	EBA048
01029-00	12cm x 5m	1 x 5	24 x 5	204-1929	EBA049

Tensopress®

High Compression Extensible Bandage (Type 3C)

Available on Drug Tariff

Tensopress delivers measured, safe, sub-bandage pressure for the treatment of venous leg ulcers by increasing blood flow velocity and improving venous return.

- A cotton viscose mix
- Central green line guides and simplifies precise and even application



Code	Size	Pack Qty	Case Qty	PIP Code	NHS SC Code
71723-01	7.5cm x 3m	1 x 10	6 x 10	043-5339	EBA075
71723-00	10cm x 3m	1 x 10	8 x 10	043-5891	EBA077

JOBST® Foam

Foam Rubber Kidney Shaped Pad

JOBST® Foam is designed to evenly distribute the pressure from compression stockings over bony prominences.

- Latex-free formulation



Code	Size	Pack Qty
78496-00	9cm x 1cm	1 x 1

JOBST® Compri2

Two-Layer Bandage System

JOBST® Compri2 is latex-free and has accessible guidance tools to aid application and promote the correct compression. Two compression levels, ensuring the right level can be selected for both venous and mixed aetiology leg ulcers.

Compression on demand

Two choices to achieve the compression needed

Type of disease	Incidence	ABPI	Compression level	
Venous	~60%	0.8 - 1.3	40 mmHg	JOBST® Compri2
Mixed aetiology	~20%	>0.5 - <0.8	20 - 30 mmHg	JOBST® Compri2 lite
Arterial	~20%	<0.5	No compression	-






1st layer

- Application around the heel becomes much easier because the bandage stretches in both dimensions
- Tearable to allow extra padding to be easily applied where necessary







2nd layer

- A simple indicator to assist achieving the correct compression level every time
- The bandage can be adjusted during application



Code	Ankle Size	PIP Code	NHS SC Code
JOBST® Compri2 - regular compression (ABPI >0.8)			
76271-00 (Latex-free)	18 - 25cm	400-6177	EBA145
76271-01 (Latex-free)	25 - 32cm	400-6185	EBA146
JOBST® Compri2 lite - light compression (ABPI 0.5 - 0.8)			
76271-02 (Latex-free)	18 - 25cm	400-6193	EBA147
76271-03 (Latex-free)	25 - 32cm	400-6201	EBA148





JOBST® Comprifore

Available on Drug Tariff

Multi-Layer Compression Bandage System for Effective Venous Leg Ulcer Therapy

The JOBST® Comprifore multi-layer bandage kit contains all the components needed for the treatment of venous leg ulcers for an ankle circumference between 18 - 25cm, providing sustained, graduated compression of 40mmHg[1].

Sustained compression

JOBST® Comprifore provides sustained compression for up to seven days for an excellent cost-benefit ratio. The bandage system delivers a high stiffness factor for effective compression and allows for easy and safe application to minimise nursing time.

Options			
JOBST® Comprifore is available as a kit for the treatment of venous leg ulcers with an ABPI of 0.8 or above	JOBST® Comprifore lite Kit provides reduced compression for mixed aetiology ulcers with an ABPI of between 0.5 and 0.8	Both kits are available in latex-free versions to reduce risk of allergies	JOBST® Comprifore components are also available individually
✓	✓	✓	✓

1 Targeted nominal compression for an average ankle size.

JOBST® Comprifore

Product	Code	Description	Contents	Case Unit	PIP Code	NHS SC Code
JOBST® Comprifore Kits						
JOBST® Comprifore	72661-04	4-Layer kit	1 Kit	8 Kits	369-3678	EBA066
JOBST® Comprifore LF	72661-05	Latex-Free 4-Layer kit	1 Kit	8 Kits	369-3694	EBA067
JOBST® Comprifore lite	72661-06	3-Layer kit	1 Kit	8 Kits	369-3702	EBA068
JOBST® Comprifore lite LF	72661-07	Latex-Free 3-Layer kit	1 Kit	8 Kits	369-3710	EBA069
Note: suitable for an ankle circumference of 18 - 25cm						

Product	JOBST® Comprifore #1	JOBST® Comprifore #2	JOBST® Comprifore #3	JOBST® Comprifore #4
JOBST® Comprifore	✓	✓	✓	✓
JOBST® Comprifore LF	✓	✓	✓ (LF)	✓ (LF)
JOBST® Comprifore lite	✓	✓	-	✓
JOBST® Comprifore lite LF	✓	✓	-	✓ (LF)

Components within the JOBST® Comprifore multi-layer kits

JOBST® Comprifore #1

Absorbent padding bandage

- Distributes pressure evenly
- Protects bony prominences
- Absorbs exudate



JOBST® Comprifore #3

Light compression bandage

- Provides the first layer of compression
- Central coloured line facilitates correct application
- Also available in a latex-free version



JOBST® Comprifore #2

Light support bandage

- Aids absorbency
- Smooths padding bandage



JOBST® Comprifore #4

Flexible cohesive compression bandage

- Provides the second layer of compression
- Keeps bandages securely in place
- No additional fixation necessary
- Also available in latex-free version



JOBST® Classic

Available on Drug Tariff

Ready-to-Wear RAL Compression Hosiery

JOBST® Classic is a range of ready-to-wear, circular-knit hosiery for the treatment of moderate lymphoedema and venous disease.

Firm feel fabric

JOBST® Classic is manufactured using a firmer feel fabric with stiff resistance, to apply a stronger level of compression. The reinforced toe area and robust elastic heel provides high durability with a comfortable fit.

Compression Classes	
RAL CCL 2 = 23-32mmHg	RAL CCL 3 = 34-46mmHg
✓	✓

Sizes and Options / Materials - Knitted or dotted silicone top band for thigh high				
Six Sizes	Regular	Petite	Open Toe	Closed Toe
✓	✓	✓	✓	✓






Knitted top band



Dotted silicone top band

Colours	
Sand	Pearl

Styles				
				
AD Knee High	AG Thigh High	AT Tights	AG-T Chap Style*	AT Men*

*not available via Drug Tariff

JOBST® Classic - CCL 2

Size	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC
	Knee High - Sand						**Thigh High - Sand**					
	RAL CCL 2 - Open Toe						**RAL CCL 2 - Open Toe**					
	Petite			**Regular**			**Dotted silicone band**					
I	77402-12	415-2401	-	77401-12	415-2278	-	77404-12	415-2468	-			
II	77402-13	415-2419	-	77401-13	415-2351	-	77404-13	415-2476	-			
III	77402-14	415-2427	-	77401-14	415-2369	-	77404-14	415-2484	-			
IV	77402-15	415-2435	-	77401-15	415-2377	-	77404-15	415-2492	-			
V	77402-16	415-2443	-	77401-16	415-2385	-	77404-16	415-2500	-			
VI	77402-17	415-2450	-	77401-17	415-2393	-	77404-17	415-2518	-			
	Knee High - Sand											
	RAL CCL 2 - Closed Toe											
	Petite			**Regular**								
I	77402-30	-	-	77401-30	-	-						
II	77402-31	-	-	77401-31	-	-						
III	77402-32	-	-	77401-32	-	-						
IV	77402-33	-	-	77401-33	-	-						
V	77402-34	-	-	77401-34	-	-						
VI	77402-35	-	-	77401-35	-	-						
	Tights - Sand						**Knee High - Pearl**					
	RAL CCL 2 - Open Toe						**RAL CCL 2 - Open Toe**					
	Petite			**Regular**			**Petite**			**Regular**		
I	79488-18	415-2286	-	79489-18	415-2344	-	77402-00	415-1635	-	77401-00	415-1692	-
II	79488-19	415-2302	-	79489-19	415-2229	-	77402-01	415-1643	-	77401-01	415-1700	-
III	79488-20	415-2294	-	79489-20	415-2237	-	77402-02	415-1650	-	77401-02	415-1718	-
IV	79488-21	415-2310	-	79489-21	415-2245	-	77402-03	415-1668	-	77401-03	415-1726	-
V	79488-22	415-2328	-	79489-22	415-2252	-	77402-04	415-1676	-	77401-04	415-2211	-
VI	79488-23	415-2336	-	79489-23	415-2260	-	77402-05	415-1684	-	77401-05	415-1734	-
	Thigh High - Pearl						**Tights - Pearl**					
	RAL CCL 2 - Open Toe						**RAL CCL 2 - Open Toe**					
	Dotted silicone band						**Petite**			**Regular**		
I	77404-00	415-1742	-				79488-00	415-1577	-	79489-00	415-1528	-
II	77404-01	415-1759	-				79488-01	415-1585	-	79489-01	415-1536	-
III	77404-02	415-1767	-				79488-02	415-1593	-	79489-02	415-1544	-
IV	77404-03	415-1775	-				79488-03	415-1601	-	79489-03	415-1551	-
V	77404-04	415-1783	-				79488-04	415-1619	-	79489-04	415-1569	-
VI	77404-05	415-1791	-				79488-05	415-1627	-	79489-05	415-1510	-

JOBST® Classic - CCL 3

Size	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC
Knee High - Sand							**Thigh High - Sand**					
	RAL CCL 3 - Open Toe						**RAL CCL 3 - Open Toe**					
	Petite			**Regular**			**Dotted silicone band**					
I	73530-12	415-2583	-	73531-12	415-2526	-	73533-12	415-2641	-			
II	73530-13	415-2591	-	73531-13	415-2534	-	73533-13	415-2658	-			
III	73530-14	415-2609	-	73531-14	415-2542	-	73533-14	415-2666	-			
IV	73530-15	415-2617	-	73531-15	415-2559	-	73533-15	415-2674	-			
V	73530-16	415-2625	-	73531-16	415-2567	-	73533-16	415-2682	-			
VI	73530-17	415-2633	-	73531-17	415-2575	-	73533-17	415-2690	-			
Knee High - Sand												
	RAL CCL 3 - Closed Toe											
	Petite			**Regular**								
I	73530-24	-	-	73531-24	-	-						
II	73530-25	-	-	73531-25	-	-						
III	73530-26	-	-	73531-26	-	-						
IV	73530-27	-	-	73531-27	-	-						
V	73530-28	-	-	73531-28	-	-						
VI	73530-29	-	-	73531-29	-	-						
Tights - Sand							**Knee High - Pearl**					
	RAL CCL 3 - Open Toe						**RAL CCL 3 - Open Toe**					
	Petite			**Regular**			**Petite**			**Regular**		
I	73537-12	415-2880	-	73538-12	415-2765	-	73530-00	415-1866	-	73531-00	415-1809	-
II	73537-13	415-2898	-	73538-13	415-2773	-	73530-01	415-1874	-	73531-01	415-1817	-
III	73537-14	415-2906	-	73538-14	415-2781	-	73530-02	415-1882	-	73531-02	415-1825	-
IV	73537-15	415-2914	-	73538-15	415-2799	-	73530-03	415-1890	-	73531-03	415-1833	-
V	73537-16	415-2922	-	73538-16	415-2807	-	73530-04	415-1908	-	73531-04	415-1841	-
VI	73537-17	415-2930	-	73538-17	415-2815	-	73530-05	415-1916	-	73531-05	415-1858	-
Thigh High - Pearl							**Tights - Pearl**					
	RAL CCL 3 - Open Toe						**RAL CCL 3 - Open Toe**					
	Dotted silicone band						**Petite**			**Regular**		
I	73533-00	415-1924	-				73537-00	415-2823	-	73538-00	415-2708	-
II	73533-01	415-1932	-				73537-01	415-2831	-	73538-01	415-2716	-
III	73533-02	415-1940	-				73537-02	415-2849	-	73538-02	415-2724	-
IV	73533-03	415-1957	-				73537-03	415-2856	-	73538-03	415-2732	-
V	73533-04	415-1965	-				73537-04	415-2864	-	73538-04	415-2740	-
VI	73533-05	415-1973	-				73537-05	415-2872	-	73538-05	415-2757	-


RAL compression classes 2 (23-32mmHg) and 3 (34-46mmHg)							
		Sizes					
		I	II	III	IV	V	VI
Circumference in cm	**cG**	43-48	45-52	49-56	53-60	56-64	60-68
	cF	39-52	41-56	44-60	47-65	50-69	53-73
	cE	30-37	33-40	35-43	37-45	39-48	41-51
	cD	27-33	29-36	32-39	34-42	36-45	38-48
	cC	28-34	30-37	33-40	35-43	37-46	39-49
	cB1	23-27	24-29	26-32	29-35	31-37	33-39
	cB	18-20	20-22	22-24	24-26	26-28	28-30
	cY	26-31	28-33	29-35	31-37	32-38	33-40
	cA	17-22	19-24	21-26	23-29	25-32	27-34

Length measurement for all sizes:

AD
ID petite: 34-39 cm
ID regular: 39-45 cm

AG
Length: 70-80 cm

AT
IK petite: 65-75 cm
IK regular: 75-85 cm



JOBST® Bella™ Lite

Ready-to-Wear Upper Limb Compression Garments

Available on Drug Tariff

JOBST® Bella™ Lite circular-knit, seamless compression garments for the upper limb are suitable for treating mild to moderate levels of lymphoedema. These comfortable garments are available as an armsleeve or a gauntlet that can be worn together or separately, and a combined armsleeve.

Comfortable and effective

The soft, fine knit construction of JOBST® Bella™ Lite makes it particularly durable, air permeable and easy to don. JOBST® Bella™ Lite sizes fit a wide range of hand and arm measurements, making it convenient and cost effective.

Compression Levels	
CCL 1 = 15-20mmHg	CCL 2 = 20-30mmHg
✓	✓

Sizes and Options - Knitted or dotted silicone band for the armsleeve	
Six Sizes	Regular
✓	✓




Knitted top band



Dotted silicone top band

Colours	
Beige	Black

Styles		
		
CG Armsleeve	AC1 Gauntlet with Thumb	AG Combined Armsleeve

JOBST® Bella™ Lite - CCL 1

	Code	Option	Size	Length	Colour	PIP	NHS SC
Armsleeve							
	75273-02	Knitted band	Small	Regular	Beige	382-1196	EGD9103
	75273-00	Knitted band	Small	Long	Beige	325-7011	EGD9106
	75274-02	Knitted band	Medium	Regular	Beige	382-1188	EGD9102
	75274-00	Knitted band	Medium	Long	Beige	325-7037	EGD9105
	75275-02	Knitted band	Large	Regular	Beige	382-1170	EGD9101
	75275-00	Knitted band	Large	Long	Beige	325-7029	EGD9104
	75273-03	Silicone band	Small	Regular	Beige	382-1238	EGD9112
	75273-01	Silicone band	Small	Long	Beige	325-6955	EGD9109
	75274-03	Silicone band	Medium	Regular	Beige	382-1220	EGD9111
	75274-01	Silicone band	Medium	Long	Beige	325-6971	EGD9108
	75275-03	Silicone band	Large	Regular	Beige	382-1212	EGD9110
	75275-01	Silicone band	Large	Long	Beige	325-6963	EGD9107
Gauntlet							
	75279-00	With thumb	Small	-	Beige	325-6823	-
	75280-00	With thumb	Medium	-	Beige	325-6831	-
	75281-00	With thumb	Large	-	Beige	325-6922	-
Combined Armsleeve							
	77678-00	Knitted band	Small	Regular	Beige	398-9571	-
	77678-03	Knitted band	Small	Long	Beige	398-9605	-
	77678-01	Knitted band	Medium	Regular	Beige	398-9589	-
	77678-04	Knitted band	Medium	Long	Beige	398-9613	-
	77678-02	Knitted band	Large	Regular	Beige	398-9597	-
	77678-05	Knitted band	Large	Long	Beige	398-9621	-
	77678-06	Silicone band	Small	Regular	Beige	398-9639	-
	77678-09	Silicone band	Small	Long	Beige	398-9662	-
	77678-07	Silicone band	Medium	Regular	Beige	398-9647	-
	77678-10	Silicone band	Medium	Long	Beige	398-9670	-
	77678-08	Silicone band	Large	Regular	Beige	398-9654	-
	77678-11	Silicone band	Large	Long	Beige	398-9688	-
	77678-00	Knitted band	Small	Regular	Beige	398-9571	-
	77678-03	Knitted band	Small	Long	Beige	398-9605	-

Armsleeve

Size	 Wrist	 Mid-Lower	 Mid-Upper
Small	14 - 20cm	18 - 23cm	22 - 40cm
Medium	16 - 22cm	23 - 29cm	27 - 45cm
Large	18 - 24cm	29 - 34cm	33 - 51cm

Length



Note: length measurement taken on inside of arm

Regular up to (‹) 43cm | **Long** over (›) 43cm

Gauntlet

Size	 Palm	 Wrist
Small	15 - 19cm	14 - 20cm
Medium	19 - 22cm	16 - 22cm
Large	22 - 26cm	18 - 24cm

JOBST® Bella™ Lite - CCL 2

	Code	Option	Size	Length	Colour	PIP	NHS SC
Armsleeve							
	75276-05	Knitted band	Small	Regular	Beige	382-1261	EGD9115
	75276-00	Knitted band	Small	Long	Beige	325-6989	EGD9118
	75277-05	Knitted band	Medium	Regular	Beige	382-1204	EGD9114
	75277-00	Knitted band	Medium	Long	Beige	325-7003	EGD9117
	75278-05	Knitted band	Large	Regular	Beige	382-1246	EGD9113
	75278-00	Knitted band	Large	Long	Beige	325-6997	EGD9116
	75276-06	Silicone band	Small	Regular	Beige	382-1287	EGD9127
	75276-01	Silicone band	Small	Long	Beige	325-6898	EGD9124
	75277-06	Silicone band	Medium	Regular	Beige	382-1279	EGD9126
	75277-01	Silicone band	Medium	Long	Beige	325-6914	EGD9123
	75278-06	Silicone band	Large	Regular	Beige	382-1253	EGD9125
	75278-01	Silicone band	Large	Long	Beige	325-6906	EGD9122
	75276-04	Knitted band	Small	Long	Black	363-7568	EGD9121
	75277-04	Knitted band	Medium	Long	Black	363-7576	EGD9120
	75278-04	Knitted band	Large	Long	Black	363-7550	EGD9119
Gauntlet							
	75282-00	With thumb	Small	-	Beige	325-6773	-
	75283-00	With thumb	Medium	-	Beige	325-6781	-
	75284-00	With thumb	Large	-	Beige	325-6799	-
	75282-03	With thumb	Small	-	Black	363-7543	-
	75283-03	With thumb	Medium	-	Black	363-7535	-
	75284-03	With thumb	Large	-	Black	363-7527	-
Combined Armsleeve							
	77680-00	Knitted band	Small	Regular	Beige	398-9696	-
	77680-03	Knitted band	Small	Long	Beige	398-9720	-
	77680-01	Knitted band	Medium	Regular	Beige	398-9704	-
	77680-04	Knitted band	Medium	Long	Beige	398-9738	-
	77680-02	Knitted band	Large	Regular	Beige	398-9712	-
	77680-05	Knitted band	Large	Long	Beige	398-9746	-
	77680-06	Silicone band	Small	Regular	Beige	398-9753	-
	77680-09	Silicone band	Small	Long	Beige	398-9787	-
	77680-07	Silicone band	Medium	Regular	Beige	398-9761	-
	77680-10	Silicone band	Medium	Long	Beige	398-9795	-
	77680-08	Silicone band	Large	Regular	Beige	398-9779	-
	77680-11	Silicone band	Large	Long	Beige	398-9803	-

Combined Armsleeve			
			
Size	**Palm**	**Mid-Lower**	**Mid-Upper**
Small	15 - 19cm	18 - 23cm	22 - 40cm
Medium	19 - 22cm	23 - 29cm	27 - 45cm
Large	22 - 26cm	29 - 34cm	33 - 51cm

Length	
	
Note: length measurement taken on inside of arm	
Regular up to (‹) 43cm	**Long** over (›) 43cm

Gauntlet		
		
Size	**Palm**	**Wrist**
Small	15 - 19cm	14 - 20cm
Medium	19 - 22cm	16 - 22cm
Large	22 - 26cm	18 - 24cm

JOBST® forMen Ambition

Available on Drug Tariff

Ready-to-Wear RAL Compression Hosiery

JOBST® forMen is a range of ready to wear, circular-knit hosiery for the treatment of mild to moderate lymphoedema and venous disease.

Stylish and effective


Inspired by a man's dress sock the timeless ribbed design, available in a variety of colours, offers better patient compliance. The balanced mix of soft yarns and high quality cotton ensures high elasticity for ease of donning and increased comfort. For the effective management of lymphatic and vascular conditions and helps prevent the recurrence of leg ulcers.

Compression Classes	
RAL CCL 1 = 18-21mmHg	RAL CCL 2 = 23-32mmHg
✓	✓


Sizes and Options / Materials			
Six Sizes	Regular	Long	Closed Toe
✓	✓	✓	✓

Colours				
Black	Dark Grey	Brown	Navy	Khaki

Style



AD Knee High



	Measurement Point	I	II	III	IV	V	VI
	RAL Compression Classes 1 (18-21mmHg) and 2 (23-32mmHg)						
circumference in cm	cD	28 - 34	31 - 38	34 - 41	37 - 46	40 - 49	42 - 52
	cC	30 - 36	32 - 40	34 - 43	38 - 47	40 - 50	43 - 53
	cB1	23.5 - 27.5	26 - 30.5	28.5 - 33	31 - 36.5	33.5 - 39	35.5 - 41.5
	cB	18 - 20	20 - 22	22 - 24	24 - 26	26 - 28	28 - 30
	cY	26 - 31	28 - 33	29 - 35	31 - 37	32 - 38	33 - 40
	cA	17 - 22	19 - 24	21 - 26	23 - 29	25 - 32	27 - 34
	Leg length measurements for all sizes in cm: regular 39 – 45cm, long 45 – 51cm						

JOBST® forMen Ambition

Size	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC
Below Knee - Black												
	RAL CCL 1 - Closed Toe						**RAL CCL 2 - Closed Toe**					
	Long			**Regular**			**Long**			**Regular**		
I	79483-24	393-0930	-	79482-24	393-1227	-	79485-24	393-0500	-	79484-24	393-0625	-
II	79483-25	393-0948	-	79482-25	393-1235	-	79485-25	393-0518	-	79484-25	393-0633	-
III	79483-26	393-0955	-	79482-26	393-1243	-	79485-26	393-0526	-	79484-26	393-0641	-
IV	79483-27	393-0963	-	79482-27	393-1250	-	79485-27	393-0534	-	79484-27	393-0658	-
V	79483-28	393-0971	-	79482-28	393-1268	-	79485-28	393-0542	-	79484-28	393-0666	-
VI	79483-29	393-0989	-	79482-29	393-1276	-	79485-29	393-0559	-	79484-29	393-0674	-
Below Knee - Dark Grey												
	RAL CCL 1 - Closed Toe						**RAL CCL 2 - Closed Toe**					
	Long			**Regular**			**Long**			**Regular**		
I	79483-06	393-0872	-	79482-06	393-1169	-	79485-06	393-0328	-	79484-06	393-0567	-
II	79483-07	393-0880	-	79482-07	393-1177	-	79485-07	393-0336	-	79484-07	393-0575	-
III	79483-08	393-0898	-	79482-08	393-1185	-	79485-08	393-0344	-	79484-08	393-0583	-
IV	79483-09	393-0906	-	79482-09	393-1193	-	79485-09	393-0351	-	79484-09	393-0591	-
V	79483-10	393-0914	-	79482-10	393-1201	-	79485-10	393-0369	-	79484-10	393-0609	-
VI	79483-11	393-0922	-	79482-11	393-1219	-	79485-11	393-0377	-	79484-11	393-0617	-
Below Knee - Brown												
	RAL CCL 1 - Closed Toe						**RAL CCL 2 - Closed Toe**					
	Long			**Regular**			**Long**			**Regular**		
I	79483-18	393-0997	-	79482-18	393-1284	-	79485-18	393-0443	-	79484-18	393-0682	-
II	79483-19	393-1003	-	79482-19	393-1292	-	79485-19	393-0450	-	79484-19	393-0690	-
III	79483-20	393-1011	-	79482-20	393-1300	-	79485-20	393-0468	-	79484-20	393-0708	-
IV	79483-21	393-1029	-	79482-21	393-1318	-	79485-21	393-0476	-	79484-21	393-0716	-
V	79483-22	393-1037	-	79482-22	393-1326	-	79485-22	393-0484	-	79484-22	393-0724	-
VI	79483-23	393-1045	-	79482-23	393-1334	-	79485-23	393-0492	-	79484-23	393-0732	-
Below Knee - Navy												
	RAL CCL 1 - Closed Toe						**RAL CCL 2 - Closed Toe**					
	Long			**Regular**			**Long**			**Regular**		
I	79483-12	393-1102	-	79482-12	393-1409	-	79485-12	393-0385	-	79484-12	393-0807	-
II	79483-13	393-1110	-	79482-13	393-1417	-	79485-13	393-0393	-	79484-13	393-0815	-
III	79483-14	393-1128	-	79482-14	393-1425	-	79485-14	393-0401	-	79484-14	393-0823	-
IV	79483-15	393-1136	-	79482-15	393-1433	-	79485-15	393-0419	-	79484-15	393-0831	-
V	79483-16	393-1144	-	79482-16	393-1441	-	79485-16	393-0427	-	79484-16	393-0849	-
VI	79483-17	393-1151	-	79482-17	393-1458	-	79485-17	393-0435	-	79484-17	393-0856	-
Below Knee - Khaki												
	RAL CCL 1 - Closed Toe						**RAL CCL 2 - Closed Toe**					
	Long			**Regular**			**Long**			**Regular**		
I	79483-00	393-0864	-	79482-00	393-1342	-	79485-00	393-0260	-	79484-00	393-0740	-
II	79483-01	393-1052	-	79482-01	393-1359	-	79485-01	393-0278	-	79484-01	393-0757	-
III	79483-02	393-1060	-	79482-02	393-1367	-	79485-02	393-0286	-	79484-02	393-0765	-
IV	79483-03	393-1078	-	79482-03	393-1375	-	79485-03	393-0294	-	79484-03	393-0773	-
V	79483-04	393-1086	-	79482-04	393-1383	-	79485-04	393-0302	-	79484-04	393-0781	-
VI	79483-05	393-1094	-	79482-05	393-1391	-	79485-05	393-0310	-	79484-05	393-0799	-

JOBST® forMen Explore

Available on Drug Tariff

Ready-to-Wear RAL Compression Hosiery

JOBST® forMen is a range of ready to wear, circular-knit hosiery for the treatment of mild to moderate lymphoedema and venous disease.

Durable and effective


Reinforced foot zones offer increased strength and durability for active wear. The high cotton content of JOBST® forMen Explore – the classic of all clothing fabrics – ensures a natural feel against the skin and reduces the potential for skin irritations. For the effective management of lymphatic and vascular conditions and helps prevent the recurrence of leg ulcers.

Compression Classes	
RAL CCL 1 = 18-21mmHg	RAL CCL 2 = 23-32mmHg
✓	✓


Sizes and Options / Materials			
Six Sizes	Regular	Long	Closed Toe
✓	✓	✓	✓

Colours			
Black	Dark grey	Navy	Khaki

Style



AD Knee High



Measurement Point		I	II	III	IV	V	VI
RAL Compression Classes 1 (18-21mmHg) and 2 (23-32mmHg)							
circumference in cm	cD	28 - 34	31 - 38	34 - 41	37 - 46	40 - 49	42 - 52
	cC	30 - 36	32 - 40	34 - 43	38 - 47	40 - 50	43 - 53
	cB1	23.5 - 27.5	26 - 30.5	28.5 - 33	31 - 36.5	33.5 - 39	35.5 - 41.5
	cB	18 - 20	20 - 22	22 - 24	24 - 26	26 - 28	28 - 30
	cY	26 - 31	28 - 33	29 - 35	31 - 37	32 - 38	33 - 40
	cA	17 - 22	19 - 24	21 - 26	23 - 29	25 - 32	27 - 34
	Leg length measurements for all sizes in cm: regular 39 – 45cm, long 45 – 51cm						

JOBST® forMen Explore

Size	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC
Below Knee - Black												
	RAL CCL 1 - Closed Toe						**RAL CCL 2 - Closed Toe**					
	Long			**Regular**			**Long**			**Regular**		
I	79509-24	393-2464	-	79508-24	393-2704	-	79511-24	393-1987	-	79510-24	393-2225	-
II	79509-25	393-2472	-	79508-25	393-2712	-	79511-25	393-1995	-	79510-25	393-2233	-
III	79509-26	393-2480	-	79508-26	393-2720	-	79511-26	393-2001	-	79510-26	393-2241	-
IV	79509-27	393-2498	-	79508-27	393-2738	-	79511-27	393-2019	-	79510-27	393-2258	-
V	79509-28	393-2506	-	79508-28	393-2746	-	79511-28	393-2027	-	79510-28	393-2266	-
VI	79509-29	393-2514	-	79508-29	393-2753	-	79511-29	393-2035	-	79510-29	393-2274	-
Below Knee - Dark Grey												
	RAL CCL 1 - Closed Toe						**RAL CCL 2 - Closed Toe**					
	Long			**Regular**			**Long**			**Regular**		
I	79509-06	393-2407	-	79508-06	393-2647	-	79511-06	393-1920	-	79510-06	393-2167	-
II	79509-07	393-2415	-	79508-07	393-2654	-	79511-07	393-1938	-	79510-07	393-2175	-
III	79509-08	393-2423	-	79508-08	393-2662	-	79511-08	393-1946	-	79510-08	393-2183	-
IV	79509-09	393-2431	-	79508-09	393-2670	-	79511-09	393-1953	-	79510-09	393-2191	-
V	79509-10	393-2449	-	79508-10	393-2688	-	79511-10	393-1961	-	79510-10	393-2209	-
VI	79509-11	393-2456	-	79508-11	393-2696	-	79511-11	393-1979	-	79510-11	393-2217	-
Below Knee - Navy												
	RAL CCL 1 - Closed Toe						**RAL CCL 2 - Closed Toe**					
	Long			**Regular**			**Long**			**Regular**		
I	79509-12	393-2589	-	79508-12	393-2829	-	79511-12	393-2100	-	79510-12	393-2340	-
II	79509-13	393-2597	-	79508-13	393-2837	-	79511-13	393-2118	-	79510-13	393-2357	-
III	79509-14	393-2605	-	79508-14	393-2845	-	79511-14	393-2126	-	79510-14	393-2365	-
IV	79509-15	393-2613	-	79508-15	393-2852	-	79511-15	393-2134	-	79510-15	393-2373	-
V	79509-16	393-2621	-	79508-16	393-2860	-	79511-16	393-2142	-	79510-16	393-2381	-
VI	79509-17	393-2639	-	79508-17	393-2878	-	79511-17	393-2159	-	79510-17	393-2399	-
Below Knee - Khaki												
	RAL CCL 1 - Closed Toe						**RAL CCL 2 - Closed Toe**					
	Long			**Regular**			**Long**			**Regular**		
I	79509-00	393-2522	-	79508-00	393-2761	-	79511-00	393-2043	-	79510-00	393-2282	-
II	79509-01	393-2530	-	79508-01	393-2779	-	79511-01	393-2050	-	79510-01	393-2290	-
III	79509-02	393-2548	-	79508-02	393-2787	-	79511-02	393-2068	-	79510-02	393-2308	-
IV	79509-03	393-2555	-	79508-03	393-2795	-	79511-03	393-2076	-	79510-03	393-2316	-
V	79509-04	393-2563	-	79508-04	393-2803	-	79511-04	393-2084	-	79510-04	393-2324	-
VI	79509-05	393-2571	-	79508-05	393-2811	-	79511-05	393-2092	-	79510-05	393-2332	-

JOBST® Opaque

Ready-to-Wear RAL Compression Hosiery

Available on Drug Tariff

JOBST® Opaque is a range of ready-to-wear, circular-knit hosiery for the treatment of mild to moderate lymphoedema and venous disease.

Soft and comfortable

JOBST® Advanced Comfort combines high-quality yarns and a superior knit structure, to make JOBST® Opaque even softer than ever and creates unbelievable wearing comfort. The Nilit® Aquarius fibres effectively wick away moisture and the double-covered inlay yarns make JOBST® Opaque incredibly easy to don.


Compression Classes	
RAL CCL 1 = 18-21mmHg	RAL CCL 2 = 23-32mmHg
✓	✓

Sizes and Options / Materials						
Six Sizes	Regular	Petite	Wider Circumference	Made-to-Measure	Open Toe	Closed Toe
✓	✓	✓	✓	✓	✓	✓


Knitted top band



Lace silicone top band



Dotted silicone top band


Sensitive top band


Colours				
Natural	Black	Caramel	Bronze	Navy

Styles


AD Knee High


AG Thigh High


AT Tights



AG-T Chap Style


AT Men


AT Maternity

JOBST® Opaque

	Measurement Point	I	II	III	IV	V	VI
		RAL Compression Classes 1 (18-21mmHg) and 2 (23-32mmHg)					
circumference in cm	cG silicone band	43 - 48	45 - 52	49 - 56	53 - 60	56 - 64	60 - 68
	cG silicone band wide	49 - 57	53 - 62	57 - 67	61 - 72	65 - 77	69 - 81
	cG	43 - 57	45 - 62	49 - 67	53 - 72	56 - 77	60 - 81
	cF	39 - 52	41 - 56	44 - 60	47 - 65	50 - 69	53 - 73
	cE	30 - 37	33 - 40	35 - 43	37 - 45	39 - 48	41 - 51
	cD	27 - 33	29 - 36	32 - 39	34 - 42	36 - 45	38 - 48
	cC	28 - 34	30 - 37	33 - 40	35 - 43	37 - 46	39 - 49
	cB1	23 - 27	24 - 29	26 - 32	29 - 35	31 - 37	33 - 39
	cB	18 - 20	20 - 22	22 - 24	24 - 26	26 - 28	28 - 30
	cY	26 - 31	28 - 33	29 - 35	31 - 37	32 - 38	33 - 40
	cA	17 - 22	19 - 24	21 - 26	23 - 29	25 - 32	27 - 34
	Length measurements for all sizes (in cm): AD: ID regular 39 - 45, ID petite 34 - 39 AG: IG regular 70 - 80, IG petite 60 - 69 AT / AG-T: IK regular 75 - 85, IK petite 65 - 75						





JOBST® Opaque - CCL 1

Size	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC
Knee High - Natural												
	RAL CCL 1 - Open Toe						**RAL CCL 1 - Closed Toe**					
	Petite			**Regular**			**Petite**			**Regular**		
I	70616-48	326-1757	-	70615-48	326-1633	-	70616-18	326-2052	-	70615-18	326-1930	-
II	70616-49	326-1765	-	70615-49	326-1641	-	70616-19	326-2060	-	70615-19	326-1948	-
III	70616-50	326-1773	-	70615-50	326-1658	-	70616-20	326-2078	-	70615-20	326-1955	-
IV	70616-51	326-1781	-	70615-51	326-1666	-	70616-21	326-2086	-	70615-21	326-1963	-
V	70616-52	326-1799	-	70615-52	326-1674	-	70616-22	326-2094	-	70615-22	326-1971	-
VI	70616-53	326-1864	-	70615-53	326-1682	-	70616-23	326-2102	-	70615-23	326-1989	-
Knee High - Natural												
	RAL CCL 1 - Open Toe, Dotted Silicone Band						**RAL CCL 1 - Closed Toe, Dotted Silicone Band**					
	Petite			**Regular**			**Petite**			**Regular**		
I	70616-84	-	-	70615-84	-	-	70616-108	-	-	70615-108	-	-
II	70616-85	-	-	70615-85	-	-	70616-109	-	-	70615-109	-	-
III	70616-86	-	-	70615-86	-	-	70616-110	-	-	70615-110	-	-
IV	70616-87	-	-	70615-87	-	-	70616-111	-	-	70615-111	-	-
V	70616-88	-	-	70615-88	-	-	70616-112	-	-	70615-112	-	-
VI	70616-89	-	-	70615-89	-	-	70616-113	-	-	70615-113	-	-
Knee High - Black												
	RAL CCL 1 - Open Toe						**RAL CCL 1 - Closed Toe**					
	Petite			**Regular**			**Petite**			**Regular**		
I	70616-42	326-1807	-	70615-42	326-1690	-	70616-24	326-1997	-	70615-24	326-1872	-
II	70616-43	326-1815	-	70615-43	326-1708	-	70616-25	326-2003	-	70615-25	326-1880	-
III	70616-44	326-1823	-	70615-44	326-1716	-	70616-26	326-2011	-	70615-26	326-1898	-
IV	70616-45	326-1831	-	70615-45	326-1724	-	70616-27	326-2029	-	70615-27	326-1906	-
V	70616-46	326-1849	-	70615-46	326-1732	-	70616-28	326-2037	-	70615-28	326-1914	-
VI	70616-47	326-1856	-	70615-47	326-1740	-	70616-29	326-2045	-	70615-29	326-1922	-
Knee High - Black												
	RAL CCL 1 - Open Toe, Dotted Silicone Band						**RAL CCL 1 - Closed Toe, Dotted Silicone Band**					
	Petite			**Regular**			**Petite**			**Regular**		
I	70616-78	-	-	70615-78	-	-	70616-114	-	-	70615-114	-	-
II	70616-79	-	-	70615-79	-	-	70616-115	-	-	70615-115	-	-
III	70616-80	-	-	70615-80	-	-	70616-116	-	-	70615-116	-	-
IV	70616-81	-	-	70615-81	-	-	70616-117	-	-	70615-117	-	-
V	70616-82	-	-	70615-82	-	-	70616-118	-	-	70615-118	-	-
VI	70616-83	-	-	70615-83	-	-	70616-119	-	-	70615-119	-	-
Knee High - Navy												
	RAL CCL 1 - Open Toe						**RAL CCL 1 - Closed Toe**					
	Petite			**Regular**			**Petite**			**Regular**		
I	70616-54	360-7256	-	70615-54	360-7314	-	70616-12	360-7132	-	70615-12	360-7199	-
II	70616-55	360-7264	-	70615-55	360-7322	-	70616-13	360-7140	-	70615-13	360-7207	-
III	70616-56	360-7272	-	70615-56	360-7330	-	70616-14	360-7157	-	70615-14	360-7215	-
IV	70616-57	360-7280	-	70615-57	360-7348	-	70616-15	360-7165	-	70615-15	360-7223	-
V	70616-58	360-7298	-	70615-58	360-7355	-	70616-16	360-7173	-	70615-16	360-7231	-
VI	70616-59	360-7306	-	70615-59	360-7363	-	70616-17	360-7181	-	70615-17	360-7249	-

JOBST® Opaque - CCL 1

Size	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC
Knee High - Navy												
	RAL CCL 1 - Open Toe, Dotted Silicone Band						**RAL CCL 1 - Closed Toe, Dotted Silicone Band**					
	Petite			**Regular**			**Petite**			**Regular**		
I	70616-72	-	-	70615-72	-	-	70616-102	-	-	70615-102	-	-
II	70616-73	-	-	70615-73	-	-	70616-103	-	-	70615-103	-	-
III	70616-74	-	-	70615-74	-	-	70616-104	-	-	70615-104	-	-
IV	70616-75	-	-	70615-75	-	-	70616-105	-	-	70615-105	-	-
V	70616-76	-	-	70615-76	-	-	70616-106	-	-	70615-106	-	-
VI	70616-77	-	-	70615-77	-	-	70616-107	-	-	70615-107	-	-
Knee High - Bronze												
	RAL CCL 1 - Open Toe						**RAL CCL 1 - Closed Toe**					
	Petite			**Regular**			**Petite**			**Regular**		
I	70616-120	380-7021	-	70615-120	380-6478	-	70616-126	380-7088	-	70615-126	380-6544	-
II	70616-121	380-7039	-	70615-121	380-6486	-	70616-127	380-7096	-	70615-127	380-6551	-
III	70616-122	380-7047	-	70615-122	380-6494	-	70616-128	380-7104	-	70615-128	380-6569	-
IV	70616-123	380-7054	-	70615-123	380-6502	-	70616-129	380-7112	-	70615-129	380-6577	-
V	70616-124	380-7062	-	70615-124	380-6510	-	70616-130	380-7120	-	70615-130	380-6585	-
VI	70616-125	380-7070	-	70615-125	380-6528	-	70616-131	380-7138	-	70615-131	380-6593	-
Knee High - Caramel												
	RAL CCL 1 - Open Toe						**RAL CCL 1 - Closed Toe**					
	Petite			**Regular**			**Petite**			**Regular**		
I	70616-36	384-0691	-	70615-36	384-0451	-	70616-00	384-0766	-	70615-00	384-0519	-
II	70616-37	384-0709	-	70615-37	384-0469	-	70616-01	384-0774	-	70615-01	384-0527	-
III	70616-38	384-0717	-	70615-38	384-0477	-	70616-02	384-0782	-	70615-02	384-0535	-
IV	70616-39	384-0725	-	70615-39	384-0485	-	70616-03	384-0790	-	70615-03	384-0543	-
V	70616-40	384-0733	-	70615-40	384-0493	-	70616-04	384-0808	-	70615-04	384-0550	-
VI	70616-41	384-0741	-	70615-41	384-0501	-	70616-05	384-0816	-	70615-05	384-0568	-
Thigh High - Natural												
	RAL CCL 1 - Open Toe, Dotted Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70618-252	363-8293	-	70618-384	363-8418	-	70617-255	363-8350	-	70617-387	363-8475	-
II	70618-253	363-8301	-	70618-385	363-8426	-	70617-256	363-8368	-	70617-388	363-8483	-
III	70618-254	363-8319	-	70618-386	363-8434	-	70617-257	363-8376	-	70617-389	363-8491	-
IV	70618-255	363-8327	-	70618-387	363-8442	-	70617-258	363-8384	-	70617-390	363-8509	-
V	70618-256	363-8335	-	70618-388	363-8459	-	70617-259	363-8392	-	70617-391	363-8517	-
VI	70618-257	363-8343	-	70618-389	363-8467	-	70617-260	363-8400	-	70617-392	363-8525	-
Thigh High - Natural												
	RAL CCL 1 - Closed Toe, Dotted Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70618-276	363-8533	-	70618-408	363-8657	-	70617-279	363-8590	-	70617-411	363-8715	-
II	70618-277	363-8541	-	70618-409	363-8665	-	70617-280	363-8608	-	70617-412	363-8723	-
III	70618-278	363-8558	-	70618-410	363-8673	-	70617-281	363-8616	-	70617-413	363-8731	-
IV	70618-279	363-8566	-	70618-411	363-8681	-	70617-282	363-8624	-	70617-414	363-8749	-
V	70618-280	363-8574	-	70618-412	363-8699	-	70617-283	363-8632	-	70617-415	363-8756	-
VI	70618-281	363-8582	-	70618-413	363-8707	-	70617-284	363-8640	-	70617-416	363-8764	-

JOBST® Opaque - CCL 1

Size	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC
Thigh High - Natural												
	RAL CCL 1 - Open Toe, Sensitive Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70618-582	326-9099	-	70618-630	326-9164	-	70617-561	326-8927	-	70617-609	326-9024	-
II	70618-583	326-9107	-	70618-631	326-9172	-	70617-562	326-8935	-	70617-610	326-9032	-
III	70618-584	326-9115	-	70618-632	326-9180	-	70617-563	326-8943	-	70617-611	326-9040	-
IV	70618-585	326-9131	-	70618-633	326-9198	-	70617-564	326-8950	-	70617-612	326-9057	-
V	70618-586	326-9149	-	70618-634	326-9206	-	70617-565	326-8968	-	70617-613	326-9065	-
VI	70618-587	326-9156	-	70618-635	326-9214	-	70617-566	326-8976	-	70617-614	326-9073	-
Thigh High - Natural												
	RAL CCL 1 - Closed Toe, Sensitive Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70618-600	357-1791	-	70618-648	357-2054	-	70617-579	357-1924	-	70617-627	357-2435	-
II	70618-601	357-1809	-	70618-649	357-2070	-	70617-580	357-1932	-	70617-628	357-2443	-
III	70618-602	357-1817	-	70618-650	357-2088	-	70617-581	357-1940	-	70617-629	357-2450	-
IV	70618-603	357-1825	-	70618-651	357-2104	-	70617-582	357-1957	-	70617-630	357-2468	-
V	70618-604	357-1833	-	70618-652	357-2112	-	70617-583	357-1965	-	70617-631	357-2476	-
VI	70618-605	357-1841	-	70618-653	357-2120	-	70617-584	357-1973	-	70617-632	357-2484	-
Thigh High - Natural												
	RAL CCL 1 - Open Toe, Lace Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70618-192	357-3896	-	70618-210	357-4316	-	70617-195	357-4258	-	70617-213	357-4373	-
II	70618-193	357-3904	-	70618-211	357-4324	-	70617-196	357-4266	-	70617-214	357-4381	-
III	70618-194	357-3912	-	70618-212	357-4332	-	70617-197	357-4274	-	70617-215	357-4399	-
IV	70618-195	357-3920	-	70618-213	357-4340	-	70617-198	357-4282	-	70617-216	357-4407	-
V	70618-196	357-3938	-	70618-214	357-4357	-	70617-199	357-4290	-	70617-217	357-4415	-
VI	70618-197	357-3946	-	70618-215	357-4365	-	70617-200	357-4308	-	70617-218	357-4423	-
Thigh High - Natural												
	RAL CCL 1 - Closed Toe, Lace Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70618-42	357-2518	-	70618-126	357-3078	-	70617-42	357-2658	-	70617-126	357-3201	-
II	70618-43	357-2526	-	70618-127	357-3094	-	70617-43	357-2666	-	70617-127	357-3219	-
III	70618-44	357-2542	-	70618-128	357-3102	-	70617-44	357-2674	-	70617-128	357-3235	-
IV	70618-45	357-2559	-	70618-129	357-3128	-	70617-45	357-2690	-	70617-129	357-3250	-
V	70618-46	357-2567	-	70618-130	357-3144	-	70617-46	357-2708	-	70617-130	357-3268	-
VI	70618-47	357-2575	-	70618-131	357-3169	-	70617-47	357-2724	-	70617-131	357-3284	-
Thigh High - Black												
	RAL CCL 1 - Open Toe, Dotted Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70618-246	363-8772	-	70618-378	363-8897	-	70617-249	363-8830	-	70617-381	363-8954	-
II	70618-247	363-8780	-	70618-379	363-8905	-	70617-250	363-8848	-	70617-382	363-8962	-
III	70618-248	363-8798	-	70618-380	363-8913	-	70617-251	363-8855	-	70617-383	363-8970	-
IV	70618-249	363-8806	-	70618-381	363-8921	-	70617-252	363-8863	-	70617-384	363-8988	-
V	70618-250	363-8814	-	70618-382	363-8939	-	70617-253	363-8871	-	70617-385	363-8996	-
VI	70618-251	363-8822	-	70618-383	363-8947	-	70617-254	363-8889	-	70617-386	363-9002	-

JOBST® Opaque - CCL 1

Size	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC
Thigh High - Black												
	RAL CCL 1 - Closed Toe, Dotted Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70618-282	363-9010	-	70618-414	363-9762	-	70617-285	363-9077	-	70617-417	363-9846	-
II	70618-283	363-9028	-	70618-415	363-9770	-	70617-286	363-9085	-	70617-418	363-9853	-
III	70618-284	363-9036	-	70618-416	363-9788	-	70617-287	363-9093	-	70617-419	363-9861	-
IV	70618-285	363-9044	-	70618-417	363-9796	-	70617-288	363-9101	-	70617-420	363-9879	-
V	70618-286	363-9051	-	70618-418	363-9804	-	70617-289	363-9119	-	70617-421	363-9887	-
VI	70618-287	363-9069	-	70618-419	363-9812	-	70617-290	363-9127	-	70617-422	363-9895	-
Thigh High - Black												
	RAL CCL 1 - Open Toe, Sensitive Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70618-576	327-5492	-	70618-624	327-5351	-	70617-555	327-5773	-	70617-603	327-5856	-
II	70618-577	327-5658	-	70618-625	327-5377	-	70617-556	327-5781	-	70617-604	327-5880	-
III	70618-578	327-5674	-	70618-626	327-5385	-	70617-557	327-5807	-	70617-605	327-5898	-
IV	70618-579	327-5708	-	70618-627	327-5393	-	70617-558	327-5823	-	70617-606	327-5914	-
V	70618-580	327-5724	-	70618-628	327-5401	-	70617-559	327-5831	-	70617-607	327-5930	-
VI	70618-581	327-5732	-	70618-629	327-5427	-	70617-560	327-5849	-	70617-608	327-5948	-
Thigh High - Black												
	RAL CCL 1 - Closed Toe, Sensitive Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70618-606	357-1866	-	70618-654	357-2146	-	70617-585	357-1981	-	70617-633	357-2237	-
II	70618-607	357-1874	-	70618-655	357-2161	-	70617-586	357-1999	-	70617-634	357-2245	-
III	70618-608	357-1882	-	70618-656	357-2179	-	70617-587	357-2005	-	70617-635	357-2260	-
IV	70618-609	357-1890	-	70618-657	357-2195	-	70617-588	357-2013	-	70617-636	357-2278	-
V	70618-610	357-1908	-	70618-658	357-2203	-	70617-589	357-2021	-	70617-637	357-2294	-
VI	70618-611	357-1916	-	70618-659	357-2211	-	70617-590	357-2039	-	70617-638	357-2310	-
Thigh High - Black												
	RAL CCL 1 - Open Toe, Lace Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70618-186	357-3557	-	70618-204	357-4084	-	70617-189	357-3953	-	70617-207	357-4027	-
II	70618-187	357-3581	-	70618-205	357-4092	-	70617-190	357-3961	-	70617-208	357-4035	-
III	70618-188	357-3599	-	70618-206	357-4100	-	70617-191	357-3979	-	70617-209	357-4043	-
IV	70618-189	357-3615	-	70618-207	357-4118	-	70617-192	357-3987	-	70617-210	357-4050	-
V	70618-190	357-3631	-	70618-208	357-4126	-	70617-193	357-3995	-	70617-211	357-4068	-
VI	70618-191	357-3649	-	70618-209	357-4134	-	70617-194	357-4019	-	70617-212	357-4076	-
Thigh High - Black												
	RAL CCL 1 - Closed Toe, Lace Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70618-54	357-2591	-	70618-138	357-3300	-	70617-54	357-2732	-	70617-138	357-3441	-
II	70618-55	357-2609	-	70618-139	357-3326	-	70617-55	357-2757	-	70617-139	357-3458	-
III	70618-56	357-2617	-	70618-140	357-3334	-	70617-56	357-2773	-	70617-140	357-3466	-
IV	70618-57	357-2625	-	70618-141	357-3342	-	70617-57	357-2799	-	70617-141	357-3474	-
V	70618-58	357-2633	-	70618-142	357-3409	-	70617-58	357-2815	-	70617-142	357-3482	-
VI	70618-59	357-2641	-	70618-143	357-3433	-	70617-59	357-2823	-	70617-143	357-3490	-

JOBST® Opaque - CCL 1

Size	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC
Thigh High - Navy												
	RAL CCL 1 - Open Toe, Dotted Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70618-240	363-9903	-	70618-372	364-0026	-	70617-243	363-9960	-	70617-375	364-0083	-
II	70618-241	363-9911	-	70618-373	364-0034	-	70617-244	363-9978	-	70617-376	364-0091	-
III	70618-242	363-9929	-	70618-374	364-0042	-	70617-245	363-9986	-	70617-377	364-0109	-
IV	70618-243	363-9937	-	70618-375	364-0059	-	70617-246	363-9994	-	70617-378	364-0117	-
V	70618-244	363-9945	-	70618-376	364-0067	-	70617-247	363-0000	-	70617-379	364-0125	-
VI	70618-245	363-9952	-	70618-377	364-0075	-	70617-248	363-0018	-	70617-380	364-0133	-
Thigh High - Navy												
	RAL CCL 1 - Closed Toe, Dotted Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70618-270	364-0141	-	70618-402	364-0265	-	70617-273	364-0208	-	70617-405	364-0323	-
II	70618-271	364-0158	-	70618-403	364-0273	-	70617-274	364-0216	-	70617-406	364-0331	-
III	70618-272	364-0166	-	70618-404	364-0281	-	70617-275	364-0224	-	70617-407	364-0349	-
IV	70618-273	364-0174	-	70618-405	364-0299	-	70617-276	364-0232	-	70617-408	364-0356	-
V	70618-274	364-0182	-	70618-406	364-0307	-	70617-277	364-0240	-	70617-409	364-0364	-
VI	70618-275	364-0190	-	70618-407	364-0315	-	70617-278	364-0257	-	70617-410	364-0372	-
Thigh High - Navy												
	RAL CCL 1 - Open Toe, Sensitive Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70618-570	360-8593	-	70618-618	360-8650	-	70617-549	360-8718	-	70617-597	360-8775	-
II	70618-571	360-8601	-	70618-619	360-8668	-	70617-550	360-8726	-	70617-598	360-8783	-
III	70618-572	360-8619	-	70618-620	360-8676	-	70617-551	360-8734	-	70617-599	360-8791	-
IV	70618-573	360-8627	-	70618-621	360-8684	-	70617-552	360-8742	-	70617-600	360-8809	-
V	70618-574	360-8635	-	70618-622	360-8692	-	70617-553	360-8759	-	70617-601	360-8817	-
VI	70618-575	360-8643	-	70618-623	360-8700	-	70617-554	360-8767	-	70617-602	360-8825	-
Thigh High - Navy												
	RAL CCL 1 - Closed Toe, Sensitive Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70618-594	360-7611	-	70618-642	360-7678	-	70617-573	360-7736	-	70617-621	360-7793	-
II	70618-595	360-7629	-	70618-643	360-7686	-	70617-574	360-7744	-	70617-622	360-7801	-
III	70618-596	360-7637	-	70618-644	360-7694	-	70617-575	360-7751	-	70617-623	360-7819	-
IV	70618-597	360-7645	-	70618-645	360-7702	-	70617-576	360-7769	-	70617-624	360-7827	-
V	70618-598	360-7652	-	70618-646	360-7710	-	70617-577	360-7777	-	70617-625	360-7835	-
VI	70618-599	360-7660	-	70618-647	360-7728	-	70617-578	360-7785	-	70617-626	360-7843	-
Thigh High - Navy												
	RAL CCL 1 - Open Toe, Lace Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70618-222	360-8833	-	70618-354	360-8890	-	70617-225	360-8957	-	70617-357	360-9013	-
II	70618-223	360-8841	-	70618-355	360-8908	-	70617-226	360-8965	-	70617-358	360-9021	-
III	70618-224	360-8858	-	70618-356	360-8916	-	70617-227	360-8973	-	70617-359	360-9039	-
IV	70618-225	360-8866	-	70618-357	360-8924	-	70617-228	360-8981	-	70617-360	360-9047	-
V	70618-226	360-8874	-	70618-358	360-8932	-	70617-229	360-8999	-	70617-361	360-9054	-
VI	70618-227	360-8882	-	70618-359	360-8940	-	70617-230	360-9005	-	70617-362	360-9062	-

JOBST® Opaque - CCL 1

Size	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC
Thigh High - Navy												
	RAL CCL 1 - Closed Toe, Lace Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70618-30	360-7850	-	70618-114	360-7918	-	70617-30	360-7975	-	70617-114	360-8031	-
II	70618-31	360-7868	-	70618-115	360-7926	-	70617-31	360-7983	-	70617-115	360-8049	-
III	70618-32	360-7876	-	70618-116	360-7934	-	70617-32	360-7991	-	70617-116	360-8056	-
IV	70618-33	360-7884	-	70618-117	360-7942	-	70617-33	360-8007	-	70617-117	360-8064	-
V	70618-34	360-7892	-	70618-118	360-7959	-	70617-34	360-8015	-	70617-118	360-8072	-
VI	70618-35	360-7900	-	70618-119	360-7967	-	70617-35	360-8023	-	70617-119	360-8080	-
Thigh High - Bronze												
	RAL CCL 1 - Open Toe, Dotted Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70618-444	383-8711	-	70618-492	383-9073	-	70617-447	383-7069	-	70617-495	383-7630	-
II	70618-445	383-8729	-	70618-493	383-9081	-	70617-448	383-7077	-	70617-496	383-7648	-
III	70618-446	383-8737	-	70618-494	383-9099	-	70617-449	383-7085	-	70617-497	383-7655	-
IV	70618-447	383-8745	-	70618-495	383-9107	-	70617-450	383-7093	-	70617-498	383-7663	-
V	70618-448	383-8752	-	70618-496	383-9115	-	70617-451	383-7101	-	70617-499	383-7671	-
VI	70618-449	383-8760	-	70618-497	383-9123	-	70617-452	383-7119	-	70617-500	383-7689	-
Thigh High - Bronze												
	RAL CCL 1 - Closed Toe, Dotted Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70618-450	383-8778	-	70618-498	383-9131	-	70617-453	383-7135	-	70617-501	383-7697	-
II	70618-451	383-8786	-	70618-499	383-9149	-	70617-454	383-7143	-	70617-502	383-7705	-
III	70618-452	383-8794	-	70618-500	383-9156	-	70617-455	383-7150	-	70617-503	383-7713	-
IV	70618-453	383-8802	-	70618-501	383-9164	-	70617-456	383-7168	-	70617-504	383-7721	-
V	70618-454	383-8810	-	70618-502	383-9172	-	70617-457	383-7176	-	70617-505	383-7739	-
VI	70618-455	383-8828	-	70618-503	383-9180	-	70617-458	383-7184	-	70617-506	383-7747	-
Thigh High - Bronze												
	RAL CCL 1 - Open Toe, Sensitive Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70618-660	383-8471	-	70618-672	383-8836	-	70617-639	380-7278	-	70617-651	383-7192	-
II	70618-661	383-8489	-	70618-673	383-8844	-	70617-640	380-7286	-	70617-652	383-7200	-
III	70618-662	383-8497	-	70618-674	383-8851	-	70617-641	380-7294	-	70617-653	383-7218	-
IV	70618-663	383-8505	-	70618-675	383-8869	-	70617-642	380-7302	-	70617-654	383-7226	-
V	70618-664	383-8513	-	70618-676	383-8877	-	70617-643	380-7310	-	70617-655	383-7234	-
VI	70618-665	383-8521	-	70618-677	383-8885	-	70617-644	380-7328	-	70617-656	383-7242	-
Thigh High - Bronze												
	RAL CCL 1 - Closed Toe, Sensitive Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70618-666	383-8539	-	70618-678	383-8893	-	70617-645	380-7336	-	70617-657	383-7259	-
II	70618-667	383-8547	-	70618-679	383-8901	-	70617-646	380-7344	-	70617-658	383-7267	-
III	70618-668	383-8554	-	70618-680	383-8919	-	70617-647	380-7351	-	70617-659	383-7275	-
IV	70618-669	383-8562	-	70618-681	383-8927	-	70617-648	380-7369	-	70617-660	383-7283	-
V	70618-670	383-8570	-	70618-682	383-8935	-	70617-649	380-7377	-	70617-661	383-7291	-
VI	70618-671	383-8588	-	70618-683	383-8943	-	70617-650	380-7385	-	70617-662	383-7309	-

JOBST® Opaque - CCL 1

Size	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC
Thigh High - Bronze												
	RAL CCL 1 - Open Toe, Lace Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70618-432	383-8596	-	70618-480	383-8950	-	70617-435	380-7393	-	70617-483	383-7333	-
II	70618-433	383-8604	-	70618-481	383-8968	-	70617-436	380-7401	-	70617-484	383-7341	-
III	70618-434	383-8612	-	70618-482	383-8976	-	70617-437	380-7419	-	70617-485	383-7358	-
IV	70618-435	383-8620	-	70618-483	383-8984	-	70617-438	380-7427	-	70617-486	383-7366	-
V	70618-436	383-8638	-	70618-484	383-8992	-	70617-439	380-7435	-	70617-487	383-7374	-
VI	70618-437	383-8646	-	70618-485	383-9008	-	70617-440	380-7443	-	70617-488	383-7382	-
Thigh High - Bronze												
	RAL CCL 1 - Closed Toe, Lace Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70618-438	383-8653	-	70618-486	383-9016	-	70617-441	380-7450	-	70617-489	383-7390	-
II	70618-439	383-8661	-	70618-487	383-9024	-	70617-442	380-7468	-	70617-490	383-7408	-
III	70618-440	383-8679	-	70618-488	383-9032	-	70617-443	380-7476	-	70617-491	383-7416	-
IV	70618-441	383-8687	-	70618-489	383-9040	-	70617-444	380-7484	-	70617-492	383-7424	-
V	70618-442	383-8695	-	70618-490	383-9057	-	70617-445	380-7492	-	70617-493	383-7432	-
VI	70618-443	383-8703	-	70618-491	383-9065	-	70617-446	380-7500	-	70617-494	383-7440	-
Thigh High - Caramel												
	RAL CCL 1 - Open Toe, Dotted Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70618-228	384-2903	-	70618-360	384-3265	-	70617-231	384-1186	-	70617-363	384-1640	-
II	70618-229	384-2911	-	70618-361	384-3273	-	70617-232	384-1194	-	70617-364	384-1657	-
III	70618-230	384-2929	-	70618-362	384-3281	-	70617-233	384-1202	-	70617-365	384-1665	-
IV	70618-231	384-2937	-	70618-363	384-3299	-	70617-234	384-1210	-	70617-366	384-1673	-
V	70618-232	384-2945	-	70618-364	384-3307	-	70617-235	384-1228	-	70617-367	384-1681	-
VI	70618-233	384-2952	-	70618-365	384-3315	-	70617-236	384-1236	-	70617-368	384-1699	-
Thigh High - Caramel												
	RAL CCL 1 - Closed Toe, Dotted Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70618-258	384-2960	-	70618-390	384-3323	-	70617-261	384-1244	-	70617-393	384-1731	-
II	70618-259	384-2978	-	70618-391	384-3331	-	70617-262	384-1251	-	70617-394	384-1749	-
III	70618-260	384-2986	-	70618-392	384-3349	-	70617-263	384-1269	-	70617-395	384-1756	-
IV	70618-261	384-2994	-	70618-393	384-3356	-	70617-264	384-1277	-	70617-396	384-1764	-
V	70618-262	384-3000	-	70618-394	384-3364	-	70617-265	384-1285	-	70617-397	384-1772	-
VI	70618-263	384-3018	-	70618-395	384-3372	-	70617-266	384-1293	-	70617-398	384-1780	-
Thigh High - Caramel												
	RAL CCL 1 - Open Toe, Sensitive Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70618-564	384-0949	-	70618-612	384-3026	-	70617-543	384-2663	-	70617-591	384-1335	-
II	70618-565	384-0956	-	70618-613	384-3034	-	70617-544	384-2671	-	70617-592	384-1343	-
III	70618-566	384-0964	-	70618-614	384-3042	-	70617-545	384-2689	-	70617-593	384-1350	-
IV	70618-567	384-0972	-	70618-615	384-3059	-	70617-546	384-2697	-	70617-594	384-1368	-
V	70618-568	384-0980	-	70618-616	384-3067	-	70617-547	384-2705	-	70617-595	384-1376	-
VI	70618-569	384-0998	-	70618-617	384-3075	-	70617-548	384-2713	-	70617-596	384-1384	-

JOBST® Opaque - CCL 1

Size	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC
Thigh High - Caramel												
	RAL CCL 1 - Closed Toe, Sensitive Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70618-588	384-2721	-	70618-636	384-3083	-	70617-567	384-1004	-	70617-615	384-1426	-
II	70618-589	384-2739	-	70618-637	384-3091	-	70617-568	384-1012	-	70617-616	384-1434	-
III	70618-590	384-2747	-	70618-638	384-3109	-	70617-569	384-1020	-	70617-617	384-1442	-
IV	70618-591	384-2754	-	70618-639	384-3117	-	70617-570	384-1038	-	70617-618	384-1459	-
V	70618-592	384-2762	-	70618-640	384-3125	-	70617-571	384-1046	-	70617-619	384-1467	-
VI	70618-593	384-2770	-	70618-641	384-3133	-	70617-572	384-1053	-	70617-620	384-1475	-
Thigh High - Caramel												
	RAL CCL 1 - Open Toe, Lace Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70618-180	384-2788	-	70618-198	384-3141	-	70617-183	384-1061	-	70617-201	384-1491	-
II	70618-181	384-2796	-	70618-199	384-3158	-	70617-184	384-1079	-	70617-202	384-1509	-
III	70618-182	384-2804	-	70618-200	384-3166	-	70617-185	384-1087	-	70617-203	384-1517	-
IV	70618-183	384-2812	-	70618-201	384-3174	-	70617-186	384-1095	-	70617-204	384-1525	-
V	70618-184	384-2820	-	70618-202	384-3182	-	70617-187	384-1103	-	70617-205	384-1533	-
VI	70618-185	384-2838	-	70618-203	384-3190	-	70617-188	384-1111	-	70617-206	384-1541	-
Thigh High - Caramel												
	RAL CCL 1 - Closed Toe, Lace Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70618-06	384-2846	-	70618-90	384-3208	-	70617-06	384-1129	-	70617-90	384-1574	-
II	70618-07	384-2853	-	70618-91	384-3216	-	70617-07	384-1137	-	70617-91	384-1582	-
III	70618-08	384-2861	-	70618-92	384-3224	-	70617-08	384-1145	-	70617-92	384-1590	-
IV	70618-09	384-2879	-	70618-93	384-3232	-	70617-09	384-1152	-	70617-93	384-1608	-
V	70618-10	384-2887	-	70618-94	384-3240	-	70617-10	384-1160	-	70617-94	384-1616	-
VI	70618-11	384-2895	-	70618-95	384-3257	-	70617-11	384-1178	-	70617-95	384-1624	-
Tights - Natural												
	RAL CCL 1 - Open Toe						**RAL CCL 1 - Closed Toe**					
	Petite			**Regular**			**Petite**			**Regular**		
I	70622-84	357-2898	-	70621-84	357-2955	-	70622-36	357-2682	-	70621-36	357-2831	-
II	70622-85	357-2906	-	70621-85	357-2963	-	70622-37	357-2716	-	70621-37	357-2849	-
III	70622-86	357-2914	-	70621-86	357-2971	-	70622-38	357-2740	-	70621-38	357-2856	-
IV	70622-87	357-2922	-	70621-87	357-2989	-	70622-39	357-2765	-	70621-39	357-2864	-
V	70622-88	357-2930	-	70621-88	357-2997	-	70622-40	357-2781	-	70621-40	357-2872	-
VI	70622-89	357-2948	-	70621-89	357-3003	-	70622-41	357-2807	-	70621-41	357-2880	-
Tights - Black												
	RAL CCL 1 - Open Toe						**RAL CCL 1 - Closed Toe**					
	Petite			**Regular**			**Petite**			**Regular**		
I	70622-78	357-3763	-	70621-78	357-3839	-	70622-48	357-3508	-	70621-48	357-3656	-
II	70622-79	357-3789	-	70621-79	357-3847	-	70622-49	357-3516	-	70621-49	357-3672	-
III	70622-80	357-3797	-	70621-80	357-3854	-	70622-50	357-3524	-	70621-50	357-3722	-
IV	70622-81	357-3805	-	70621-81	357-3862	-	70622-51	357-3540	-	70621-51	357-3730	-
V	70622-82	357-3813	-	70621-82	357-3870	-	70622-52	357-3573	-	70621-52	357-3748	-
VI	70622-83	357-3821	-	70621-83	357-3888	-	70622-53	357-3607	-	70621-53	357-3755	-

JOBST® Opaque - CCL 1

Size	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC
Tights - Navy												
	RAL CCL 1 - Open Toe						**RAL CCL 1 - Closed Toe**					
	Petite			**Regular**			**Petite**			**Regular**		
I	70622-108	361-1597	-	70621-108	361-1654	-	70622-24	361-1472	-	70621-24	361-1530	-
II	70622-109	361-1605	-	70621-109	361-1662	-	70622-25	361-1480	-	70621-25	361-1548	-
III	70622-110	361-1613	-	70621-110	361-1670	-	70622-26	361-1498	-	70621-26	361-1555	-
IV	70622-111	361-1621	-	70621-111	361-1688	-	70622-27	361-1506	-	70621-27	361-1563	-
V	70622-112	361-1639	-	70621-112	361-1696	-	70622-28	361-1514	-	70621-28	361-1571	-
VI	70622-113	361-1647	-	70621-113	361-1704	-	70622-29	361-1522	-	70621-29	361-1589	-
Tights - Bronze												
	RAL CCL 1 - Open Toe						**RAL CCL 1 - Closed Toe**					
	Petite			**Regular**			**Petite**			**Regular**		
I	70622-240	384-0337	-	70621-240	383-9933	-	70622-246	384-0394	-	70621-246	383-9990	-
II	70622-241	384-0345	-	70621-241	383-9941	-	70622-247	384-0402	-	70621-247	384-0006	-
III	70622-242	384-0352	-	70621-242	383-9958	-	70622-248	384-0410	-	70621-248	384-0014	-
IV	70622-243	384-0360	-	70621-243	383-9966	-	70622-249	384-0428	-	70621-249	384-0022	-
V	70622-244	384-0378	-	70621-244	383-9974	-	70622-250	384-0436	-	70621-250	384-0030	-
VI	70622-245	384-0386	-	70621-245	383-9982	-	70622-251	384-0444	-	70621-251	384-0048	-
Tights - Caramel												
	RAL CCL 1 - Open Toe						**RAL CCL 1 - Closed Toe**					
	Petite			**Regular**			**Petite**			**Regular**		
I	70622-72	384-4610	-	70621-72	384-4289	-	70622-00	384-4701	-	70621-00	384-4354	-
II	70622-73	384-4628	-	70621-73	384-4297	-	70622-01	384-4719	-	70621-01	384-4362	-
III	70622-74	384-4636	-	70621-74	384-4305	-	70622-02	384-4727	-	70621-02	384-4370	-
IV	70622-75	384-4644	-	70621-75	384-4313	-	70622-03	384-4735	-	70621-03	384-4388	-
V	70622-76	384-4651	-	70621-76	384-4321	-	70622-04	384-4743	-	70621-04	384-4396	-
VI	70622-77	384-4669	-	70621-77	384-4339	-	70622-05	384-4750	-	70621-05	384-4404	-
Tights for Men - Natural												
	RAL CCL 1 - Open Toe						**RAL CCL 1 - Closed Toe**					
	Petite			**Regular**			**Petite**			**Regular**		
I	70622-204	-	-	70621-204	-	-	70622-234	-	-	70621-234	-	-
II	70622-205	-	-	70621-205	-	-	70622-235	-	-	70621-235	-	-
III	70622-206	-	-	70621-206	-	-	70622-236	-	-	70621-236	-	-
IV	70622-207	-	-	70621-207	-	-	70622-237	-	-	70621-237	-	-
V	70622-208	-	-	70621-208	-	-	70622-238	-	-	70621-238	-	-
VI	70622-209	-	-	70621-209	-	-	70622-239	-	-	70621-239	-	-
Tights for Men - Black												
	RAL CCL 1 - Open Toe						**RAL CCL 1 - Closed Toe**					
	Petite			**Regular**			**Petite**			**Regular**		
I	70622-198	-	-	70621-198	-	-	70622-228	-	-	70621-228	-	-
II	70622-199	-	-	70621-199	-	-	70622-229	-	-	70621-229	-	-
III	70622-200	-	-	70621-200	-	-	70622-230	-	-	70621-230	-	-
IV	70622-201	-	-	70621-201	-	-	70622-231	-	-	70621-231	-	-
V	70622-202	-	-	70621-202	-	-	70622-232	-	-	70621-232	-	-
VI	70622-203	-	-	70621-203	-	-	70622-233	-	-	70621-233	-	-

JOBST® Opaque - CCL 1

Size	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC
Chap Style - Natural												
	RAL CCL 1 - Open Toe											
	Left, Petite			**Left, Regular**			**Right, Petite**			**Right, Regular**		
I	70620-54	-	-	70619-54	-	-	70620-24	-	-	70619-24	-	-
II	70620-55	-	-	70619-55	-	-	70620-25	-	-	70619-25	-	-
III	70620-56	-	-	70619-56	-	-	70620-26	-	-	70619-26	-	-
IV	70620-57	-	-	70619-57	-	-	70620-27	-	-	70619-27	-	-
V	70620-58	-	-	70619-58	-	-	70620-28	-	-	70619-28	-	-
VI	70620-59	-	-	70619-59	-	-	70620-29	-	-	70619-29	-	-
Chap Style - Natural												
	RAL CCL 1 - Closed Toe											
	Left, Petite			**Left, Regular**			**Right, Petite**			**Right, Regular**		
I	70620-144	-	-	70619-144	-	-	70620-114	-	-	70619-114	-	-
II	70620-145	-	-	70619-145	-	-	70620-115	-	-	70619-115	-	-
III	70620-146	-	-	70619-146	-	-	70620-116	-	-	70619-116	-	-
IV	70620-147	-	-	70619-147	-	-	70620-117	-	-	70619-117	-	-
V	70620-148	-	-	70619-148	-	-	70620-118	-	-	70619-118	-	-
VI	70620-149	-	-	70619-149	-	-	70620-119	-	-	70619-119	-	-
Chap Style - Natural												
	RAL CCL 1 - Open Toe						**RAL CCL 1 - Closed Toe**					
	Pair, Petite			**Pair, Regular**			**Pair, Petite**			**Pair, Regular**		
I	70620-84	-	-	70619-84	-	-	70620-174	-	-	70619-174	-	-
II	70620-85	-	-	70619-85	-	-	70620-175	-	-	70619-175	-	-
III	70620-86	-	-	70619-86	-	-	70620-176	-	-	70619-176	-	-
IV	70620-87	-	-	70619-87	-	-	70620-177	-	-	70619-177	-	-
V	70620-88	-	-	70619-88	-	-	70620-178	-	-	70619-178	-	-
VI	70620-89	-	-	70619-89	-	-	70620-179	-	-	70619-179	-	-



JOBST® Opaque - CCL 2

Size	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC
Knee High - Natural												
	RAL CCL 2 - Open Toe						**RAL CCL 2 - Closed Toe**					
	Petite			**Regular**			**Petite**			**Regular**		
I	70626-24	326-2441	-	70625-24	326-2326	-	70626-30	326-3746	-	70625-30	326-3621	-
II	70626-25	326-2458	-	70625-25	326-2334	-	70626-31	326-3753	-	70625-31	326-3639	-
III	70626-26	326-2466	-	70625-26	326-2342	-	70626-32	326-3761	-	70625-32	326-3647	-
IV	70626-27	326-2474	-	70625-27	326-2359	-	70626-33	326-3779	-	70625-33	326-3654	-
V	70626-28	326-2482	-	70625-28	326-2367	-	70626-34	326-3787	-	70625-34	326-3662	-
VI	70626-29	326-2490	-	70625-29	326-2375	-	70626-35	326-3795	-	70625-35	326-3670	-
Knee High - Natural												
	RAL CCL 2 - Open Toe, Dotted Silicone Band						**RAL CCL 2 - Closed Toe, Dotted Silicone Band**					
	Petite			**Regular**			**Petite**			**Regular**		
I	70626-66	-	-	70625-66	-	EGD7740	70626-78	-	-	70625-78	-	-
II	70626-67	-	-	70625-67	-	EGD7741	70626-79	-	-	70625-79	-	-
III	70626-68	-	-	70625-68	-	EGD7742	70626-80	-	-	70625-80	-	-
IV	70626-69	-	-	70625-69	-	EGD7743	70626-81	-	-	70625-81	-	-
V	70626-70	-	-	70625-70	-	EGD7744	70626-82	-	-	70625-82	-	-
VI	70626-71	-	-	70625-71	-	EGD7745	70626-83	-	-	70625-83	-	-
Knee High - Black												
	RAL CCL 2 - Open Toe						**RAL CCL 2 - Closed Toe**					
	Petite			**Regular**			**Petite**			**Regular**		
I	70626-48	326-2383	-	70625-48	326-2268	-	70626-36	326-3688	-	70625-36	326-3563	-
II	70626-49	326-2391	-	70625-49	326-2276	-	70626-37	326-3696	-	70625-37	326-3571	-
III	70626-50	326-2409	-	70625-50	326-2284	-	70626-38	326-3704	-	70625-38	326-3589	-
IV	70626-51	326-2417	-	70625-51	326-2292	-	70626-39	326-3712	-	70625-39	326-3597	-
V	70626-52	326-2425	-	70625-52	326-2300	-	70626-40	326-3720	-	70625-40	326-3605	-
VI	70626-53	326-2433	-	70625-53	326-2318	-	70626-41	326-3738	-	70625-41	326-3613	-
Knee High - Black												
	RAL CCL 2 - Open Toe, Dotted Silicone Band						**RAL CCL 2 - Closed Toe, Dotted Silicone Band**					
	Petite			**Regular**			**Petite**			**Regular**		
I	70626-60	-	-	70625-60	-	-	70626-84	-	-	70625-84	-	-
II	70626-61	-	-	70625-61	-	-	70626-85	-	-	70625-85	-	-
III	70626-62	-	-	70625-62	-	-	70626-86	-	-	70625-86	-	-
IV	70626-63	-	-	70625-63	-	-	70626-87	-	-	70625-87	-	-
V	70626-64	-	-	70625-64	-	-	70626-88	-	-	70625-88	-	-
VI	70626-65	-	-	70625-65	-	-	70626-89	-	-	70625-89	-	-
Knee High - Navy												
	RAL CCL 2 - Open Toe						**RAL CCL 2 - Closed Toe**					
	Petite			**Regular**			**Petite**			**Regular**		
I	70626-90	360-7496	-	70625-90	360-7553	-	70626-18	360-7371	-	70625-18	360-7439	-
II	70626-91	360-7504	-	70625-91	360-7561	-	70626-19	360-7389	-	70625-19	360-7447	-
III	70626-92	360-7512	-	70625-92	360-7579	-	70626-20	360-7397	-	70625-20	360-7454	-
IV	70626-93	360-7520	-	70625-93	360-7587	-	70626-21	360-7405	-	70625-21	360-7462	-
V	70626-94	360-7538	-	70625-94	360-7595	-	70626-22	360-7413	-	70625-22	360-7470	-
VI	70626-95	360-7546	-	70625-95	360-7603	-	70626-23	360-7421	-	70625-23	360-7488	-

JOBST® Opaque - CCL 2

Size	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC
Knee High - Navy												
	RAL CCL 2 - Open Toe, Dotted Silicone Band						**RAL CCL 2 - Closed Toe, Dotted Silicone Band**					
	Petite			**Regular**			**Petite**			**Regular**		
I	70626-102	-	-	70625-102	-	-	70626-114	-	-	70625-114	-	-
II	70626-103	-	-	70625-103	-	-	70626-115	-	-	70625-115	-	-
III	70626-104	-	-	70625-104	-	-	70626-116	-	-	70625-116	-	-
IV	70626-105	-	-	70625-105	-	-	70626-117	-	-	70625-117	-	-
V	70626-106	-	-	70625-106	-	-	70626-118	-	-	70625-118	-	-
VI	70626-107	-	-	70625-107	-	-	70626-119	-	-	70625-119	-	-
Knee High - Bronze												
	RAL CCL 2 - Open Toe						**RAL CCL 2 - Closed Toe**					
	Petite			**Regular**			**Petite**			**Regular**		
I	70626-120	380-7146	-	70625-120	380-6601	-	70626-126	380-7203	-	70625-126	380-6841	-
II	70626-121	380-7153	-	70625-121	380-6619	-	70626-127	380-7211	-	70625-127	380-6858	-
III	70626-122	380-7161	-	70625-122	380-6627	-	70626-128	380-7229	-	70625-128	380-6866	-
IV	70626-123	380-7179	-	70625-123	380-6635	-	70626-129	380-7237	-	70625-129	380-6874	-
V	70626-124	380-7187	-	70625-124	380-6643	-	70626-130	380-7245	-	70625-130	380-6882	-
VI	70626-125	380-7195	-	70625-125	380-6650	-	70626-131	380-7252	-	70625-131	380-6890	-
Knee High - Caramel												
	RAL CCL 2 - Open Toe						**RAL CCL 2 - Closed Toe**					
	Petite			**Regular**			**Petite**			**Regular**		
I	70626-00	384-0824	-	70625-00	384-0584	-	70626-06	384-0881	-	70625-06	384-0634	-
II	70626-01	384-0832	-	70625-01	384-0576	-	70626-07	384-0899	-	70625-07	384-0642	-
III	70626-02	384-0840	-	70625-02	384-0592	-	70626-08	384-0907	-	70625-08	384-0659	-
IV	70626-03	384-0857	-	70625-03	384-0600	-	70626-09	384-0915	-	70625-09	384-0667	-
V	70626-04	384-0865	-	70625-04	384-0618	-	70626-10	384-0923	-	70625-10	384-0675	-
VI	70626-05	384-0873	-	70625-05	384-0626	-	70626-11	384-0931	-	70625-11	384-0683	-
Thigh High - Natural												
	RAL CCL 2 - Open Toe, Dotted Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70628-126	326-8877	-	70628-270	326-3142	-	70627-126	326-2599	-	70627-270	326-3027	-
II	70628-127	326-8885	-	70628-271	326-3159	-	70627-127	326-2607	-	70627-271	326-3035	-
III	70628-128	326-8893	-	70628-272	326-3167	-	70627-128	326-2615	-	70627-272	326-3043	-
IV	70628-129	326-8901	-	70628-273	326-3175	-	70627-129	326-2623	-	70627-273	326-3050	-
V	70628-130	326-8919	-	70628-274	326-3183	-	70627-130	326-2631	-	70627-274	326-3068	-
VI	70628-131	326-8869	-	70628-275	326-3191	-	70627-131	326-2649	-	70627-275	326-3076	-
Thigh High - Natural												
	RAL CCL 2 - Closed Toe, Dotted Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70628-150	357-5057	-	70628-294	357-8150	-	70627-150	357-5115	-	70627-294	357-8218	-
II	70628-151	357-5065	-	70628-295	357-8168	-	70627-151	357-5123	-	70627-295	357-8226	-
III	70628-152	357-5073	-	70628-296	357-8176	-	70627-152	357-5131	-	70627-296	357-8234	-
IV	70628-153	357-5081	-	70628-297	357-8184	-	70627-153	357-5149	-	70627-297	357-8242	-
V	70628-154	357-5099	-	70628-298	357-8192	-	70627-154	357-5156	-	70627-298	357-8259	-
VI	70628-155	357-5107	-	70628-299	357-8200	-	70627-155	357-5164	-	70627-299	357-8267	-

JOBST® Opaque - CCL 2

Size	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC
Thigh High - Natural												
	RAL CCL 2 - Open Toe, Sensitive Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70628-582	358-3499	-	70628-630	358-3614	-	70627-558	358-3739	-	70627-606	358-3861	-
II	70628-583	358-3507	-	70628-631	358-3622	-	70627-559	358-3747	-	70627-607	358-3879	-
III	70628-584	358-3515	-	70628-632	358-3630	-	70627-560	358-3754	-	70627-608	358-3887	-
IV	70628-585	358-3523	-	70628-633	358-3648	-	70627-561	358-3762	-	70627-609	358-3895	-
V	70628-586	358-3531	-	70628-634	358-3655	-	70627-562	358-3770	-	70627-610	358-3911	-
VI	70628-587	358-3549	-	70628-635	358-3663	-	70627-563	358-3788	-	70627-611	358-3929	-
Thigh High - Natural												
	RAL CCL 2 - Closed Toe, Sensitive Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70628-600	357-4431	-	70628-648	357-4936	-	70627-576	357-4613	-	70627-624	357-8036	-
II	70628-601	357-4449	-	70628-649	357-4944	-	70627-577	357-4621	-	70627-625	357-8044	-
III	70628-602	357-4456	-	70628-650	357-4951	-	70627-578	357-4639	-	70627-626	357-8051	-
IV	70628-603	357-4464	-	70628-651	357-4969	-	70627-579	357-4647	-	70627-627	357-8069	-
V	70628-604	357-4472	-	70628-652	357-4977	-	70627-580	357-4654	-	70627-628	357-8077	-
VI	70628-605	357-4480	-	70628-653	357-4985	-	70627-581	357-4662	-	70627-629	357-8085	-
Thigh High - Natural												
	RAL CCL 2 - Open Toe, Lace Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70628-54	357-5875	-	70628-222	357-7236	-	70627-54	357-5172	-	70627-222	357-6626	-
II	70628-55	357-5883	-	70628-223	357-7244	-	70627-55	357-5180	-	70627-223	357-6634	-
III	70628-56	357-5891	-	70628-224	357-7251	-	70627-56	357-5198	-	70627-224	357-6642	-
IV	70628-57	357-5917	-	70628-225	357-7269	-	70627-57	357-5206	-	70627-225	357-6659	-
V	70628-58	357-5925	-	70628-226	357-7277	-	70627-58	357-5214	-	70627-226	357-6675	-
VI	70628-59	357-5933	-	70628-227	357-7285	-	70627-59	357-5222	-	70627-227	357-6683	-
Thigh High - Natural												
	RAL CCL 2 - Closed Toe, Lace Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70628-66	357-5958	-	70628-234	357-6733	-	70627-66	357-6436	-	70627-234	357-7178	-
II	70628-67	357-5966	-	70628-235	357-6741	-	70627-67	357-6444	-	70627-235	357-7186	-
III	70628-68	357-5982	-	70628-236	357-6766	-	70627-68	357-6451	-	70627-236	357-7194	-
IV	70628-69	357-5990	-	70628-237	357-6774	-	70627-69	357-6469	-	70627-237	357-7202	-
V	70628-70	357-6006	-	70628-238	357-6782	-	70627-70	357-6477	-	70627-238	357-7210	-
VI	70628-71	357-6014	-	70628-239	357-6808	-	70627-71	357-6485	-	70627-239	357-7228	-
Thigh High - Black												
	RAL CCL 2 - Open Toe, Dotted Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70628-330	326-8786	-	70628-348	326-3084	-	70627-330	326-2532	-	70627-348	326-2656	-
II	70628-331	326-8794	-	70628-349	326-3092	-	70627-331	326-2540	-	70627-349	326-2664	-
III	70628-332	326-8810	-	70628-350	326-3100	-	70627-332	326-2557	-	70627-350	326-2672	-
IV	70628-333	326-8828	-	70628-351	326-3118	-	70627-333	326-2565	-	70627-351	326-2680	-
V	70628-334	326-8836	-	70628-352	326-3126	-	70627-334	326-2573	-	70627-352	326-2698	-
VI	70628-335	326-8851	-	70628-353	326-3134	-	70627-335	326-2581	-	70627-353	326-2706	-

JOBST® Opaque - CCL 2

Size	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC
Thigh High - Black												
	RAL CCL 2 - Closed Toe, Dotted Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70628-156	357-4670	-	70628-300	357-8093	-	70627-156	357-4738	-	70627-300	357-7970	-
II	70628-157	357-4688	-	70628-301	357-8101	-	70627-157	357-4746	-	70627-301	357-7988	-
III	70628-158	357-4696	-	70628-302	357-8119	-	70627-158	357-4753	-	70627-302	357-7996	-
IV	70628-159	357-4704	-	70628-303	357-8127	-	70627-159	357-4761	-	70627-303	357-8002	-
V	70628-160	357-4712	-	70628-304	357-8135	-	70627-160	357-4779	-	70627-304	357-8010	-
VI	70628-161	357-4720	-	70628-305	357-8143	-	70627-161	357-4787	-	70627-305	357-8028	-
Thigh High - Black												
	RAL CCL 2 - Open Toe, Sensitive Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70628-576	358-3556	-	70628-624	358-3671	-	70627-552	358-3796	-	70627-600	358-3937	-
II	70628-577	358-3564	-	70628-625	358-3689	-	70627-553	358-3804	-	70627-601	358-3945	-
III	70628-578	358-3572	-	70628-626	358-3697	-	70627-554	358-3812	-	70627-602	358-3952	-
IV	70628-579	358-3580	-	70628-627	358-3705	-	70627-555	358-3820	-	70627-603	358-3960	-
V	70628-580	358-3598	-	70628-628	358-3713	-	70627-556	358-3838	-	70627-604	358-3978	-
VI	70628-581	358-3606	-	70628-629	358-3721	-	70627-557	358-3846	-	70627-605	358-3986	-
Thigh High - Black												
	RAL CCL 2 - Closed Toe, Sensitive Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70628-606	357-4498	-	70628-654	357-4993	-	70627-582	357-4555	-	70627-630	326-3027	-
II	70628-607	357-4506	-	70628-655	357-5008	-	70627-583	357-4563	-	70627-631	326-3035	-
III	70628-608	357-4514	-	70628-656	357-5016	-	70627-584	357-4571	-	70627-632	326-3043	-
IV	70628-609	357-4522	-	70628-657	357-5024	-	70627-585	357-4589	-	70627-633	326-3050	-
V	70628-610	357-4530	-	70628-658	357-5032	-	70627-586	357-4597	-	70627-634	326-3068	-
VI	70628-611	357-4548	-	70628-659	357-5040	-	70627-587	357-4605	-	70627-635	326-3076	-
Thigh High - Black												
	RAL CCL 2 - Open Toe, Lace Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70628-317	357-5800	-	70628-342	357-6493	-	70627-318	357-5230	-	70627-342	357-6568	-
II	70628-319	357-5818	-	70628-343	357-6501	-	70627-319	357-5248	-	70627-343	357-6576	-
III	70628-320	357-5826	-	70628-344	357-6527	-	70627-320	357-5255	-	70627-344	357-6584	-
IV	70628-321	357-5842	-	70628-345	357-6535	-	70627-321	357-5263	-	70627-345	357-6592	-
V	70628-322	357-5859	-	70628-346	357-6543	-	70627-322	357-5271	-	70627-346	357-6600	-
VI	70628-323	357-5867	-	70628-347	357-6550	-	70627-323	357-5289	-	70627-347	357-6618	-
Thigh High - Black												
	RAL CCL 2 - Closed Toe, Lace Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70628-78	357-6022	-	70628-246	357-6824	-	70627-78	357-6352	-	70627-246	357-7350	-
II	70628-79	357-6030	-	70628-247	357-6832	-	70627-79	357-6360	-	70627-247	357-7368	-
III	70628-80	357-6048	-	70628-248	357-6840	-	70627-80	357-6378	-	70627-248	357-7376	-
IV	70628-81	357-6063	-	70628-249	357-6865	-	70627-81	357-6394	-	70627-249	357-7384	-
V	70628-82	357-6089	-	70628-250	357-6873	-	70627-82	357-6402	-	70627-250	357-7392	-
VI	70628-83	357-6097	-	70628-251	357-6881	-	70627-83	357-6410	-	70627-251	357-7400	-

JOBST® Opaque - CCL 2

Size	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC
Thigh High - Navy												
	RAL CCL 2 - Open Toe, Dotted Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70628-378	361-1233	-	70628-414	361-1290	-	70627-372	361-1357	-	70627-414	361-1415	-
II	70628-379	361-1241	-	70628-415	361-1308	-	70627-373	361-1365	-	70627-415	361-1423	-
III	70628-380	361-1258	-	70628-416	361-1316	-	70627-374	361-1373	-	70627-416	361-1431	-
IV	70628-381	361-1266	-	70628-417	361-1324	-	70627-375	361-1381	-	70627-417	361-1449	-
V	70628-382	361-1274	-	70628-418	361-1332	-	70627-376	361-1399	-	70627-418	361-1456	-
VI	70628-383	361-1282	-	70628-419	361-1340	-	70627-377	361-1407	-	70627-419	361-1464	-
Thigh High - Navy												
	RAL CCL 2 - Closed Toe, Dotted Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70628-144	360-9898	-	70628-288	360-9971	-	70627-144	361-0037	-	70627-288	361-0094	-
II	70628-145	360-9914	-	70628-289	360-9989	-	70627-145	361-0045	-	70627-289	361-0102	-
III	70628-146	360-9922	-	70628-290	360-9997	-	70627-146	361-0052	-	70627-290	361-0110	-
IV	70628-147	360-9930	-	70628-291	361-0003	-	70627-147	361-0060	-	70627-291	361-0128	-
V	70628-148	360-9955	-	70628-292	361-0011	-	70627-148	361-0078	-	70627-292	361-0136	-
VI	70628-149	360-9963	-	70628-293	361-0029	-	70627-149	361-0086	-	70627-293	361-0144	-
Thigh High - Navy												
	RAL CCL 2 - Open Toe, Sensitive Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70628-570	361-0151	-	70628-618	361-0219	-	70627-546	361-0276	-	70627-594	361-0334	-
II	70628-571	361-0169	-	70628-619	361-0227	-	70627-547	361-0284	-	70627-595	361-0342	-
III	70628-572	361-0177	-	70628-620	361-0235	-	70627-548	361-0292	-	70627-596	361-0359	-
IV	70628-573	361-0185	-	70628-621	361-0243	-	70627-549	361-0300	-	70627-597	361-0367	-
V	70628-574	361-0193	-	70628-622	361-0250	-	70627-550	361-0318	-	70627-598	361-0375	-
VI	70628-575	361-0201	-	70628-623	361-0268	-	70627-551	361-0326	-	70627-599	361-0383	-
Thigh High - Navy												
	RAL CCL 2 - Closed Toe, Sensitive Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70628-594	360-9344	-	70628-642	360-9476	-	70627-570	360-9534	-	70627-618	360-9591	-
II	70628-595	360-9351	-	70628-643	360-9484	-	70627-571	360-9542	-	70627-619	360-9609	-
III	70628-596	360-9369	-	70628-644	360-9492	-	70627-572	360-9559	-	70627-620	360-9617	-
IV	70628-597	360-9377	-	70628-645	360-9500	-	70627-573	360-9567	-	70627-621	360-9625	-
V	70628-598	360-9385	-	70628-646	360-9518	-	70627-574	360-9575	-	70627-622	360-9633	-
VI	70628-599	360-9393	-	70628-647	360-9526	-	70627-575	360-9583	-	70627-623	360-9641	-
Thigh High - Navy												
	RAL CCL 2 - Open Toe, Lace Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70628-372	361-0391	-	70628-408	361-0458	-	70627-366	361-0516	-	70627-408	361-0573	-
II	70628-373	361-0409	-	70628-409	361-0466	-	70627-367	361-0524	-	70627-409	361-0581	-
III	70628-374	361-0417	-	70628-410	361-0474	-	70627-368	361-0532	-	70627-410	361-0599	-
IV	70628-375	361-0425	-	70628-411	361-0482	-	70627-369	361-0540	-	70627-411	361-0607	-
V	70628-376	361-0433	-	70628-412	361-0490	-	70627-370	361-0557	-	70627-412	361-0615	-
VI	70628-377	361-0441	-	70628-413	361-0508	-	70627-371	361-0565	-	70627-413	361-0623	-

JOBST® Opaque - CCL 2

Size	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC
Thigh High - Navy												
	RAL CCL 2 - Closed Toe, Lace Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70628-42	360-9658	-	70628-210	360-9716	-	70627-42	360-9773	-	70627-210	360-9831	-
II	70628-43	360-9666	-	70628-211	360-9724	-	70627-43	360-9781	-	70627-211	360-9849	-
III	70628-44	360-9674	-	70628-212	360-9732	-	70627-44	360-9799	-	70627-212	360-9856	-
IV	70628-45	360-9682	-	70628-213	360-9740	-	70627-45	360-9807	-	70627-213	360-9864	-
V	70628-46	360-9690	-	70628-214	360-9757	-	70627-46	360-9815	-	70627-214	360-9872	-
VI	70628-47	360-9708	-	70628-215	360-9765	-	70627-47	360-9823	-	70627-215	360-9880	-
Thigh High - Bronze												
	RAL CCL 2 - Open Toe, Dotted Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70628-444	383-9446	-	70628-492	383-9818	-	70627-444	383-7994	-	70627-492	383-8356	-
II	70628-445	383-9453	-	70628-493	383-9826	-	70627-445	383-8000	-	70627-493	383-8364	-
III	70628-446	383-9461	-	70628-494	383-9834	-	70627-446	383-8018	-	70627-494	383-8372	-
IV	70628-447	383-9479	-	70628-495	383-9842	-	70627-447	383-8026	-	70627-495	383-8380	-
V	70628-448	383-9487	-	70628-496	383-9859	-	70627-448	383-8034	-	70627-496	383-8398	-
VI	70628-449	383-9495	-	70628-497	383-9867	-	70627-449	383-8042	-	70627-497	383-8406	-
Thigh High - Bronze												
	RAL CCL 2 - Closed Toe, Dotted Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70628-450	383-9511	-	70628-498	383-9875	-	70627-450	383-8059	-	70627-498	383-8414	-
II	70628-451	383-9529	-	70628-499	383-9883	-	70627-451	383-8067	-	70627-499	383-8422	-
III	70628-452	383-9537	-	70628-500	383-9891	-	70627-452	383-8075	-	70627-500	383-8430	-
IV	70628-453	383-9545	-	70628-501	383-9909	-	70627-453	383-8083	-	70627-501	383-8448	-
V	70628-454	383-9552	-	70628-502	383-9917	-	70627-454	383-8091	-	70627-502	383-8455	-
VI	70628-455	383-9560	-	70628-503	383-9925	-	70627-455	383-8109	-	70627-503	383-8463	-
Thigh High - Bronze												
	RAL CCL 2 - Open Toe, Sensitive Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70628-660	383-9198	-	70628-672	383-9578	-	70627-636	383-7754	-	70627-648	383-8117	-
II	70628-661	383-9198	-	70628-673	383-9586	-	70627-637	383-7762	-	70627-649	383-8125	-
III	70628-662	383-9214	-	70628-674	383-9594	-	70627-638	383-7770	-	70627-650	383-8133	-
IV	70628-663	383-9222	-	70628-675	383-9602	-	70627-639	383-7788	-	70627-651	383-8141	-
V	70628-664	383-9230	-	70628-676	383-9610	-	70627-640	383-7796	-	70627-652	383-8158	-
VI	70628-665	383-9248	-	70628-677	383-9628	-	70627-641	383-7804	-	70627-653	383-8166	-
Thigh High - Bronze												
	RAL CCL 2 - Closed Toe, Sensitive Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70628-666	383-9263	-	70628-678	383-9636	-	70627-642	383-7812	-	70627-654	383-8174	-
II	70628-667	383-9271	-	70628-679	383-9644	-	70627-643	383-7820	-	70627-655	383-8182	-
III	70628-668	383-9289	-	70628-680	383-9651	-	70627-644	383-7838	-	70627-656	383-8190	-
IV	70628-669	383-9297	-	70628-681	383-9669	-	70627-645	383-7846	-	70627-657	383-8208	-
V	70628-670	383-9305	-	70628-682	383-9677	-	70627-646	383-7853	-	70627-658	383-8216	-
VI	70628-671	383-9313	-	70628-683	383-9685	-	70627-647	383-7861	-	70627-659	383-8224	-

JOBST® Opaque - CCL 2

Size	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC
Thigh High - Bronze												
	RAL CCL 2 - Open Toe, Lace Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70628-432	383-9321	-	70628-480	383-9693	-	70627-432	383-7879	-	70627-480	383-8232	-
II	70628-433	383-9339	-	70628-481	383-9701	-	70627-433	383-7887	-	70627-481	383-8240	-
III	70628-434	383-9347	-	70628-482	383-9719	-	70627-434	383-7895	-	70627-482	383-8257	-
IV	70628-435	383-9354	-	70628-483	383-9727	-	70627-435	383-7903	-	70627-483	383-8265	-
V	70628-436	383-9362	-	70628-484	383-9735	-	70627-436	383-7911	-	70627-484	383-8273	-
VI	70628-437	383-9370	-	70628-485	383-9743	-	70627-437	383-7929	-	70627-485	383-8281	-
Thigh High - Bronze												
	RAL CCL 2 - Closed Toe, Lace Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70628-438	383-9388	-	70628-486	383-9750	-	70627-438	383-7937	-	70627-486	383-8299	-
II	70628-439	383-9396	-	70628-487	383-9768	-	70627-439	383-7945	-	70627-487	383-8307	-
III	70628-440	383-9404	-	70628-488	383-9776	-	70627-440	383-7952	-	70627-488	383-8315	-
IV	70628-441	383-9412	-	70628-489	383-9784	-	70627-441	383-7960	-	70627-489	383-8323	-
V	70628-442	383-9420	-	70628-490	383-9792	-	70627-442	383-7978	-	70627-490	383-8331	-
VI	70628-443	383-9438	-	70628-491	383-9800	-	70627-443	383-7986	-	70627-491	383-8349	-
Thigh High - Caramel												
	RAL CCL 2 - Open Toe, Dotted Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70628-120	384-3638	-	70628-264	384-4156	-	70627-120	384-2135	-	70627-264	384-2549	-
II	70628-121	384-3646	-	70628-265	384-4164	-	70627-121	384-2143	-	70627-265	384-2556	-
III	70628-122	384-3653	-	70628-266	384-4172	-	70627-122	384-2150	-	70627-266	384-2564	-
IV	70628-123	384-3687	-	70628-267	384-4180	-	70627-123	384-2168	-	70627-267	384-2572	-
V	70628-124	384-3661	-	70628-268	384-4198	-	70627-124	384-2176	-	70627-268	384-2580	-
VI	70628-125	384-3679	-	70628-269	384-4206	-	70627-125	384-2184	-	70627-269	384-2598	-
Thigh High - Caramel												
	RAL CCL 2 - Closed Toe, Dotted Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70628-132	384-3703	-	70628-276	384-4214	-	70627-132	384-2192	-	70627-276	384-2606	-
II	70628-133	384-3711	-	70628-277	384-4222	-	70627-133	384-2200	-	70627-277	384-2614	-
III	70628-134	384-3729	-	70628-278	384-4230	-	70627-134	384-2218	-	70627-278	384-2622	-
IV	70628-135	384-3737	-	70628-279	384-4248	-	70627-135	384-2226	-	70627-279	384-2630	-
V	70628-136	384-3745	-	70628-280	384-4255	-	70627-136	384-2234	-	70627-280	384-2648	-
VI	70628-137	384-3752	-	70628-281	384-4263	-	70627-137	384-2242	-	70627-281	384-2655	-
Thigh High - Caramel												
	RAL CCL 2 - Open Toe, Sensitive Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70628-564	384-3380	-	70628-612	384-3802	-	70627-540	384-1822	-	70627-588	384-2259	-
II	70628-565	384-3398	-	70628-613	384-3810	-	70627-541	384-1830	-	70627-589	384-2267	-
III	70628-566	384-3406	-	70628-614	384-3828	-	70627-542	384-1848	-	70627-590	384-2275	-
IV	70628-567	384-3414	-	70628-615	384-3836	-	70627-543	384-1855	-	70627-591	384-2283	-
V	70628-568	384-3422	-	70628-616	384-3844	-	70627-544	384-1863	-	70627-592	384-2291	-
VI	70628-569	384-3430	-	70628-617	384-3851	-	70627-545	384-1871	-	70627-593	384-2309	-

JOBST® Opaque - CCL 2

Size	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC
Thigh High - Caramel												
	RAL CCL 2 - Closed Toe, Sensitive Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70628-588	384-1905	-	70628-636	384-3893	-	70627-564	384-3448	-	70627-612	384-2341	-
II	70628-589	384-1913	-	70628-637	384-3901	-	70627-565	384-3455	-	70627-613	384-2358	-
III	70628-590	384-1921	-	70628-638	384-3919	-	70627-566	384-3463	-	70627-614	384-2366	-
IV	70628-591	384-1939	-	70628-639	384-3927	-	70627-567	384-3471	-	70627-615	384-2374	-
V	70628-592	384-1947	-	70628-640	384-3935	-	70627-568	384-3489	-	70627-616	384-2382	-
VI	70628-593	384-1954	-	70628-641	384-3943	-	70627-569	384-3497	-	70627-617	384-2390	-
Thigh High - Caramel												
	RAL CCL 2 - Open Toe, Lace Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70628-06	384-3505	-	70628-174	384-3992	-	70627-06	384-1962	-	70627-174	384-2408	-
II	70628-07	384-3513	-	70628-175	384-4008	-	70627-07	384-1970	-	70627-175	384-2416	-
III	70628-08	384-3521	-	70628-176	384-4016	-	70627-08	384-1988	-	70627-177	384-2424	-
IV	70628-09	384-3539	-	70628-177	384-4024	-	70627-09	384-1996	-	70627-176	384-2432	-
V	70628-10	384-3547	-	70628-178	384-4032	-	70627-10	384-2002	-	70627-178	384-2440	-
VI	70628-11	384-3554	-	70628-179	384-4040	-	70627-11	384-2010	-	70627-179	384-2457	-
Thigh High - Caramel												
	RAL CCL 2 - Closed Toe, Lace Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	70628-18	384-3570	-	70628-186	384-4099	-	70627-18	384-2044	-	70627-186	384-2473	-
II	70628-19	384-3588	-	70628-187	384-4107	-	70627-19	384-2051	-	70627-187	384-2481	-
III	70628-20	384-3596	-	70628-188	384-4115	-	70627-20	384-2069	-	70627-188	384-2499	-
IV	70628-21	384-3604	-	70628-189	384-4123	-	70627-21	384-2077	-	70627-189	384-2507	-
V	70628-22	384-3612	-	70628-190	384-4131	-	70627-22	384-2085	-	70627-190	384-2515	-
VI	70628-23	384-3620	-	70628-191	384-4149	-	70627-23	384-2093	-	70627-191	384-2523	-
Tights - Natural												
	RAL CCL 2 - Open Toe						**RAL CCL 2 - Closed Toe**					
	Petite			**Regular**			**Petite**			**Regular**		
I	70628-582	357-3193	-	70631-48	357-3359	-	70632-54	357-3011	-	70631-54	357-3086	-
II	70628-583	357-3227	-	70631-49	357-3367	-	70632-55	357-3029	-	70631-55	357-3110	-
III	70628-584	357-3243	-	70631-50	357-3375	-	70632-56	357-3037	-	70631-56	357-3136	-
IV	70628-585	357-3276	-	70631-51	357-3383	-	70632-57	357-3045	-	70631-57	357-3151	-
V	70628-586	357-3292	-	70631-52	357-3417	-	70632-58	357-3052	-	70631-58	357-3177	-
VI	70628-587	357-3318	-	70631-53	357-3425	-	70632-59	357-3060	-	70631-59	357-3185	-
Tights - Black												
	RAL CCL 2 - Open Toe						**RAL CCL 2 - Closed Toe**					
	Petite			**Regular**			**Petite**			**Regular**		
I	70632-108	357-2344	-	70631-108	357-2401	-	70632-72	357-2047	-	70631-72	357-2229	-
II	70632-109	357-2351	-	70631-109	357-2419	-	70632-73	357-2062	-	70631-73	357-2252	-
III	70632-110	357-2369	-	70631-110	357-2427	-	70632-74	357-2096	-	70631-74	357-2286	-
IV	70632-111	357-2377	-	70631-111	357-2492	-	70632-75	357-2138	-	70631-75	357-2302	-
V	70632-112	357-2385	-	70631-112	357-2500	-	70632-76	357-2153	-	70631-76	357-2328	-
VI	70632-113	357-2393	-	70631-113	357-2534	-	70632-77	357-2187	-	70631-77	357-2336	-

JOBST® Opaque - CCL 2

Size	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC
Tights - Navy												
	RAL CCL 2 - Open Toe						**RAL CCL 2 - Closed Toe**					
	Petite			**Regular**			**Petite**			**Regular**		
I	70632-126	361-1837	-	70631-150	361-1894	-	70632-36	361-1712	-	70631-36	361-1779	-
II	70632-127	361-1845	-	70631-151	361-1902	-	70632-37	361-1720	-	70631-37	361-1787	-
III	70632-128	361-1852	-	70631-152	361-1910	-	70632-38	361-1738	-	70631-38	361-1795	-
IV	70632-129	361-1860	-	70631-153	361-1928	-	70632-39	361-1746	-	70631-39	361-1803	-
V	70632-130	361-1878	-	70631-154	361-1936	-	70632-40	361-1753	-	70631-40	361-1811	-
VI	70632-131	361-1886	-	70631-155	361-1944	-	70632-41	361-1761	-	70631-41	361-1829	-
Tights - Bronze												
	RAL CCL 2 - Open Toe						**RAL CCL 2 - Closed Toe**					
	Petite			**Regular**			**Petite**			**Regular**		
I	70632-240	384-0204	-	70631-240	384-0055	-	70632-246	384-0279	-	70631-246	384-0113	-
II	70632-241	384-0212	-	70631-241	384-0063	-	70632-247	384-0287	-	70631-247	384-0121	-
III	70632-242	384-0220	-	70631-242	384-0071	-	70632-248	384-0295	-	70631-248	384-0139	-
IV	70632-243	384-0238	-	70631-243	384-0089	-	70632-249	384-0303	-	70631-249	384-0147	-
V	70632-244	384-0246	-	70631-244	384-0097	-	70632-250	384-0311	-	70631-250	384-0154	-
VI	70632-245	384-0253	-	70631-245	384-0105	-	70632-251	384-0329	-	70631-251	384-0162	-
Tights - Caramel												
	RAL CCL 2 - Open Toe						**RAL CCL 2 - Closed Toe**					
	Petite			**Regular**			**Petite**			**Regular**		
I	70632-00	384-4768	-	70631-00	384-4412	-	70632-06	384-4842	-	70631-06	384-4537	-
II	70632-01	384-4776	-	70631-01	384-4420	-	70632-07	384-4859	-	70631-07	384-4545	-
III	70632-02	384-4784	-	70631-02	384-4438	-	70632-08	384-4867	-	70631-08	384-4552	-
IV	70632-03	384-4792	-	70631-03	384-4446	-	70632-09	384-4875	-	70631-09	384-4560	-
V	70632-04	384-4800	-	70631-04	384-4453	-	70632-10	384-4883	-	70631-10	384-4578	-
VI	70632-05	384-4818	-	70631-05	384-4461	-	70632-11	384-4891	-	70631-11	384-4586	-
Tights for Men - Natural												
	RAL CCL 2 - Open Toe						**RAL CCL 2 - Closed Toe**					
	Petite			**Regular**			**Petite**			**Regular**		
I	70632-102	-	-	70631-102	-	-	70632-234	-	-	70631-234	-	-
II	70632-103	-	-	70631-103	-	-	70632-235	-	-	70631-235	-	-
III	70632-104	-	-	70631-104	-	-	70632-236	-	-	70631-236	-	-
IV	70632-105	-	-	70631-105	-	-	70632-237	-	-	70631-237	-	-
V	70632-106	-	-	70631-106	-	-	70632-238	-	-	70631-238	-	-
VI	70632-107	-	-	70631-107	-	-	70632-239	-	-	70631-239	-	-
Tights for Men - Black												
	RAL CCL 2 - Open Toe						**RAL CCL 2 - Closed Toe**					
	Petite			**Regular**			**Petite**			**Regular**		
I	70632-114	-	-	70631-114	-	-	70632-228	-	-	70631-228	-	-
II	70632-115	-	-	70631-1 15	-	-	70632-229	-	-	70631-229	-	-
III	70632-116	-	-	70631-116	-	-	70632-230	-	-	70631-230	-	-
IV	70632-117	-	-	70631-117	-	-	70632-231	-	-	70631-231	-	-
V	70632-118	-	-	70631-118	-	-	70632-232	-	-	70631-232	-	-
VI	70632-119	-	-	70631-119	-	-	70632-233	-	-	70631-233	-	-

JOBST® Opaque - CCL 2

Size	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC
	Chap Style - Natural											
	RAL CCL 2 - Open Toe											
	Left, Petite			**Left, Regular**			**Right, Petite**			**Right, Regular**		
I	70630-42	-	-	70629-42	-	-	70630-36	-	-	70629-36	-	-
II	70630-43	-	-	70629-43	-	-	70630-37	-	-	70629-37	-	-
III	70630-44	-	-	70629-44	-	-	70630-38	-	-	70629-38	-	-
IV	70630-45	-	-	70629-45	-	-	70630-39	-	-	70629-39	-	-
V	70630-46	-	-	70629-46	-	-	70630-40	-	-	70629-40	-	-
VI	70630-47	-	-	70629-47	-	-	70630-41	-	-	70629-41	-	-
	Chap Style - Natural											
	RAL CCL 2 - Closed Toe											
	Left, Petite			**Left, Regular**			**Right, Petite**			**Right, Regular**		
I	70630-60	-	-	70629-60	-	-	70630-54	-	-	70629-54	-	-
II	70630-61	-	-	70629-61	-	-	70630-55	-	-	70629-55	-	-
III	70630-62	-	-	70629-62	-	-	70630-56	-	-	70629-56	-	-
IV	70630-63	-	-	70629-63	-	-	70630-57	-	-	70629-57	-	-
V	70630-64	-	-	70629-64	-	-	70630-58	-	-	70629-58	-	-
VI	70630-65	-	-	70629-65	-	-	70630-59	-	-	70629-59	-	-
	Chap Style - Natural											
	RAL CCL 2 - Open Toe						**RAL CCL 2 - Closed Toe**					
	Pair, Petite			**Pair, Regular**			**Pair, Petite**			**Pair, Regular**		
I	70630-48	-	-	70629-48	-	-	70630-66	-	-	70629-66	-	-
II	70630-49	-	-	70629-49	-	-	70630-67	-	-	70629-67	-	-
III	70630-50	-	-	70629-50	-	-	70630-68	-	-	70629-68	-	-
IV	70630-51	-	-	70629-51	-	-	70630-69	-	-	70629-69	-	-
V	70630-52	-	-	70629-52	-	-	70630-70	-	-	70629-70	-	-
VI	70630-53	-	-	70629-53	-	-	70630-71	-	-	70629-71	-	-



JOBST® UltraSheer

Available on Drug Tariff

Ready-to-Wear RAL Compression Hosiery

JOBST® Ultrasheer is a range of ready-to-wear, circular-knit hosiery for the treatment of mild to moderate lymphoedema and venous disease.

Elegant and stylish

JOBST® UltraSheer provides outstanding sheerness and a silky look for an elegant and stylish appearance that encourages better patient compliance. The fine yarns deliver superior breathability and moisture management making the garments more comfortable to wear. With six fashionable colours for increased patient choice, JOBST® UltraSheer is suitable for any wardrobe.

Compression Classes	
RAL CCL 1 = 18-21mmHg	RAL CCL 2 = 23-32mmHg
✓	✓

Sizes and Options / Materials					
Six Sizes	Regular	Petite	Wider Circumference	Open Toe	Closed Toe
✓	✓	✓	✓	✓	✓



Lace silicone top band





Sensitive top band




Decorative dotted

Colours					
Natural	Caramel	Bronze	Espresso	Anthracite	Black

Styles			
			
AD Knee High	AG Thigh High	Tights	Maternity Tights

JOBST® UltraSheer

	Measurement Point	I	II	III	IV	V	VI
		RAL Compression Classes 1 (18-21mmHg) and 2 (23-32mmHg)					
circumference in cm	cG silicone band	43 - 48	45 - 52	49 - 56	53 - 60	56 - 64	60 - 68
	cG silicone band wide	49 - 57	53 - 62	57 - 67	61 - 72	65 - 77	69 - 81
	cG	43 - 57	45 - 62	49 - 67	53 - 72	56 - 77	60 - 81
	cF	39 - 52	41 - 56	44 - 60	47 - 65	50 - 69	53 - 73
	cE	30 - 37	33 - 40	35 - 43	37 - 45	39 - 48	41 - 51
	cD	27 - 33	29 - 36	32 - 39	34 - 42	36 - 45	38 - 48
	cC	28 - 34	30 - 37	33 - 40	35 - 43	37 - 46	39 - 49
	cB1	23 - 27	24 - 29	26 - 32	29 - 35	31 - 37	33 - 39
	cB	18 - 20	20 - 22	22 - 24	24 - 26	26 - 28	28 - 30
	cY	26 - 31	28 - 33	29 - 35	31 - 37	32 - 38	33 - 40
	cA	17 - 22	19 - 24	21 - 26	23 - 29	25 - 32	27 - 34
	Length measurements for all sizes (in cm): AD: ID regular 39 - 45, ID petite 34 - 39 AG: IG regular 70 - 80, IG petite 60 - 69 AT / AG-T: IK regular 75 - 85, IK petite 65 - 75						





JOBST® UltraSheer - CCL 1

Size	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC
Knee High - Anthracite							**Knee High - Black**					
	RAL CCL 1 - Closed Toe						**RAL CCL 1 - Closed Toe**					
	Petite			**Regular**			**Petite**			**Regular**		
I	73576-30	389-1546	-	73575-30	389-1538	-	73576-24	389-2098	-	73575-24	389-1744	EGD4973
II	73576-31	389-2049	-	73575-31	389-1694	-	73576-25	389-2106	-	73575-25	389-1751	EGD4974
III	73576-32	389-2056	-	73575-32	389-1702	-	73576-26	389-2114	-	73575-26	389-1769	EGD4891
IV	73576-33	389-2064	-	73575-33	389-1710	-	73576-27	389-2122	-	73575-27	389-1777	EGD4936
V	73576-34	389-2072	-	73575-34	389-1728	-	73576-28	389-2130	-	73575-28	389-1785	EGD4937
VI	73576-35	389-2080	-	73575-35	389-1736	-	73576-29	389-2148	-	73575-29	389-1793	EGD4938
Knee High - Bronze							**Knee High - Espresso**					
	RAL CCL 1 - Closed Toe						**RAL CCL 1 - Closed Toe**					
	Petite			**Regular**			**Petite**			**Regular**		
I	73576-96	389-2155	-	73575-84	389-1801	-	73576-12	389-2213	-	73575-12	389-1868	-
II	73576-97	389-2163	-	73575-85	389-1819	-	73576-13	389-2221	-	73575-13	389-1876	-
III	73576-98	389-2171	-	73575-86	389-1827	-	73576-14	389-2239	-	73575-14	389-1884	-
IV	73576-99	389-2189	-	73575-87	389-1835	-	73576-15	389-2247	-	73575-15	389-1892	-
V	73576-100	389-2197	-	73575-88	389-1843	-	73576-16	389-2254	-	73575-16	389-1900	-
VI	73576-101	389-2205	-	73575-89	389-1850	-	73576-17	389-2262	-	73575-17	389-1918	-
Knee High - Caramel							**Knee High - Natural**					
	RAL CCL 1 - Closed Toe						**RAL CCL 1 - Closed Toe**					
	Petite			**Regular**			**Petite**			**Regular**		
I	73576-00	389-2270	-	73575-00	389-1926	-	73576-18	389-2338	-	73575-18	389-1983	-
II	73576-01	389-2288	-	73575-01	389-1934	-	73576-19	389-2346	-	73575-19	389-1991	-
III	73576-02	389-2296	-	73575-02	389-1942	-	73576-20	389-2353	-	73575-20	389-2007	-
IV	73576-03	389-2304	-	73575-03	389-1959	-	73576-21	389-2361	-	73575-21	389-2015	-
V	73576-04	389-2312	-	73575-04	389-1967	-	73576-22	389-2379	-	73575-22	389-2023	-
VI	73576-05	389-2320	-	73575-05	389-1975	-	73576-23	389-2387	-	73575-23	389-2031	-
Thigh High - Anthracite												
	RAL CCL 1 - Closed Toe, Sensitive Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	73578-602	389-1637	-	73578-644	389-1652	-	73577-594	389-1579	-	73577-636	389-1595	-
II	73578-603	390-4463	-	73578-645	389-5463	-	73577-595	389-3997	-	73577-637	389-4946	-
III	73578-604	390-4471	-	73578-646	389-5471	-	73577-596	389-4003	-	73577-638	389-4953	-
IV	73578-605	390-4489	-	73578-647	389-5489	-	73577-597	389-4011	-	73577-639	389-4961	-
V	73578-606	390-4497	-	73578-648	389-5497	-	73577-598	389-4029	-	73577-640	389-4979	-
VI	73578-607	390-4505	-	73578-649	389-5505	-	73577-599	389-4037	-	73577-641	389-4987	-
Thigh High - Black												
	RAL CCL 1 - Open Toe, Sensitive Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	73578-662	389-1645	-	73578-668	389-1660	-	73577-654	389-1587	-	73577-660	389-1603	-
II	73578-663	390-4810	-	73578-669	389-5810	-	73577-655	389-4383	-	73577-661	389-5299	-
III	73578-664	390-4828	-	73578-670	389-5828	-	73577-656	389-4391	-	73577-662	389-5307	-
IV	73578-665	390-4836	-	73578-671	389-5836	-	73577-657	389-4417	-	73577-663	389-5315	-
V	73578-666	390-4844	-	73578-672	389-5844	-	73577-658	389-4425	-	73577-664	389-5323	-
VI	73578-667	390-4851	-	73578-673	389-5851	-	73577-659	389-4433	-	73577-665	389-5331	-

JOBST® UltraSheer - CCL 1

Size	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC
Thigh High - Black												
	RAL CCL 1 - Closed Toe, Sensitive Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	73578-596	390-4513	-	73578-638	389-5513	-	73577-588	389-4045	-	73577-630	389-4995	-
II	73578-597	390-4521	-	73578-639	389-5521	-	73577-589	389-4052	-	73577-631	389-5000	-
III	73578-598	390-4539	-	73578-640	389-5539	-	73577-590	389-4060	-	73577-632	389-5018	-
IV	73578-599	390-4547	-	73578-641	389-5547	-	73577-591	389-4078	-	73577-633	389-5026	-
V	73578-600	390-4554	-	73578-642	389-5554	-	73577-592	389-4086	-	73577-634	389-5034	-
VI	73578-601	390-4562	-	73578-643	389-5562	-	73577-593	389-4094	-	73577-635	389-5042	-
Thigh High - Bronze												
	RAL CCL 1 - Closed Toe, Sensitive Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	73578-674	390-4570	-	73578-680	389-5570	-	73577-666	389-4102	-	73577-672	389-5059	-
II	73578-675	390-4588	-	73578-681	389-5588	-	73577-667	389-4110	-	73577-673	389-5067	-
III	73578-676	390-4596	-	73578-682	389-5596	-	73577-668	389-4128	-	73577-674	389-5075	-
IV	73578-677	390-4604	-	73578-683	389-5604	-	73577-669	389-4136	-	73577-675	389-5083	-
V	73578-678	390-4612	-	73578-684	389-5612	-	73577-670	389-4144	-	73577-676	389-5091	-
VI	73578-679	390-4620	-	73578-685	389-5620	-	73577-671	389-4151	-	73577-677	389-5109	-
Thigh High - Espresso												
	RAL CCL 1 - Closed Toe, Sensitive Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	73578-584	390-4638	-	73578-626	389-5638	-	73577-576	389-4185	-	73577-618	389-5117	-
II	73578-585	390-4646	-	73578-627	389-5646	-	73577-577	389-4193	-	73577-619	389-5125	-
III	73578-586	390-4653	-	73578-628	389-5653	-	73577-578	389-4201	-	73577-620	389-5133	-
IV	73578-587	390-4661	-	73578-629	389-5661	-	73577-579	389-4219	-	73577-621	389-5141	-
V	73578-588	390-4679	-	73578-630	389-5679	-	73577-580	389-4227	-	73577-622	389-5158	-
VI	73578-589	390-4687	-	73578-631	389-5687	-	73577-581	389-4235	-	73577-623	389-5166	-
Thigh High - Caramel												
	RAL CCL 1 - Open Toe, Sensitive Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	73578-608	390-4869	-	73578-650	389-5869	-	73577-600	389-4441	-	73577-642	389-5349	-
II	73578-609	390-4877	-	73578-651	389-5877	-	73577-601	389-4458	-	73577-643	389-5356	-
III	73578-610	390-4885	-	73578-652	389-5885	-	73577-602	389-4466	-	73577-644	389-5364	-
IV	73578-611	390-4893	-	73578-653	389-5893	-	73577-603	389-4474	-	73577-645	389-5372	-
V	73578-612	390-4901	-	73578-654	389-5901	-	73577-604	389-4482	-	73577-646	389-5380	-
VI	73578-613	390-4919	-	73578-655	389-5919	-	73577-605	389-4490	-	73577-647	389-5398	-
Thigh High - Caramel												
	RAL CCL 1 - Closed Toe, Sensitive Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	73578-578	390-4695	-	73578-620	389-5695	-	73577-570	389-4243	-	73577-612	389-5174	-
II	73578-579	390-4703	-	73578-621	389-5703	-	73577-571	389-4250	-	73577-613	389-5182	-
III	73578-580	390-4711	-	73578-622	389-5711	-	73577-572	389-4268	-	73577-614	389-5190	-
IV	73578-581	390-4729	-	73578-623	389-5729	-	73577-573	389-4276	-	73577-615	389-5208	-
V	73578-582	390-4737	-	73578-624	389-5737	-	73577-574	389-4284	-	73577-616	389-5216	-
VI	73578-583	390-4745	-	73578-625	389-5745	-	73577-575	389-4292	-	73577-617	389-5224	-

JOBST® UltraSheer - CCL 1

Size	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC
Thigh High - Natural												
	RAL CCL 1 - Open Toe, Sensitive Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	73578-614	390-4927	-	73578-656	389-5927	-	73577-606	389-4516	-	73577-648	389-5406	-
II	73578-615	390-4935	-	73578-657	389-5935	-	73577-607	389-4524	-	73577-649	389-5414	-
III	73578-616	390-4943	-	73578-658	389-5943	-	73577-608	389-4532	-	73577-650	389-5422	-
IV	73578-617	390-4950	-	73578-659	389-5950	-	73577-609	389-4540	-	73577-651	389-5430	-
V	73578-618	390-4968	-	73578-660	389-5968	-	73577-610	389-4557	-	73577-652	389-5448	-
VI	73578-619	390-4976	-	73578-661	389-5976	-	73577-611	389-4565	-	73577-653	389-5455	-
Thigh High - Natural												
	RAL CCL 1 - Closed Toe, Sensitive Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	73578-590	390-4752	-	73578-632	389-5752	-	73577-582	389-4326	-	73577-624	389-5232	-
II	73578-591	390-4760	-	73578-633	389-5760	-	73577-583	389-4334	-	73577-625	389-5240	-
III	73578-592	390-4778	-	73578-634	389-5778	-	73577-584	389-4342	-	73577-626	389-5257	-
IV	73578-593	390-4786	-	73578-635	389-5786	-	73577-585	389-4359	-	73577-627	389-5265	-
V	73578-594	390-4794	-	73578-636	389-5794	-	73577-586	389-4367	-	73577-628	389-5273	-
VI	73578-595	390-4802	-	73578-637	389-5802	-	73577-587	389-4375	-	73577-629	389-5281	-
Thigh High - Anthracite												
	RAL CCL 1 - Closed Toe, Lace Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	73578-66	389-1611	-	73578-164	389-1629	-	73577-66	389-1553	-	73577-156	389-1561	-
II	73578-67	389-9721	-	73578-165	390-0248	-	73577-67	389-2478	-	73577-157	389-2866	-
III	73578-68	389-9739	-	73578-166	390-0255	-	73577-68	389-2486	-	73577-158	389-2874	-
IV	73578-69	389-9747	-	73578-167	390-0263	-	73577-69	389-2494	-	73577-159	389-2882	-
V	73578-70	389-9754	-	73578-168	390-0271	-	73577-70	389-2502	-	73577-160	389-2890	-
VI	73578-71	389-9762	-	73578-169	390-0289	-	73577-71	389-2510	-	73577-161	389-2908	-
Thigh High - Black												
	RAL CCL 1 - Closed Toe, Lace Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	73578-54	389-9770	-	73578-152	390-0297	-	73577-54	389-2528	-	73577-144	389-2916	-
II	73578-55	389-9788	-	73578-153	390-0305	-	73577-55	389-2536	-	73577-145	389-2924	-
III	73578-56	389-9796	-	73578-154	390-0313	-	73577-56	389-2544	-	73577-146	389-2932	-
IV	73578-57	389-9804	-	73578-155	390-0321	-	73577-57	389-2551	-	73577-147	389-2940	-
V	73578-58	389-9812	-	73578-156	390-0339	-	73577-58	389-2569	-	73577-148	389-2957	-
VI	73578-59	389-9820	-	73578-157	390-0347	-	73577-59	389-2577	-	73577-149	389-2965	-
Thigh High - Bronze												
	RAL CCL 1 - Closed Toe, Lace Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	73578-308	389-9838	-	73578-320	390-0354	-	73577-300	389-2585	-	73577-312	389-2999	-
II	73578-309	389-9846	-	73578-321	390-0362	-	73577-301	389-2593	-	73577-313	389-3005	-
III	73578-310	389-9853	-	73578-322	390-0370	-	73577-302	389-2601	-	73577-314	389-3013	-
IV	73578-311	389-9861	-	73578-323	390-0388	-	73577-303	389-2619	-	73577-315	389-3021	-
V	73578-312	389-9879	-	73578-324	390-0396	-	73577-304	389-2627	-	73577-316	389-3039	-
VI	73578-313	389-9887	-	73578-325	390-0404	-	73577-305	389-2635	-	73577-317	389-3047	-

JOBST® UltraSheer - CCL 1

Size	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC
Thigh High - Espresso												
	RAL CCL 1 - Closed Toe, Lace Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	73578-30	389-9895	-	73578-128	390-0412	-	73577-30	389-2650	-	73577-120	389-3096	-
II	73578-31	389-9903	-	73578-129	390-0420	-	73577-31	389-2668	-	73577-121	389-3104	-
III	73578-32	389-9911	-	73578-130	390-0438	-	73577-32	389-2676	-	73577-122	389-3112	-
IV	73578-33	389-9929	-	73578-131	390-0446	-	73577-33	389-2684	-	73577-123	389-3120	-
V	73578-34	389-9937	-	73578-132	390-0453	-	73577-34	389-2692	-	73577-124	389-3138	-
VI	73578-35	389-9945	-	73578-133	390-0461	-	73577-35	389-2700	-	73577-125	389-3146	-
Thigh High - Caramel												
	RAL CCL 1 - Closed Toe, Lace Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	73578-06	389-9952	-	73578-104	390-0479	-	73577-06	389-2726	-	73577-96	389-3153	-
II	73578-07	389-9960	-	73578-105	390-0487	-	73577-07	389-2734	-	73577-97	389-3161	-
III	73578-08	389-9978	-	73578-106	390-0495	-	73577-08	389-2742	-	73577-98	389-3179	-
IV	73578-09	389-9986	-	73578-107	390-0503	-	73577-09	389-2759	-	73577-99	389-3187	-
V	73578-10	389-9994	-	73578-108	390-0511	-	73577-10	389-2767	-	73577-100	389-3195	-
VI	73578-11	390-0008	-	73578-109	390-0529	-	73577-11	389-2775	-	73577-101	389-3203	-
Thigh High - Natural												
	RAL CCL 1 - Closed Toe, Lace Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	73578-42	390-0016	-	73578-140	390-0537	-	73577-42	389-2809	-	73577-132	389-3930	-
II	73578-43	390-0024	-	73578-141	390-0545	-	73577-43	389-2817	-	73577-133	389-3948	-
III	73578-44	390-0032	-	73578-142	390-0552	-	73577-44	389-2825	-	73577-134	389-3955	-
IV	73578-45	390-0040	-	73578-143	390-0560	-	73577-45	389-2833	-	73577-135	389-3963	-
V	73578-46	390-0057	-	73578-144	390-0578	-	73577-46	389-2841	-	73577-136	389-3971	-
VI	73578-47	390-0065	-	73578-145	390-0586	-	73577-47	389-2858	-	73577-137	389-3989	-
Tights - Anthracite							**Tights - Black**					
	RAL CCL 1 - Closed Toe						**RAL CCL 1 - Closed Toe**					
	Petite			**Regular**			**Petite**			**Regular**		
I	73582-60	389-1686	-	73581-60	389-1678	-	73582-48	389-9234	-	73581-48	389-8822	-
II	73582-61	389-9184	-	73581-61	389-8772	-	73582-49	389-9242	-	73581-49	389-8830	-
III	73582-62	389-9192	-	73581-62	389-8780	-	73582-50	389-9259	-	73581-50	389-8848	-
IV	73582-63	389-9200	-	73581-63	389-8798	-	73582-51	389-9267	-	73581-51	389-8855	-
V	73582-64	389-9218	-	73581-64	389-8806	-	73582-52	389-9275	-	73581-52	389-8863	-
VI	73582-65	389-9226	-	73581-65	389-8814	-	73582-53	389-9283	-	73581-53	389-8871	-
Tights - Bronze							**Tights - Espresso**					
	RAL CCL 1 - Closed Toe						**RAL CCL 1 - Closed Toe**					
	Petite			**Regular**			**Petite**			**Regular**		
I	73582-216	389-9291	-	73581-204	389-8889	-	73582-24	389-9358	-	73581-24	389-8947	-
II	73582-217	389-9309	-	73581-205	389-8897	-	73582-25	389-9366	-	73581-25	389-8954	-
III	73582-218	389-9317	-	73581-206	389-8905	-	73582-26	389-9374	-	73581-26	389-8962	-
IV	73582-219	389-9325	-	73581-207	389-8913	-	73582-27	389-9382	-	73581-27	389-8970	-
V	73582-220	389-9333	-	73581-208	389-8921	-	73582-28	389-9390	-	73581-28	389-8988	-
VI	73582-221	389-9341	-	73581-209	389-8939	-	73582-29	389-9408	-	73581-29	389-8996	-

JOBST® UltraSheer - CCL 1

Size	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC
	Tights - Caramel						**Tights - Natural**					
	RAL CCL 1 - Closed Toe						**RAL CCL 1 - Closed Toe**					
	Petite			**Regular**			**Petite**			**Regular**		
I	73582-00	389-9416	-	73581-00	389-9002	-	73582-36	389-9473	-	73581-36	389-9069	-
II	73582-01	389-9424	-	73581-01	389-9010	-	73582-37	389-9481	-	73581-37	389-9077	-
III	73582-02	389-9432	-	73581-02	389-9028	-	73582-38	389-9499	-	73581-38	389-9085	-
IV	73582-03	389-9440	-	73581-03	389-9036	-	73582-39	389-9507	-	73581-39	389-9093	-
V	73582-04	389-9457	-	73581-04	389-9044	-	73582-40	389-9515	-	73581-40	389-9101	-
VI	73582-05	389-9465	-	73581-05	389-9051	-	73582-41	389-9523	-	73581-41	389-9119	-
	Maternity Tights - Anthracite						**Maternity Tights - Black**					
	RAL CCL 1 - Closed Toe						**RAL CCL 1 - Closed Toe**					
	Petite			**Regular**			**Petite**			**Regular**		
I	73582-258	-	-	73581-246	-	-	73582-252	-	-	73581-240	-	-
II	73582-259	-	-	73581-247	-	-	73582-253	-	-	73581-241	-	-
III	73582-260	-	-	73581-248	-	-	73582-254	-	-	73581-242	-	-
IV	73582-261	-	-	73581-249	-	-	73582-255	-	-	73581-243	-	-
V	73582-262	-	-	73581-250	-	-	73582-256	-	-	73581-244	-	-
VI	73582-263	-	-	73581-251	-	-	73582-257	-	-	73581-245	-	-
	Maternity Tights - Bronze						**Maternity Tights - Espresso**					
	RAL CCL 1 - Closed Toe						**RAL CCL 1 - Closed Toe**					
	Petite			**Regular**			**Petite**			**Regular**		
I	73582-288	-	-	73581-276	-	-	73582-240	-	-	73581-228	-	-
II	73582-289	-	-	73581-277	-	-	73582-241	-	-	73581-229	-	-
III	73582-290	-	-	73581-278	-	-	73582-242	-	-	73581-230	-	-
IV	73582-291	-	-	73581-279	-	-	73582-243	-	-	73581-231	-	-
V	73582-292	-	-	73581-280	-	-	73582-244	-	-	73581-232	-	-
VI	73582-293	-	-	73581-281	-	-	73582-245	-	-	73581-233	-	-
	Maternity Tights - Caramel						**Maternity Tights - Natural**					
	RAL CCL 1 - Closed Toe						**RAL CCL 1 - Closed Toe**					
	Petite			**Regular**			**Petite**			**Regular**		
I	73582-234	-	-	73581-222	-	-	73582-246	-	-	73581-234	-	-
II	73582-235	-	-	73581-223	-	-	73582-247	-	-	73581-235	-	-
III	73582-236	-	-	73581-224	-	-	73582-248	-	-	73581-236	-	-
IV	73582-237	-	-	73581-225	-	-	73582-249	-	-	73581-237	-	-
V	73582-238	-	-	73581-226	-	-	73582-250	-	-	73581-238	-	-
VI	73582-239	-	-	73581-227	-	-	73582-251	-	-	73581-239	-	-

JOBST® UltraSheer - CCL 2

Size	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC
Knee High - Anthracite							**Knee High - Bronze**					
	RAL CCL 2 - Closed Toe						**RAL CCL 2 - Closed Toe**					
	Petite			**Regular**			**Petite**			**Regular**		
I	73584-42	389-3575	-	73583-42	389-3211	-	73584-115	389-3690	-	73583-102	389-3336	-
II	73584-43	389-3583	-	73583-43	389-3229	-	73584-116	389-3708	-	73583-103	389-3344	-
III	73584-44	389-3591	-	73583-44	389-3237	-	73584-117	389-3716	-	73583-104	389-3351	-
IV	73584-45	389-3609	-	73583-45	389-3245	-	73584-118	389-3724	-	73583-105	389-3369	-
V	73584-46	389-3617	-	73583-46	389-3252	-	73584-119	389-3732	-	73583-106	389-3377	-
VI	73584-47	389-3625	-	73583-47	389-3260	-	73584-120	389-3740	-	73583-107	389-3385	-
Knee High - Black												
	RAL CCL 2 - Open Toe						**RAL CCL 2 - Closed Toe**					
	Petite			**Regular**			**Petite**			**Regular**		
I	73584-52	389-4763	-	73583-48	389-4581	-	73584-36	389-3633	-	73583-36	389-3278	-
II	73584-48	389-4771	-	73583-49	389-4599	-	73584-37	389-3641	-	73583-37	389-3286	-
III	73584-49	389-4789	-	73583-50	389-4607	-	73584-38	389-3658	-	73583-38	389-3294	-
IV	73584-50	389-4797	-	73583-51	389-4615	-	73584-39	389-3666	-	73583-39	389-3302	-
V	73584-51	389-4805	-	73583-52	389-4623	-	73584-40	389-3674	-	73583-40	389-3310	-
VI	73584-53	389-4813	-	73583-53	389-4631	-	73584-41	389-3682	-	73583-41	389-3328	-
Knee High - Espresso												
	RAL CCL 2 - Closed Toe											
	Petite			**Regular**								
I	73584-18	389-3757	-	73583-18	389-3393	-						
II	73584-19	389-3765	-	73583-19	389-3401	-						
III	73584-20	389-3773	-	73583-20	389-3419	-						
IV	73584-21	389-3781	-	73583-21	389-3427	-						
V	73584-22	389-3799	-	73583-22	389-3435	-						
VI	73584-23	389-3807	-	73583-23	389-3443	-						
Knee High - Caramel												
	RAL CCL 2 - Open Toe						**RAL CCL 2 - Closed Toe**					
	Petite			**Regular**			**Petite**			**Regular**		
I	73584-00	389-4821	-	73583-00	389-4649	-	73584-06	389-3815	-	73583-06	389-3450	-
II	73584-01	389-4839	-	73583-01	389-4656	-	73584-07	389-3823	-	73583-07	389-3468	-
III	73584-02	389-4847	-	73583-02	389-4664	-	73584-08	389-3831	-	73583-08	389-3476	-
IV	73584-03	389-4854	-	73583-03	389-4672	-	73584-09	389-3849	-	73583-09	389-3484	-
V	73584-04	389-4862	-	73583-04	389-4680	-	73584-10	389-3856	-	73583-10	389-3492	-
VI	73584-05	389-4870	-	73583-05	389-4698	-	73584-11	389-3864	-	73583-11	389-3500	-
Knee High - Natural												
	RAL CCL 2 - Open Toe						**RAL CCL 2 - Closed Toe**					
	Petite			**Regular**			**Petite**			**Regular**		
I	73584-24	389-4888	-	73583-24	389-4706	EGD4918	73584-30	389-3872	-	73583-30	389-3518	-
II	73584-25	389-4896	-	73583-25	389-4714	EGD4919	73584-31	389-3880	-	73583-31	389-3526	-
III	73584-26	389-4904	-	73583-26	389-4722	EGD4939	73584-32	389-3898	-	73583-32	389-3534	-
IV	73584-27	389-4912	-	73583-27	389-4730	EGD4940	73584-33	389-3906	-	73583-33	389-3542	-
V	73584-28	389-4920	-	73583-28	389-4748	EGD4920	73584-34	389-3914	-	73583-34	389-3559	-
VI	73584-29	389-4938	-	73583-29	389-4755	EGD4977	73584-35	389-3922	-	73583-35	389-3567	-

JOBST® UltraSheer - CCL 2

Size	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC
Thigh High - Anthracite												
	RAL CCL 2 - Closed Toe, Decorative Dotted Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	73586-252	390-7292	-	73586-300	390-7656	-	73585-244	390-6211	-	73585-292	390-6757	-
II	73586-253	390-7300	-	73586-301	390-7664	-	73585-245	390-6229	-	73585-293	390-6765	-
III	73586-254	390-7318	-	73586-302	390-7672	-	73585-246	390-6237	-	73585-294	390-6773	-
IV	73586-255	390-7326	-	73586-303	390-7680	-	73585-247	390-6245	-	73585-295	390-6781	-
V	73586-256	390-7334	-	73586-304	390-7698	-	73585-248	390-6252	-	73585-296	390-6799	-
VI	73586-257	390-7342	-	73586-305	390-7706	-	73585-249	390-6260	-	73585-297	390-6807	-
Thigh High - Black												
	RAL CCL 2 - Open Toe, Decorative Dotted Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	73586-326	390-8019	-	73586-332	390-8191	-	73585-318	390-6575	-	73585-324	390-7110	-
II	73586-327	390-8027	-	73586-333	390-8209	-	73585-319	390-6583	-	73585-325	390-7128	-
III	73586-328	390-8035	-	73586-334	390-8217	-	73585-320	390-6591	-	73585-326	390-7136	-
IV	73586-329	390-8043	-	73586-335	390-8225	-	73585-321	390-6609	-	73585-327	390-7144	-
V	73586-330	390-8050	-	73586-336	390-8233	-	73585-322	390-6617	-	73585-328	390-7151	-
VI	73586-331	390-8068	-	73586-337	390-8241	-	73585-323	390-6625	-	73585-329	390-7169	-
Thigh High - Black												
	RAL CCL 2 - Closed Toe, Decorative Dotted Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	73586-246	390-7359	-	73586-294	390-7714	-	73585-238	390-6278	-	73585-286	390-6815	-
II	73586-247	390-7367	-	73586-295	390-7722	-	73585-239	390-6286	-	73585-287	390-6823	-
III	73586-248	390-7375	-	73586-296	390-7730	-	73585-240	390-6294	-	73585-288	390-6831	-
IV	73586-249	390-7383	-	73586-297	390-7748	-	73585-241	390-6302	-	73585-289	390-6849	-
V	73586-250	390-7391	-	73586-298	390-7755	-	73585-242	390-6310	-	73585-290	390-6856	-
VI	73586-251	390-7409	-	73586-299	390-7763	-	73585-243	390-6328	-	73585-291	390-6864	-
Thigh High - Bronze												
	RAL CCL 2 - Closed Toe, Decorative Dotted Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	73586-481	390-7417	-	73586-487	390-7771	-	73585-473	390-6336	-	73585-479	390-6872	-
II	73586-482	390-7425	-	73586-488	390-7789	-	73585-474	390-6344	-	73585-480	390-6880	-
III	73586-483	390-7433	-	73586-489	390-7797	-	73585-475	390-6351	-	73585-481	390-6898	-
IV	73586-484	390-7441	-	73586-490	390-7805	-	73585-476	390-6369	-	73585-482	390-6906	-
V	73586-485	390-7458	-	73586-491	390-7813	-	73585-477	390-6377	-	73585-483	390-6914	-
VI	73586-486	390-7466	-	73586-492	390-7821	-	73585-478	390-6385	-	73585-484	390-6922	-
Thigh High - Espresso												
	RAL CCL 2 - Closed Toe, Decorative Dotted Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	73586-234	390-7474	-	73586-282	390-7839	-	73585-226	390-6393	-	73585-274	390-6930	-
II	73586-235	390-7482	-	73586-283	390-7847	-	73585-227	390-6401	-	73585-275	390-6948	-
III	73586-236	390-7490	-	73586-284	390-7854	-	73585-228	390-6419	-	73585-276	390-6955	-
IV	73586-237	390-7508	-	73586-285	390-7862	-	73585-229	390-6427	-	73585-277	390-6963	-
V	73586-238	390-7516	-	73586-286	390-7870	-	73585-230	390-6435	-	73585-278	390-6971	-
VI	73586-239	390-7524	-	73586-287	390-7888	-	73585-231	390-6443	-	73585-279	390-6989	-

JOBST® UltraSheer - CCL 2

Size	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC
Thigh High - Caramel												
	RAL CCL 2 - Open Toe, Decorative Dotted Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	73586-210	390-8076	-	73586-258	390-8258	-	73585-202	390-6633	-	73585-250	390-7177	-
II	73586-211	390-8084	-	73586-259	390-8266	-	73585-203	390-6641	-	73585-251	390-7185	-
III	73586-212	390-8092	-	73586-260	390-8274	-	73585-204	390-6658	-	73585-252	390-7193	-
IV	73586-213	390-8100	-	73586-261	390-8282	-	73585-205	390-6666	-	73585-253	390-7201	-
V	73586-214	390-8118	-	73586-262	390-8290	-	73585-206	390-6674	-	73585-254	390-7219	-
VI	73586-215	390-8126	-	73586-263	390-8308	-	73585-207	390-6682	-	73585-255	390-7227	-
Thigh High - Caramel												
	RAL CCL 2 - Closed Toe, Decorative Dotted Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	73586-222	390-7532	-	73586-270	390-7896	-	73585-214	390-6450	-	73585-262	390-6997	-
II	73586-223	390-7540	-	73586-271	390-7904	-	73585-215	390-6468	-	73585-263	390-7003	-
III	73586-224	390-7557	-	73586-272	390-7912	-	73585-216	390-6476	-	73585-264	390-7011	-
IV	73586-225	390-7565	-	73586-273	390-7920	-	73585-217	390-6484	-	73585-265	390-7029	-
V	73586-226	390-7573	-	73586-274	390-7938	-	73585-218	390-6492	-	73585-266	390-7037	-
VI	73586-227	390-7581	-	73586-275	390-7946	-	73585-219	390-6500	-	73585-267	390-7045	-
Thigh High - Natural												
	RAL CCL 2 - Open Toe, Decorative Dotted Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	73586-216	390-8134	-	73586-264	390-8316	-	73585-208	390-6690	-	73585-256	390-7235	-
II	73586-217	390-8142	-	73586-265	390-8324	-	73585-209	390-6708	-	73585-257	390-7243	-
III	73586-218	390-8159	-	73586-266	390-8332	-	73585-210	390-6716	-	73585-258	390-7250	-
IV	73586-219	390-8167	-	73586-267	390-8340	-	73585-211	390-6724	-	73585-259	390-7268	-
V	73586-220	390-8175	-	73586-268	390-8357	-	73585-212	390-6732	-	73585-260	390-7276	-
VI	73586-221	390-8183	-	73586-269	390-8365	-	73585-213	390-6740	-	73585-261	390-7284	-
Thigh High - Natural												
	RAL CCL 2 - Closed Toe, Decorative Dotted Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	73586-240	390-7599	-	73586-288	390-7953	-	73585-232	390-6518	-	73585-280	390-7052	-
II	73586-241	390-7607	-	73586-289	390-7961	-	73585-233	390-6526	-	73585-281	390-7060	-
III	73586-242	390-7615	-	73586-290	390-7979	-	73585-234	390-6534	-	73585-282	390-7078	-
IV	73586-243	390-7623	-	73586-291	390-7987	-	73585-235	390-6542	-	73585-283	390-7086	-
V	73586-244	390-7631	-	73586-292	390-7995	-	73585-236	390-6559	-	73585-284	390-7094	-
VI	73586-245	390-7649	-	73586-293	390-8001	-	73585-237	390-6567	-	73585-285	390-7102	-
Thigh High - Anthracite												
	RAL CCL 2 - Closed Toe, Sensitive Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	73586-637	390-4984	-	73586-679	390-5494	-	73585-629	390-3366	-	73585-671	390-3721	-
II	73586-638	390-4992	-	73586-680	390-5502	-	73585-630	390-3374	-	73585-672	390-3739	-
III	73586-639	390-5007	-	73586-681	390-5510	-	73585-631	390-3382	-	73585-673	390-3747	-
IV	73586-640	390-5015	-	73586-682	390-5528	-	73585-632	390-3390	-	73585-674	390-3754	-
V	73586-641	390-5023	-	73586-683	390-5536	-	73585-633	390-3408	-	73585-675	390-3762	-
VI	73586-642	390-5031	-	73586-684	390-5544	-	73585-634	390-3416	-	73585-676	390-3770	-

JOBST® UltraSheer - CCL 2

Size	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC
Thigh High - Black												
	RAL CCL 2 - Open Toe, Sensitive Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	73586-685	390-5858	-	73586-691	390-6039	-	73585-677	390-4083	-	73585-683	390-4273	-
II	73586-686	390-5866	-	73586-692	390-6047	-	73585-678	390-4091	-	73585-684	390-4281	-
III	73586-687	390-5874	-	73586-693	390-6054	-	73585-679	390-4109	-	73585-685	390-4299	-
IV	73586-688	390-5882	-	73586-694	390-6062	-	73585-680	390-4117	-	73585-686	390-4307	-
V	73586-689	390-5890	-	73586-695	390-6070	-	73585-681	390-4125	-	73585-687	390-4315	-
VI	73586-690	390-5908	-	73586-696	390-6088	-	73585-682	390-4133	-	73585-688	390-4323	-
Thigh High - Black												
	RAL CCL 2 - Closed Toe, Sensitive Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	73586-631	390-5049	-	73586-673	390-5551	-	73585-623	390-3424	-	73585-665	390-3788	-
II	73586-632	390-5056	-	73586-674	390-5569	-	73585-624	390-3432	-	73585-666	390-3796	-
III	73586-633	390-5064	-	73586-675	390-5577	-	73585-625	390-3440	-	73585-667	390-3804	-
IV	73586-634	390-5072	-	73586-676	390-5585	-	73585-626	390-3457	-	73585-668	390-3812	-
V	73586-635	390-5080	-	73586-677	390-5593	-	73585-627	390-3465	-	73585-669	390-3820	-
VI	73586-636	390-5098	-	73586-678	390-5601	-	73585-628	390-3473	-	73585-670	390-3838	-
Thigh High - Bronze												
	RAL CCL 2 - Closed Toe, Sensitive Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	73586-702	390-5106	-	73586-708	390-5619	-	73585-694	390-3481	-	73585-700	390-3846	-
II	73586-703	390-5114	-	73586-709	390-5627	-	73585-695	390-3499	-	73585-701	390-3853	-
III	73586-704	390-5122	-	73586-710	390-5635	-	73585-696	390-3507	-	73585-702	390-3861	-
IV	73586-705	390-5130	-	73586-711	390-5643	-	73585-697	390-3515	-	73585-703	390-3879	-
V	73586-706	390-5148	-	73586-712	390-5650	-	73585-698	390-3523	-	73585-704	390-3887	-
VI	73586-707	390-5155	-	73586-713	390-5668	-	73585-699	390-3531	-	73585-705	390-3895	-
Thigh High - Espresso												
	RAL CCL 2 - Closed Toe, Sensitive Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	73586-619	390-5163	-	73586-661	390-5676	-	73585-611	390-3549	-	73585-653	390-3903	-
II	73586-620	390-5171	-	73586-662	390-5684	-	73585-612	390-3556	-	73585-654	390-3911	-
III	73586-621	390-5189	-	73586-663	390-5692	-	73585-613	390-3564	-	73585-655	390-3929	-
IV	73586-622	390-5197	-	73586-664	390-5700	-	73585-614	390-3572	-	73585-656	390-3937	-
V	73586-623	390-5205	-	73586-665	390-5718	-	73585-615	390-3580	-	73585-657	390-3945	-
VI	73586-624	390-5213	-	73586-666	390-5726	-	73585-616	390-3598	-	73585-658	390-3952	-
Thigh High - Caramel												
	RAL CCL 2 - Open Toe, Sensitive Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	73586-601	390-5916	-	73586-643	390-6096	-	73585-593	390-4141	-	73585-635	390-4331	-
II	73586-602	390-5924	-	73586-644	390-6104	-	73585-594	390-4158	-	73585-636	390-4349	-
III	73586-603	390-5932	-	73586-645	390-6112	-	73585-595	390-4166	-	73585-637	390-4356	-
IV	73586-604	390-5940	-	73586-646	390-6120	-	73585-596	390-4174	-	73585-638	390-4364	-
V	73586-605	390-5957	-	73586-647	390-6138	-	73585-597	390-4182	-	73585-639	390-4372	-
VI	73586-606	390-5965	-	73586-648	390-6146	-	73585-598	390-4190	-	73585-640	390-4380	-

JOBST® UltraSheer - CCL 2

Size	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC
Thigh High - Caramel												
	RAL CCL 2 - Closed Toe, Sensitive Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	73586-613	390-5221	-	73586-655	390-5734	-	73585-605	390-3606	-	73585-647	390-3960	-
II	73586-614	390-5239	-	73586-656	390-5742	-	73585-606	390-3614	-	73585-648	390-3978	-
III	73586-615	390-5247	-	73586-657	390-5759	-	73585-607	390-3622	-	73585-649	390-3986	-
IV	73586-616	390-5254	-	73586-658	390-5767	-	73585-608	390-3630	-	73585-650	390-3994	-
V	73586-617	390-5262	-	73586-659	390-5775	-	73585-609	390-3648	-	73585-651	390-4000	-
VI	73586-618	390-5270	-	73586-660	390-5783	-	73585-610	390-3655	-	73585-652	390-4018	-
Thigh High - Natural												
	RAL CCL 2 - Open Toe, Sensitive Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	73586-607	390-5973	-	73586-649	390-6153	-	73585-599	390-4208	-	73585-641	390-4398	-
II	73586-608	390-5981	-	73586-650	390-6161	-	73585-600	390-4216	-	73585-642	390-4406	-
III	73586-609	390-5999	-	73586-651	390-6179	-	73585-601	390-4224	-	73585-643	390-4414	-
IV	73586-610	390-6005	-	73586-652	390-6187	-	73585-602	390-4232	-	73585-644	390-4422	-
V	73586-611	390-6013	-	73586-653	390-6195	-	73585-603	390-4240	-	73585-645	390-4430	-
VI	73586-612	390-6021	-	73586-654	390-6203	-	73585-604	390-4257	-	73585-646	390-4448	-
Thigh High - Natural												
	RAL CCL 2 - Open Toe, Decorative Dotted Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	73586-1542	390-5288	-	73586-1548	390-5791	-	73585-1534	390-3663	-	73585-1540	390-4026	-
II	73586-1543	390-5296	-	73586-1549	390-5809	-	73585-1535	390-3671	-	73585-1541	390-4034	-
III	73586-1544	390-5304	-	73586-1550	390-5817	-	73585-1536	390-3689	-	73585-1542	390-4042	-
IV	73586-1545	390-5312	-	73586-1551	390-5825	-	73585-1537	390-3697	-	73585-1543	390-4059	-
V	73586-1546	390-5320	-	73586-1552	390-5833	-	73585-1538	390-3705	-	73585-1544	390-4067	-
VI	73586-1547	390-5338	-	73586-1553	390-5841	-	73585-1539	390-3713	-	73585-1545	390-4075	-
Thigh High - Anthracite												
	RAL CCL 2 - Closed Toe, Lace Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	73586-90	390-2020	-	73586-204	390-2384	-	73585-90	390-0610	-	73585-196	390-1402	-
II	73586-91	390-2038	-	73586-205	390-2392	-	73585-91	390-0628	-	73585-197	390-1410	-
III	73586-92	390-2046	-	73586-206	390-2400	-	73585-92	390-0636	-	73585-198	390-1428	-
IV	73586-93	390-2053	-	73586-207	390-2418	-	73585-93	390-0644	-	73585-199	390-1436	-
V	73586-94	390-2061	-	73586-208	390-2426	-	73585-94	390-0651	-	73585-200	390-1444	-
VI	73586-95	390-2079	-	73586-209	390-2434	-	73585-95	390-0669	-	73585-201	390-1451	-
Thigh High - Black												
	RAL CCL 2 - Open Toe, Lace Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	73586-308	390-2988	-	73586-320	390-3168	-	73585-300	390-1022	-	73585-312	390-1204	-
II	73586-309	390-2996	-	73586-321	390-3176	-	73585-301	390-1030	-	73585-313	390-1212	-
III	73586-310	390-3002	-	73586-322	390-3184	-	73585-302	390-1048	-	73585-314	390-1220	-
IV	73586-311	390-3010	-	73586-323	390-3192	-	73585-303	390-1055	-	73585-315	390-1238	-
V	73586-312	390-3028	-	73586-324	390-3200	-	73585-304	390-1063	-	73585-316	390-1246	-
VI	73586-313	390-3036	-	73586-325	390-3218	-	73585-305	390-1071	-	73585-317	390-1253	-

JOBST® UltraSheer - CCL 2

Size	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC
Thigh High - Black												
	RAL CCL 2 - Closed Toe, Lace Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	73586-78	390-2087	-	73586-192	390-2442	-	73585-78	390-0677	-	73585-184	390-1469	-
II	73586-79	390-2095	-	73586-193	390-2459	-	73585-79	390-0685	-	73585-185	390-1477	-
III	73586-80	390-2103	-	73586-194	390-2467	-	73585-80	390-0693	-	73585-186	390-1485	-
IV	73586-81	390-2111	-	73586-195	390-2475	-	73585-81	390-0701	-	73585-187	390-1493	-
V	73586-82	390-2129	-	73586-196	390-2483	-	73585-82	390-0719	-	73585-188	390-1501	-
VI	73586-83	390-2137	-	73586-197	390-2491	-	73585-83	390-0727	-	73585-189	390-1519	-
Thigh High - Bronze												
	RAL CCL 2 - Closed Toe, Lace Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	73586-463	390-2145	-	73586-475	390-2509	-	73585-455	390-0735	-	73585-467	390-1667	-
II	73586-464	390-2152	-	73586-476	390-2517	-	73585-456	390-0743	-	73585-468	390-1675	-
III	73586-465	390-2160	-	73586-477	390-2525	-	73585-457	390-0750	-	73585-469	390-1683	-
IV	73586-466	390-2178	-	73586-478	390-2533	-	73585-458	390-0768	-	73585-470	390-1691	-
V	73586-467	390-2186	-	73586-479	390-2541	-	73585-459	390-0776	-	73585-471	390-1709	-
VI	73586-468	390-2194	-	73586-480	390-2558	-	73585-460	390-0784	-	73585-472	390-1717	-
Thigh High - Espresso												
	RAL CCL 2 - Closed Toe, Lace Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	73586-42	390-2202	-	73586-156	390-2566	-	73585-42	390-0792	-	73585-148	390-1725	-
II	73586-43	390-2210	-	73586-157	390-2574	-	73585-43	390-0800	-	73585-149	390-1733	-
III	73586-44	390-2228	-	73586-158	390-2582	-	73585-44	390-0818	-	73585-150	390-1741	-
IV	73586-45	390-2236	-	73586-159	390-2590	-	73585-45	390-0826	-	73585-151	390-1758	-
V	73586-46	390-2244	-	73586-160	390-2608	-	73585-46	390-0834	-	73585-152	390-1766	-
VI	73586-47	390-2251	-	73586-161	390-2616	-	73585-47	390-0842	-	73585-153	390-1774	-
Thigh High - Caramel												
	RAL CCL 2 - Open Toe, Lace Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	73586-06	390-3044	-	73586-120	390-3226	-	73585-06	390-1089	-	73585-112	390-1261	-
II	73586-07	390-3051	-	73586-121	390-3234	-	73585-07	390-1097	-	73585-113	390-1279	-
III	73586-08	390-3069	-	73586-122	390-3242	-	73585-08	390-1105	-	73585-114	390-1287	-
IV	73586-09	390-3077	-	73586-123	390-3259	-	73585-09	390-1113	-	73585-115	390-1295	-
V	73586-10	390-3085	-	73586-124	390-3267	-	73585-10	390-1121	-	73585-116	390-1303	-
VI	73586-11	390-3093	-	73586-125	390-3275	-	73585-11	390-1139	-	73585-117	390-1311	-
Thigh High - Caramel												
	RAL CCL 2 - Closed Toe, Lace Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	73586-18	390-2269	-	73586-132	390-2624	-	73585-18	390-0859	-	73585-124	390-1782	-
II	73586-19	390-2277	-	73586-133	390-2632	-	73585-19	390-0867	-	73585-125	390-1790	-
III	73586-20	390-2285	-	73586-134	390-2640	-	73585-20	390-0875	-	73585-126	390-1808	-
IV	73586-21	390-2293	-	73586-135	390-2657	-	73585-21	390-0883	-	73585-127	390-1816	-
V	73586-22	390-2301	-	73586-136	390-2665	-	73585-22	390-0891	-	73585-128	390-1824	-
VI	73586-23	390-2319	-	73586-137	390-2673	-	73585-23	390-0909	-	73585-129	390-1832	-

JOBST® UltraSheer - CCL 2

Size	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC
Thigh High - Natural												
	RAL CCL 2 - Open Toe, Lace Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	-	390-3101	EGD7931	73586-168	390-3283	-	73585-54	390-1147	-	73585-160	390-1329	-
II	-	390-3119	EGD7932	73586-169	390-3291	-	73585-55	390-1154	-	73585-161	390-1337	-
III	-	390-3127	EGD4885	73586-170	390-3309	-	73585-56	390-1162	-	73585-162	390-1345	-
IV	-	390-3135	EGD7933	73586-171	390-3317	-	73585-57	390-1170	-	73585-163	390-1352	-
V	-	390-3143	EGD7934	73586-172	390-3325	-	73585-58	390-1188	-	73585-164	390-1360	-
VI	-	390-3150	EGD7935	73586-173	390-3333	-	73585-59	390-1196	-	73585-165	390-1378	-
Thigh High - Natural												
	RAL CCL 2 - Closed Toe, Lace Silicone Band											
	Petite			**Petite, wider circumference**			**Regular**			**Regular, wider circumference**		
I	73586-66	390-2327	-	73586-180	390-2681	-	73585-66	390-0917	-	73585-172	390-1840	-
II	73586-67	390-2335	-	73586-181	390-2699	-	73585-67	390-0925	-	73585-173	390-1857	-
III	73586-68	390-2343	-	73586-182	390-2707	-	73585-68	390-0933	-	73585-174	390-1865	-
IV	73586-69	390-2350	-	73586-183	390-2715	-	73585-69	390-0941	-	73585-175	390-1873	-
V	73586-70	390-2368	-	73586-184	390-2723	-	73585-70	390-0958	-	73585-176	390-1881	-
VI	73586-71	390-2376	-	73586-185	390-2731	-	73585-71	390-0966	-	73585-177	390-1899	-
Tights - Anthracite							**Tights - Bronze**					
	RAL CCL 2 - Closed Toe						**RAL CCL 2 - Closed Toe**					
	Petite			**Regular**			**Petite**			**Regular**		
I	73590-84	389-6586	-	73589-84	389-6040	-	73590-318	389-6701	-	73589-306	389-6222	-
II	73590-85	389-6594	-	73589-85	389-6057	-	73590-319	389-6719	-	73589-307	389-6230	-
III	73590-86	389-6602	-	73589-86	389-6065	-	73590-320	389-6727	-	73589-308	389-6248	-
IV	73590-87	389-6610	-	73589-87	389-6073	-	73590-321	389-6735	-	73589-309	389-6255	-
V	73590-88	389-6628	-	73589-88	389-6081	-	73590-322	389-6743	-	73589-310	389-6263	-
VI	73590-89	389-6636	-	73589-89	389-6099	-	73590-323	389-6750	-	73589-311	389-6271	-
Tights - Black												
	RAL CCL 2 - Open Toe						**RAL CCL 2 - Closed Toe**					
	Petite			**Regular**			**Petite**			**Regular**		
I	73590-96	389-6941	-	73589-96	389-6164	-	73590-72	389-6644	-	73589-72	389-6107	-
II	73590-97	389-6958	-	73589-97	389-6172	-	73590-73	389-6651	-	73589-73	389-6115	-
III	73590-98	389-6966	-	73589-98	389-6180	-	73590-74	389-6669	-	73589-74	389-6123	-
IV	73590-99	389-6974	-	73589-99	389-6198	-	73590-75	389-6677	-	73589-75	389-6131	-
V	73590-100	389-6982	-	73589-100	389-6206	-	73590-76	389-6685	-	73589-76	389-6149	-
VI	73590-101	389-6990	-	73589-101	389-6214	-	73590-77	389-6693	-	73589-77	389-6156	-
Tights - Espresso												
	RAL CCL 2 - Closed Toe											
	Petite			**Regular**								
I	73590-36	389-6768	-	73589-36	389-6289	-						
II	73590-37	389-6776	-	73589-37	389-6297	-						
III	73590-38	389-6784	-	73589-38	389-6305	-						
IV	73590-39	389-6792	-	73589-39	389-6313	-						
V	73590-40	389-6800	-	73589-40	389-6321	-						
VI	73590-41	389-6818	-	73589-41	389-6339	-						

JOBST® UltraSheer - CCL 2

Size	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC
Tights - Caramel												
	RAL CCL 2 - Open Toe						**RAL CCL 2 - Closed Toe**					
	Petite			**Regular**			**Petite**			**Regular**		
I	73590-00	389-7006	-	73589-00	389-6461	-	73590-06	389-6826	-	73589-06	389-6347	-
II	73590-01	389-7014	-	73589-01	389-6479	-	73590-07	389-6834	-	73589-07	389-6354	-
III	73590-02	389-7022	-	73589-02	389-6487	-	73590-08	389-6842	-	73589-08	389-6362	-
IV	73590-03	389-7030	-	73589-03	389-6495	-	73590-09	389-6859	-	73589-09	389-6370	-
V	73590-04	389-7048	-	73589-04	389-6503	-	73590-10	389-6867	-	73589-10	389-6388	-
VI	73590-05	389-7055	-	73589-05	389-6511	-	73590-11	389-6875	-	73589-11	389-6396	-
Tights - Natural												
	RAL CCL 2 - Open Toe						**RAL CCL 2 - Closed Toe**					
	Petite			**Regular**			**Petite**			**Regular**		
I	73590-48	389-7063	-	73589-48	389-6529	-	73590-54	389-6883	-	73589-54	389-6404	-
II	73590-49	389-7071	-	73589-49	389-6537	-	73590-55	389-6891	-	73589-55	389-6412	-
III	73590-50	389-7089	-	73589-50	389-6545	-	73590-56	389-6909	-	73589-56	389-6420	-
IV	73590-51	389-7097	-	73589-51	389-6552	-	73590-57	389-6917	-	73589-57	389-6438	-
V	73590-52	389-7105	-	73589-52	389-6560	-	73590-58	389-6925	-	73589-58	389-6446	-
VI	73590-53	389-7113	-	73589-53	389-6578	-	73590-59	389-6933	-	73589-59	389-6453	-
Maternity Tights - Anthracite							**Maternity Tights - Bronze**					
	RAL CCL 2 - Closed Toe						**RAL CCL 2 - Closed Toe**					
	Petite			**Regular**			**Petite**			**Regular**		
I	73590-450	-	-	73589-438	-	-	73590-486	-	-	73589-474	-	-
II	73590-451	-	-	73589-439	-	-	73590-487	-	-	73589-475	-	-
III	73590-452	-	-	73589-440	-	-	73590-488	-	-	73589-476	-	-
IV	73590-453	-	-	73589-441	-	-	73590-489	-	-	73589-477	-	-
V	73590-454	-	-	73589-442	-	-	73590-490	-	-	73589-478	-	-
VI	73590-455	-	-	73589-443	-	-	73590-491	-	-	73589-479	-	-
Maternity Tights - Black												
	RAL CCL 2 - Open Toe						**RAL CCL 2 - Closed Toe**					
	Petite			**Regular**			**Petite**			**Regular**		
I	73590-456	-	-	73589-444	-	-	73590-444	-	-	73589-432	-	-
II	73590-457	-	-	73589-445	-	-	73590-445	-	-	73589-433	-	-
III	73590-458	-	-	73589-446	-	-	73590-446	-	-	73589-434	-	-
IV	73590-459	-	-	73589-447	-	-	73590-447	-	-	73589-435	-	-
V	73590-460	-	-	73589-448	-	-	73590-448	-	-	73589-436	-	-
VI	73590-461	-	-	73589-449	-	-	73590-449	-	-	73589-437	-	-
Maternity Tights - Espresso												
	RAL CCL 2 - Closed Toe											
	Petite			**Regular**								
I	73590-426	-	-	73589-414	-	-						
II	73590-427	-	-	73589-415	-	-						
III	73590-428	-	-	73589-416	-	-						
IV	73590-429	-	-	73589-417	-	-						
V	73590-430	-	-	73589-418	-	-						
VI	73590-431	-	-	73589-419	-	-						

JOBST® UltraSheer - CCL 2

Size	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC	Code	PIP	NHS SC
Maternity Tights - Caramel												
	RAL CCL 2 - Open Toe						**RAL CCL 2 - Closed Toe**					
	Petite			**Regular**			**Petite**			**Regular**		
I	73590-414	-	-	73589-402	-	-	73590-420	-	-	73589-408	-	-
II	73590-415	-	-	73589-403	-	-	73590-421	-	-	73589-409	-	-
III	73590-416	-	-	73589-404	-	-	73590-422	-	-	73589-410	-	-
IV	73590-417	-	-	73589-405	-	-	73590-423	-	-	73589-411	-	-
V	73590-418	-	-	73589-406	-	-	73590-424	-	-	73589-412	-	-
VI	73590-419	-	-	73589-407	-	-	73590-425	-	-	73589-413	-	-
Maternity Tights - Natural												
	RAL CCL 2 - Open Toe						**RAL CCL 2 - Closed Toe**					
	Petite			**Regular**			**Petite**			**Regular**		
I	73590-432	-	-	73589-420	-	-	73590-438	-	-	73589-426	-	-
II	73590-433	-	-	73589-421	-	-	73590-439	-	-	73589-427	-	-
III	73590-434	-	-	73589-422	-	-	73590-440	-	-	73589-428	-	-
IV	73590-435	-	-	73589-423	-	-	73590-441	-	-	73589-429	-	-
V	73590-436	-	-	73589-424	-	-	73590-442	-	-	73589-430	-	-
VI	73590-437	-	-	73589-425	-	-	73590-443	-	-	73589-431	-	-



JOBST® UlcerCare

Ready-to-Wear and Custom-Fit

Available on Drug Tariff

Innovative two-in-one compression system for the progressive management of venous leg ulcers.

Nominal compression of 40mmHg* at the ankle



JOBST® UlcerCare consists of two components, a medical stocking and a compression liner. Together these provide a nominal compression of 40mmHg at the ankle. JOBST® UlcerCare delivers a high stiffness factor to provide effective compression for actively healing leg ulcers, whist continuing to encourage the supported self-management of patients. JOBST® UlcerCare kit and JOBST® UlcerCare Liner comply with the Lower Limb Reccomendations from the National Wound Care Strategy Program.[1]

Compression
40mmHg* at the ankle

Ready-to-Wear - Sizes, Options, Materials			
Seven Sizes	Choice of Zipper	Left or Right Zipper	Zipper Pull Tag
✓	✓	✓	✓


Custom-Fit - Sizes, Options, Materials				
Custom Size	Straight Open Toe	Slant Open Toe	Lateral Zipper Position	Medial Zipper Position
✓	✓	✓	✓	✓

Colours - Liner	Colours - Outer Stocking	
White	Beige	Black

Styles	
	
White Liner (Contains Silk)	Outer Stocking - With or Without Zipper

*Targeted nominal compression for an average ankle size.

Size	S	M	L	XL	2XL	3XL	4XL
Ankle	18 – 21 cm	21 – 24 cm	24 – 27 cm	27 – 30 cm	30 – 33 cm	33 – 36 cm	36 – 39 cm
Calf	29 – 35 cm	34 – 41 cm	40 – 47 cm	46 – 54 cm	48 – 56 cm	50 – 58 cm	52 – 60 cm



1 National Wound Care Strategy Programme (2020), Lower Limb: Recommendations for Clinical Care, NWCSP.

JOBST® UlcerCare

Code	Size	PIP	Contents Stocking/ Liner	NHS SC	Code	Size	PIP	Contents Stocking/ Liner	NHS SC
Ready-to-Wear Stocking and Liner									
Black					**Beige**				
73630-41	S	364-4556	1/2	-	73630-21	S	364-4622	1/2	-
73630-42	M	364-4549	1/2	-	73630-22	M	364-4614	1/2	-
73630-43	L	364-4531	1/2	-	73630-23	L	364-4606	1/2	-
73630-44	XL	364-4523	1/2	-	73630-24	XL	364-4598	1/2	-
73630-45	2XL	364-4515	1/2	-	73630-25	2XL	364-4580	1/2	-
73630-46	3XL	364-4507	1/2	-	73630-26	3XL	364-4572	1/2	-
73630-47	4XL	364-4499	1/2	-	73630-27	4XL	364-4564	1/2	-
Ready-to-Wear Zippered Stocking and Liner									
Zipper to be placed on side opposite to the ulceration									
Zipper Left - Black					**Zipper Left - Beige**				
73631-41	S	364-1263	1/2	-	73631-21	S	364-1271	1/2	-
73631-42	M	364-1248	1/2	-	73631-22	M	364-1255	1/2	-
73631-43	L	364-1222	1/2	-	73631-23	L	364-1230	1/2	-
73631-44	XL	364-1206	1/2	-	73631-24	XL	364-1214	1/2	-
73631-45	2XL	364-1172	1/2	-	73631-25	2XL	364-1180	1/2	-
73631-46	3XL	364-1156	1/2	-	73631-26	3XL	364-1164	1/2	-
73631-47	4XL	364-1131	1/2	-	73631-27	4XL	364-1149	1/2	-
Zipper Right - Black					**Zipper Right - Beige**				
73631-71	S	364-1404	1/2	-	73631-61	S	364-1412	1/2	-
73631-72	M	364-1388	1/2	-	73631-62	M	364-1396	1/2	-
73631-73	L	364-1362	1/2	-	73631-63	L	364-1370	1/2	-
73631-74	XL	364-1347	1/2	-	73631-64	XL	364-1354	1/2	-
73631-75	2XL	364-1321	1/2	-	73631-65	2XL	364-1339	1/2	-
73631-76	3XL	364-1305	1/2	-	73631-66	3XL	364-1313	1/2	-
73631-77	4XL	364-1289	1/2	-	73631-67	4XL	364-1297	1/2	-
Ready-to-Wear									
Replacement Liners - White									
73632-21	S	364-7666	3	-	73632-25	2XL	364-1123	3	-
73632-22	M	364-7674	3	-	73632-26	3XL	364-1115	3	-
73632-23	L	364-7682	3	-	73632-27	4XL	364-1107	3	-
73632-24	XL	364-7708	3	-					

Custom-Fit	
Black or Caramel	
LZ-07-04	Medical stocking
LL-07-04	Pack of two liners



JOBST® FarrowWrap®

Ready-to-Wear Wrap Compression System

Available on Drug Tariff

JOBST® FarrowWrap® garments are a range of adjustable wrap compression systems that are available in a variety of styles, fabrics, sizes and compression classes. Made with durable materials, they are designed to help reduce limb swelling and maintain volume reductions, as well as promote healing of venous leg ulcers. They are designed to provide the same strong support as multi-layer bandaging but in a simplified, intuitive design.

Compression Levels	
Lite = 20-30mmHg	**Strong and Classic** = 30-40mmHg
✓	✓



Sizes and Options / Materials				
Five Sizes	Regular	Tall / Long	Short	Made-to-Measure*
✓	✓	✓	✓	✓


Classic


Strong


Lite

Colours	
Tan	Black**

Styles					
					
Arm	Gaunlet	Foot	Leg	Thigh	Liners*

*Not available on FP10 / GP10
**Armpiece only

JOBST® FarrowWrap® LITE

Upper Limb, Ready-to-Wear Wrap Compression System

Compression: **20-30mmHg**

Mild to moderate lymphoedema and chronic oedema.

A double-laminated fabric with built in liner is especially suitable for palliative and paediatric patients. The fabric is both durable and breathable.





Code	Option	Size	Length	Colour	PIP Code
Armpiece					
76650-11	Left	Small	Short	Tan	392-6953
76650-05	Left	Small	Regular	Tan	392-7019
76650-35	Left	Small	Long	Tan	392-7076
76650-13	Left	Medium	Short	Tan	392-6961
76650-07	Left	Medium	Regular	Tan	392-7027
76650-01	Left	Medium	Long	Tan	392-7084
76650-15	Left	Large	Short	Tan	392-6979
76650-09	Left	Large	Regular	Tan	392-7035
76650-03	Left	Large	Long	Tan	392-7092
76650-29	Left	Small	Short	Black	392-6987
76650-23	Left	Small	Regular	Black	392-7043
76650-17	Left	Small	Long	Black	392-7100
76650-31	Left	Medium	Short	Black	392-6995
76650-25	Left	Medium	Regular	Black	392-7050
76650-19	Left	Medium	Long	Black	392-7118
76650-33	Left	Large	Short	Black	392-7001
76650-27	Left	Large	Regular	Black	392-7068
76650-21	Left	Large	Long	Black	392-7126
76650-12	Right	Small	Short	Tan	392-6771
76650-06	Right	Small	Regular	Tan	392-6839
76650-00	Right	Small	Long	Tan	392-6896
76650-14	Right	Medium	Short	Tan	392-6789
76650-08	Right	Medium	Regular	Tan	392-6847
76650-02	Right	Medium	Long	Tan	392-6904
76650-16	Right	Large	Short	Tan	392-6797
76650-10	Right	Large	Regular	Tan	392-6854
76650-04	Right	Large	Long	Tan	392-6912
76650-30	Right	Small	Short	Black	392-6805
76650-24	Right	Small	Regular	Black	392-6862
76650-18	Right	Small	Long	Black	392-6920
76650-32	Right	Medium	Short	Black	392-6813
76650-26	Right	Medium	Regular	Black	392-6870
76650-20	Right	Medium	Long	Black	392-6938
76650-34	Right	Large	Short	Black	392-6821
76650-28	Right	Large	Regular	Black	392-6888
76650-22	Right	Large	Long	Black	392-6946
Gauntlet					
76653-00	Ambidextrous	Extra small	-	Tan	404-5647
76653-01	Ambidextrous	Small	-	Tan	404-5654
76653-02	Ambidextrous	Medium	-	Tan	404-5662
76653-03	Ambidextrous	Large	-	Tan	404-5670

Note: if your patient does not fit within the RTW sizing chart for LITE or requires a STRONG or CLASSIC armpiece, these are available as a custom-fit variant. Further details available on request.

JOBST® FarrowWrap®

Ready-to-wear are pre-made pieces to fit the majority of patient sizes. Most variants are fully reimbursable on FP10 / GP10. The sizing charts below can be used to determine if your patient fits within a ready-to-wear size. If not, a custom-fit (CF) option is available. CF means pieces are manufactured to the limb measurements of each patient to provide a better fit for those limb shapes that do not fit within the RTW size range.

Please note: custom-fit is not available on FP10 / GP10.


JOBST® FarrowWrap® Armpiece



	Wrist	Elbow	Axilla		Length
Small	14 - 18cm	20 - 27cm	22 - 31cm	Short	40 - 43cm
Medium	16 - 21cm	25 - 34cm	29 - 39cm	Regular	44 - 47cm
Large	19 - 25cm	30 - 40cm	32 - 45cm	Long	48 - 51cm

NOTE: Measure length from the outside of the wrist to the axilla (G) with a slight bend in arm


JOBST® FarrowWrap® Gauntlet



Measurement Points	X Small	Small	Medium	Large
a - Palm at base of fingers	< 18cm	18 - 20cm	20 - 22cm	> 22cm
b - Palm at fold of thumb	< 19cm	19 - 21cm	21 - 23cm	> 24cm


Accessories

Foam Squares Pack*



Code	Size
76675-01	Small - 8 x 8cm
76675-02	Large - 10 x 10cm

GarmentGrip Roll*



Code	Size
76672-00	14 x 80cm
76672-01	14 x 120cm
76672-02	14 x 550cm
76672-03	20 x 120cm
76672-04	30 x 150cm
76672-05	10 x 60cm

*Not available on FP10

JOBST® FarrowWrap® LITE

Lower Limb, Ready-to-Wear Wrap Compression System

Compression: **20-30mmHg**

Mild to moderate lymphoedema / chronic oedema and chronic venous insufficiency.

A double-laminated fabric with built in liner is especially suitable for prevention of venous leg ulcer recurrence. The fabric is both durable and breathable.





Code	Size	Length	Colour	NHS SC Code	PIP Code
Footpiece					
76652-05	Extra small	Regular	Tan	EGD8743	360-6332
76652-00	Extra small	Long	Tan	EGD8738	360-6381
76652-06	Small	Regular	Tan	EGD8744	360-6340
76652-01	Small	Long	Tan	EGD8739	360-6399
76652-07	Medium	Regular	Tan	EGD8745	360-6357
76652-02	Medium	Long	Tan	EGD8740	360-6407
76652-08	Large	Regular	Tan	EGD8746	360-6365
76652-03	Large	Long	Tan	EGD8741	360-6415
76652-09	Extra large	Regular	Tan	EGD8747	360-6373
76652-04	Extra large	Long	Tan	EGD8742	360-6423
Legpiece					
76654-00	Extra small	Regular	Tan	EGD8748	360-6233
76654-05	Extra small	Tall	Tan	EGD8753	360-6282
76654-01	Small	Regular	Tan	EGD8749	360-6241
76654-06	Small	Tall	Tan	EGD8754	360-6290
76654-02	Medium	Regular	Tan	EGD8750	360-6258
76654-07	Medium	Tall	Tan	EGD8755	360-6308
76654-03	Large	Regular	Tan	EGD8751	360-6266
76654-08	Large	Tall	Tan	EGD8756	360-6316
76654-04	Extra large	Regular	Tan	EGD8752	360-6274
76654-09	Extra large	Tall	Tan	EGD8757	360-6324
Kneepiece (seperate)*					
73367-00	Extra small	-	Tan	-	-
73367-01	Small	-	Tan	-	-
73367-02	Medium	-	Tan	-	-
73367-03	Large	-	Tan	-	-
73367-04	Extra large	-	Tan	-	-
Thighpiece (including kneepiece)					
76655-05	Extra small	Short	Tan	EGD8763	360-6084
76655-00	Extra small	Regular	Tan	EGD8758	360-6134
76655-10	Extra small	Tall	Tan	EGD8768	360-6183
76655-06	Small	Short	Tan	EGD8764	360-6092
76655-01	Small	Regular	Tan	EGD8759	360-6142
76655-11	Small	Tall	Tan	EGD8769	360-6191
76655-07	Medium	Short	Tan	EGD8765	360-6100
76655-02	Medium	Regular	Tan	EGD8760	360-6159
76655-12	Medium	Tall	Tan	EGD8770	360-6209
76655-08	Large	Short	Tan	EGD8766	360-6118
76655-03	Large	Regular	Tan	EGD8761	360-6167
76655-13	Large	Tall	Tan	EGD8771	360-6217
76655-09	Extra large	Short	Tan	EGD8767	360-6126
76655-04	Extra large	Regular	Tan	EGD8762	360-6175
76655-14	Extra large	Tall	Tan	EGD8772	360-6225

Note: if your patient does not fit within the RTW sizing chart, LITE is also available as a custom-fit variant. Further details available on request. *Not available on FP10 / GP10.

JOBST® FarrowWrap® Strong

Lower Limb, Ready-to-Wear Wrap Compression System

Compression: **30-40mmHg**

Moderate to severe lymphoedema / chronic oedema with some limb shape distortion, chronic venous insufficiency and venous leg ulceration.

Premium fabric with strong compression and a soft inner layer suitable for patients with skin folds, some shape distortion and moderate to severe oedema.





Code	Size	Length	Colour	NHS SC Code	PIP Code
Footpiece					
76646-04	Extra small	Regular	Tan	EGD8803	381-9976
76646-09	Extra small	Long	Tan	EGD8808	381-9927
76646-05	Small	Regular	Tan	EGD8804	382-0008
76646-00	Small	Long	Tan	EGD8799	381-9950
76646-06	Medium	Regular	Tan	EGD8834	381-9992
76646-01	Medium	Long	Tan	EGD8800	381-9943
76646-07	Large	Regular	Tan	EGD8806	381-9984
76646-02	Large	Long	Tan	EGD8801	381-9935
76646-08	Extra large	Regular	Tan	EGD8807	381-9968
76646-03	Extra large	Long	Tan	EGD8802	381-9919
Legpiece					
76648-00	Extra small	Regular	Tan	EGD8809	382-0073
76648-05	Extra small	Tall	Tan	EGD8814	382-0024
76648-01	Small	Regular	Tan	EGD8810	382-0107
76648-06	Small	Tall	Tan	EGD8815	382-0057
76648-02	Medium	Regular	Tan	-	382-0099
76648-07	Medium	Tall	Tan	EGD8816	382-0040
76648-03	Large	Regular	Tan	EGD8812	382-0081
76648-08	Large	Tall	Tan	EGD8817	382-0032
76648-04	Extra large	Regular	Tan	EGD8813	382-0065
76648-09	Extra large	Tall	Tan	EGD8818	382-0016
Kneepiece (seperate)*					
73368-00	Extra small	-	Tan	-	-
73368-01	Small	-	Tan	-	-
73368-02	Medium	-	Tan	-	-
73368-03	Large	-	Tan	-	-
73368-04	Extra large	-	Tan	-	-
73368-05	Extra extra large	-	Tan	-	-
Thighpiece (including kneepiece)					
76649-05	Extra small	Short	Tan	EGD8728	382-0222
76649-00	Extra small	Regular	Tan	EGD8819	382-0172
76649-10	Extra small	Tall	Tan	EGD8733	382-0123
76649-06	Small	Short	Tan	EGD8729	382-0255
76649-01	Small	Regular	Tan	EGD8820	382-0206
76649-11	Small	Tall	Tan	EGD8734	382-0156
76649-07	Medium	Short	Tan	EGD8730	382-0248
76649-02	Medium	Regular	Tan	EGD8836	382-0198
76649-12	Medium	Tall	Tan	EGD8735	382-0149
76649-08	Large	Short	Tan	EGD8731	382-0230
76649-03	Large	Regular	Tan	EGD8822	382-0180
76649-13	Large	Tall	Tan	EGD8736	382-0131
76649-09	Extra large	Short	Tan	EGD8732	382-0214
76649-04	Extra large	Regular	Tan	EGD8727	382-0164
76649-14	Extra large	Tall	Tan	EGD8737	382-0115

Note: if your patient does not fit within the RTW sizing chart, STRONG is also available as a custom-fit variant. Further details available on request. *Not available on FP10 / GP10.

JOBST® FarrowWrap® Classic

Lower Limb, Ready-to-Wear Wrap Compression System

Compression: **30-40mmHg**

Moderate to severe lymphoedema / chronic oedema, chronic venous insufficiency and venous leg ulceration.

The original material that has been clinically proven over many years of use on thousands of patients worldwide. The rigid nature of the material offers performance beyond the compression range making it excellent for patients with deep skin folds and moderate to severe oedema.





Code	Size	Length	Colour	NHS SC Code	PIP Code
Footpiece					
76656-05	Extra small	Regular	Tan	EGD8860	360-6688
76656-00	Extra small	Long	Tan	EGD8773	360-6738
76656-06	Small	Regular	Tan	EGD8861	360-6696
76656-01	Small	Long	Tan	EGD8774	360-6746
76656-07	Medium	Regular	Tan	EGD8862	360-6704
76656-02	Medium	Long	Tan	EGD8857	360-6753
76656-08	Large	Regular	Tan	EGD8863	360-6712
76656-03	Large	Long	Tan	EGD8858	360-6761
76656-09	Extra large	Regular	Tan	EGD8864	360-6720
76656-04	Extra large	Long	Tan	EGD8859	360-6779
Legpiece					
76657-00	Extra small	Regular	Tan	EGD8865	360-6589
76657-05	Extra small	Tall	Tan	EGD8870	360-6639
76657-01	Small	Regular	Tan	EGD8866	360-6597
76657-06	Small	Tall	Tan	EGD8871	360-6647
76657-02	Medium	Regular	Tan	EGD8867	360-6605
76657-07	Medium	Tall	Tan	EGD8872	360-6654
76657-03	Large	Regular	Tan	EGD8868	360-6613
76657-08	Large	Tall	Tan	EGD8873	360-6662
76657-04	Extra large	Regular	Tan	EGD8869	360-6621
76657-09	Extra large	Tall	Tan	EGD8874	360-6670
Thighpiece (including kneepiece)					
76658-06	Extra small	Short	Tan	EGD8881	360-6431
76658-00	Extra small	Regular	Tan	EGD8875	360-6480
76658-12	Extra small	Tall	Tan	EGD8887	360-6530
76658-07	Small	Short	Tan	EGD8882	360-6449
76658-01	Small	Regular	Tan	EGD8876	360-6498
76658-13	Small	Tall	Tan	EGD8888	360-6548
76658-08	Medium	Short	Tan	EGD8883	360-6456
76658-02	Medium	Regular	Tan	EGD8877	360-6506
76658-14	Medium	Tall	Tan	EGD8889	360-6555
76658-09	Large	Short	Tan	EGD8884	360-6464
76658-03	Large	Regular	Tan	EGD8878	360-6514
76658-15	Large	Tall	Tan	EGD8890	360-6563
76658-10	Extra large	Short	Tan	EGD8885	360-6472
76658-04	Extra large	Regular	Tan	EGD8879	360-6522
76658-16	Extra large	Tall	Tan	EGD8891	360-6571

Note: if your patient does not fit within the RTW sizing chart, CLASSIC is also available as a custom-fit variant. Further details available on request.

JOBST® FarrowWrap®

Ready-to-Wear Liner*

A pair of liners is supplied with each JOBST® FarrowWrap® piece with the exception of ready-to-wear upper limb armpiece^. It provides a smooth base and skin protection prior to applying the outer wrap piece. In addition, daily washing of the outer wrap piece is not recommended due to the nature of the material. However, a fresh liner can be worn daily as the liner can be washed daily. Additional liners are available to order separately but they are not available on Drug Tariff.

- Beige - open toe knee high liner
- TG® Soft - soft and stretchy open toe liner

Liner	Code	NHS SC Code	Piece	Size
Beige closed toe	76666-21	EGD8826	Leg AD	Small (15-58cm)
Beige closed toe	76666-22	EGD8827	Leg AD	Large (52-70cm)
TG Soft Open toe	76669-00	EGD8828	Leg AD	Small (<40cm)
TG Soft Open toe	76669-02	EGD8829	Leg AD	Medium (40-70cm)
TG Soft Open toe	76669-04	EGD8830	Leg AD	Large (70-120cm)
TG Soft Open toe	76670-02	EGD8832	Thigh AG	Medium (40-70cm)
TG Soft Open toe	76670-04	EGD8833	Thigh AG	Large (70-120cm)

JOBST® FarrowHybrid®

Ready-to-Wear Liner With Compression in the Foot Only

For the treatment of mild to moderate oedema of the foot. JOBST® FarrowHybrid® provides even compression around the foot and ankle whilst remaining non-compressive above the ankle.

Code	Size	Colour	NHS SC Code	PIP Code
ADI (20-30mmHg)				
76664-00	Small (shoe size 2-5)	Taupe	EGD8893	408-1626
76664-01	Medium (shoe size 6-8)	Taupe	EGD8894	408-1634
76664-02	Large (shoe size 9-11)	Taupe	EGD8895	408-1642
76664-03	Extra large (shoe size 12-13)	Taupe	EGD8896	408-1659
76664-12	Small - wide (shoe size 2-5)	Taupe	EGD8905	408-1667
76664-13	Medium - wide (shoe size 6-8)	Taupe	EGD8823	408-1675
76664-14	Large - wide (shoe size 9-11)	Taupe	EGD8824	408-2384
76664-15	Extra large - wide (shoe size 12-13)	Taupe	EGD8825	408-1683
76664-16	Small - wide (shoe size 2-5)	Black	-	408-6765
76664-17	Medium - wide (shoe size 6-8)	Black	-	408-6757
76664-18	Large - wide (shoe size 9-11)	Black	-	408-6740
76664-19	Extra large - wide (shoe size 12-13)	Black	-	408-6732

Do not apply JOBST® FarrowHybrid® if calf circumference exceeds 60cm.
Consider a JOBST® FarrowWrap® footpiece in this circumstance.

^RTW upper limb has a liner built into the LITE armpiece
*Not available on FP10 / GP10

Ready-to-Wear Lower Limb Measuring Guide

For Ready-to-Wear

3 x Measurements required for legpiece
3 x Measurements required for thighpiece
2 x Measurements required for footpiece

For Custom-Fit


All measurements required for custom-fit.

Note:

For A-D and D-E, measure posteriorly and follow leg contour.

Note:

For E-G measure medially and follow contours.



	Extra Small	Small	Medium	Large	Extra Large
Footpiece					
Mid-foot	22 - 24cm	25 - 27cm	28 - 30cm	31 - 34cm	35 - 40cm
Regular	15 - 17cm	16 - 19cm	18 - 20cm	19 - 22cm	21 - 24cm
Tall	18 - 20cm	20 - 22cm	21 - 23cm	23 - 25cm	25 - 27cm
Legpiece					
Widest calf	36 - 43cm	42 - 50cm	48 - 58cm	53 - 63cm	58 - 68cm
Ankle	21 - 25cm	25 - 30cm	30 - 36cm	36 - 42cm	42 - 50cm
Regular	35 - 39cm	37 - 41cm	39 - 43cm	41 - 45cm	42 - 46cm
Tall	40 - 44cm	42 - 46cm	44 - 48cm	46 - 50cm	47 - 51cm
Kneepiece (seperate)					
Above knee	47 - 52cm	53 - 59cm	60 - 66cm	67 - 73cm	74 - 80cm
Widest calf	36 - 43cm	42 - 50cm	48 - 58cm	53 - 63cm	58 - 68cm
Thighpiece (including kneepiece)					
Groin	60 - 70cm	65 - 75cm	70 - 80cm	75 - 90cm	85 - 100cm
Above knee	47 - 52cm	53 - 59cm	60 - 66cm	67 - 73cm	74 - 80cm
Short	16 - 20cm	16 - 20cm	16 - 20cm	16 - 20cm	16 - 20cm
Regular	21 - 26cm	21 - 26cm	21 - 26cm	21 - 26cm	21 - 26cm
Tall	27 - 31cm	27 - 31cm	27 - 31cm	27 - 31cm	27 - 31cm
Liners*					
Beige closed toe	-	15 - 58cm		52 - 70cm	-
TG® Soft open toe	-	<40cm	40 - 70cm	70 - 120cm	-
FarrowHybrid® shoe size	-	2 - 5	6 - 8	9 - 11	12 - 13

*Supplied in pairs
Notes for JOBST® FarrowHybrid®: Compression in the foot. Standard foot = 20-27cm. Wide foot = 28-35cm.
JOBST® FarrowHybrid® maximum leg circumference 60cm

JOBST® FarrowWrap® 4000

Available on Drug Tariff

JOBST® FarrowWrap® 4000 is a four band, interlocking wrap compression system that provides 30-40mmHg, for the treatment of mild to moderate oedema, with or without a venous leg ulcer present.


JOBST® FarrowWrap® 4000 comes as a kit including 1 x JOBST® FarrowHybrid® 20-30mmHg compression sock and 1 x JOBST® FarrowWrap® 4000 outer wrap.

The JOBST® FarrowHybrid® provides compression from foot to ankle only, to address oedema in the foot.

Compression Level
CCL 3 = 30-40mmHg
✓


Sizes and Options		
Four Sizes	Regular	Tall
✓	✓	✓

Colours	
Tan	Black

Styles

Leg



3 simple measurements



Measuring Information

	Extra Small	Small	Medium	Large
JOBST® FarrowWrap® 4000 Legpiece				
cB Ankle	17 - 23cm	19 - 28cm	25 - 37cm	29 - 43cm
cC Calf	25 - 38cm	32 - 48cm	37 - 56cm	43 - 67cm
B-D Regular	30 - 33cm	30 - 33cm	30 - 33cm	30 - 33cm
B-D Tall	34 - 37cm	34 - 37cm	34 - 37cm	34 - 37cm

Ordering Information

Code	Size	Length	Colour	PIP Code
JOBST® FarrowWrap® 4000 Legpiece				
76661-00	Extra small	Regular	Tan	408-1543
76661-08	Extra small	Regular	Black	408-1584
76661-04	Extra small	Tall	Tan	408-1501
76661-12	Extra small	Tall	Black	408-1469
76661-01	Small	Regular	Tan	408-1550
76661-09	Small	Regular	Black	408-1592
76661-05	Small	Tall	Tan	408-1519
76661-13	Small	Tall	Black	408-1477
76661-02	Medium	Regular	Tan	408-1568
76661-10	Medium	Regular	Black	408-1600
76661-06	Medium	Tall	Tan	408-1527
76661-14	Medium	Tall	Black	408-1485
76661-03	Large	Regular	Tan	408-1576
76661-11	Large	Regular	Black	408-1618
76661-07	Large	Tall	Tan	408-1535
76661-15	Large	Tall	Black	408-1493

†A1 measurement only required for JOBST® FarrowHybrid®

JOBST® Farrow® Toe Cap

The JOBST® Farrow® toe cap is a ready-to-wear garment that provides compression for patients with mild to moderate oedema in the toes and foot. The comfortable fabric and ultra thin seams ensure high wearing comfort.

JOBST® Farrow® toe cap can be worn on either foot by turning the garment inside out.

This garment is a fantastic addition to the JOBST® FarrowWrap® range and can be worn together with the JOBST® FarrowWrap® footpiece.

Compression Levels	
CCL 1 = 15-20mmHg	CCL 2 = 20-30mmHg
✓	✓


Sizes and Options	
Four Sizes	Toes Can be Trimmed to The Correct Length
✓	✓

Colours	
Tan	Black

Measuring Information

Size	Extra Small	Small	Medium	Large
cA	20 - 22cm	22 - 24cm	24 - 26cm	26 - 28cm





Ordering Information

Colour	Compression Level	X Small	Small	Medium	Large
Tan	15-20mmHg	71027-00	71027-01	71027-02	71027-03
	20-30mmHg	72947-00	72947-01	72947-02	72947-03
Black	15-20mmHg	71027-04	71027-05	71027-06	71027-07
	20-30mmHg	72947-04	72947-05	72947-06	72947-07

JOBST® Farrow® Glove

The JOBST® Farrow® glove is a ready-to-wear garment that provides gentle compression for patients with mild to moderate oedema. The comfortable fabric and ultra thin seams ensure high wearing comfort.

The JOBST® Farrow® glove is available in a broad range of sizes, and can be worn on either hand by turning the garment inside out.

This garment is a fantastic addition to the JOBST® FarrowWrap® range and can be worn together with the JOBST® FarrowWrap® gauntlet.


Compression Levels	
CCL 1 = 15-20mmHg	CCL 2 = 20-30mmHg
✓	✓

Sizes and Options	
Eight Sizes	Fingers Can be Trimmed to The Correct Length
✓	✓

Colours	
Tan	Black

Measuring Information

Size	A	B
	Circumference (cm)	
1	16.6 - 18.0cm	15.6 - 16.7cm
2	16.6 - 18.0cm	16.7 - 18.1cm
3	18.0 - 19.5cm	15.9 - 17.0cm
4	18.0 - 19.5cm	17.0 - 18.5cm
5	18.0 - 19.5cm	18.5 - 20.4cm
6	19.5 - 21.5cm	17.2 - 18.7cm
7	19.5 - 21.5cm	18.7 - 20.6cm
8	21.5 - 23.1cm	18.3 - 20.1cm





Ordering Information

Colour	Compression Level	Size 1	Size 2	Size 3	Size 4	Size 5	Size 6	Size 7	Size 8
Tan	15-20mmHg	75463-00	75463-01	75463-02	75463-03	75463-04	75463-05	75463-06	75463-07
	20-30mmHg	75924-00	75924-01	75924-02	75924-03	75924-04	75924-05	75924-06	75924-07
Black	15-20mmHg	75463-10	75463-11	75463-12	75463-13	75463-14	75463-15	75463-16	75463-17
	20-30mmHg	75924-10	75924-11	75924-12	75924-13	75924-14	75924-15	75924-16	75924-17

JOBST® Elvarex® Custom-Fit

How to create a custom-fit garment for your patient

Each patient you measure for a custom-fit garment will have different requirements and preferences and these can be accommodated quickly and easily with the JOBST® Elvarex® custom-fit range. The following information will guide you through the process of 'building' a patient's custom-fit garment. Along with the Essity compression therapy price list, it will enable you to calculate the price of each individual garment you prescribe and, if required, select the appropriate Drug Tariff codes that are essential for a community prescription.

Select product

Firstly, select the product you require - JOBST® Elvarex®, JOBST® Elvarex® Plus, JOBST® Elvarex® Soft or JOBST® Elvarex® Soft Seamless.

Select basic style and compression class

Then, select the basic garment style required from the choices available under each product, for example 'Knee High'. For each garment style, you will find a list of codes relating to the compression class. This list of codes relates to the basic garment in the standard colour (eg: knee high with slant open toe in beige).

Select options

For each basic garment style, there is also a list of options. The options and codes relate to additional garment features that can be added to the basic garment you have selected, in order to provide the individually specified finished garment. The options include T-heel, slipform, zipper, silicone band, closed toe, colour, etc.

Example

Your patient requires a pair of JOBST® Elvarex® custom-fit tights - tights are made up of 2 x thigh high stockings and 1 x body bandage. The requirement is for CCL 1 on the right leg and CCL 2 on the left. Also, the patient prefers closed toes and needs a T-heel on the left leg. The garment would be created as follows:

Description	Code
Thigh high stocking CCL 1	L1-01-06
Thigh high stocking CCL 2	L2-02-06
Body bandage CCL 1	L1-01-08
Closed toes	L-A001 x 2
T-heel	L-A010

Chap and tights styles

For chap and tights styles, you will see an extra column headed "Add". These codes must be included in the order to achieve the final 'chap style' or 'tights' garment. For example, to order a chap style garment, you must select the basic garment in the required compression class (thigh high stocking A-G) and add a waist attachment with garter belt. You can also choose to add other options to customise the final garment design.

Drug Tariff codes / prescription

We have a wide range of styles and options available in the JOBST® Elvarex® range. The majority of these are available on Drug Tariff and so are reimbursable against a prescription. It is important that you refer to the latest copy of the Essity compression therapy Drug Tariff price list to ensure that the styles and options you require are all available on Drug Tariff. The Drug Tariff codes need to be included on your prescription.

Schema number

Each JOBST® Elvarex® garment order is provided with a schema number. This is a unique identification number for each garment manufactured and the number can be found on the order sheet inside the packaging. By requesting a re-order quoting this schema number, you will receive an exact copy of the original garment. The pharmacist still needs the prescription form quoting the 'schema' number and the order codes.

Further assistance

Please do not hesitate to contact our Customer Service team if you require any further information on 0845 122 3600.

JOBST® Elvarex® Custom-Fit

Available on Drug Tariff

Custom-fit, flat-knit, RAL compression garments

For over 70 years, JOBST® Elvarex® has been a recognised standard of therapy providing reliable flat-knitted, graduated compression garments for the treatment of lymphoedema and chronic oedema. Through our high degree of innovation and manufacturing expertise we are now the UK's leading flat-knit compression garment provider. JOBST® Elvarex® is a time-proven classic and seen as the absolute specialist amongst flat-knitted compression garments.

Unique stitch pattern

JOBST® Elvarex® is knitted on straight rows of needles and only becomes a garment when sewn together. The number of stitches can be varied to achieve virtually any desired shape or accomomodate any level of limb shape distortion. A high-working pressure to resting pressure ratio guarantees effective mobilisation of oedema and improves lymph drainage. The open, wide-mesh structure of the unique stitch pattern provides high air permeability and breathability. The coarse-knitted loops ensure the garment lies flat over skin folds, will not tourniquet ensuring a reduced risk of trauma and skin damage. The complex yarn formula is very robust and provides high resistance, especially important for more stubborn oedema and rebound oedema.

Compression Classes - Upper Limb		
RAL CCL 1 = 15-21mmHg	RAL CCL 2 = 23-32mmHg	RAL CCL 2 Forte = 23-32mmHg
✓	✓	✓

Compression Classes - Lower Limb					
RAL CCL 1 18-21mmHg	RAL CCL 2 23-32mmHg	RAL CCL 3 34-46mmHg	RAL CCL 3 Forte 34-46mmHg	RAL CCL 4 49-70mmHg	RAL CCL 4 Super 60-90mmHg
✓	✓	✓	✓	✓	✓

Sizes	
Custom-Fit	Available as custom-fit garments for a precise anatomical fit for each individual patient. Majority of patients can be catered for regardless of limb size or shape.
✓	

Options
There are over 13,000 variations available with the majority being available on Drug Tariff, including the most popular options, dotted silicone band, slipform, T-heel and zipper. For full details, please see the 'Options' section alongside the basic garment style required.


Colours - Upper Limb					
Beige	Black	Cranberry	Macchiato	Ruby red	Dark blue

Colours - Lower Limb					
Beige	Black	Cranberry	Dark brown	Dark blue	Stone
Grey	Henna	Denim	Graphite	Aubergine	

Styles
Available in 14 basic styles, with all conceivable variants, to accommodate a broad therapeutic range for lymphoedema management of the lower and upper extremities.

JOBST® Elvarex® Custom-Fit


Knee High



Basic style - slant open toe in beige		
CCL 1 - 18-21mmHg	L1-01-04	1
CCL 2 - 23-32mmHg	L2-02-04	1
CCL 3 - 34-46mmHg	L3-03-04	1
CCL 3 Forte - 34-46mmHg	L3-04-04	1
CCL 4 - 49-70mmHg	L4-05-04	1
CCL 4 Super - 60-90mmHg	L5-06-04	1

Options		
Closed toe	L-A001	1
2 Ankle pads (profile)	L-A002	1
Zipper	L-A003	1
Silicone band	L-A004	1
Non-standard colour	L-A008	1
T-heel	L-A010	1
SoftFit	L-A016	1
Silk pocket		


Knee High Footless



Basic style - in beige		
CCL 1 - 18-21mmHg	L1-01-24	1
CCL 2 - 23-32mmHg	L2-02-24	1
CCL 3 - 34-46mmHg	L3-03-24	1
CCL 3 Forte - 34-46mmHg	L3-04-24	1
CCL 4 - 49-70mmHg	L4-05-24	1
CCL 4 Super - 60-90mmHg	L5-06-24	1

Options *Not available on Drug Tariff		
Zipper	L-A003	1
Silicone band	L-A004	1
Non-standard colour	L-A008	1
SoftFit	L-A016	1
Silk pocket		


Mid Thigh*



Basic style - slant open toe in beige		
CCL 1 - 18-21mmHg	L1-01-20	1
CCL 2 - 23-32mmHg	L2-02-20	1
CCL 3 - 34-46mmHg	L3-03-20	1
CCL 3 Forte - 34-46mmHg	L3-04-20	1
CCL 4 - 49-70mmHg	L4-05-20	1
CCL 4 Super - 60-90mmHg	L5-06-20	1

Options *Not available on Drug Tariff		
Closed toe	L-A001	1
2 Ankle pads (profile)	L-A002	1
Zipper	L-A003	1
Silicone band	L-A004	1
Non-standard colour	L-A008	1
T-Heel	L-A010	1
Knee functional zone	L-A015	1
No volant		
Silk pocket		


Thigh High



Basic style - slant open toe in beige		
CCL 1 - 18-21mmHg	L1-01-06	1
CCL 2 - 23-32mmHg	L2-02-06	1
CCL 3 - 34-46mmHg	L3-03-06	1
CCL 3 Forte - 34-46mmHg	L3-04-06	1
CCL 4 - 49-70mmHg	L4-05-06	1
CCL 4 Super - 60-90mmHg	L5-06-06	1

Options *Not available on Drug Tariff		
Closed toe	L-A001	1
2 Ankle pads (profile)	L-A002	1
Zipper	L-A003	1
Silicone band	L-A004	1
Non-standard colour	L-A008	1
Slipform	L-A009	1
T-heel	L-A010	1
Top functional zone*	L-A014	1
Knee functional zone	L-A015	1
Straight top		
No volant		
Silk pocket		

Thigh High Footless



Basic style - in beige		
CCL 1 - 18-21mmHg	L1-01-16	1
CCL 2 - 23-32mmHg	L2-02-16	1
CCL 3 - 34-46mmHg	L3-03-16	1
CCL 3 Forte - 34-46mmHg	L4-04-16	1
CCL 4 - 49-70mmHg	L5-05-16	1
CCL 4 Super - 60-90mmHg	L6-06-16	1

A pair of footless tights consists of 2 x thigh footless and 1 x body bandage. The body bandage codes can be found in the tights section overleaf.

Options *Not available on Drug Tariff		
Zipper	L-A003	1
Silicone band	L-A004	1
Non-standard colour	L-A008	1
Slipform	L-A009	1
Top functional zone*	L-A014	1
Knee functional zone	L-A015	1
No volant		
Silk pocket		

JOBST® Elvarex® Custom-Fit

Tights




Basic style - slant open toe in beige		
CCL 1 - 18-21mmHg	L1-01-06	1
CCL 2 - 23-32mmHg	L2-02-06	1
CCL 3 - 34-46mmHg	L3-03-06	1
CCL 3 Forte - 34-46mmHg	L3-04-06	1
CCL 4 - 49-70mmHg	L4-05-06	1
CCL 4 Super - 60-90mmHg	L5-06-06	1
Add body bandage		
CCL 1 - 18-21mmHg	L1-01-08	
CCL 2 - 23-32mmHg	L2-02-08	
CCL 3 - 34-46mmHg	L3-03-08	
CCL 3 Forte - 34-46mmHg	L3-04-08	
CCL 4 - 49-70mmHg	L4-05-08	
CCL 4 Super - 60-90mmHg	L5-06-08	

A pair of tights consists of 2 legs and a body bandage sewn together. One-legged tights can also be constructed by adding a leg extension to a body bandage (one-legged garment can offer stability to the unaffected limb). Compression classes for each leg and the body bandage can be mixed on the same garment.

Options *Not available on Drug Tariff		
Closed toe	L-A001	1
2 Ankle pads (profile)	L-A002	1
Zipper	L-A003	1
Silicone band	L-A004	1
Leg extension	L-A007	1
Non-standard colour	L-A008	1
Slipform	L-A009	1
T-heel	L-A010	1
Fly for men	L-A012	1
Adjustable waistband	L-A013	1
Top functional zone*	L-A014	1
Knee functional zone	L-A015	1
Reinforced gusset		
Straight top		
No volant		
Silk pocket		
Open pubis		
Crotch for men		

Capri Leg*



Basic style - in beige		
CCL 1 - 18-21mmHg	L1-01-14	1
CCL 2 - 23-32mmHg	L2-02-14	1
CCL 3 - 34-46mmHg	L3-03-14	1
CCL 3 Forte - 34-46mmHg	L3-04-14	1
CCL 4 - 49-70mmHg	L4-05-14	1
CCL 4 Super - 60-90mmHg	L5-06-14	1

Options *Not available on Drug Tariff		
Zipper	L-A003	1
Non-standard colour	L-A008	1
Slipform	L-A009	1
Top functional zone*	L-A014	1
Knee functional zone	L-A015	1
Silk pocket		

A pair of Capri tights for CCL3Forte to 4Super can be ordered by combining 2 Capri legs and a body bandage sewn together. The body bandage codes can be found in the tights section above. A Capri chap style garment consists of a Capri leg (or 2 x Capri legs for a pair) with a waist attachment with garter belt for all compression classes. The waist attachment with garter belt code can be found in the Chap style garment section below.


Capri Tights (CCL 1-3 only)



Basic style - in beige		
CCL 1 - 18-21mmHg	L1-01-28	1
CCL 2 - 23-32mmHg	L2-02-28	1
CCL 3 - 34-46mmHg	L3-03-28	1

Options		
Non-standard colour	L-A008	1
Slipform	L-A009	1
Adjustable waistband	L-A013	1

Chap Style



Basic style - slant open toe in beige		
CCL 1 - 18-21mmHg	L1-01-06	1
CCL 2 - 23-32mmHg	L2-02-06	1
CCL 3 - 34-46mmHg	L3-03-06	1
CCL 3 Forte - 34-46mmHg	L3-04-06	1
CCL 4 - 49-70mmHg	L4-05-06	1
CCL 4 Super - 60-90mmHg	L5-06-06	1
Add waist attachment with garter belt		
	L-A005	1

Chap style garments consist of one thigh high stocking and a waist attachment with garter belt or two thigh high stockings for a pair of chap style garments

Options		
Closed toe	L-A001	1
2 Ankle pads (profile)	L-A002	1
Zipper	L-A003	1
Silicone band	L-A004	1
Non-standard colour	L-A008	1
Slipform	L-A009	1
T-heel	L-A010	1
Adjustable waistband	L-A013	1
Knee functional zone	L-A015	1
Straight top		
No volant		
Silk pocket		
Full compression in waist attachment		

Foot Cap with Toes




Basic style - with toes in beige†		
CCL 1 - 18-21mmHg	L1-01-02	1
CCL 2 - 23-32mmHg	L2-02-02	1
CCL 3 - 34-46mmHg	L3-03-02	1

Options		
Non-standard colour	L-A008	1

†1.5cm slanted little toe opening as standard

JOBST® Elvarex® Custom-Fit

Bermuda



Basic style - body bandage and two leg extensions in beige		
CCL 1 - 18-21mmHg	L1-01-26	1
CCL 2 - 23-32mmHg	L2-02-26	1
CCL 3 - 34-46mmHg	L3-03-26	1
CCL 3 Forte - 34-46mmHg	L3-04-26	1
CCL 4 - 49-70mmHg	L4-05-26	1
CCL 4 Super - 60-90mmHg	L5-06-26	1

Options *Not available on Drug Tariff		
Zipper	L-A003	1
Non-standard colour	L-A008	1
Slipform	L-A009	1
Fly for men	L-A012	1
Adjustable waistband	L-A013	1
Top functional zone*	L-A014	1
Reinforced gusset		
Open pubis		
Crotch for men		
Silk pocket		

Glove / Gauntlet to Wrist



Basic style - no fingers / thumb in beige		
CCL 1 - 15-21mmHg	U1-01-01	1
CCL 2 - 23-32mmHg	U2-02-01	1
CCL 2 Forte - 23-32mmHg	U3-03-01	1

Options		
Finger - CCL 1 - 15-21mmHg	U1-A009	1
Finger - CCL 2 - 23-32mmHg	U2-A011	1
Finger - CCL 2 Forte - 23-32mmHg	U3-A013	1
Zipper	U-A001	1
Non-standard colour	U-A007	1

Glove / Gauntlet to Elbow



Basic style - no fingers / thumb in beige		
CCL 1 - 15-21mmHg	U1-01-03	1
CCL 2 - 23-32mmHg	U2-02-03	1
CCL 2 Forte - 23-32mmHg	U3-03-03	1

Options		
Finger - CCL 1 - 15-21mmHg	U1-A009	1
Finger - CCL 2 - 23-32mmHg	U2-A011	1
Finger - CCL 2 Forte - 23-32mmHg	U3-A013	1
Zipper	U-A001	1
Silicone band	U-A002	1
Non-standard colour	U-A007	1
Bias top		

Armsleeve to Axilla with Glove / Gauntlet



Basic style - no fingers / thumb in beige		
CCL 1 - 15-21mmHg	U1-01-05	1
CCL 2 - 23-32mmHg	U2-02-05	1
CCL 2 Forte - 23-32mmHg	U3-03-05	1

Options		
Finger - CCL 1 - 15-21mmHg	U1-A009	1
Finger - CCL 2 - 23-32mmHg	U2-A011	1
Finger - CCL 2 Forte - 23-32mmHg	U3-A013	1
Zipper	U-A001	1
Silicone band	U-A002	1
Shoulder cap	U-A003	1
Non-standard colour	U-A007	1
SoftFit	U-A016	1
Full compression shoulder cap		
Straight top		

Armsleeve from Wrist to Axilla



Basic style - in beige		
CCL 1 - 15-21mmHg	U1-01-07	1
CCL 2 - 23-32mmHg	U2-02-07	1
CCL 2 Forte - 23-32mmHg	U3-03-07	1

Options		
Zipper	U-A001	1
Silicone band	U-A002	1
Shoulder cap	U-A003	1
Non-standard colour	U-A007	1
Elbow functional zone	U-A008	1
SoftFit	U-A016	1
Full compression shoulder cap		
Straight top		

JOBST® Elvarex® Plus

Custom-Fit, Flat-Knit, RAL Compression Garments

Available on Drug Tariff

JOBST® Elvarex® Plus is part of a robust family of products which through innovative product development, continues the JOBST® Elvarex® success story in the management of lymphoedema and chronic oedema.

A new generation in flat-knit garments

JOBST® Elvarex® Plus has succeeded in adapting the inner strength of the JOBST® Elvarex® material to the skin's sensory perception and the complex movements of the hands and feet. With the unique 3D knitting technology, the custom-fit, flat-knit, seamless compression garments reproduce each body shape accurately and fit like a second skin.


Compression Classes - Upper Limb	
RAL CCL 1 = 15-21mmHg	RAL CCL 2 = 23-32mmHg
✓	✓

Compression Classes - Lower Limb	
RAL CCL 1 = 18-21mmHg	RAL CCL 2 = 23-32mmHg
✓	✓

Sizes and Options			
Custom-Fit	Foot Cap With Four or Five Toes	Glove / Gaunlet to Elbow With / Without Bias Cut	Individual Right / Left Garments
✓	✓	✓	✓

Colours - Upper Limb					
Beige	Black	Cranberry	Macchiato	Ruby red	Dark blue

Colours - Lower Limb					
Beige	Black	Cranberry	Dark brown	Dark blue	Grey

Styles			
			
AA Foot Cap	AC1 Glove With Thumb and Fingers to Wrist	AC1 Gauntlet With Thumb to Wrist	AE Glove or Gauntlet With Thumb to Elbow

JOBST® Elvarex® Plus

Foot Cap with Toes



Compression Class	Code	Pieces
Basic style – foot cap with toes in beige†		
CCL 1 – 18-21mmHg	L1-12-02	1
CCL 2 – 23-32mmHg	L2-13-02	1

Options	Code	Pieces
Options		
Non-standard colour	L-A008P	1

†1.5cm slanted little toe opening as standard


Glove / Gauntlet to Wrist



Compression Class	Code	Pieces
Basic style – no fingers / thumb in beige		
CCL 1 – 15-21mmHg	U1-12-01	1
CCL 2 – 23-32mmHg	U2-13-01	1

Options	Code	Pieces
Options		
Finger CCL 1 – 15-21mmHg	U1-A016	1
Finger CCL 2 – 23-32mmHg	U2-A017	1
Non-standard colour	U-A007P	1

Glove / Gauntlet to Wrist



Compression Class	Code	Pieces
Basic style – no fingers / thumb in beige		
CCL 1 – 15-21mmHg	U1-12-03	1
CCL 2 – 23-32mmHg	U2-13-03	1

Options	Code	Pieces
Options		
Finger CCL 1 – 15-21mmHg	U1-A016	1
Finger CCL 2 – 23-32mmHg	U2-A017	1
Non-standard colour	U-A007P	1
Bias top		



Seamless

3D knitting

JOBST® Elvarex® material

JOBST® Elvarex® Soft

Available on Drug Tariff

Custom-Fit, Flat-Knit, RAL Compression Garments

Upper and lower limb custom-fit, flat-knit garments for the treatment of mild to moderate lymphoedema and specifically designed for those patients with dry, sensitive or fragile skin. JOBST® Elvarex® Soft provides choice and allows patients to live with oedema the comfortable way.

Soft and skin friendly

The unique knitted fabric contains a soft textured yarn which creates an extremely smooth fabric surface and makes the inside of the seam particularly soft. This combination of skin-protecting softness and excellent breathability makes the garments particularly well suited for patients with dry, sensitive or fragile skin. There are a number of specifically designed material enhancements, such as the knitted edge without an edging seam and the comfort dotted silicone band, which all combine to ensure that the garment provides the optimum fit and comfort.

Compression Classes - Upper Limb	
RAL CCL 1 = 15-21mmHg	RAL CCL 2 = 23-32mmHg
✓	✓

Compression Classes - Lower Limb		
RAL CCL 1 = 18-21mmHg	RAL CCL 2 = 23-32mmHg	RAL CCL 3 = 34-46mmHg
✓	✓	✓

Sizes and Options	
Custom-Fit	Available as custom-fit garments for a precise anatomical fit for each individual patient. Majority of patients can be catered for regardless of limb size or shape.
✓	


Colours - Upper Limb				
Beige	Black	Cranberry	Dark brown	Macchiato
Grey	Ruby red	Pine green	Sunflower yellow	

Colours - Lower Limb				
Beige	Black	Cranberry	Dark brown	Dark blue
Grey	Ruby red	Pine green	Sunflower yellow	

Styles
Available in eight basic styles, with many options to accommodate a wide range of clinical and patient needs.

JOBST® Elvarex® Soft


Knee High



Basic style - slant open toe in beige		
CCL 1 – 18-21mmHg	L1-10-04	1
CCL 2 – 23-32mmHg	L2-08-04	1
CCL 3 – 34-46mmHg	L3-09-04	1

Options	Code	Pieces
Closed toe	L-A001S	1
Silicone band	L-A004S	1
Non-standard colour	L-A008S	1
T-heel	L-A010S	1
SoftFit	L-A016S	1

Thigh High



Basic style - slant open toe in beige		
CCL 1 – 18-21mmHg	L1-10-06	1
CCL 2 – 23-32mmHg	L2-08-06	1
CCL 3 – 34-46mmHg	L3-09-06	1

Options	Code	Pieces
Closed toe	L-A001S	1
Silicone band	L-A004S	1
Slipform†	L-A011S	1
Non-standard colour	L-A008S	1
T-heel	L-A010S	1

†Please note thigh high, tights, Capri tights and Bermuda style garments are manufactured with slipform.

Tights




Basic style - slant open toe in beige		
CCL 1 – 18-21mmHg	L1-10-06	1
CCL 2 – 23-32mmHg	L2-08-06	1
CCL 3 – 34-46mmHg	L3-09-06	1
Add body bandage		
CCL 1 – 18-21mmHg	L1-10-08	1
CCL 2 – 23-32mmHg	L2-08-08	1
CCL 3 – 34-46mmHg	L3-08-08	1

Options	Code	Pieces
Closed toe	L-A001S	1
Slipform†	L-A011S	1
Non-standard colour	L-A008S	1
T-heel	L-A010S	1
Adjustable waistband	L-A013S	1

†Please note thigh high, tights, Capri tights and Bermuda style garments are manufactured with slipform.

Thigh High Footless (Footless Tights)*



Basic style - footless thigh high in beige		
CCL 1 - 18-21mmHg	L1-10-16	1
CCL 2 - 23-32mmHg	L2-08-16	1
CCL 3 - 34-46mmHg	L3-09-16	1
Add body bandage		
CCL 1 - 18-21mmHg	L1-10-08	1
CCL 2 - 23-32mmHg	L2-08-08	1
CCL 3 - 34-46mmHg	L3-08-08	1

Options	Code	Pieces
Non-standard colour	L-A008S	1
Slipform†	L-A011S	1
Adjustable waistband	L-A013S	1

*Footless tights consist of 2 x thigh high footless and 1 x body bandage sewn together.
†Please note thigh high, tights and Capri tights are manufactured with slipform. The Drug Tariff code must be included on your prescription.

Capri Leg (Capri Tights)




Basic style - Capri leg in beige		
CCL 1 - 18-21mmHg	L1-10-14	1
CCL 2 - 23-32mmHg	L2-08-14	1
CCL 3 - 34-46mmHg	L3-09-14	1

Options	Code	Pieces
Non-standard colour	L-A008S	1
Slipform†	L-A011S	1
Adjustable waistband	L-A013S	1

†Please note thigh high, tights, Capri tights and Bermuda style garments are manufactured with slipform.

JOBST® Elvarex® Soft

Chap Style



Basic style - foot cap with toes in beige†			Options		
CCL 1 - 18-21mmHg	L1-10-06	1	Closed toe	L-A001S	1
CCL 2 - 23-32mmHg	L2-08-06	1	Non-standard colour	L-A008S	1
CCL 3 - 34-46mmHg	L3-09-06	1	T-heel	L-A010S	1
			Full compression in waist band		
Add waist attachment with garter belt					
	L-A005S	1			

Bermuda



Basic style - body bandage and two leg extensions in beige			Options		
CCL 1 - 15-21mmHg	L1-10-26	1	Non-standard colour	U-A008S	1
CCL 2 - 23-32mmHg	L2-08-26	1	Slipform†	L-A011S	1
CCL 3 - 34-46mmHg	L3-09-26	1			

†Please note thigh high, tights, Capri tights and Bermuda style garments are manufactured with slipform.

Armsleeve from Wrist to Axilla



Basic style - with bias cut in beige			Options		
CCL 1 - 15-21mmHg	U1-10-07	1	Silicone band	U-A002S	1
CCL 2 - 23-32mmHg	U2-08-07	1	Non-standard colour	U-A007S	1
			SoftFit	U-A018S	1



JOBST® Elvarex® Soft Seamless

Available on Drug Tariff

Custom-Fit, Flat-Knit, RAL Compression Garments

A technological innovation from Essity. The world's first seamless flat-knit compression glove and foot cap.

Seamless flat-knit design

Constructed entirely without a seam, these unique custom-fit, flat-knit garments provide unrivalled patient comfort during wear. The perfect anatomical fit guarantees even, uninterrupted and consistent compression for optimum oedema management, even in critical areas, such as the fingers and the ball of the foot. The real advantage of this innovation is that the garments are especially comfortable to wear.


Compression Classes - Upper Limb	
RAL CCL 1 = 15-21mmHg	RAL CCL 2 = 23-32mmHg
✓	✓

Compression Classes - Lower Limb	
RAL CCL 1 = 18-21mmHg	RAL CCL 2 = 23-32mmHg
✓	✓

Sizes and Options	
Custom-Fit	Individual Right / Left Garments
✓	✓

Colours - Upper Limb						
Beige	Black	Cranberry	Dark brown	Macchiato	Grey	Ruby red

Colours - Lower Limb						
Beige	Black	Cranberry	Dark brown	Dark blue	Grey	Ruby red

Styles			
			
AA Foot Cap	AC1 Glove With Thumb and Fingers to Wrist	AC1 Gauntlet With Thumb to Wrist	AE Glove or Gauntlet With Thumb to Elbow

JOBST® Elvarex® Soft Seamless

	Compression Class	Code	Pieces	Options	Code	Pieces
Foot Cap with Toes						
	**Basic style - foot cap with toes in beige†**			**Options**		
	CCL 1 - 18-21mmHg	L1-11-02	1	Non-standard colour	L-A008S	1
	CCL 2 - 23-32mmHg	L2-07-02	1			
				†1.5cm slanted little toe opening as standard		
Glove / Gauntlet to Wrist						
	**Basic style - no fingers / thumb in beige**			**Options**		
	CCL 1 - 15-21mmHg	U1-11-01	1	Finger CCL 1 - 15-21mmHg	U1-A015	1
	CCL 2 - 23-32mmHg	AC1-WF-WT	1	Finger CCL 2 - 23-32mmHg	U2-A014	1
				Non-standard colour	U-A007S	1
Glove / Gauntlet to Wrist						
	**Basic style - no fingers / thumb in beige**			**Options**		
	CCL 1 - 15-21mmHg	U1-11-03	1	Finger CCL 1 - 15-21mmHg	U1-A015	1
	CCL 2 - 23-32mmHg	U2-07-03	1	Finger CCL 2 - 23-32mmHg	U2-A014	1
				Non-standard colour	U-A007S	1



JOBST® Custom-Fit

Styles and Options - Upper Limb

Popular styles



AC1
Glove



AE
Glove to elbow with bias top



CG1
Armsleeve with bias top



AG1
Combined armsleeve and glove / gaunlet



CH
Armsleeve with shoulder cap (and bra loop)

Armsleeve options



Silicone dotted band, on top



Silicone dotted band, inside



3/4 silicone dotted band, inside



Shoulder cap with bra loop



Shoulder cap with body strap



Elbow functional zone



Armsleeve with silk pocket at bend of elbow



Bias top

Glove / gauntlet options



Gaunlet without fingers, with thumb



Glove with fingers, with thumb



Zipper



Silk pocket (to allow padding to be inserted)

JOBST® Custom-Fit

Styles and Options - Lower Limb

Popular styles



AD
Knee high



AG
Thigh high
(with slip form)



AG-T
Chap style garment (thigh high with waist attachment). Also available as a pair



AG-HT
One-legged tights. Thigh high with body bandage and leg extension



AT
Tights (2 x thigh high with 1 x body bandage)



Bermuda style



Capri tights

Top band options



Straight cut



Slipform (French cut / outside higher)



Silicone dotted band; inside



Silicone dotted band; on top



3/4 silicone dotted band; inside



2.5cm silicone dotted band, on top

JOBST® Custom-Fit

Styles and Options - Lower Limb

Silicone dotted band examples



5cm silicone dotted pieces sewn at G into the area of the slipform



Silicone dotted strips sewn at G into the area of the slipform (silicone strips are available in lengths of 10, 12, 15, 20 and 25cm)



10cm silicone vertical strips stitched at the back on top of the seam and at the front, from measuring point G downwards

Body bandage options



Knitted top of the body bandage. Full compression in body bandage



Adjustable waistband



Crotch for men



Fly for men



Top functional zone

SlipForm examples



Body bandage



Thigh high with silicone dotted band, inside



Thigh high with silicone dotted band, on top

JOBST® Custom-Fit

Styles and Options - Lower Limb

Foot options



Straight foot



Slant foot



Closed toe



Open toe



T-Heel



Ankle pad (profile)



Foot cap with toes



Silk pocket at ankle flex (Y measuring point)

Other lower limb options



Zipper



Silk pocket at back of knee



Leg extension



Waist attachment with garter belt for Chap style garments



Knee functional zone

JOBST® Confidence

Available on Drug Tariff

Custom-Fit, Flat-Knit, RAL Compression Garments

Introducing a new era of comfort with the next-generation of flat-knit technology for mild to moderate oedema, lipoedema and venous disease.


Inspired by patient needs and developed with next-generation technology

JOBST® Confidence helps patients feel unrestricted in body and mind, with a one-of-a-kind fit garment that is as unique as they are. The new generation of flat-knit technology used in JOBST® Confidence allows for a precise fit to the real anatomical body shapes of each of your individual patients. This personally refined fit not only provides comfort, but also supports your patients with a greater degree of movement


Compression Classes - Lower Limb		
RAL CCL 1 = 18-21mmHg	RAL CCL 2 = 23-32mmHg	RAL CCL 3 = 34-46mmHg
✓	✓	✓

Compression Classes - Upper Limb	
RAL CCL 1 = 15-21mmHg	RAL CCL 2 = 23-32mmHg
✓	✓

Sizes and Options						
Knee functional zone	Ankle functional zone	Closed toe	Slant toe	Halux Valgus	Elbow Functional Zone	Extra Angled Elbow
✓	✓	✓	✓	✓	✓	✓





Deco line




Initials

Colours					
Beige	Black	Jeans heather	Red heather	Anthracite heather	Caramel

Styles		
		
Knee high	Thigh high	Armsleeve C-G

JOBST® Confidence


Knee High



Basic style	Code	Pieces
CCL 1 - 18-21mmHg	L1-16-04	1
CCL 2 - 23-32mmHg	L2-17-04	1
CCL 3 - 34-46mmHg	L3-18-04	1

Options	Code	Pieces
Closed toe	L-A001C	1
SoftFit	L-A016C	1
Slipform	L-A011C	1
Ankle Functional zone	L-A010C	1

Thigh High



Basic style	Code	Pieces
CCL 1 - 18-21mmHg	L1-16-06	1
CCL 2 - 23-32mmHg	L2-17-06	1
CCL 3 - 34-46mmHg	L3-18-06	1

Options	Code	Pieces
Closed toe	L-A001C	1
Slipform	L-A011C	1
Ankle Functional zone	L-A010C	1
Knee functional zone	L-A015C	1

Armsleeve



Basic style	Code	Pieces
CCL 1 - 15-21mmHg	U1-16-07	1
CCL 2 - 23-32mmHg	U2-17-07	1

Options	Code	Pieces
Elbow Functional Zone	U-A008C	1
Silicone Band	U-A002C	1
SoftFit	U-A018C	1



JOBST® Relax

Slimline, Lightweight and Flexible

JOBST® Relax is designed to complement lymphoedema management during the day and offers an alternative to night-time compression solutions, such as compression bandages. JOBST® Relax is a custom-made, compression garment that provides the oedema control patients need, both at night-time and during periods of relaxation, combined with the look and comfort that patients want. It is intended for mild to severe lymphoedema patients in the maintenance phase, who are looking for a garment that will better suit their lifestyle.

JOBST® Relax can be measured in the same way as JOBST® Elvarex®.

Compression Levels	
Upper Limb	
CCL 1 = 15-20mmHg	
✓	
Lower Limb	
CCL 1 = 15-20mmHg	CCL 2 = 20-30mmHg
✓	✓

Sizes and Options / Materials
Zipper
✓

Colours	
Beige	Rose

Styles			
			
Armsleeve with gauntlet	Armsleeve	Knee High	Thigh High

More than just a good night's sleep

Slimline, lightweight and flexible design makes JOBST® Relax ideal for relaxing before bedtime



Anatomically shaped for a comfortable fit



Freedom of movement to accommodate the patient's lifestyle



Coolmax® and Lycra® fibres provide an improved fit and drying times



Easy to fold and fits comfortably into an overnight bag



Benefits of compression at night

Night-time compression can provide comfort and a reassuring sense of being in control. It might offer benefits which can even be felt the next day. Patients have reported benefits of reduced swelling, improved pain management and sleep when night-time compression was used, compared to an increase in swelling documented when night-time compression was not used[1].

Ordering information

Style / Option	Code	Style / Option	Code
JOBST® Relax			
Knee High Basic style - open toe			
CCL 1 - 15-20mmHg	L1-14-04	Zipper	L-A003R
CCL 2 - 20-30mmHg	L2-15-04		
Thigh High Basic style - open toe			
CCL 1 - 15-20mmHg	L1-14-06	Zipper	L-A003R
CCL 2 - 20-30mmHg	L2-15-06		
Armsleeve to Axilla with Gauntlet Basic style with bias top and thumb			
CCL 1 - 15-20mmHg	U1-14-05	Zipper	U-A001R
Armsleeve from Wrist to Axilla Basic style with bias top			
CCL 1 - 15-20mmHg	U1-14-07	Zipper	U-A001R

1 Hibbard, J.H. and Greene, J. (2013) 'What the evidence show about patient activation: Better health outcomes and care experiences; fewer data on costs', Health Affairs, 32 (2), pp. 207 – 214.

JOBST® Surgical Vest

Post-Operative Compression Garment



The JOBST® Surgical Vest is designed in luxuriously soft fabric, to enhance patient comfort and satisfaction. Patient comfort is of paramount importance, particularly during post-operative recovery.

Firm and supportive

The JOBST® Surgical Vest has a firm and supportive construction that supports breast and chest tissue after surgery, helping to manage post-operative oedema and minimising haematoma and seroma formation. In addition, the durable, nylon-zipper front-fastening enables easy donning and removal. The Velcro® brand hook-and-loop fastener system allows for easy adjustments of the garment.

Sizes			
Five Sizes	Both Cups	Without Cups	Left or Right
✓	✓	✓	✓

Colours
White

Ordering information

		Code	Code	Code
With Cups				
Size	**Chest and cup Sizes**	**Both cups**	**Left cup only**	**Right cup only**
1	30C, 32B, 32C, 34A, 34B, 36A	46111-01801-01	46111-01819-01	46111-01806-01
2	34C, 34D, 36B, 36C, 38A, 38B	46111-01802-01	46111-01820-01	46111-01807-01
3	34DD, 36D, 36DD, 38C, 38D, 40C	46111-01803-01	46111-01844-01	46111-01808-01
4	38DD, 40D, 40DD, 42C, 42D, 44C	46111-01804-01	46111-01845-01	46111-01809-01
5	42DD, 44D, 44DD, 46C, 46D, 48C	46111-01805-01	46111-01846-01	46111-01814-01
Without Cups				
Size	**Chest circumference^**			
1	79-89 cm (31"- 35")	46111-01891-01	-	-
2	89-99 cm (35 1/8"- 39")	46111-01892-01	-	-
3	99-109 cm (39 1/8"- 43")	46111-01893-01	-	-
4	109-119 cm (43 1/8"- 47")	46111-01894-01	-	-
5	120-130 cm (47 1/8"- 51")	46111-01895-01	-	-

^under-arm measurement

Bellisse®

Comfortable, Supportive Treatment for Breast Oedema

Bellisse® is an innovative treatment approach for chest wall, breast oedema and / or fibrosis, following breast cancer or surgery. Do not confuse this revolutionary "compression garment" with a mastectomy or post-operative bra. This is a medical compression garment that helps to reduce discomfort and provides gentle, even compression around the torso, facilitating drainage of excess lymphatic fluid from the breast and chest wall area. It is available in a large range of sizes and available in two colours, buff and black.

Bellisse® can be worn day and night, is fully adjustable, and provides gentle compression with or without the addition of JOBST® JoViPads. A number of anatomically shaped pads are available and can be used in combination with Bellisse®.

Sizes	
Eleven Band Sizes	Three Cup Sizes
✓	✓

Colours	
Black	Buff



Bellisse®

Measuring Information

It is recommended that measurements are taken when the oedema is at its most reduced state. Often this is found first thing in the morning.

Measurement Points


1. Measure in cm the circumference at Point 1, located directly under the bust (where a regular bra sits)

2. Move 5 - 7cm down towards the rib cage from Point 1 and measure circumference at Point 2

3. To determine the band size, take the larger of the two measurements from Point 1 or Point 2

4. Measure the circumference at Point 3 at the fullest part of the bust

5. To determine the cup size, subtract the larger of the two measurements from Point 1 or Point 2 from the bust measurement (Point 3). Use the difference in the chart below to select the cup size



Cup Size	A/B	C/D	DD/E
	5cm or less	5 - 10cm	10 - 15cm

If more than 15cm, go up one one band size

Band measurement (cm)	75 - 80	81 - 85	86 - 90	91 - 95	96 - 100	101 - 105	106 - 110	111 - 115	116 - 120	121 - 125	126 - 130
Bellisse Band Size	30	32	34	36	38	40	42	44	46	48	50

Ordering Information

Bellisse Size	Black	Buff	Bellisse Size	Black	Buff
30 A/B	N/A	76125-29	42 A/B	76125-15	76125-46
30 C/D	N/A	76125-30	42 C/D	76125-16	76125-47
32 A/B	76125-00	76125-31	42 DD/E	76125-17	76125-48
32 C/D	76125-01	76125-32	44 A/B	76125-18	76125-49
32 DD/E	76125-02	76125-33	44 C/D	76125-19	76125-50
34 A/B	76125-03	76125-34	44 DD/E	76125-20	76125-51
34 C/D	76125-04	76125-35	46 A/B	76125-21	76125-52
34 DD/E	76125-05	76125-36	46 C/D	76125-22	76125-53
36 A/B	76125-06	76125-37	46 DD/E	76125-23	76125-54
36 C/D	76125-07	76125-38	48 A/B	76125-24	76125-55
36 DD/E	76125-08	76125-39	48 C/D	76125-25	76125-56
38 A/B	76125-09	76125-40	48 DD/E	76125-26	76125-57
38 C/D	76125-10	76125-41	50 C/D	76125-27	76125-58
38 DD/E	76125-11	76125-42	50 DD/E	76125-28	76125-59
40 A/B	76125-12	76125-43			
40 C/D	76125-13	76125-44			
40 DD/E	76125-14	76125-45			

JOBST® JoViPads

Working in Combination

Available in different shapes, JOBST® JoViPads address particular areas of the breast and chest wall affected by pitting oedema and / or fibrosis. The pads provide gentle and prolonged pressure to soften congested tissue and encourage evacuation of lymph fluid to functional pathways. Designed for use under Bellisse®.



Cleavage Pad

	Black	Buff
One size	75928-00	75928-01



Lateral Pad

	Black	Buff
Small	76166-00	76166-03
Medium	76166-01	76166-04
Large	76166-02	76166-05



Mini Axilla Pad

	Black	Buff
Small	75926-00	75926-03
Medium	75926-01	75926-04
Large	75926-02	75926-05
Extra large	75926-07	75926-06



Inframammary Pad

	Black	Buff
Small	75927-00	75927-04
Medium	75927-01	75927-05
Large	75927-02	75927-06
Extra large	75927-03	75927-07



Post-Lumpectomy Pad

	Black	Buff
A/B - Small	75931-00	75931-04
C - Medium	75931-01	75931-05
D/DD - Large	75931-02	75931-06
E - Extra large	75931-03	75931-07



Unilateral Post-Mastectomy Pad

	Black	Buff
Small	75932-00	75932-04
Medium	75932-01	75932-05
Large	75932-02	75932-06
Extra large	75932-03	75932-07



Chest-Wall Pocket Pad

	Black	Buff
Small	75934-00	75934-04
Medium	75934-01	75934-05
Large	75934-02	75934-06
Extra large	75934-03	75934-07



Serratus Anterior Pad

	Black	Buff
Small	75929-00	75929-04
Medium	75929-01	75929-05
Large	75929-02	75929-06
Extra large	75929-03	75929-07



Axilla Pad

	Black	Buff
30-38cm - S/M	75925-00	75925-02
38-46cm - L/XL	75925-01	75925-03



Drain Pocket Pad

	Black	Buff
One size	76126-00	76126-01



Double Mastectomy Pad

	Black	Buff
Small	75933-00	75933-03
Medium	75933-01	75933-04
Large	N/A	75933-05



Bellisse Extender

	Black	Buff
One size	76127-00	76127-01

JOBST® JoViPads are available in both black and buff varients and a number of sizes ranging from small to extra large.

JOBST® JoViPads

The control of lymphoedema is a challenge for both the patient and health care professional. JOBST® JoViPads provide additional compression to the oedema and can help soften fibrotic tissues.

The irregular foam chip padding soften fibrotic tissues by alternating high and low pressure, creating a massage effect to help reduce oedema as well as addressing and preventing fibrosis. They are very useful in treatment of breast, chest wall, genital and head and neck lymphoedema.

Head, Neck and Eye Pads

Chin Strap Standard



Size	Code
Small	76133-00
Medium	76133-01
Large	76133-02
Extra large	76133-03

Chin Strap Extended



Size	Code
Small	76134-00
Medium	76134-01
Large	76134-02
Extra large	76134-03

Neck Pad Universal

Covers neck and back ear



Size	Code
Small / Medium	77610-00
Large / Extra large	77610-01

Neck Pad Peri-Auricular

Covers the neck front and back ear



Size	Code
Small / Medium	76407-00
Large / Extra large	76407-01

Eye Pad

With band




Size	Size	Code
Small / Medium	53 - 58.5cm	76135-00
Large / Extra large	58.5 - 66cm	76135-01

Measuring Information

Facial and Neck Garment Sizes

Size	Head circumference	Neck circumference
Small	53 - 56cm	28 - 33cm
Medium	57 - 58cm	34 - 38cm
Large	59 - 61cm	39 - 43cm
Extra large	62 - 66cm	44 - 48.5cm



JOBST® JoViPads are available in both black and buff varients and a number of sizes ranging from small to extra large.

JOBST® JoViPads

Genital Pads

Female Genital Pads - Suitable for the treatment of labia oedema and / or fibrosis. The width of the pad at the front has been designed to encompass any hysterectomy scar at the bikini line by completely covering it.

Male Genital Pads - Suitable for the treatment of scrotum oedema and / or fibrosis.

Female

Size	Size*	Code
Small	86 - 99cm	75936-01
Medium	100 - 110cm	75936-02
Large	111 - 120cm	75936-03
Extra large	121 - 132cm	75936-04

*Hip circumference

Male

Size	Size*	Code
Small	28 - 35cm	75937-00
Medium	36 - 41cm	75937-01
Large	42 - 47cm	75937-02
Extra large	48 - 52cm	75937-03

*Size as measured in scrotal circumference



Front



Back



Front



Back

Multi-Purpose Pads

Multi-purpose pads are slim line and can be applied in a variety of ways, to almost any area of the body, where oedema and / or fibrosis is a problem. The pads are ivory coloured, foam-filled pads made from organic materials cotton and spandex.

Multi-purpose pads are most effective under compression bandages, wrap compression systems or compression garments.

Round



Size	Code
One size - diameter 12.5cm	76411-00

Rectangle



Size	Code
Small / Medium - 10 X 20cm	76417-00
Large / Extra large - 10 x 30cm	76417-01

Palm



Size	Code
One size	76416-00

Dorsum



Size	Code
One size	76413-00

Oval



Size	Code
One size - 20 x 12.5cm	76412-0

Kidney



Size	Code
Small - 7.5 x 17.5cm	76414-00
Medium - 12.5 x 25cm	76414-01
Large - 17.5 x 32.5cm	76414-03

JOBST® JoViPads are available in both black and buff varients and a number of sizes ranging from small to extra large.

JOBST® JoViPads

Multi-Purpose PitPaks

PitPaks have been specially developed for the treatment of persistent, long term fibrosis or lipodematosclerosis, and are available in a variety of shapes making them anatomically suitable for various body parts.

PitPaks are harder than foam pads so it is advised that the patient's tolerance should be taken into account when determining the amount of desired compression

PitPaks Cilinders

Complete filling

Size	Circumference	Length	Code
Small	18 - 23cm	15cm	76426-00
Medium	23 - 30.5cm	16.5 cm	76426-01
Large	28 - 38cm	20.5cm	76426-02
Extra large	35.5 - 43cm	23cm	76426-03



Half-filling

Size	Circumference	Length	Code
Small	18 -25.5cm	19cm	76426-00
Medium	20.5 - 30.5cm	21.5cm	76426-01



PitPaks

Round



Size	Code
One size - 15cm	76418-01

Rectangle



Size	Code
One size - 15 x 22.5cm	76422-00

Inframammary



Size	Code
Small - 10 x 21.5cm	76425-00
Medium - 11.5 x 22.5cm	76425-01
Large - 12.5 x 24cm	76425-02
Extra large - 14 x 25cm	76425-03

Dorsum



Size	Code
One size - 10 x 12.5cm	76424-00

Oval



Size	Code
One size - 14 x 21.5cm	76419-01

Half Circle



Size	Code
Small - 15cm	76421-00
Medium - 20cm	76421-01

Tapard



Size	Code
One size - 5.5 x 17.5cm	76423-00

Crescent



Size	Code
One size - 5 x 12.5cm	76421-00

Terms and Conditions

1. FORMATION OF CONTRACT

1.1. The following conditions (the "Conditions") apply to the contract between the Seller and its customer (the "Buyer") for the supply of goods ("Goods") (the "Contract") to the exclusion of any other terms that the Buyer seeks to impose or incorporate, or which are implied by trade, custom, practice or course of dealing.

1.2. The Buyer acknowledges, that, in entering into the Contract, it is acting in the course of its business, trade or profession and that it is not acting as a "consumer" for the purposes of, and as defined in, any applicable laws including (without limitation) the Sale of Goods Act, 1979 (as amended) the Supply of Goods and Services Act, 1982 (as amended), the Unfair Contract Terms Act ,1977 (as amended), the Consumer Contracts Regulations, 2013, or otherwise.

1.3. All quotations, offers and tenders are made and all orders are accepted subject to the Conditions. Except as otherwise provided in the Conditions, all other terms, conditions or warranties express or implied are excluded from the Contract unless expressly accepted in writing by the Seller. In particular, the Seller is not be responsible for ensuring that the Goods are suitable for the Buyer's purposes.

1.4. The "Seller" shall mean BSN medical Ltd or any holding company, subsidiary or subsidiary of a holding company of BSN medical Ltd (for the purposes of this clause, the terms holding company and subsidiary shall be as defined is section 1159 of the Companies Act 2006).

1.5. The Buyer's purchase order or written acceptance of the Seller's quotation constitutes an offer by the Buyer to purchase Goods in accordance with these Conditions (the "Order").

1.6. The Order shall only be deemed to be accepted when the Seller issues written acceptance of the Order ("Despatch Note") at which point and on which date the Contract shall come into existence.

1.7. Any quotation given by the Seller shall not constitute an offer, shall only be valid for a maximum period of 30 days from date of issue, and may be withdrawn by the Seller immediately at any time by written or oral notice to the Buyer.

1.8. All references to INCOTERMS shall be to the INCOTERMS 2010.

2. PRICES

2.1. The price of the Goods shall be the price set out in the Despatch Note, or, if no price is quoted, either the price set out in the Seller's published price list in force as at the date of delivery or the price agreed between the parties in writing from time to time.

2.2. There is a minimum order value of £100.00 ("Minimum Order Value") for each Order. The Seller may accept orders below the Minimum Order Value but any such orders accepted by the Seller will be subject to an administration charge of up to, but not exceeding, £25.00. For the avoidance of doubt, the Minimum Order Value does not apply to orders of Compression Garments as set out in Condition 24.3.

2.3. In respect of delivery (or arrangement of delivery) of the Goods:
2.3.1. if delivery is not included in the price, the Buyer shall pay all transportation and insurance costs incurred by the Seller in making or arranging such delivery; and
2.3.2.the Buyer shall be liable to the Seller for further carriage costs and any demurrage costs incurred by the Seller if vehicles are unduly delayed at the place of delivery.

2.4. Unless otherwise agreed in writing, all prices are quoted exclusive of VAT.

2.5. The Seller reserves the right to increase the price of the Goods, by giving notice to the Buyer at any time before delivery to reflect an increase in the cost of the Goods to the Seller in excess of 10%, that is due to:
2.5.1. any factor beyond the control of the Seller (including foreign exchange fluctuations, increases in taxes and duties, the action of any government or authority, and increases in labour, materials and other manufacturing costs);
2.5.2.any request by the Buyer to change the delivery date(s), quantities or types of Goods ordered, or the specification of the Goods; or
2.5.3.any delay caused by any instructions of the Buyer in respect of the Goods or failure of the Buyer to give the Seller adequate or accurate information or instructions in respect of the Goods.

2.6. The Seller reserves the right to make an additional charge to cover any increase in transport costs occurring before the date of delivery.

3. PAYMENT

3.1. Unless otherwise agreed by the Seller in writing, payment shall be due and payable on the 20th day of the month following the date of invoice and shall be made without any set-off, counterclaim or deduction whatsoever. Where any dispute arises, the Buyer shall not oppose and shall consent to the immediate enforcement of any judgement entered against it regardless of any cross-claim, set-off or counterclaim.

3.2. The Seller shall be entitled to submit its invoice with its delivery advice note or at any time afterwards except that where delivery has been postponed at the request of, or by the default of, the Buyer, the Seller may submit its invoice at any time after the Goods are ready for delivery or would have been ready in the ordinary course but for the Buyer's request or default.

3.3. Where Goods are delivered by instalments the Seller may invoice each instalment separately and the Buyer shall pay such invoices in accordance with these Conditions.

3.4. No disputes arising under the Contract nor delays (other than delays acknowledged by the Seller in writing), shall interfere with prompt payment in full by the Buyer.

3.5. If the Buyer shall default in payment whether under this Contract or any other contract with the Seller, the Seller shall be entitled, without prejudice to any other right or remedy, to do all or any of the following:-
3.5.1. to suspend any or all further deliveries under the Contract and under any other contract between the Seller and the Buyer, without notice;
3.5.2.to charge interest on any amount outstanding at the rate of 3% per annum above the base rate of Royal Bank of Scotland plc, such interest being charged as a separate, continuing obligation not merging with any judgment, together with all debt recovery costs;
3.5.3.to serve notice on the Buyer requiring immediate payment for all Goods supplied or to be supplied by the Seller under this Contract and all other contracts with the Buyer whether or not payment is otherwise due;
3.5.4.by its servants or agents, enter the Buyer's premises and recover and dispose of the Goods, and the Buyer shall not be entitled to make any claim against the Seller in respect of such entry or disposal; and/or
3.5.5.to sue for the price of the Goods and any other sums then due or rendered due under the terms of any contract, even though title may not have passed to the Buyer or the Goods have been recovered, provided that the Seller shall give credit for the cost (excluding all profit) of any Goods not supplied or recovered from the Buyer as certified by the Seller, such certificate to be conclusive proof of the cost of such Goods.

3.6. Unless agreed in writing by the Seller, the Buyer shall have no right to set-off, withhold or deduct from monies due to the Seller whether in respect of any sums due from the Seller, any sums claimed to be due from the Seller or otherwise.

4. DELIVERY

4.1. The Seller shall deliver the Goods to the location set out in the Order or such other location as the parties may agree ("Delivery Location").

4.2. Delivery of the Goods within the European Union shall be delivered DAP unless the Seller has specifically agreed otherwise in writing.

4.3. Delivery of the Goods anywhere outside the European Union shall be EXW at the Seller's premises, unless the Seller specifically agrees in writing to deliver the Goods to the Buyer, in which case delivery of the Goods shall be DAP at the Buyer's premises.

4.4. Time for delivery is given as accurately as possible but is not guaranteed. The Buyer shall have no right to damages or to cancel the Contract for failure for any cause to meet any delivery time stated nor shall the Buyer be entitled to make, or to purport to make, time for delivery of the essence of the Contract. The Seller shall be entitled to deliver within 30 days of any delivery time stated.

4.5. The date for delivery shall in every case be dependent upon prompt receipt of all necessary information, final instructions or approvals from the Buyer in connection with any act required to enable the Seller to despatch the Goods.

4.6. If any order for Goods is to be delivered by several instalments to the Buyer, each such instalment shall be treated as a separate and identifiable contract and the rights of either party under each contract shall be construed accordingly, save only that the Seller may suspend delivery whilst payment is overdue in respect of any previous instalment.

4.7. If the Seller fails to make delivery or makes defective delivery of any instalment, such failure or defective delivery shall not affect the obligations of the parties under the Contract in respect of the other instalments.

4.8. The Seller will endeavour to comply with reasonable requests by the Buyer for postponement of delivery of the Goods but shall be under no obligation to do so. Where delivery is postponed otherwise than due to default by the Seller then, without prejudice to all other rights and remedies available to the Seller, the Buyer shall pay all costs and expenses, including a reasonable charge for storage and transportation.

4.9. The Seller is under no obligation to accept Goods which are returned to them due to an error on the part of the Buyer (including but not limited to over ordering of Goods or holding an excess of stock) ("Returned Goods"). Where the Seller agrees in writing to accept Returned Goods the Buyer shall be liable to pay the full price for the Goods, less the following, to the Seller:
4.9.1. Goods returned within 30 days of delivery – 90% of the price will be refunded;
4.9.2.Goods returned between 30 and 90 days after delivery – 80% of the price will be refunded;
4.9.3.Goods returned between 90 and 180 days after delivery – 60% of the price will be refunded;
4.9.4.Goods returned between 180 and 270 days after delivery – 40% of the price will be refunded;
4.9.5.Goods returned between 270 and 365 days after delivery – 20% of the price will be refunded; and
4.9.6.for Goods returned 365 days or more after delivery there will be no refund.
For the purposes of this Condition 4.9, the Goods will be considered 'returned' when the Seller is in physical possession of the Returned Goods.

4.10.If the Buyer refuses to take delivery of any Goods with no valid reason (whether such reason is valid is to be determined at the absolute discretion of the Seller) ("Refused Goods") the Buyer is liable to pay to the Seller the full delivery costs for the Refused Goods plus a handling charge of 10% of the price of the Refused Goods or £15, whichever is the higher.

4.11. Subject to Conditions 4.9 and 4.10, the Seller will only accept Returned Goods and Refused Goods if:
4.11.1. the Goods are complete units of sale; and
4.11.2.the Goods are in the same condition as at the date of delivery; and
4.11.3.the Goods are saleable.

4.12.The Seller will organise for the collection of all Goods returned in accordance with Conditions 4.9, 4.10 and 4.11.

5. PACKAGING

5.1. Packaging supplied by the Seller, unless otherwise expressly agreed in writing, is intended to provide adequate protection in normal conditions of transit of expected duration.

5.2. The Buyer shall be responsible for ensuring the Goods are kept in adequate storage conditions once delivered, including but not limited to protecting the Goods against the effects of humidity and temperature.

5.3. The Seller shall be entitled to deliver the total weight, area or volume of any Contract in containers of such varying weights and sizes as it may decide at its absolute discretion.

5.4. Where packages, pallets, drums, beams or other containers ("Packaging") are stated to be returnable, they shall be returned in good order and condition at the Buyer's expense to the destination named on the Seller's invoice and the Buyer shall notify the Seller of the date of despatch. If any such Packaging was charged for on delivery, credit will be given when it is received at the named destination in good order and condition. If the Buyer fails within a reasonable period to return any such Packaging in good order and condition, the Buyer shall, if not so charged on delivery, be charged for the Packaging at the Seller's standard rate prevailing at the time.

6. RISK AND TITLE

6.1. Risk shall pass to the Buyer (so that the Buyer is then responsible for all loss or deterioration of the Goods or for any damage occurring):-
6.1.1. if the Seller is responsible for delivering or arranging for delivery of the Goods, at the time when the Goods arrive at the Delivery Location; or
6.1.2. in all other circumstances, on the completion of loading of the Goods at the premises of the Seller.

6.2. Title to the Goods shall only pass to the Buyer if the Buyer has paid to the Seller all sums (including any default interest) due from it to the Seller under the Contract and under all other contracts between the Seller and the Buyer (including any sums due under contracts made after the Contract) whether or not the same are immediately payable.

6.3. Until title to the Goods has passed to the Buyer, the Buyer shall:
6.3.1. store the Goods separately from all other goods held by the Buyer so that they remain readily identifiable as the Seller's property;
6.3.2.not remove, deface or obscure any identifying mark or packaging on or relating to the Goods;
6.3.3.maintain the Goods in satisfactory condition and keep them insured against all risks for their full price on the Seller's behalf from the date of delivery;
6.3.4. notify the Seller immediately if it becomes subject to any of the events listed in Condition 14.1; and
6.3.5. give the Seller such information relating to the Goods as the Seller may require from time to time.

6.4. Subject to Condition 6.5, the Buyer may resell or use the Goods in the ordinary course of its business (but not otherwise) before the Seller receives payment for the Goods. However, if the Buyer resells the Goods before that time:
6.4.1. it does so as principal and not as the Seller's agent; and
6.4.2. title to the Goods shall pass from the Seller to the Buyer immediately before the time at which resale by the Buyer occurs.

6.5. If before title to the Goods passes to the Buyer the Buyer becomes subject to any of the events listed in Condition 14.1 then, without limiting any other right or remedy the Seller may have:
6.5.1. the Buyer's right to resell Goods or use them in the ordinary course of its business ceases immediately; and
6.5.2. the Seller may at any time:
(a) require the Buyer to deliver up all Goods in its possession which have not been resold, or irrevocably incorporated into another product; and
(b) if the Customer fails to do so promptly, enter any premises of the Buyer or of any third party where the Goods are stored in order to recover them.

6.5. The Seller may recover Goods in respect of which title has not passed to the Buyer at any time and the Buyer irrevocably licences the Seller, its officers, employees and agents to enter upon any premises of the Buyer, with or without vehicles, for the purpose either of satisfying itself that Condition 6.3 is being complied with by the Buyer or of recovering any Goods in respect of which title has not passed to the Buyer.

7. MANUFACTURER GUARANTEES

7.1. Some of the Goods sold by the Seller come with a manufacturer's guarantee. The Buyer should refer to the manufacturer's guarantee provided with the Goods for details of the applicable terms and conditions.

8. VARIATIONS

8.1. No variations to the Goods required by the Buyer shall be binding on the Seller unless agreed by the Seller in writing, including, as part of such agreement, any change to the contract price arising in consequence of the variation to the Goods.

8.2. Following agreement of any such variation, any dates quoted for delivery shall be extended accordingly.

9. SAMPLE AND TESTING

9.1. The Seller reserves the right to alter the dimensions or composition of the Goods supplied to conform to applicable standards or laws.

9.2. The Seller shall not alter the specification of any Goods to those of its current specification for such Goods from time to time ordered by the Buyer, without prior agreement of the Buyer, provided that the Seller may alter raw materials or the provider of raw materials such that the product remains in compliance with the original specification.

9.3. The Buyer is responsible for testing and examining the Goods upon delivery to ensure that:-
9.3.1. the Goods are suitable and fit for any purpose for which the Buyer intends to use said Goods; and
9.3.2. the Goods delivered do not contain any contamination likely to cause damage to the Buyer's property or to have an adverse effect on the Buyer's processes or products.

9.4. The Buyer shall indemnify the Seller against any loss arising or damage suffered by the Seller pursuant to any failure by the Buyer to :-
9.4.1. comply with Condition 9.3; or
9.4.2. notify the Seller pursuant to Condition 11.2 of any findings pursuant to its testing and examination of the Goods under Condition 9.3.

9.5. Any samples, drawings, descriptive matter or advertising issued by the Seller and any descriptions or illustrations of the Goods contained in the Seller's catalogues, brochures, on the Seller's Website, elsewhere online or in any other digital format, are issued or published for the sole purpose of giving an approximate idea of the Goods described in them. They shall not form part of the Contract or have any contractual force.

10. USE OF THE GOODS

10.1. All recommendations and advice given by or on behalf of the Seller to the Buyer as to the methods of storing, applying or using the Goods, the purposes to which the Goods may be applied and the suitability of using the Goods in any manufacturing process or in conjunction with any other materials, are given without liability on the part of the Seller its servants or agents and where the Goods are intended for use in connection with any drugs or other Goods or substance whatever, or for the packaging or securing of goods or articles of any kind the Buyer shall be responsible for ensuring that such articles will not be, or not liable to be adversely affected by contact with or proximity to the goods supplied and all liability to the Buyer or any third party in respect of any damage so caused is hereby expressly excluded.

10.2. Any instructions given by the Seller is not intended to substitute for advice given by a medical practitioner, pharmacist or other licensed healthcare professional. Information and statements about Goods are not intended to be used to diagnose, treat, cure or prevent any disease or health condition. The Seller accepts no liability for inaccuracies or misstatements about Goods by manufacturers or other third parties.

11. DEFECTS AND DELIVERY DISCREPANCIES

11.1. The Seller warrants that the Goods supplied shall comply in all material respects with its standard specification for those Goods from time to time.

11.2. The Buyer shall notify the Seller of:-
11.2.1. any discrepancy between the delivered quantity of Goods and the stated quantity set out on the delivery documentation (a "Discrepancy"); and/or
11.2.2. any reason (other than a change in the Buyer's requirements) whereby the Goods delivered are not in accordance with the Seller's standard specification for the Goods (a "Defect"),
and where the Discrepancy or Defect is or should be apparent on delivery, the Buyer shall note the Discrepancy or Defect on the Seller's delivery documents.

11.3. The Buyer shall only be entitled to claim (and then subject to Condition 12) for Defects which are or would be apparent on visual inspection at the time of delivery if:-
11.3.1. the Buyer inspects and undertakes its testing of the Goods within 3 (three) working days (and in any event before use of any part of the Goods) following the date of their arrival at the Delivery Location;
11.3.2. the Seller is notified in writing within 10 (ten) working days after the inspection or testing under Condition 11.3.1 of the particular Defect; and
11.3.3. the Seller is given a reasonable opportunity to inspect the Goods and investigate any complaint.

11.4. The Buyer shall only be entitled to claim (and then subject to Condition 12) for Discrepancies which are or would be apparent on visual inspection at the time of delivery if the Seller is notified in writing within 3 (three) working days after the time of delivery.

11.5. The Buyer shall only be entitled to claim (and then subject to Condition 12) for Defects and Discrepancies which are not apparent on visual inspection at the time of delivery if:-
11.5.1. the Seller is notified in writing:
(a) within three (3) days of the Buyer discovering the Defect or Discrepancy; and
(b) the Buyer discovers the Defect or Discrepancy within 3 (three) months following the date of delivery of the Goods; and
11.5.2. the Buyer makes no subsequent use of the Goods or alteration or interference to or with the Goods before the Seller has a reasonable opportunity to inspect the Goods and investigate any complaint.

11.6. The Buyer shall indemnify the Seller against any loss arising or damage suffered through failure to inform the Seller within 5 (five) working days of the estimated date of delivery notified by the Seller in the event of non-delivery of a whole consignment.

11.7. Where the Buyer makes notification in accordance with Conditions 11.3, 11.4 or 11.5, of a Defect or Discrepancy which is objectively justified, the Seller shall (at the Seller's option):-
11.7.1. remove the delivered Goods and provide replacement Goods; or
11.7.2. remove the delivered Goods and reimburse the Buyer (whether by way of credit note or otherwise) the price paid for those Goods.

12. EXCLUSION AND LIMITATION OF LIABILITY

12.1. Subject to Conditions 12.2 to 12.3, the Seller's aggregate liability to the Buyer in respect of all claims arising out of or in connection with any Contract (including without limitation as a result of breach of contract, negligence, or any other tort, under statute or otherwise) will be limited to the greater of a sum equal to 500% (five hundred per cent) of the aggregate amount paid or payable under the Contract or £80,000 (eighty thousand pounds).

12.2. Subject to condition 12.3, the Seller shall under no circumstances whatever be liable to the Customer, whether in contract, tort (including negligence), breach of statutory duty, or otherwise, for:
12.2.1. loss of profits (whether direct or indirect);
12.2.2. loss of sales or business (whether direct or indirect);
12.2.3. loss of agreements or contracts (whether direct or indirect);
12.2.4. loss of anticipated savings (whether direct or indirect);
12.2.5. loss of or damage to goodwill (whether direct or indirect);
12.2.6. loss of use or corruption of software, data or information; and/or
12.2.7. any indirect or consequential loss.

12.3. Nothing in these Conditions shall limit or exclude the Seller's liability for:
12.3.1. death or personal injury caused by its negligence, or the negligence of its employees, agents or subcontractors;
12.3.2. fraud or fraudulent misrepresentation;
12.3.3. breach of the terms implied by section 12 of the Sale of Goods Act 1979 (title and quiet possession);
12.3.4. any breach of the Contract that results from the wilful act or wilful omission of the Seller or its employees, agents or subcontractors;
12.3.5. defective products under the Consumer Protection Act 1987; or
12.3.6. any other matter in respect of which liability cannot by applicable law be limited.

12.4. The terms implied by sections 13 to 15 of the Sale of Goods Act 1979 and the terms implied by sections 3 to 5 of the Supply of Goods and Services Act 1982 are, to the fullest extent permitted by law, excluded from the Contract.

12.5. This clause 12 shall survive termination of the Contract.

13. CONFIDENTIAL INFORMATION AND INTELLECTUAL PROPERTY RIGHTS

13.1. All drawings, documents, confidential records, computer software and other information supplied by the Seller, whether produced by itself or a third party, are supplied on the express understanding that copyright is reserved to the Seller (or the third party) and that the Buyer will not, without the written consent of the Seller, give away, loan, exhibit or sell any drawings, documents, records, software or other information or extracts from them or copies of them or use them in any way except in connection with the Goods in respect of which they are issued.

13.2. All claims for alleged infringement of patents, trade marks, registered designs, design right or copyright received by the Buyer relating to the Goods must be notified immediately to the Seller. If requested by the Seller, the Seller shall be entitled to have conduct of any proceedings relating to any such claim in such manner as the Seller thinks fit and the Buyer will provide to the Seller such reasonable assistance as the Seller may request. The cost of any such proceedings will be borne by the Seller.

13.3. If any allegation shall be made against the Buyer to the effect that the supply of such of the Goods as are manufactured by the Seller infringes the intellectual property rights of any third party or the Seller has reason to believe that any such allegation is likely to be made, the Seller may at its option and expense modify or replace the Goods so as to avoid the infringement (but without adversely effecting the overall performance of the Goods), or obtain for the benefit of the Buyer the right to continue to use the Goods, or repurchase the Goods at the contract price as reduced by a reasonable provision for depreciation. If the Seller pursues any of such options, the Buyer will have no rights or remedies against the Seller arising directly or indirectly out of the alleged infringement.

13.4. To the extent that the Goods are to be manufactured in accordance with a specification supplied by the Buyer, the Buyer shall indemnify the Seller against all liabilities, costs, expenses, damages and losses (including any direct, indirect or consequential losses, loss of profit, loss of reputation and all interest, penalties and legal and other professional costs and expenses) suffered or incurred by the Seller in connection with any claim made against the Seller for actual or alleged infringement of a third party's intellectual property rights arising out of or in connection with the Seller's use of the specification. This clause 13.4 shall survive termination of the Contract.

14. INSOLVENCY

14.1. The Seller may terminate the Contract with immediate effect by giving the Buyer written notice at any time if any encumbrancer takes possession of or a receiver, administrative receiver or similar officer is appointed over any of the property or assets of the Buyer or if the Buyer makes any voluntary arrangement with its creditors or becomes subject to an administration order or has an administrator appointed or goes into liquidation or has a resolution for its winding-up passed (except for the purpose of amalgamation or reconstruction not involving insolvency where the resulting entity agrees to be bound by or assumes the obligations imposed on the Buyer) or anything analogous to any of these events under the law of any jurisdiction occurs in relation to the Buyer or if the Buyer ceases or threatens to cease to carry on business.

14.2. The termination of the Contract pursuant to Condition 14.1 will not affect any rights and/or liabilities of the Buyer which have accrued before termination.

15. FORCE MAJEURE

15.1. The Seller shall be under no liability for any delay or failure to perform any of its obligations under the Contract arising from any cause or causes beyond its reasonable control including without limitation any of the following: act of God, war, riots, civil commotions, strikes, lock-outs, trade disputes, fires, breakdowns, failure of a utility service, interruptions of transport, governmental action or restriction, malicious damage, shortages of labour or materials or breakdown of machinery, delay in delivery by Seller's suppliers or any other cause whatsoever (whether or not similar to the foregoing) ("Force Majeure Event").

15.1. If the Force Majeure Event prevents the Seller from providing supply of Goods under the Contract for more than 1 month, either party shall have the right to terminate the effected part of the Contract immediately by giving written notice to the other party.

16. CONSUMER PROTECTION ACT 1987 (THE "ACT")
In circumstances in which the Seller supplies Goods to the Buyer for incorporation with, or use ancillary to, any composite or other products to be produced, manufactured, processed or supplied by the Buyer or a third party then the Buyer shall indemnify the Seller against all liabilities, costs, expenses, damages and losses (including any direct, indirect or consequential losses, loss of profit, loss of reputation and all interest, penalties and legal and other professional costs and expenses) that the Seller incurs, or has to bear, if any claim or claims are made against the Seller, pursuant to the Act.

17. NOTICE

17.1. Any notice or other communication given under or in connection with the Contract shall be in writing and shall be delivered personally, sent by pre-paid first class post or other next working day delivery service or e-mail.

17.2. A notice or other communication shall be deemed to have been received: if delivered personally, when left at the Seller's registered office; if sent by pre-paid first class post or other next working day delivery service, at 9.00 am on the second Business Day after posting or if sent by e-mail, one Business Day after transmission.

17.3. In proving the service of any notice, it will be sufficient to prove, in the case of a letter, that such letter was properly addressed, stamped and placed in the post and, in the case of an e-mail, that such e-mail was sent to the specified e-mail address of the addressee.

17.4. The provisions of this clause shall not apply to the service of any proceedings or other documents in any legal action.

18. THE SELLER'S RIGHT TO VARY THESE TERMS
18.1. The Seller amends these Conditions from time to time.
18.2. The Conditions in force at the time of the Buyer's Order will apply to the Contract.
19. ASSIGNMENT
19.1. The Seller may at any time assign, transfer, mortgage, charge, subcontract or deal in any other manner with all or any of its rights under the Contract and may subcontract or delegate in any manner any or all of its obligations under the Contract to any third party.
19.2. The Buyer shall not, without the prior written consent of the Seller, sub-license, assign or transfer in any way any of its rights, liabilities and/or obligations under the Contract on a temporary or permanent basis to any third party.
20. SEVERABILITY
If any provision or part-provision of the Contract is or becomes invalid, illegal or unenforceable, it shall be deemed modified to the minimum extent necessary to make it valid, legal and enforceable. If such modification is not possible, the relevant provision or part-provision shall be deemed deleted. Any modification to or deletion of a provision or part-provision under this Condition shall not affect the validity and enforceability of the rest of the Contract.
21. WAIVER
21.1. A failure by the Seller to exercise, or a delay in exercising, any right or remedy under the Contract shall not constitute a waiver of the right or remedy or a waiver of any other rights or remedies which the Seller may otherwise have and no single or partial exercise of any right or remedy under the Contract shall prevent any further exercise of the right or remedy or the exercise of any other right or remedy.
21.2. Any waiver by the Seller of a breach of any of the terms of the Contract or of any default under the Contract shall not be deemed a waiver by the Seller of any subsequent breach or default and shall not affect the other terms of the Contract.
22. THIRD PARTY RIGHTS
No person who is not a party to the Contract is intended to receive a benefit under, or be entitled to enforce, the Contract pursuant to the Contracts (Rights of Third Parties) Act 1999. Any right or remedy of a third party which exists or is available apart from the Act is not affected.
23. LAW AND JURISDICTION
The parties agree that any disputes arising or in any way connected with the subject matter of the Contract (whether of a contractual or tortious nature or otherwise) shall be subject to the laws of England and in the case of proceedings issued against the Seller shall be subject to the jurisdiction of the English courts only. Schedule 1 to the Uniform Law on International Sales Act 1967 shall not apply to the Contract.
24. SPECIAL CONDITIONS RELATING TO THE SUPPLY OF 'ready to wear' COMPRESSION GARMENTS
24.1. The provisions of this Condition 24 apply to 'Ready to Wear' compression garments ("Compression Garments") only. For the avoidance of doubt, the Conditions and this Condition 24 do not apply to 'Made to Measure' compression garments.
24.2. In the event of a conflict between this Condition 24 and the remainder of the Contract, the provisions in this Condition 24 will prevail.
24.3. The Minimum Order Value set out in Condition 2.1 shall not apply to Compression Garments.
April 2015

Notes


Notes

Notes

Notes

The Lymph Drainage System

Simplified diagram of the most important anatomical areas



Conception: Peter D. Asmussen Graphic design: Robert Darroll © JOBST Emmerich 2002

JOBST GmbH
Beiersdorfstr.1 · 46446
Emmerich

Information



Customer Service

Our experienced UK Customer Service team use their expertise to process orders efficiently and answer product queries.

- Call: 0345 122 3600
- Direct fax: 0345 122 3450
- Email: compression.uk@jobst.com

PATH

Professional Accredited Training for Healthcare Professionals

Training and Education

PATH is a gold standard education platform with a diverse selection of therapy areas and learning styles, to help you provide the best possible patient care.

- Contact your local Compression Therapy Specialist or scan the QR code for more details

JOBST Delivered

JOBST® Delivered

JOBST Delivered provides a prompt and hassle-free delivery of JOBST garments and can dispense the garment to any UK address. This can be the patient's home, clinic or hospice – and it's delivered free of charge! For more information and a WELCOME PACK - including new envelopes:

- Email: jobstdelivered@jobst.com
- Call: 0808 196 6250
- Visit: www.jobst.co.uk



Compression Therapy Specialist

Your local Essity Compression Therapy Specialist is your first point of contact for our entire range of compression therapy products and as such is able to offer you advice, support and detailed product information.

- Call your local Compression Therapy Specialist directly or call the Concierge Service on 01482 670177 to schedule a visit

essity

JOBST®, an Essity brand

Essity
T/A BSN medical Limited
PO Box 258 • Willerby • Hull • HU10 6WT

www.jobst.co.uk
Tel: 01482 670100 • Fax: 01482 670111
E-mail: compression.uk@jobst.com

/JOBSTUK
@JOBSTUK