
Niti-S Biliary Covered Stent [Giobor]

for malignant biliary strictures



Features

- **Partially covered design**

(1) The long covered portion (70%)
(2) The bare portion (30%)

- **Anti-migration flared end**

- **Radiopaque markers:** Three (3) at both ends & Two (2) at the end of the covered part

TaeWoong MEDICAL

Ordering Information

Endoscopic				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
BS0806FP	8	6	8.5	180
BS0808FP		8		
BS0810FP		10		
BS1006FP	10	6	8.5	180
BS1008FP		8		
BS1010FP		10		

Aquilant
A HC2I Company
csuk@hc21.group
01256 306506
healthcare21.eu

ⓒ 2022, Taewoong Medical Co., Ltd. All rights reserved. TW-MBB-CG (Rev.00) Printed in Korea. May 04. 2022